





TABLE OF CONTENTS

GENERAL	4
Introduction	4
Spare parts and adaptations	4
Tightening Allen screws	4
Torque	5
Checks	5
Identifying and alleviating malfunctions	5
FRAME	6
Changing Frame	6
SEAT	7
Fitting seat cover to frame	7
Front seat height (SHv)	7
Rear seat height (SHh)	8
Checking the folding unit for play and adjusting	9
BACKREST	10
Tension adjustable backrest	10
Backrest height (RH)	13
Backrest angle (RW)	14
Push handles / backrest telescopes	14
Stabilisation bar	16
FOOTRESTS	17
Footplate	17
Lower leg length (UL)	17
Footrest mounted in high position	18
SIDE PARTS	19
Clothes-guard assembly	19
Mud guard assembly	20
Siderest assembly and adjustment	21
Armrest assembly and adjustment	21
FRONT WHEELS	22
Replacing a front wheel	22
Replacing a front wheel fork	22
Front wheels, repositioning to the frame	23
Front wheel supporter, adjusting the steering error angle	23
Front wheel supporter, adjusting the trail angle	23
REAR WHEELS	24
Rear wheels, repositioning	24
Wheel camber, adapter sleeves	24
Adjusting the removable axle	24
Distance sleeves	25
Adapterplate	25
BRAKES	26
Fitting / adjusting the parking brake	26
OPTIONS & ACCESSORIES	27
Antitipper	27
Tipper aid	28
Transit wheels	28
Cane holder	29
Fitting the pelvic belt	29



GENERAL

Introduction

This service manual is part of the instructions and contains the technical information for servicing, configuring and repairing a Kuschall® wheelchair.



WARNING!

Danger of accident and severe injuries.

If the wheelchair is improperly set it can cause accidents and severe injuries.

► **Changes to the wheelchair may only be carried out by the dealer.**


To guarantee the required safety and reliability, all wheelchairs must be comprehensively checked once a year.

In part, assembly and adjustment require extensive experience. For this reason, the following assembly instructions have been split into three categories:

Requirement	Symbol
Easy – technical understanding required	●○○
Intermediate – specialist knowledge required	●●○
Difficult – specialist wheelchair assembly knowledge and experience required	●●●

The required tools and their respective sizes are listed above each instruction. The instructions include information on the torques with which the respective screw connections must be tightened. Adhering to the given torques requires the use of a torque spanner.

Tools	Symbol
Allen key	● 2x3, 4, 5, 6
Phillips screwdriver	✕ 2
Straddle spanner	—C 19
Socket spanner/ring spanner	⬡ 8, 10, 14, 22



Spare parts and adaptations


All spare parts can be purchased from Kuschall®'s Customer Services. An electronic spare parts catalog is available by logging onto www.kueschall.com. Only original spare parts may be used. Installing additional adaptations to a Kuschall® wheelchair requires the prior written approval of Kuschall® AG.

Tightening Allen screws

Allen keys are not designed for greater forces. When tightening or loosening an Allen screw, it is therefore advisable to apply force to the nut to prevent the hexagon socket from being damaged.

Tightening and loosening


- Turn the nut with a socket spanner (only use a straddle spanner if there is insufficient space) and merely hold the screw tight with the Allen key.






Tightening and loosening if there is no nut

If an Allen screw is directly screwed into a screw thread, the screw must be tightened using an Allen key.

 Ensure that the Allen key is of good quality and not worn.



Torque

■ All screw connections must be tightened with the torques specified in the following instructions.

Checks

Visual check

■ Check the entire frame for cracks, especially the areas around joints and welded seams.

Checking the screw connections

■ Check all bolts with the torques specified in the instructions regularly, and adjust if required.



CAUTION!

Several screw connections have been secured with safety seals. If these are opened, they must be secured again using new safety seals. High-power and low power adhesives are available. For torque entries notice shall be made whether an adhesive and which adhesive needs to be used.

Identifying and alleviating malfunctions

Malfunction	Possible cause	Measure
The wheelchair will not move in a straight line	Incorrect tire pressure in a rear wheel	Correct tire pressure
	One or more spokes broken	Replace defective spoke(s)
	Spokes unevenly tensioned	Tighten excessively loose spokes
	Dirty or damaged wheel bearings	Clean or replace bearings
	Bearing block of castor fork is not vertical	Align bearing block vertically
	Front wheels not set to the same height	Position the front wheels in such a way that they touch the ground at the same time
	Rear wheels not parallel or axes not aligned	Adjust the prestress load on the scissor mechanism and/or the trail
The wheelchair tips backwards too easily	Rear wheels have been fitted too far forward	Fit rear wheels further back
	Backrest angle too great	Reduce backrest angle
	Seat angle too great	Mount the adapter plate lower on the side profile Mount the smaller castor fork
The brakes engage poorly or asymmetrically	Incorrect tire pressure in one or both rear wheels	Correct tire pressure
	Brake setting incorrect	Correct brake setting
Roll resistance is too great	Insufficient tire pressure in the rear wheels	Correct tire pressure
	Rear wheels are not parallel	Ensure that the rear wheels are parallel (4° wheel camber)
The front wheels wobble when moving fast	Insufficient tension in the front wheel bearings block	Lightly tighten the nut in the bearings block axle
	Front wheel is worn flat	Replace front wheel
The front wheel is stiff or stuck	Dirty or damaged bearings	Clean or replace the bearings
The wheelchair is very difficult to unfold	The backrestcover is too tight	Loosen the backrestbands a little
Handling seems imprecise	The scissor mechanism is not closed properly	If required, remove dirt from scissor mechanism
	The scissor mechanism is misaligned	Realign the scissor mechanism



FRAME

The Champion 08 frame is available in aluminium, titan or carbon. Aluminium or titan frames are available with frame angles of 75° and 90°, carbon frames with a frame angle of 87°

Changing Frame

Fitting frame tubes

- Fit both frame tubes into the side supporters^① using bolts ^② and tighten lightly.



When assembling the titan frames, make sure an aluminium spacer sleeve ^④ is placed between frame tube ^③ and side supporter.

Drilling frame mounting holes

Difficulty: ●●●

Tools: drill, drill bit: Ø 5,9 mm

- Before drilling the frame mounting holes, fit and adjust the following components:

- Front wheels, → Chapter Front wheels;
 - Front wheels, repositioning to the frame
 - Front wheel supporter, adjusting the steering error angle
 - Front wheel supporter, adjusting the trail angle
- Footrests, → Chapter Footrests;
 - Footrest, height adjustment
- Adapter sleeves, Distance sleeves, Adapterplates,
 - Chapter Rear Wheels;
 - Adjustment of rear wheel parallelism
 - Adapterplate, adjustment to folding unit or assembly of new adapterplate
 - Adapterplate spacer assembly
- Seat cover, → Chapter Seat;
 - Fitting seat cover to frame



Make sure that all settings on the wheelchair are correct. To guarantee that the wheelchair runs straight, make sure that the frame is parallel and in plumb.



- Drill the frame mounting holes ^⑤ carefully on both sides through the side supporters and the seat cover from the inside out to Ø 5.9 mm.

Attaching the side supporters, frame tubes and seat cover



The purpose of the threaded connection through the frame mounting holes is to provide additional fixing. The frame is thus easier to guide and is better secured. Any misalignment of the toe angle of the front wheels is counteracted.

- Attach the side supporters, frame tubes, seat cover and optionally, the holders for the mudguard or siderests using bolts through the frame mounting holes ^⑥ on both sides and tighten.






SEAT

Fitting seat cover to frame


- Fit seat cover to side supports on both sides using bolts ①.
















Front seat height (SHv)

The following possibilities are available to adjust the front seat height (SHv):

- Replace the front wheel with a larger or smaller one,
→ Chap. Front wheels; Replacing a front wheel.
- Replace the front fork with a larger or smaller one,
→ Chap. Front wheels; Replacing a front wheel fork.

 Adjusting the front seat height changes the seat angle. It may be necessary to adjust the rear seat height correspondingly.

 In the case of a 4° wheel camber, it must be ensured that the rear wheels are parallel after changing the front seat height. If required, they must be readjusted,
→ Chap. Rear wheels; Adjustment of rear wheel parallelism.

Front seat height (SHv) per frame size, front fork and front wheel												
SHv [mm]	Aluminium / Titan frame 75° (active)											
	(short frame)						(long frame)					
	 5"	 5"	 6"	 6"	 8"	 8"	 5"	 5"	 6"	 6"	 8"	 8"
450		●5"	●4"	●3"			●3"					
460			●5"	●4"	●3"		●4"	●3"				
470			●5"	●4"	●3"		●4"	●3"				
480			●6"	●5"	●4"	●3"	●4"	●3"				
490				●6"	●5"	●4"	●5"	●4"	●3"			
500					●6"	●5"		●5"	●4"	●3"		
510						●6"		●6"	●5"	●4"	●3"	
520						●6"			●6"	●5"	●4"	
530									●6"	●5"	●4"	
540										●6"	●5"	



Front seat height (SHv) per frame size, front fork and front wheel												
Aluminium / Titan frame 90° (dynamic)												
SHv [mm]	(short frame)						(long frame)					
	5"	5"	6"	6"	8"	8"	5"	5"	6"	6"	8"	8"
450		4"	3"									
460		5"	4"	3"			3"					
470			5"	4"	3"		4"	3"				
480			5"	4"	3"			4"	3"			
490				5"	4"	3"		4"	3"			
500					5"	4"		5"	4"	3"		
510						5"			5"	4"	3"	
520										5"	4"	3"
530										5"	4"	3"
540											5"	4"

Front seat height (SHv) per frame size, front fork and front wheel						
Carbon frame 87°						
UL 400 – 500						
SHv [mm]	5"	5"	6"	6"	8"	8"
450	3"					
460	4"	3"				
480		4"	3"			
490		5"	4"	3"		
500			5"	4"	3"	
510				5"	4"	3"
520					5"	4"
530					5"	4"
540						5"

Rear seat height (SHh)


The following possibilities are available to adjust the rear seat height (SHh):

- Adjustment of the adapter plate to the folding unit, → Chap. Rear wheels; Adapterplate, adjustment to folding unit or assembly of new adapterplate.
- Replacement of the rear wheel with a larger or smaller one.



- i** In the case of a 4° wheel camber, it must be ensured that the rear wheels are parallel after changing the rear seat height. If required, they must be readjusted.
 → Chap. Rear wheels; Adjustment of rear wheel parallelism.

Rear seat height (SHh) by position in the adapter plate and rear wheel size			
SHh [mm]	22" wheel	24" wheel	26" wheel
390	3	1	
400	4	2	
410	4	2	
420	5	3	1
430	6	4	2
440	7	5	3
450		6	4
460		6	4
470		7	5
480			6
490			7



Assumption: Seat angle = 60 mm

Checking the folding unit for play and adjusting

The bolt of the scissor mechanism must be checked regularly for play and tightened, if necessary.

- !** Both scissor levers must lie full on top of one another.
 The scissor levers must not have any play.
 The scissor mechanism must open and close easily.

Difficulty: ●●○

Tools: ● 6 ○ 14

- Check the bolt ① for play.

If the scissor mechanism is too easy to move or the scissor levers have play:

- Loosen the bolt and remove.
- Check the bolt and and replace, if necessary.

- !** **WARNING!** Danger of severe injury to the user as a result of using the wrong bolt.
 When replacing the bolt ①, ensure that only original parts (spare part kit 1540893) are used.

- Insert bolt ① with washer ② through the scissor lever and screw onto the nut ③ ensuring there is no play.

i To obtain a play-free bolt connection, tighten bolt ① firmly and then loosen by a quarter turn.


- Hold bolt ① firmly with an Allen key and tighten counter nut ④.

If the scissor mechanism is too difficult to move:

- Loosen counter nut ④ and bolt ① slightly.
- Tighten counter nut ④ again.

Check:

Always ensure that the screws ③ are tightened in such a way that they are seated on the scissor mechanism when it is open. Please note that the screws ③ are secured with adhesive and this must be replaced as appropriate.



④ → 13 Nm



BACKREST

Tension adjustable backrest

Tension adjustable backrest with push handles standard


RH [mm]	Push-handle	Inter-mediate tube	Velcro straps without stabilisation bar		Velcro straps with stabilisation bar		
	①	②	③	④	③	④	
300	S	S	L	M			S
315	S	S	L	M			S
330	S	S	L	S+M	not		S
345	S	L	L	S+M	possible		M
360	S	L	L	S+M			M
375	S	L	L	L			M
390	S	L	L	L			M
405	L	L	L	L	S	S+L	L
420	L	L	L	S+L	S	S+L	L
435	L	L	L	S+L	S	S+L	L
450	L	L	L	S+L	S	M+L	L
465	L	L	L	S+L	S	M+L	L


Tension adjustable backrest with foldable push handles

RH [mm]	teles-copic tube	Inter-mediate tube	Velcro straps without stabilisation bar		Velcro straps with stabilisation bar		
	①	②	③	④			
300	S	S	L	M			S
315	S	S	L	M			S
330	S	L	L	S+M			S
345	S	L	L	S+M			M
360	S	L	L	S+M			M
375	S	L	L	L	not		M
390	L	L	L	L	possible		M
405	L	L	L	L			L
420	L	L	L	S+L			L
435	L	L	L	S+L			L
450	L	L	L	S+L			L
465	L	L	L	S+L			L





Tension adjustable backrest without push handles

RH [mm]	teles- copic tube	Inter- mediate tube	Velcro straps without stabilisation bar		Velcro straps with stabilisation bar		
	①	②	③	④	③	④	
300	S	S	L	M			S
315	S	S	L	M			S
330	L	S	L	S+M	not		S
345	L	S	L	S+M	possible		M
360	L	S	L	S+M			M
375	L	L	L	L			M
390	L	L	L	L			M
405	L	L	L	L	S	S+L	L
420	L	L	L	S+L	S	S+L	L
435	L	L	L	S+L	S	M+L	L
450	L	L	L	S+L	S	M+L	L
465	L	L	L	S+L	S	M+L	L




Tension adjustable backrest with rearsset push handles height adjustable


RH [mm]	teles- copic tube	Inter- mediate tube	Velcro straps without stabilisation bar		Velcro straps with stabilisation bar		
	①	②	③	④	③	④	
300	S	S	S	M			S
315	S	S	S	M			S
330	S	L	L	M	not		S
345	S	L	L	M	possible		M
360	S	L	L	S+M			M
375	S	L	S	L			M
390	L	L	S	L			M
405	L	L	S	S+L	S	2S+M	L
420	L	L	L	L	S	2S+M	L
435	L	L	L	L	S	S+L	L
450	L	L	L	S+L	S	S+L	L
465	L	L	L	S+L	S	S+L	L





Tension adjustable backrest with push handles, height adjustable

RH [mm]	Push-handle	Backrest-tube	Velcro straps without stabilisation bar		Velcro straps with stabilisation bar		
	①	②	③	④	③	④	
300	S	different for every backrest height	L	M			S
315	S		L	M			S
330	M		L	S+M	not		S
345	M		L	S+M	possible		M
360	M		L	S+M			M
375	M		L	L			M
390	M		L	L			M
405	L		L	L	S	S+L	L
420	L		L	S+L	S	S+L	L
435	L		L	S+L	S	S+L	L
450	L		L	S+L	S	M+L	L
465	L		L	S+L	S	M+L	L





Backrest height (RH)


In order to adjust the backrest height (RH), the backrest tubes have to be fixed in a different position or they have to be exchanged.

Standard backrest cover, height adjustment

Difficulty: ●○○

Tools: ●3 8

- Push the backrest cover so that bolt and nut on the backrest tube are revealed.
- Remove bolt and nut ①.
- Adjust the backrest tube ② to the required height then insert the screw ① into the closest hole and tighten it again.
- Carry out the same setting on both sides.
- Reposition the backrest cover correctly.





Tension adjustable back, height adjustment

Difficulty: ●○○

Tools: ●3 8

- Remove the backrest cover ①.
- Slide the backrest straps either up or down to locate the fixing bolt ②.
- Remove bolt and nut ②.
- Adjust the backrest tube to the required height then insert the screw ② into the closest hole and tighten it again.
- Carry out the same setting on both sides.
- Reposition the backrest cover.

 In the case of major changes of the backrest height (RH), bands may have to be inserted or removed and a larger/smaller backrest cover may be necessary.





Backrest angle (RW)

Backrest, angle adjustment


Difficulty: ●●○

Tools: ●4 □10

- Remove the lower bolt ① from the side supporter ② and move the backrest ③ to the required position. The spacer ④ remains fixed in the latch bolt ⑤.
- Insert bolt and nut ① into the closest hole and tighten it again.
- Carry out the same setting on both sides.

Function control:

The backrest must fold easily and the ratchet bolt must sit tightly against the side supporter.




Push handles / backrest telescopes

Push handle / telescopic tube replacement

Difficulty: ●●○

Tools: ●3 □8 ×2

- Remove the backrest cover ①.
- Remove the screws ③ holding the uppermost backrest band ② (or standard backrest cover) to the push handles ⑤ (or telescopic tubes).
- Slide the backrest straps ④ (or standard backrest cover) either up or down to locate the fixing bolt ⑥.
- Remove bolts and nuts on both sides ⑥.
- Remove push handles ⑤ (or telescopic tubes).
- Slide new push handle through backrest band ② and fix it on the backrest with bolt ⑥.
- Fix the upmost backrest band ② (or standard backrest cover) with screws ③.
- Carry out the same setting on both sides.
- Reposition the backrest cover correctly.






Height adjustable integrated push handles

Difficulty: ●●○

Tools: ●3 ○8 ×2

- Fit the backrest tube ① corresponding to the desired backrest height (RH) on both sides using the bolt ②.
- Slide the backrest bands ③ corresponding to the desired backrest height (RH) and, if required, the end band ④ (or the standard backrest cover) over the backrest tube ①.
- Slide the push handle ⑥ into the backrest tube ① on both sides and secure it at the desired height using the clamp bolt ⑦.
- Secure the end band ④ (or the standard backrest cover) to the backrest tube ① on both sides using the screws ⑤.



Height adjustable rear set push handles

These push handles can only be used in combination with adjustable backrests, not with standard backrests.


Difficulty: ●●○

Tools: ●3 ○8 ×2

- Fit the intermediate tube ① on both sides using the bolts ②.
- Slide the backrest bands ③ corresponding to the desired backrest height (RH) onto the intermediate tube ①.
- Fit the telescopic tube ④ corresponding to the desired backrest height (RH) on both sides using the bolt ⑤.
- Slide the holder ⑥ and the sleeve ⑦ onto the telescopic tube ④ and secure it using the screws ⑧.

For the minimum backrest height, the holder of the rear set push handle must be fitted to the intermediate tube ①. In this case, the sleeve ⑦ is not required.

- Slide the end band ⑨ onto the telescopic tube ④ and secure it on both sides using the screws ⑩.



Replacing foldable push handles

Difficulty: ●○○

Tools: ●4

- Remove the old foldable push handle and the old threaded sleeve.
- Slide the new threaded sleeve ① onto the telescopic tube ②.




For safety reasons, it is important that a new threaded sleeve is used and that the old sleeve is not left in, and reused for, the telescopic tube.

- Slide the new foldable push handle ③ onto the telescopic tube ②.
- Fix the foldable push handle with the screw ④.



The telescopic tube must be replaced as well, if no foldable push handle has been fitted and a foldable push handle is to be newly fitted.





Replacing the handle

An adhesive (e.g. hair spray) is used in these instructions. When applied to the handle, this substance works as a lubricant and as an adhesive once dry.




After drying, the adhesive used must be able to resist a pull-off force of 750 N. If in doubt, contact Kuschall® AG.

Difficulty: ●●○

- Remove the old handle.
- Remove any residue (residual adhesive, grease, dust) from the push handle tube.
- Apply a thin layer of hair spray all over the surface of the push handle tube onto which the handle is to be slid.
- Apply a thin layer of hair spray to the inside of the handle.
- Slide the new handle onto the push handle tube.
- Move the handle into the correct position (grooves facing upwards).



If a long handle has been fitted and this is to be replaced with a short one, the push handle tube must be shortened by 35 mm. The push handle tube must be replaced when switching from a short to a long handle.



Stabilisation bar


Stabilisation bar assembly

Difficulty: ●●○

Tools: ● 3 8 x 2

With a back height RH405 or larger it is possible to assemble a stabilisation bar to increase the rigidity of the backrest handles.

- Remove backrest cover ①, push handles ② and 10 cm backrest band (or end band, if no push handles are assembled).
- Assemble a 5 cm backrest band ③ or end band with the screws ⑩ to the push handles ②.
- Assemble the push handles ② with bolts and nuts ⑫.
- Attach the clamps ④ together with the RH socket ⑥ and the LH socket ⑧ with the screws ⑦ below the backrest band ③ to the push handles ②.
- Replace the push handle/backrest band/socket assembly.
- Press pin ⑤ and slide the stabilisation bar ⑩ into the RH socket ⑥ then swing the stabilisation bar upwards, press pin ⑨ and click the stabilisation bar into the LH socket ⑧.



⑦ → 4 Nm

⑫ → 7 Nm

⑬ → 4 Nm



FOOTRESTS

Lower leg length (UL)

To change the lower leg length, the footrest can be fixed at a higher or lower position, → Chap. Footrests, Footrest, height adjustment.

Short lower leg lengths (UL) can be set using a high-mounted footrest, → Chap. Footrests, Footrest mounted in high position.

Aluminium frame: UL220 – UL310 high mounted footrest
 UL320 – UL390 standard footrest, short frame
 UL400 – UL500 standard footrest, long frame

Titan frame: UL220 – UL340 high mounted footrest
 UL360 – UL390 standard footrest, short frame
 UL400 – UL500 standard footrest, long frame

Carbon frame: UL300 – UL340 high mounted footrest
 UL400 – UL500 standard footrest

Footrest, height adjustment


Difficulty: ●○○

Tools: ●4 8

- Remove on both sides bolt and nut ①, which fix the telescopic tube ② to the frame ③.
- Extend the footrest telescope ② to the required length, then insert the bolts ① into the closest holes.
- Carry out the same setting on both sides.
- Tighten the locking bolts ① on both sides.

Function control

Check that the footrest is firmly attached but that it can fold easily.



Footplate

Footplate replacement


Difficulty: ●○○

Tools: ●4 8

- Remove the screws ① and the screws ②.
- Remove the footplate ③ and replace it by the new one.
- Attach the footplate ③ with the screws ① and the screws ②.

Function control

Check that the footrest is firmly attached but that it can fold easily.






Fitting the footplate set to the front

The footplate is fitted set back as standard.

Difficulty: ●●○

Tools: ●4, 5

- Remove the screws ① and the screws ②.
- Remove the screws and washers ③ on both sides.
- Rotate the footplate mounting ④ by 180°. The elongated side of the footplate mounting ④ is now pointing forwards.
- Fit the footplate mounting ④ again using the screws and washers ③.
- Secure the footplate using the screws ① and the screws ②.
- Set the desired footplate angle and tighten the screws ③, → Chap. Footrests, Angle adjustment of the footplate.



① → 7 Nm

② → 7 Nm

③ → 7 Nm

Function control


Check that the footrest is firmly attached but that it can fold easily.

Angle adjustment of the footplate

Difficulty: ●○○

Tools: ●5

- Slightly loosen the bolts ① under the footplate which secure the telescopic tube ② to the footplate mounting ③.
- Set the footplate to the desired angle.
- Tighten the bolt ① on both sides.



① → 7 Nm

Function control


Check that the footplate can fold easily.

Footrest mounted in high position

Difficulty: ●●○

Tools: ●4, 5, 10

- 1 Remove the footrest from the frame.
- Detach the telescopic tube and footplate fixture ③ from the footrest ① by removing the screw ②.
- 2 Attach the mounting bracket ④ to the frame using the screws and nuts ⑤ in the holes which are the closest to the frame.
- 3 Attach the footrest ① at the desired height on the mounting brackets ④ with the screws ⑥ and the conical spacers ⑦.
- Fit caps ⑧ to the ends of the frame.
- Carry out the same setting on both sides.



⑤ → 7 Nm

⑥ → 7 Nm



SIDE PARTS

Clothes-guard assembly

Difficulty: ●●○

Tools: ●4 ⬡10

- 1 Attach the seat module fixation piece ① to the seat module with bolt and nut ② and screw ③ to the frame.
- Align the clothes guard ④ to the rear wheel and note the best hole on the clothes guard to fix it to the backrest hinge.
- 2 Remove backrest hinge screw ⑤ then re-assemble together with the clothes guard ④.
- 2 Remove the bolt ⑤ from the backrest joint.


! IMPORTANT!

If a pelvic belt is fitted, add the washer ⑥.

- Fit the clothes-guard ④ using the bolt present (and any additional washer) to the backrest joint.

Function control:

Fold and unfold the backrest (this must be done easily). Make sure the clothes-guards are positioned slightly higher than the rear wheels. Check that the clothes-guards flip up easily.



② → 7 Nm

③ → 7 Nm

⑤ → 7 Nm



Mud guard assembly


Difficulty: ●●○

Tools: ●3, 4 ◻10 ×2

If replacing a clothes guard with a mudguard, first remove the clothes guard and the bracket attached to the seat module.

- 1 Attach the bracket to the seat module with the screws ① and remount the rear wheels.
 - 2 Slightly loosen the three screws ② on the adjustment plate and move it along the mounting bracket until the mudguard is in position.
- The position of the mudguard can also be adjusted: To do this, loosen the bolts ③, adjust the mudguard position as required and retighten the bolts ③.
 - Tighten the screws ②.
 - Carry out the setting on both sides.
 - By tightening or loosening the screws ④ you can adjust how easily the armrest can be pulled out or pushed in.

! The distance between the mudguard and the wheel must be either < 8 mm or > 25 mm to prevent fingers from becoming caught between the wheel and the mudguard.



Mudguard size		Hinterrad 22"				Hinterrad 24"				Hinterrad 25"				Hinterrad 26"			
SHh [mm]	S XL L	1	2	3	4	1	2	3	4	1	2	3	4	1	2	3	4
		390	L	L	L	L	XL	XL	XL	XL	-	-	-	-	-	-	-
400	S	S	S	S	L	L	L	L	-	-	-	-	-	-	-	-	-
410	S	S	S	S	L	L	L	L	XL	XL	XL	XL	-	-	-	-	-
420	S	S	S	S	L	L	L	L	XL	XL	XL	XL	XL	XL	XL	XL	XL
430	S	S	S	S	L	L	L	L	XL	XL	L	L	XL	XL	XL	XL	XL
440	S	S	S	S	L	L	S	S	L	L	L	L	XL	XL	XL	XL	XL
450	S	S	S	S	S	S	S	S	S	S	S	S	XL	XL	XL	XL	XL
460	-	-	-	-	S	S	S	S	S	S	S	S	XL	XL	XL	XL	XL
470	-	-	-	-	S	S	S	S	S	S	S	S	L	L	L	L	L
480	-	-	-	-	-	-	-	-	S	S	S	S	L	L	L	L	L
490	-	-	-	-	-	-	-	-	S	S	S	S	S	S	S	S	S






Siderest assembly and adjustment

Difficulty: ●●○

Tools: ●3, 4 \square 10 \times 2

- 1 Attach the bracket to the seat module with the bolts and nuts ①.
- 2 Slightly loosen the three screws ② on the height adjustment bracket then move the siderest to the desired height.
- Tighten the screws ③.
- By tightening or loosening the screws ③, you can adjust how easily the armrest can be pulled out or pushed in.

 A 4° wheel camber requires that an additional set of distance plates be fitted between the adapter plate and the side profile, → Chap. Rear wheels, Adapterplate spacer assembly.



① → 7 Nm


② → 4 Nm

Armrest assembly and adjustment

Difficulty: ●●○

Tools: ●4 \square 10


- 1 Attach the vertical adaptation support ① with the bolts ②.
- 2 Attach the armrest holder ③ with the bolts and nuts ④ to the vertical adaptation support ①.
- 3 Adjust the height of the armrest by changing the position of the screw ⑤.



② → 7 Nm

④ → 7 Nm

⑤ → 4 Nm






FRONT WHEELS

Replacing a front wheel

Difficulty: ●○○

Tools: ●2x3

- 1 Remove the screw ① with disk on one side. Remove the wheel axle ②.
- 2 Remove the front wheel ③.
- 3 Place the sleeves ④ between the new front wheel ③ and the fork.
- 3 Slide the axle ② through the fork, sleeves ④ and the front wheel ③ and fix the axle using the screw ①. Here, use the new screw supplied with the wheel as this screw comes with a thread-locking device, → Chap. Front wheels, Front wheel adjustment, vertical adjustment.



Function control:


The wheel may not wobble, but must rotate easily.

Replacing a front wheel fork

Difficulty: ●●○

Tools: ⬡10

- Remove the sealing cap ① of the bearings block by inserting two screwdrivers into the grooves and flipping it off.
- Remove the nut ② with the washer ⑤.
- Remove the front wheel fork ③.
- Check the ball bearings ④ and replace them if necessary.
- Insert the new front wheel fork with the washer ⑤ and the nut ② and tighten the nut.
- Carry out the function check (see below).
- Replace the sealing cap ①.
- Check the position of the front wheel supporter, → Chap. Front wheels, Front wheel adjustment, vertical adjustment.




Function control:

Tip the wheelchair backwards by 90° so that it is lying on the backrest and the rear wheels. Turn the fork upwards (position A) and let it tip downwards.

The fork has been correctly adjusted if it easily turns to beyond the bottommost point and remains there (position B).

If the fork turns back to the lower position (position C), it has not been sufficiently tightened. There is a risk that the front wheels will start to wobble at high speeds.





Front wheels, repositioning to the frame


The front forks can be mounted in 3 different positions on the front frame. If the distance between footplate and front wheel is too small or if the footrest is mounted in high position, only the two rearmost positions can be used. With the carbon frame and the titan frame 75° having only two pairs of holes, just the rearmost position can be used.

Difficulty: ●●●


Tools: ●5 ⬡10

- Remove screws with washers ① and the nuts ②.
- Place the front wheel supporter ③ in the required position on the frame ④ then replace the screws and washers ① and the nuts ②.
- Check the position of the front wheel supporter, → Chap. Front wheels, Front wheel supporter, vertical adjustment.

The bearing block can be fitted with the holes to the front or with the holes to the back. Fitting with the holes to the back only makes sense if there is enough space between the castor wheel and the footplate.



① → 13 Nm
② → 7 Nm




Front wheel supporter, adjusting the steering error angle

Difficulty: ●●●

Tools: ●5 ⬡10

- 1 Remove nut ② and slightly loosen screw ①.
 - 2 Place a spirit level against the bearing block ③ and turn the disc ④ with an allen key so that the bearing block is perpendicular to the floor.
- Tighten screw ①.
 - Replace and tighten nut ②.



① → 13 Nm
② → 7 Nm

Front wheel supporter, adjusting the trail angle

Difficulty: ●●●

Tools: ●5 ⬡10

The trail angle of the front wheels is set and fixed when assembling the frame, → Chapter Frame, Changing frame.





REAR WHEELS

Rear wheels, repositioning


i The chair will be more stable, the further rearwards that the rear wheels are fitted.

There are 2 different adapter plates available covering 4 possible rear wheel positions. If necessary replace the adapter plate, → Chap. Rear wheels, Adapterplate, adjustment to folding unit or assembly of new adapterplate.


Difficulty: ●○○

Tools: ⬡22 —C19

- Remove the nut ① from the adapter sleeve ②.
- Insert the adapter sleeve ② in the required hole of the adapter plate ③.
- Tighten the nut ① again.
- Carry out the same setting on both sides.



① → 35 Nm



Wheel camber, adapter sleeves

1° adapter sleeves ① or 4° adapter sleeves ② can be used. The wheel camber changes accordingly.



If 4° adapter sleeves are used, the rear wheels must subsequently be adjusted until parallel.

Adjustment of rear wheel parallelism

Difficulty: ●●●



Tools: ⬡22 —C 19

- Measure the distance between the front and rear wheels at axle level.

If this distance is not the same ($x \neq y$):

- Use the open-end spanner to rotate and adjust the adapter sleeves on both sides so that the distance between the front and rear wheels at axle level is identical ($x = y$).

i If a wheelchair has been retrofitted with side parts and/or additional accessories, e.g. an antitipper, the short 1° adapter sleeve must be replaced with the long 1° adapter sleeve.




Adjusting the removable axle

Difficulty: ●●○

Tools: —C 11, 19

- Remove the rear wheel.
- Hold the end of the removable axle ② using the straddle spanner.
- Adjust the length L of the removable axle by turning the nut ①. The length is correctly adjusted if the removable axle engages correctly when fixing the wheel and wheel has just minimal clearance.

i The wheels must be exchanged (left to right side and vice versa) after adjusting both removable axles. The adjustment must now be checked or carried out again to ensure the wheels can be switched.





Distance sleeves


If additional side parts or another type of rear wheel is fitted or if the rear seat height is adjusted, the distance between the rear wheels may have to be increased by fitting additional distance sleeves.

Distance sleeve assembly

Difficulty: ●○○

Tools: ⬡22 — C 19

- Remove rear wheel.
- Loosen and remove adapter sleeve ②, nut and washer ④.
- Slide the distance sleeve ① onto the thread of the adapter sleeve ②.
- Secure the adapter sleeve to the adapter plate again using the nut and washer ④.
- In the case of a 4° wheel camber, the rear wheels must be subsequently adjusted until parallel, → Chap. Rear wheels, Adjustment of rear wheel parallelism.



② → 35 Nm


Adapterplate

Adapterplate, adjustment to folding unit or assembly of new adapterplate

Difficulty: ●●○

Tools: ●5 ⬡10


- 1 Remove the 4 bolts ① from the adapter plate ②.
 - 2 Adjust the adapter plate to the required height on the scissor axle ③. Tighten the bolts ① firmly again.
- Carry out the same setting on both sides.
 - Check the position of the front wheel supporter, → Chap. Front wheels, Front wheel supporter, vertical adjustment.



① → 13 Nm

Function control:

Check that the rear seat height is correct.

-  Note that a rear seat height adjustment changes the seat angle. If necessary the backrest angle can be adjusted accordingly.


Adapterplate spacer assembly


If a side rest is fitted with a 4° wheel camber, an additional set of distance plates must be fitted.

Difficulty: ●●○

Tools: ●5 ⬡10

- Remove the rear wheel.
- Loosen and remove the bolts ③ and remove the adapter plate ②.
- Insert both parts of the distance plate set ① between the adapter plate ② and the side profile.
- Secure the adapter plate ② and the distance plate set ① using the distance plate set bolts ③.

-  Please note that the screws of the distance plate set must be used because they are longer than the original screws.



③ → 13 Nm



BRAKES

Fitting / adjusting the parking brake

Difficulty: ●●○

Tools: ●5



Following each positioning, the rear wheel parking brakes (e.g. when changing the wheel chamber) must be readjusted.



The parking brake function is only guaranteed if the tire has the corresponding air pressure.

■ Check that the rear wheels have sufficient air.

1 Loosely fix the clamping piece 2 to the frame with the screw 1.


2 Slide the brake 3 into the correct position and tighten the screw 1.



When the brake is on, the brake shoe 4 must press approx. 4 mm into the tire.



Furthermore, please note that very little force is required for activating and deactivating the brake. If necessary, a brake lever extension can be fitted.



Visual control

Check that the parking brake is correctly positioned. The brake is correctly adjusted if the brake shoe does not press more than 4 mm into the tire when the brake is on.

Function control

Place the loaded wheelchair on a ramp with a 7° slope with the parking brake on. The wheelchair must not move. Carry out this check with the wheelchair both facing down the ramp and facing up the ramp.



OPTIONS & ACCESSORIES


Antitipper

Antitipper assembly

Difficulty: ●●○

Tools: ●6 ○22 —C 19

- 1 Remove the nut ④.
- 2 Mount the antitipper ① with the adapter sleeve ③, the nut ④, the sleeve ⑤ and the screw ⑥ to the adapter plate ②.

 If an antitipper is newly fitted, the short adapter sleeve ③ must be replaced with a longer one in certain circumstances (only applies to 1° wheel camber).


Height adjustment

- 2 Press the spring pin ⑧. Slide the inner part of the antitipper ⑦ to the required position until the spring pin ⑧ locates in the correct hole.

Function control:



The distance between the antitipper and the ground must be 40 - 60 mm. It must be easy to fold up the antitipper.

Tip the wheelchair backwards using the antitipper until the axle is perpendicular to the antitipper's point of contact with the ground. In this position, the distance between the rear wheel and the ground must be at least 50 mm.



④ → 35 Nm

⑥ → 13 Nm





Tipper aid

Tipper aid assembly and adjustment

Difficulty: ●●○

Tools: ●5

- 1 Attach the mounting block ① onto the folding unit using the 2 bolts ②.
- Attach the horizontal adaptation support ③ to the mounting block ① using the 2 bolts ④.




If a cane holder is mounted on the other side of the wheelchair, then an additional mounting block is needed to lower the tipper aid to prevent them clashing together when the wheelchair is folded.

- 2 Slide the tipper aid ⑤ into the horizontal adaptation support and fasten it with the spring clip ⑥.

Function control:

The tipper aid may not touch the floor if the wheelchair is tipped backwards.

- If the tipper aid touches the floor remove it and change the position of the spring clip ⑥. This is done by removing the tube end cap then pushing the spring clip ⑥ one hole further into the tube using a screwdriver.




Transit wheels

Transit wheels assembly and adjustment

Difficulty: ●●○

Tools: ●4 ○10

- 1 Attach the vertical adaptation support ① to the chair with the 4 bolts ②.
- 2 Place the transit wheel ③ against the vertical adapter support ① and attach it with the 2 bolts ④.
- Adjust transit wheels ③ to the correct height with the spring clip ⑤.





Cane holder

Cane holder assembly and adjustment

Difficulty: ●●○

Tools: ●5 ×2

1 Attach the mounting block ① onto the wheelchair using the 2 bolts ②.

■ Attach the horizontal adaptation support ③ to the mounting block ① using the 2 bolts ④.

i If a tipper aid is mounted on the other side of the wheelchair, then an additional mounting block is needed to lower the cane holder to prevent them clashing together when the wheelchair is folded.


2 Slide the cane holder ⑤ into the horizontal adaptation support ③ and fasten it with the spring clip ⑥.

■ Remove the screw ⑧ from the backrest cover then replace it together with the strap ⑦.

Function control:

The cane holder must not touch the floor when the wheelchair is tipped backwards.

■ If the cane holder touches the floor remove it and change the position of the spring clip ⑥. This is done by removing the tube end cap then pushing the spring clip ⑥ one hole further into the tube using a screwdriver.



Fitting the pelvic belt


Difficulty: ●○○

Tools: Plastic hammer, ●5, ⑩

■ Fit the steel strap ① with the bolt ③ on both sides of the side profile ②, pressing the upper part of the steel strap inwards, e.g. using a plastic hammer.

! Ensure that the webbing of the pelvic belt is not twisted during assembly and the locking mechanism shows towards the front.

! Ensure that the seat cover does not get caught between the steel strap and the side supporters during assembly.





Küschall AG
Benkenstrasse 260
CH-4108 Witterswil
kueschall@invacare.com
www.kuschall.com

Service Manual
Küschall Champion
ENGLISH | 2013-12

küschall[®] distributors

Belgium & Luxemburg: Invacare nv • Autobaan 22 • B-8210 Loppem

Tel: (32) (0)50 83 10 10 • Fax: (32) (0)50 83 10 11 • belgium@invacare.com

Danmark: Invacare A/S • Sdr. Ringvej 37 • DK-2605 Brøndby

Tel: (45) (0)36 90 00 00 • Fax: (45) (0)36 90 00 01 • denmark@invacare.com

Deutschland: Invacare GmbH • Alemannenstraße 10 • D-88316 Isny

Tel: (49) (0)75 62 7 00 0 • Fax: (49) (0)75 62 7 00 66 • kontakt@invacare.com

European Distributor Organisation: Invacare • Kleiststraße 49 • D-32457 Porta Westfalica

Tel: (49) (0)57 31 754 540 • Fax: (49) (0)57 31 754 541 • edo@invacare.com

España: Invacare SA • c/Areny s/n • Poligon Industrial de Celrà • E-17460 Celrà (Girona)

Tel: (34) (0)972 49 32 00 • Fax: (34) (0)972 49 32 20 • contactsp@invacare.com

France: Invacare Poirier SAS • Route de St Roch • F-37230 Fondettes

Tel: (33) (0)2 47 62 64 66 • Fax: (33) (0)2 47 42 12 24 • contactfr@invacare.com

Ireland: Invacare Ireland Ltd • Unit 5 Seatown Business Campus • Seatown Road • Swords • County Dublin – Ireland

Tel : (353) 1 810 7084 • Fax: (353) 1 810 7085 • ireland@invacare.com

Italia: Invacare Mecc San s.r.l. • Via dei Pini 62 • I-36016 Thiene (VI)

Tel: (39) 0445 38 00 59 • Fax: (39) 0445 38 00 34 • italia@invacare.com

Nederland: Invacare BV • Celsiusstraat 46 • NL-6716 BZ Ede

Tel: (31) (0)318 695 757 • Fax: (31) (0)318 695 758 • nederland@invacare.com • csede@invacare.com

Norge: Invacare AS • Grensesvingen 9 • Postboks 6230 • Etterstad • N-0603 Oslo

Tel: (47) (0)22 57 95 00 • Fax: (47) (0)22 57 95 01 • norway@invacare.com • island@invacare.com

Österreich: Invacare Austria GmbH • Herzog Odilostrasse 101 • A-5310 Mondsee

Tel.: (43) 6232 5535 0 • Fax.: (43) 6232 5535 4 • info@invacare-austria.com

Portugal: Invacare Lda • Rua Estrada Velha • 949 • P-4465-784 Leça do Balio

Tel: (351) (0)225 1059 46/47 • Fax: (351) (0)225 1057 39 • portugal@invacare.com

Sverige & Suomi: Invacare AB • Fagerstagatan 9 • S-163 91 Spånga

Tel: (46) (0)8 761 70 90 • Fax: (46) (0)8 761 81 08 • sweden@invacare.com • finland@invacare.com

Switzerland: Invacare AG • Benkenstrasse 260 • CH-4108 Witterswil

Tel.: (41) (0)61 487 70 80 • Fax.: (41) (0)61 487 70 81 • switzerland@invacare.com

United Kingdom: Invacare Limited • Pencoed Technology Park,

Pencoed, Bridgend CF35 5AQ • Switchboard Tel: (44) (0)1656 776 200, Fax: (44) (0)1656 776 201 •

Customer services Tel: (44) (0) 1656 776 222 • Fax: (44) (0) 1656 776 220

Australia: Invacare Australia Pty Ltd • ABN 45 074 676 378, PO Box 5002, 1 Lenton Place, North Rocks, NSW 2151, Australia •

Freephone: 1800 069 042, Fax: 02 8839 5353 • E-mail: sales@invacare.com.au • Web: www.invacare.com.au

New Zealand: Invacare New Zealand • PO Box 62–124, 4 Westfield Place, Mt. Wellington, Auckland, New Zealand •

Freephone: 8000 468 222, Freefax: 0800 807 788 • E-mail: sales@invacare.co.nz • Web: www.invacare.co.nz