# Invacare® Slings

|    |                                             |                                |    |
|----|---------------------------------------------|--------------------------------|----|
| EN | <b>Amputee</b>                              | User Manual . . . . .          | 1  |
| DA | <b>Amputationssejl</b>                      | Brugsanvisning . . . . .       | 3  |
| DE | <b>Patientengurt für Amputierte</b>         | Gebrauchsanweisung . . . . .   | 4  |
| ES | <b>Eslinga Amputados</b>                    | Manual del usuario . . . . .   | 6  |
| FR | <b>Sangle amputée</b>                       | Manuel d'utilisation . . . . . | 7  |
| IT | <b>Amputee</b>                              | Manuale d'uso . . . . .        | 8  |
| NL | <b>Liftmat Amputatie</b>                    | Gebruiksaanwijzing . . . . .   | 10 |
| NO | <b>Amputasjons løfteseil medium rygg</b>    | Bruksanvisning . . . . .       | 11 |
| PT | <b>Amputee</b>                              | Manual de utilização . . . . . | 12 |
| SV | <b>Amputee, Lyftsele amp mellanhög rygg</b> | Bruksanvisning . . . . .       | 14 |


**EN**

## I General


### I.1 Figures

- Figur
- Figur
- Figur
- Abbildung
- Afbeelding
- Figures
- Figura
- Figura
- Figura


1


2


3

4

5

This manual **MUST** be given to the user of the product.  
BEFORE using this product, read this manual and save for future reference.

**Yes, you can.®**

## 1.2 Symbols

|                                                                                    |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                    |
|------------------------------------------------------------------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
|    | Signal words are used in this manual and apply to hazards or unsafe practices which could result in personal injury or property damage. See the information below for definitions of the signal words.<br><b>WARNING!</b> Warning indicates a potentially hazardous situation which, if not avoided, could result in death or serious injury.<br><b>CAUTION!</b> Caution indicates a potentially hazardous situation which, if not avoided, may result in property damage or minor injury or both. |
|   | Gives useful tips, recommendations and information for efficient, trouble-free use.                                                                                                                                                                                                                                                                                                                                                                                                                |
|   | <b>IMPORTANT</b> indicates a hazardous situation that could result in damage to property if it is not avoided.                                                                                                                                                                                                                                                                                                                                                                                     |
|    | This product complies with Directive 93/42/EEC concerning medical devices.<br>The launch date of this product is stated in the CE declaration of conformity.                                                                                                                                                                                                                                                                                                                                       |
|    | Safe Working Load                                                                                                                                                                                                                                                                                                                                                                                                                                                                                  |
|    | Consult the user manual.                                                                                                                                                                                                                                                                                                                                                                                                                                                                           |
|    | DO NOT exceed 95°C washing temperature.                                                                                                                                                                                                                                                                                                                                                                                                                                                            |
|   | DO NOT bleach.                                                                                                                                                                                                                                                                                                                                                                                                                                                                                     |
|  | Tumble Dry, Normal, Low Heat                                                                                                                                                                                                                                                                                                                                                                                                                                                                       |
|  | DO NOT iron.                                                                                                                                                                                                                                                                                                                                                                                                                                                                                       |

## 1.3 Intended Use

The Invacare Amputee sling is designed to provide maximum support and security for either unilateral or bilateral amputees. For safe use of this sling, the client must have adequate head control. The sling allows easy access for personal hygiene and is easily applied in both sitting and lying positions.

## 1.4 Service Life

All Invacare slings are guaranteed for one year. Expected service life of the product is 1-5 years. Service life time varies depending on the fabric, frequency of use, laundering practices, and weight carried.

## 2 Safety

### 2.1 General Guidelines

#### **WARNING!**

- DO NOT use this product or any available optional equipment without first completely reading and understanding these instructions and any additional instructional material such as user manuals, service manuals or instruction sheets supplied with this product or optional equipment. If you are unable to understand the warnings, cautions or instructions, contact a healthcare professional, dealer or technical personnel before attempting to use this equipment - otherwise, injury or damage may occur.
- Please check your sling on delivery. If the sling is not in perfect condition, please contact your supplier.

## 2.2 Safe Working Load

|           | <b>Solid</b> |  |
|-----------|--------------|-------------------------------------------------------------------------------------|
| <b>S</b>  | 1485367      | 200 kg                                                                              |
| <b>M</b>  | 1485368      | 200 kg                                                                              |
| <b>L</b>  | 1485369      | 200 kg                                                                              |
| <b>XL</b> | 1485370      | 250 kg                                                                              |

## 3 Usage

### 3.1 Sling application

#### **WARNING!**

- Slings should be checked regularly for wear and damage (in accordance with each country's regulations). NEVER use a damaged sling.
- If you are in any doubt, please do not attempt to use the sling and seek the advice of a professional or your supplier.
- Invacare recommends that a risk assessment is always carried out by a professional prior to issuing lifting equipment. Always choose the sling design and size according to the client's weight, size and physical ability.
- Ensure that the client is not prone to spasm and there has been a risk assessment completed to address the possibility of forward pitching.

### Sling application — Sitting in chair or wheelchair

1. Lean the client forward ensuring that he is well supported to reduce risk of falling and to promote feelings of security.
2. Slide the sling down behind the back of the client in a central position with the aperture apex aligned with the base of the spine.
3. Sit the client back against the sling. **Figure 1.**
4. Pass the leg support towards the front of the chair to each side of the thigh. **Figure 2.**
5. Then place each leg strap under both thighs — one after each other.
6. Place the strap of the outer leg support in the strap on the side of the sling. **Figure 3.**

### Sling application — Lying in bed or floor

If the client is able to sit up from the bed or floor, this will ease the application.

1. Lean the client forward ensuring that they are well supported to reduce risk of falling and to promote feelings of security.
2. Slide sling down behind the back of the client in a central position with the aperture apex aligned with the base of the spine.
3. Sit the client back against the sling. **Figure 4.**

Alternatively, if the client does not have sufficient trunk control, roll the sling to mid-way point and roll the client to one side to position the rolled up sling under the lifted side. Ensure the center of the sling (green mark) is placed centrally along the spine. Roll the client to the other side, unroll the sling and roll the client back into position on the sling. **Figure 5.**

4. One at a time, bend the knees, lift the legs and apply the sling under the legs.
5. Both leg straps must be placed under each leg. **Figure 3.**

## 3.2 Attach loops to the spreader bar or lifting strap

### WARNING!

- Use only CE marked spreader bars and hoists, preferably with a rated Safe Working Load equal to or higher than that of the sling. If the SWL is lower than that of the sling, the SWL of the spreader bar or hoist is applicable.
- Always check that all the straps are properly attached to the spreader bar or lifting strap before attempting to lift the client.

1. Attach the loops to the spreader bar or lifting straps of the lifting equipment.

Attach the shoulder- and head loops before the leg hoops.

The sling can be used with spreader bars with 2 attachment points.

For security and comfort reasons, the size (width) of the spreader bar must correspond with the size of the sling.

The position of the client in the sling can be adjusted by using the different adjustment straps. Shorter loops at the shoulders will produce a more vertical lift, which will assist with positioning into a chair or wheelchair. By lowering the loops at the shoulders, a more reclined position can be achieved which is more suitable for a chair to bed transfer.

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© 2013 Invacare® Corporation. All rights reserved. Reproduction, duplication or modification in whole or in part is prohibited without prior written permission from Invacare. Trademarks are identified by ™ and ®. All trademarks are owned by or licensed to Invacare Corporation or its subsidiaries unless otherwise noted.
#### Invacare Corporation