



INVACARE® TOP END®  
**ELIMINATOR™ OSR**  
RACING WHEELCHAIR



The **Invacare® Top End® Eliminator™ OSR Racing Wheelchair** is designed for speed. The super-stiff roll-formed horizontal maintube, combined with a relaxed headtube angle for exceptional stability and aerodynamics, make it the choice of champion racers like Ernst van Dyk, Krige Schabort, Kurt Fearnley, Chantal Petitclerc and Edith Hunkeler. So, if you have your sights set on winning, the Eliminator OSR racing wheelchair will help you get to the finish line first.

## Features

6061 T6 oval tube custom aluminum frame  
Integral wrap-around fenders  
4 cage designs: V, Open-V, U or I-cage  
Kneeling, amputee or traditional upholstery  
700C or 26" high performance rear wheels  
Caliper hand brake  
Aluminum side panels  
Straight axle tube with rear wheel alignment  
Aluminum fork  
Precision track compensator system

20" front wheel with tubular tire  
Rubber-coated 3/8" diameter aluminum rims  
Aero bars  
Powder coat paint  
250 lb. weight capacity



1. I-cage
2. Open-V cage shown with cut-out axle option
3. V cage
4. U cage shown with knee plate
5. Oversized rake for maximum speed and cornering
6. Precision track compensator system



## Options

Carbon fiber rear wheels  
Corima™ carbon fiber front wheel  
Alignment gauge  
Wheel bag  
Racing suit  
Solid seat/kneel plate  
Fixed footrest

Cushion  
Cordless speedometer  
Helmet  
18" front wheel

**Invacare Corporation**  
[www.invacare.com](http://www.invacare.com)

## USA

One Invacare Way  
Elyria, Ohio  
44035-4190  
(800) 333-6900

## Canada

570 Matheson Blvd. E.,  
Unit 8  
Mississauga, Ontario  
L4Z 4G4 Canada  
(800) 668-5324

## Top End

[www.topendwheelchair.com](http://www.topendwheelchair.com)  
Sports and Recreation  
Products  
4501 63rd Circle North  
Pinellas Park, Florida  
33781-5914  
(727) 522-8677  
(800) 532-8677  
Fax: (727) 522-1007

©2013 Invacare Corporation.  
All rights reserved.  
Trademarks are identified  
by the symbols ™ and ®. All trademarks are  
owned by or licensed to  
Invacare Corporation unless  
otherwise noted. New York  
City Marathon is a registered  
trademark of New York Road  
Runners Club, Inc. Corima  
is a registered trademark  
of Corima, S.A. Zipp is  
a registered trademark  
of Compositel, Inc.  
Specifications are subject to  
change without notification.  
Form No. 01-192  
Rev 02/13  
130228

