



*küschall*<sup>®</sup> Ultra-Light

Active wheelchair



# Ultra-Light

## Content

- Frame.....2
- Backrest.....3
- Seat upholstery.....8
- Footrest.....9
- Side guards.....12
- Cushions and options.....16
- Colours.....17
- Front wheels.....18
- Rear wheels.....21
- Brakes.....27
- Safety options.....29
- Accessories.....30
- Küschall lifestyle collection.....32

## Frame

---

DDT0203      Active frame, 70°



---

DDT0212      Dynamic frame, 80°



## Backrest

### Backrest upholstery

DDT0700

Lightly upholstered backrest



DDT0701

Velcro backrest, tension adjustable



DDT0735

Aero backrest, tension adjustable



**Push handles**

DDT0750

Without push handles



DDT0702  
DDT0751

Push handles fixed,  
short or long



DDT0703

Push handles, height adjustable,  
rear set



## Ultra-Light - Backrest

---

DDT0733

Push handles, height adjustable,  
integrated



---

DDT0704

Push handles, foldable



## Back tube

DDT0705

Lumbal back



DDT0734

Foldable back



DDT0707

Back angle adjustable

DDT0718

Possible back angles:


DDT0709

$-8^\circ / -4^\circ / 0^\circ / +4^\circ / +8^\circ / +12^\circ$

DDT0719

DDT0711

DDT0712



## Backrest options

DDT1491

Backrest stabilization bar, swing away and removable






## Seat upholstery


DDT0611

Standard seat upholstery




DDT0613

Seat upholstery with central Velcro® band  
(soft band)



DDT0614

Seat upholstery with lateral Velcro® bands  
(rough band)



## Footrest

### Footrest

DDT1562  
DDT1563

Two pieces footrest,  
set to the front or to the back




DDT1564  
DDT1565

Two pieces footrest angle adjustable,  
set to the front or to the back



DDT1566  
DDT1567

One piece footrest flip-up,  
angle adjustable,  
fixed on one side when folded,  
set to the back or to the front



## Ultra-Light - Footrest

DDT1568  
DDT1569

One piece footrest flip-up,  
angle adjustable,  
fixed on both sides when folded,  
set to the back or to the front



DDT1556

Two pieces footrest high mounted, plastic



**Foot and leg support**

DDT1533

Calf-strap



DDT1530

Heel loops



DDT1536

Instep loops




## Side guards

### Clothes guards / Mud guards


DDT1400

Clothes guards, fixed



DDT1402


Mud guards, fixed



DDT1403

Mud guards, removable

*(This option adds 2 cm to the total width of the wheelchair)*




## Ultra-Light - Side guards

---

DDT1405

Mud guards carbon, fixed




---

DDT1408

Mud guards carbon, removable

*(This option adds 2 cm to the total width of the wheelchair)*



## Siderests

DDT1416  
DDT1421

Siderest height adjustable (with tools),  
removable, short or long


*(This option adds 2 cm to the total width of  
the wheelchair)*



DDT1430  
DDT1431

Siderest flip-back height adjustable,  
removable, short or long

*(This option adds 2 cm to the total width of  
the wheelchair)*



DDT1440  
DDT1441

Siderest with desk function, flip-back,  
height adjustable, removable,  
short or long

*(This option adds 2 cm to the total width of  
the wheelchair)*



## Ultra-Light - Side guards

---

DDT1411

Tubular armrest, height adjustable,  
removable, swing away, comfort padded





## Cushions and options

---

DDT0660

Flo-tech Lite cushion



DDT0661

Flo-tech Lite Visco cushion












DDT0606

Seat bag, click in and removable



## Colours

### Frame colours

Glossy black		Black bronze	
Royal blue		Royal red	
Mystic purple		Pearly sapphire	
Night shimmer		Cyber green	
Silver grey			

### Color-Power-Set (CPS)

Techno-anthracite		Lemon	
B-blue		Comet red	
Apricot			

*CPS-Set includes: Starec & High Performance rear wheel hubs with quick release protection, 22" and 24" Starec rear wheel rims, adapter plates, forks, castor housing when option is selected.*

## Front wheels

<p>DDT0800 3" Supersport</p>		<p>DDT0824 4" Skater wheel</p>	
<p>Starec wheel, available in 4" and 5" and in various CPS colours</p>		<p>DDT0811 5" Softroll</p>	
<p>DDT0810 5" Light wheel</p>		<p>DDT0820 6" Light wheel</p>	
<p>DDT0801 4" Everyday wheel</p>		<p>Pneumatic wheel, available in 6" and 7"</p>	

## Front wheel attachment

DDT0867 Standard (flow system)



DDT0868 Bearing blocks  
*(bringing the front wheels 2 cm outwards)*



## Front wheel extension

DDT0869      Front wheel extension



## Quick release fork

DDT0870      Quick release fork



# Rear wheels

## Rear wheel fixation and position


DDT0250      Rear wheel standard fixation



Rear wheel positions



DDT0249      Rear wheel extension



Rear wheel positions



DDT0909      Quad or Tetra quick release axle



## Rear wheel type

DDT0915

22", 24", 26" Kuschall Starec wheel



DDT0908  
DDT0918

24", 25" High Performance wheel



DDT0913

24" Mountain Bike wheel



Ultra-Light – Rear wheels

---

DDT0917

24", 25" Spinergy wheel



---

DDT0904

24" Spider wheel





## Rear wheels tyres

DDT1001  
Light wheel 7 bar  
22", 24", 26"



DDT1002  
Profile tyre 7 bar  
22", 24"



DDT1003  
Sports hall  
22", 24", 26"



DDT1004  
Solid tyre, profile  
24"



DDT1005  
DDT1006  
Solid tyre (KIK Type),  
grey or black,  
24"



DDT1010  
Marathon Evolution Plus  
24", 25", 26"



**Handrims**

DDT1102  
DDT1103  
Supergripp  
22", 24", 26"



DDT1104  
DDT1105  
Steel  
24", 26"



DDT1006  
DDT1007  
Titanium  
24"



DDT1009  
Maxgripp  
24"



DDT1110  
Tetra  
24", 26"



DDT1113  
DDT1114  
Light Alu  
24", 25"



## Spoke guards

DDT1202  
24" küsch



DDT1204  
24" Boa-Boa



DDT1207  
24" «Original» Carbon



DDT1210  
24" Swiss




DDT1211  
24" Pirates



# Brakes


DDT1300

Push brake, lever bent




DDT1301

Push brake, lever straight



DDT1306

Pull brake




## Ultra-Light – Brakes

---

DDT1302


Performance brake



---

DDT1303

Active brake



## Safety options

---

DDT1702  
DDT1703  
DDT1704


Antitipper, height adjustable, active and  
inactive position, removable




## Accessories

---

DDT1706  
DDT1707 Cane-holder,  
left or right



DDT1708/1709 Tipper-Aid,  
left, right or both sides



DDT1716 AirMan



## Ultra-Light – Accessories

---

DDT1728

Wheel bag





## Kuschall lifestyle collection

---

1485632  
1468279  
1451626  
1451627

Kuschall T-Shirt  
(available in S, M, L, XL, XXL)



1513357  
1513358  
1513359  
1513360  
1513361

Kuschall Jacket  
(available in S, M, L, XL, XXL)



1518297  
1518298  
1518299  
1518300  
1518301

Kuschall Sweat-shirt  
(available in S, M, L, XL, XXL)



## Ultra-Light – Kuschall lifestyle collection

---

1486608  
1486609  
1486610  
1486611  
1486612

Kuschall Body warmer  
(available in S, M, L, XL, XXL)



---

1486579                      Kuschall lanyard



---

1486578                      Kuschall Pen



## Ultra-Light – Kuschall lifestyle collection

1486577

Kuschall Carbon Pen Set



1486614

Kuschall Umbrella



1486584

Mentor Bag

1486585

(available in English, French, German)

1486583





ENGLISH

12/2008