

POST 3

Invacare® XLT  
Swing



Kunde/HMS:  
Kontaktperson:  
Adresse:  
  
Telefon:  
Mailadresse:  
Bestillingnr./SF nr:



☐ please make your selection with a cross in the box ◆ Standard feature ○ Free of Charge (NC9999) non conformance with option 9999

LPF\_NO\_XLT Swing 160101

→ Max Customer weight = 135 Kg

Code	Description	Non conformance	
BRA0010	Invacare XLT Swing		◆
BRA8000	Special (Full description and drawing by fax or E-mail)		☐
<b>PROPULSION TYPE</b>			
BRA0020	Standard wheel hub		◆
BRA0030	Assistant manoeuvred drumbrakes	0400;0695;0696;0697;0699;0710;0860;0870;0880;0930;1560;1570;1571;1580;1585;1590;1600;1660;1760;1990;2000	☐
BRA0040	One arm drive, LH	0055;0310;0320;0400;1540;1560;1570;1571;1580;1585;1590;1600;1610;1640;1660;1670;1680;1690;1700;1710;1730;1740;1760;1770;1820;1850;1871;1881;1940;1950;1990;2000	☐
BRA0050	One arm drive, RH	0055;0310;0320;0400;1540;1560;1570;1571;1580;1585;1590;1600;1610;1640;1660;1670;1680;1690;1700;1710;1730;1740;1760;1770;1810;1840;1871;1881;1940;1950;1990;2000	☐
<b>SEAT DIMENSIONS</b>			
<b>SEAT WIDTH</b>			
BRA0055	330 mm	0150;0438;0439;0440;0695;0696;0697;0699;1000;1020;1040;1060;1100;1110;1130;1360;1370;1810;1820;1900;1910;1920;1990;2000	☐
BRA0060	355 mm	0438;0439;0440;0856;1020;1060;1100;1110;1130;1360;1370;1810;1820;1990;2000	☐
BRA0070	380 mm	0856	☐
BRA0080	405 mm	0856	☐
BRA0090	430 mm	0856	◆
BRA0100	455 mm	0856	☐
BRA0110	480 mm	0856	☐
BRA0120	505 mm	0852	☐
<b>SEAT DEPTH</b>			
BRA0130	400 mm (adjustable from seat depth 360 to 400 mm)	1800	☐
BRA0140	450 mm (adjustable from seat depth 410 to 450 mm)		◆
BRA0150	500 mm (adjustable from seat depth 460 to 500 mm)		☐
<b>SEAT HEIGHT</b>			
<b>FRONT SEAT TO FLOOR HEIGHT = SHf</b>			
BRA0160	380 mm (Available as a Special- Launch date to be confirmed later.)	0320;0330;0340;0350;0360;0370;0380;0390;0400;0490;0500;0510;0520;0530;0540;0550;0560;0570;0580;0590;0600;0610;0620;0630;0640;1430;1440;1450;1470;1485;1490;1500;1510;1520	☐
BRA0170	390 mm	0330;0340;0350;0360;0370;0380;0390;0400;0500;0510;0520;0530;0540;0550;0560;0570;0580;0590;0600;0610;0620;0630;0640;1450;1460;1470;1485;1490;1500;1510;1520	☐
BRA0180	400 mm	0340;0350;0360;0370;0380;0390;0400;0510;0520;0530;0540;0550;0560;0570;0580;0590;0600;0610;0620;0630;0640;1450;1460;1470;1485;1490;1500;1510;1520	☐
BRA0190	410 mm	0360;0370;0380;0390;0400;0520;0530;0540;0550;0560;0570;0580;0590;0600;0610;0620;0630;0640;1440;1460;1470;1485;1490;1500;1510;1520	☐
BRA0200	420 mm	0360;0370;0380;0390;0400;0430;0530;0540;0550;0560;0570;0580;0590;0600;0610;0620;0630;0640;1440;1460;1470;1485;1490;1500;1510;1520	☐
BRA0210	430 mm	0370;0380;0390;0400;0540;0550;0560;0570;0580;0590;0600;0610;0620;0630;0640;1470;1485;1490;1500;1510;1520	☐
BRA0220	440 mm	0380;0390;0400;0550;0560;0570;0580;0590;0600;0610;0620;0630;0640;1485;1490;1500;1510;1520	☐
BRA0230	450 mm	0390;0400;0560;0570;0580;0590;0600;0610;0620;0630;0640;1440;1460;1470;1485;1490;1500;1510;1520	◆
BRA0240	460 mm	0310;0570;0580;0590;0600;0610;0620;0630;0640;1440;1450;1460;1485;1490;1500;1510;1520	☐
BRA0250	470 mm	0310;0320;0580;0590;0600;0610;0620;0630;0640;1510;1520	☐
BRA0260	480 mm	0310;0320;0590;0600;0610;0620;0630;0640;	☐
BRA0270	490 mm	0310;0320;0330;0600;0610;0620;0630;0640;1440;1460;1485;1490;1500	☐
BRA0280	500 mm	0310;0320;0330;0610;0620;0630;0640;1440;1450	☐
BRA0290	510 mm	0310;0320;0330;0340;0350;0620;0630;0640;1440;1450	☐
BRA0300	520 mm	0310;0320;0330;0340;0350;0360;1440;1450;1460;1470	☐
<b>REAR SEAT TO FLOOR HEIGHT =SHr, SHr ≤ SHf</b>			
BRA0310	380 mm	1550;1560;1570;1571;1580;1585;1590;1600	☐
BRA0320	390 mm	1550;1560;1570;1580;1585;1590;1600	☐
BRA0330	400 mm	1600	☐
BRA0340	410 mm		☐
BRA0350	420 mm		☐
BRA0360	430 mm		☐
BRA0370	440 mm		◆
BRA0380	450 mm	1540	☐
BRA0390	460 mm	1540	☐
BRA0400	470 mm	1540;1550;1560;1570;1580;1585;1590	☐

**SHf / Castor Size**

3" : 390-480 mm  
4" : 410-490 mm  
5" : 380, 430-510 mm  
5,5" : 440-510 mm  
6" : 450-520 mm  
7" : 480-520 mm

**SHr / Rear wheel size**

22" : 380-440 mm  
24" : 400-460 mm  
25" : 400-460 mm  
26" : 410-470 mm



<b>KNEE TO HEEL LENGTH, without seat cushion, UL ≤ SHf-60 mm</b>		
BRA0438	260 mm (Only with hanger fixed 90°)	1000;1010;1040;1050
BRA0439	270 mm (Only with hanger fixed 90°)	1000;1010;1040;1050
BRA0440	280 mm (Only with hanger fixed 90°)	1000;1010;1040;1050
BRA0450	290 mm	
BRA0460	300 mm	
BRA0470	310 mm	
BRA0480	320 mm	
BRA0490	330 mm	
BRA0500	340 mm	
BRA0510	350 mm	
BRA0520	360 mm	
BRA0530	370 mm (Only with Hanger angle adjustable or Hanger fixed 80°)	1020;1060
BRA0540	380 mm (Only with Hanger angle adjustable or Hanger fixed 80°)	1020;1060
BRA0550	390 mm (Only with Hanger angle adjustable or Hanger fixed 80°)	1020;1060
BRA0560	400 mm	◆
BRA0570	410 mm	
BRA0580	420 mm	
BRA0590	430 mm	
BRA0600	440 mm	
BRA0610	450 mm	
BRA0620	460 mm	
BRA0630	470 mm (Only with Hanger angle adjustable)	1010;1020;1050;1060
BRA0640	480 mm (Only with Hanger angle adjustable)	1010;1020;1050;1060

**BACKREST**

<b>BACKREST TYPE</b>		
BRA0680	Velcro tension adjustable, fixed height, infinitely variable angle adjustable, foldable and click in backrest (bent backrest tubes)	0750;0760;0770;0780;0790;0800;0810;0820;0830;0840;0850;0852;0854;0856;1311;1312;1313;1314;1315;1316;1317;1320;1330;1335;1336
BRA0685	Velcro tension adjustable, fixed height, infinitely variable angle adjustable, foldable and click in backrest (straight backrest tubes)	0750;0760;0770;0780;0790;0800;0810;0820;0830;0840;0850;0852;0854;0856;1311;1312;1313;1314;1315;1316;1317;1320;1330;1335;1336
BRA0690	Velcro tension adjustable, height adjustable, infinitely variable angle adjustable, foldable and click in backrest (bent backrest tubes)	0710;0720;0730;0740;0852;0854;0856;1311;1312;1313;1314;1315;1316;1317
BRA0693	Velcro tension adjustable, height adjustable, infinitely variable angle adjustable, foldable and click in backrest (straight backrest tubes)	0710;0720;0730;0740;0852;0854;0856;1311;1312;1313;1314;1315;1316;1317
BRA0695	Narrow back, velcro tension adjustable, fixed height, infinitely variable angle adjustable, foldable and click in backrest	0710;0720;0730;0750;0760;0770;0780;0790;0800;0810;0820;0830;0840;0850;0852;0854;0856;0880;0890;0900;0910;0920;1311;1312;1313;1314;1315;1316;1317;1320;1330;1335;1336
BRA0696	MatrX PB Elite std, incl hardware	0710;0720;0730;0740;0750;0760;0770;0780;0790;0800;0810;0820;0830;0840;0850;0870;0880;0890;0900;0910;0920;1310;1320;1330;1335;1336
BRA0697	MatrX PB Elite deep, incl hardware	0710;0720;0730;0740;0750;0760;0770;0780;0790;0800;0810;0820;0830;0840;0850;0870;0880;0890;0900;0910;0920;1310;1312;1313;1314;1315;1316;1317;1320;1330;1335;1336
BRA0699	Without Backrest, prepared for MatrX	0710;0720;0730;0740;0750;0760;0770;0780;0790;0800;0810;0820;0830;0840;0850;0852;0854;0856;0870;0880;0890;0900;0910;0920;1310;1311;1312;1313;1314;1315;1316;1317;1320;1330;1335;1336

**BACKREST HEIGHT**  
PLEASE NOTE THAT THE XLT CANNOT BE USED AS A SEAT DURING TRANSPORTATION IF THE BACK HEIGHT IS LOWER THAN 400 MM )

<b>FIXED BACKREST</b>		
BRA0710	200 mm	0880;0890;0900;0910;0920
BRA0720	250 mm	
BRA0730	300 mm	
BRA0740	350 mm	

<b>HEIGHT ADJUSTABLE BACKREST</b>		
BRA0750	300 mm (Height adjustable from 300 to 320 mm by 20 mm steps)	□
BRA0760	320 mm (Height adjustable from 300 to 320 mm by 20 mm steps)	□
BRA0770	350 mm (Height adjustable from 350 to 390 mm by 20 mm steps)	□
BRA0780	370 mm (Height adjustable from 350 to 390 mm by 20 mm steps)	□
BRA0790	390 mm (Height adjustable from 350 to 390 mm by 20 mm steps)	□
BRA0800	400 mm (Height adjustable from 400 to 440 mm by 20 mm steps)	◆
BRA0810	420 mm (Height adjustable from 400 to 440 mm by 20 mm steps)	□
BRA0820	440 mm (Height adjustable from 400 to 440 mm by 20 mm steps)	□
BRA0830	450 mm (Height adjustable from 450 to 490 mm by 20 mm steps)	□
BRA0840	470 mm (Height adjustable from 450 to 490 mm by 20 mm steps)	□
BRA0850	490 mm (Height adjustable from 450 to 490 mm by 20 mm steps)	□

<b>MATR X BACKREST</b>		
BRA0852	Height 310 mm	□
BRA0854	Height 410 mm	□
BRA0856	Height 460 mm	□

<b>PUSH HANDLES</b>		
BRA0860	Fixed push handles	◆
BRA0870	Foldable push handles	□
BRA0880	Height adjustable push handles, external, short grip	□
BRA0890	Height adjustable push handles, external, long grip for drum brakes	□
BRA0900	Height adjustable push handles, external, including neckrest, long grip for drum brakes	□
BRA0910	Height adjustable push handles, external, including headrest, long grip for drum brakes	□
BRA0920	Height adjustable push handles, external, including headrest with cheek support, long grip for drum brakes	□
BRA0930	Without push handles	□

**SEAT**

<b>SEAT TYPE</b>		
BRA0940	Velcro tension- and depth adjustable, infinitely variable on 40 mm.	◆

<b>SEAT CUSHIONS</b>		
BRA0960	Standard seat cushion	□
BRA0990	Without seat cushion	◆

**LEGRESTS / FOOTRESTS**

<b>LEGRESTS / HANGERS LH TYPE</b>		
BRA1000	Hanger angle adjustable - Alu2003 including calf support, LH (Only in black)	1080;1100
BRA1010	Hanger fixed 80°, swing in swing out - alu 2003, LH	◆
BRA1020	Hanger fixed 90°, swing in swing out - alu 2003, LH	1510;1520
BRA1030	Without legrest, LH	1100;1340;1350

<b>LEGRESTS / HANGERS RH TYPE</b>		
BRA1040	Hanger angle adjustable - Alu2003 including calf support, RH (Only in black)	1080;1100
BRA1050	Hanger fixed 80°, swing in swing out - alu 2003, RH	◆
BRA1060	Hanger fixed 90°, swing in swing out - alu 2003, RH	1510;1520
BRA1070	Without legrest, RH	1100;1340;1350

<b>LEGRESTS / HANGER COLOUR</b>		
BRA1080	As the chassi colour	◆
BRA1090	Black	□

<b>FOOTREST ONE PIECE incl Calfstrap with Hanger 80 deg &amp; 90 deg</b>		
BRA1100	Footrest one piece, flip-up, depth- and angle adjustable	1110;1120;1130;1140;1340;1360;1370

<b>FOOTRESTS LH TYPE</b>		
BRA1110	Footrest lower part, fixed, LH	◆
BRA1120	Footrest lower part, depth- and angle adjustable, LH	□

<b>FOOTRESTS RH TYPE</b>		
BRA1130	Footrest lower part, fixed, RH	◆
BRA1140	Footrest lower part, depth- and angle adjustable, RH	□

<b>Backrest height</b>	
Fixed	Height adjustable
200 mm	300-320 mm (20 mm step)
250 mm	350-390 mm (20 mm steps)
300 mm	400-440 mm (20 mm steps)
350 mm	450-490 mm (20 mm steps)



ARMRESTS/ SIDE GUARDS		
<b>ARMRESTS/ SIDE GUARD LH TYPE</b>		
BRA1150	Side guard, removable, LH	1900;1910;1920
BRA1160	Mudguard, removable, LH	1320;1335;1740;1900;1910;1920
BRA1170	Armrest "Rio" height adjustable, short pad (280 mm), LH	1900;1910;1920
BRA1180	Armrest "Rio" height adjustable, long pad (380 mm), LH	
BRA1190	Armrest Flip-up, height adjustable, short pad (280 mm), LH	1312;1314;1316;1320;1335;1860;1900;1910;1920
BRA1200	Armrest Flip-up, height adjustable, long pad (380 mm), LH	1312;1314;1316;1320;1335;1860
BRA1210	Armrest Hemi, LH	1900;1910;1920
BRA1220	Without armrest, LH	1860;1900;1910;1920
<b>ARMRESTS/ SIDE GUARD RH TYPE</b>		
BRA1230	Side guard, removable, RH	1900;1910;1920
BRA1240	Mudguard, removable, RH	1330;1336;1740;1900;1910;1920
BRA1250	Armrest "Rio" height adjustable, short pad (280 mm), RH	1900;1910;1920
BRA1260	Armrest "Rio" height adjustable, long pad (380 mm), RH	
BRA1270	Armrest Flip-up, height adjustable, short pad (280 mm), RH	1313;1315;1317;1330;1336;1860;1900;1910;1920
BRA1280	Armrest Flip-up, height adjustable, long pad (380 mm), RH	1313;1315;1317;1330;1336;1860
BRA1290	Armrest Hemi, RH	1900;1910;1920
BRA1300	Without armrest, RH	1860;1900;1910;1920
<b>SEATING OPTIONS</b>		
<b>BACKREST OPTIONS</b>		
BRA1310	Wedges to tension adjustable backrest	
BRA1311	Headrest MatrX, 250 mm soft pad	
BRA1312	Trunk support MatrX, fixed, offset, std, LH	
BRA1313	Trunk support MatrX, fixed, offset, std, RH	
BRA1314	Trunk support MatrX, swing away, large, LH	
BRA1315	Trunk support MatrX, swing away, large, RH	
BRA1316	Trunk support MatrX, swing away, offset, large, LH	
BRA1317	Trunk support MatrX, swing away, offset, large, RH	
BRA1320	Trunk support, LH	
BRA1330	Trunk support, RH	
BRA1335	Trunk support multi adjustable swingaway (LH)	
BRA1336	Trunk support multi adjustable swingaway (RH)	
<b>LEGREST / FOOTREST OPTIONS</b>		
BRA1340	Heel loops	
BRA1350	Calf strap	
BRA1360	Legrest amputee, LH	
BRA1370	Legrest amputee, RH	
<b>FRAME</b>		
<b>FRAME COLOURS</b>		
BRA1380	Anthracite Black	
BRA1390	Grey	
BRA1400	Blue	
BRA1410	Red	
BRA1420	Orange	
<b>CASTOR ATTACHMENT TYPES</b>		
BRA1430	Housing with angle adjustment	
<b>FRONT CASTORS</b>		
BRA1440	75x24 hard (3")	
BRA1450	100x27 medium hard (4")	
BRA1460	125x27 medium hard (5")	
BRA1470	140x40 soft (5,5")	
BRA1485	150x30 medium hard(6")	
BRA1490	150x30 soft (6")	
BRA1500	150x32 pneumatic (6")	
BRA1510	180x40 soft (7")	
BRA1520	180x40 pneumatic (7")	
<b>FORKS</b>		
BRA1530	Fixed steps	
<b>REAR WHEELS SIZE/TYPE</b>		
BRA1540	22" Standard	1660;1680;1700;1710;1990;2000
BRA1550	24" Standard	
BRA1560	24" Sporty, low profile (24x1"), aluminium handrims (Complete wheel)	1610;1620;1630;1640;1660;1680;1690;1700;1710;1770;1990;2000
BRA1570	24" Sporty, low profile (24x1"), black marathon evolution, aluminium handrims	1610;1620;1630;1640;1680;1690;1700;1710;1770;1990;2000
BRA1571	24" low profile (24x1"), stainless handrims (Complete wheel.)	1610;1620;1630;1640;1660;1670;1690;1700;1710;1770;1990;2000
BRA1580	24" Spider wheel, low profile (24x1"), aluminium handrims (Complete wheel.)	1610;1620;1630;1640;1660;1680;1690;1700;1710;1770;1990;2000
BRA1585	24" Spinerger, low profile (24x1"), black marathon evolution, titanium handrims (Complete wheel)	1610;1620;1630;1640;1670;1680;1690;1700;1710;1770;1990;2000
BRA1590	25" Low profile (25x1"), aluminium handrims (Complete wheel)	1620;1630;1640;1660;1680;1690;1700;1710;1770;1990;2000
BRA1600	26" Low profile (26x1"), aluminium handrims (Complete wheel)	1620;1630;1640;1660;1680;1690;1700;1710;1770;1990;2000
<b>TYRE TYPE, only for 22" and 24" standard rear wheels.</b>		
BRA1610	Tyres pneumatic non marking	
BRA1620	Tyres pneumatic profiled	
BRA1630	Tyres solid	
BRA1640	Tyres punctureproof	
BRA1660	Tyres black marathon evolution (Only for 24" standard rear wheel and 24" Sporty)	
<b>HANDRIM TYPE, only for 22" and 24" standard rear wheels.</b>		
BRA1670	Handrim aluminium	
BRA1680	Handrim stainless	
BRA1690	Handrim plastic coated	
BRA1700	Handrim foam coated	
BRA1710	Handrim black grip	
<b>REAR WHEEL FIXATION</b>		
BRA1730	Standard rear wheel fixation	
BRA1740	Rear wheel extension (behind backrest frame)	1760;1800;1870;1880;1990;2000
<b>CAMBER ANGLE</b>		
BRA1750	0° camber	
BRA1760	3° camber	1990;2000
<b>WHEEL OPTIONS</b>		
BRA1770	Spokeguard Frost (Only for 22" and 24" standard rear wheels)	
<b>SAFETY</b>		
<b>BRAKE TYPE</b>		
BRA1780	Standard brake	
BRA1800	Küschall active brake	1840;1850;1990;2000
BRA1810	One arm brake including brake lever extension, LH	1840;1850
BRA1820	One arm brake including brake lever extension, RH	1840;1850

SHf / Castor Size	
3"	: 390-480 mm
4"	: 410-490 mm
5"	: 380, 430-510 mm
5,5"	: 440-510 mm
6"	: 450-520 mm
7"	: 480-520 mm

SHr / Rear wheel size	
22"	: 380-440 mm
24"	: 400-460 mm
25"	: 400-460 mm
26"	: 410-470 mm

SAFETY OPTIONS			
BRA1830	Pelvic belt		<input type="checkbox"/>
BRA1840	Brake lever extension, fold away, LH		<input type="checkbox"/>
BRA1850	Brake lever extension, fold away, RH		<input type="checkbox"/>
BRA1860	Auto lock for armrest		<input type="checkbox"/>
BRA1870	Anti-tip devices, swing away, LH <i>(For safety reasons we recommend the use of two anti-tip devices.)</i>	1940	<input checked="" type="checkbox"/>
BRA1871	Anti-tip devices, extended, LH <i>(For safety reasons we recommend the use of two anti-tip devices.)</i>	1960	<input type="checkbox"/>
BRA1880	Anti-tip devices, swing away, RH <i>(For safety reasons we recommend the use of two anti-tip devices.)</i>	1950	<input checked="" type="checkbox"/>
BRA1881	Anti-tip devices, extended, RH <i>(For safety reasons we recommend the use of two anti-tip devices.)</i>	1970	<input type="checkbox"/>
BRA1890	Reflectors kit (wheels and chassis)		<input type="checkbox"/>
ACCESSORIES			
BRA1900	Table tray including mounting kit <i>(Only with long armrest pad)</i>	1920	<input type="checkbox"/>
BRA1901	Table tray cushion		<input type="checkbox"/>
BRA1910	Locking device for table tray	1920	<input type="checkbox"/>
BRA1920	Mounting kit for table (if ordered separately)		<input type="checkbox"/>
BRA1930	Pump		<input type="checkbox"/>
BRA1940	Cane holder, LH		<input type="checkbox"/>
BRA1950	Cane holder, RH		<input type="checkbox"/>
BRA1960	Step tube / Tipper aid, LH		<input type="checkbox"/>
BRA1970	Step tube / Tipper aid, RH		<input type="checkbox"/>
BRA1980	Tool kit		<input type="checkbox"/>
BRA1990	Mounting bracket E-Motion		<input type="checkbox"/>
BRA2000	Mounting bracket E-Fix		<input type="checkbox"/>
SPECIAL			
FREIGHT & PRICE			
PACKAGING			
FREIGHT PRICE			
TOTAL PRICE			