



INVACARE® TOP END®
PRO™-2
ALL SPORT WHEELCHAIR



The **Invacare® Top End® Pro™-2 All Sport Wheelchair** is all about sports. If you want to get into the game, this is the chair for you. Just name it: basketball, tennis, softball, table tennis, badmitton, floor hockey or any other sport that is played on a court and the Top End Pro-2 Wheelchair can get you started. The Pro-2 All Sport Wheelchair is offered in three seat widths with lots of adjustments which is perfect for clubs and aspiring athletes. Get in the game.

Features

Aluminum frame: choice of 3 seat widths

Adjustable seat depth 13"-17"

Screw on seat upholstery

Adjustable front and rear seat heights

Adjustable clamp-on tubular recessed footrest

Adjustable center of gravity

Removable single swivel anti-tip bar

Camber: 15°

25" High performance wheels with high pressure clincher tires

Quick release axles

Non-folding adjustable back: angle and height

Adjustable tension back upholstery

Aluminum side guards

Cushion

Positioning strap

Rear frame protectors

250 lb. weight capacity

Model #: P2AS



1. Welded bumper/recessed footrest

2. Adjustable seat depth: 13" - 17"

3. Removable single swivel anti-tip bar

4. Adjustable seat height: front and rear

5. 5K Freewheel option



Invacare Corporation
www.invacare.com

USA

One Invacare Way
Elyria, Ohio
44035-4190
(800) 333-6900

Canada

570 Matheson Blvd. E.,
Unit 8
Mississauga, Ontario
L4Z 4G4 Canada
(800) 668-5324

Top End

www.topendwheelchair.com
Sports and Recreation
Products
4501 63rd Circle North
Pinellas Park, Florida
33781-5914
(727) 522-8677
(800) 532-8677
Fax: (727) 522-1007

©2014

Invacare Corporation.

All rights reserved.

Trademarks are identified

by the symbols™ and®.

All trademarks are owned

by or licensed to Invacare

Corporation unless

otherwise noted.

Specifications are subject

to change without

notification.

Form No. 13-378

Rev. 12/14

141158

Options

Basketball carrying bag

Click straps

Metal O-Ring positioning belt

Extra seat cushion

The 5K Freewheel option is great for short road races too, so this chair truly can be used for ALL sports

Note: Some assembly is required.

