

2018-01-01

POST 2

Invacare® XLT
75/90Kunde/HMS:
Kontaktperson:
Adresse:

Telefon:
Mailadresse:
Bestillingnr./SF nr: please make your selection with a cross in the box Standard feature Free of Charge (NC9999) non conformance with option 9999

Gjeldende fra 1. juli 2018

→ Max Customer weight = 135 Kg

Code	Description	Non conformance	
BPA0010	Invacare XLT		◆
BPA8000	Special (Full description and drawing by fax or E-mail)		□
PROPULSION TYPE			
BPA0020	Standard wheel hub		◆
BPA0030	Assistant manoeuvred drumbrakes	0400;0665;0666;0667;0669;0680;0830;0840;0850; 0900;1390;1400;1410;1415;1420;1430;1490;1570; 1590; 1820;1830	□
BPA0040	One arm drive, LH	0300;0310;0400;1370;1390;1400;1410;1415; 1420;1430;1440;1470;1490;1500;1510;1520;1530; 1540;1560;1570;1590;1600;1650;1670;1701;1711; 1770; 1780;1820;1830	□
BPA0050	One arm drive, RH	0300;0310;0400;1370;1390;1400;1410;1415; 1420;1430;1440;1470;1490;1500;1510;1520;1530; 1540;1560;1570;1590;1600;1640;1680;1701;1711; 1770; 1780;1820;1830	□
FRAME TYPE			
BPA0060	Active frame, knee angle 75°		◆
BPA0070	Dynamic frame, knee angle 90°	0280;0290;0620;0630;0640;1300;1315;1320; 1330;1340;1350	□
SEAT DIMENSIONS			
SEAT WIDTH			
BPA0080	355 mm	0990;1640;1650;1820;1830	□
BPA0090	380 mm		□
BPA0100	405 mm		□
BPA0110	430 mm		□
BPA0120	455 mm		□
BPA0130	480 mm		□
SEAT DEPTH			
BPA0140	400 mm (adjustable from seat depth 360 to 400 mm)		□
BPA0150	450 mm (adjustable from seat depth 410 to 450 mm)		□
BPA0160	500 mm (adjustable from seat depth 460 to 500 mm)		□
SEAT HEIGHT			
FRONT SEAT TO FLOOR HEIGHT = SHf			
BPA0170	400 mm	0330;0340;0350;0360;0370;0380;0390;0400;0530; 0540;0550;0560;0570;0580;0590;0600;0610;0620; 0630;0640	□
BPA0180	410 mm	0340;0350;0360;0370;0380;0390;0400;0540;0550; 0560;0570;0580;0590;0600;0610;0620;0630;0640;	□
BPA0190	420 mm	0350;0360;0370;0380;0390;0400;0550;0560;0570; 0580;0590;0600;0610;0620;0630;0640	□
BPA0200	430 mm	0360;0370;0380;0390;0400;0560;0570;0580;0590; 0600;0610;0620;0630;0640	□
BPA0210	440 mm	0370;0380;0390;0400;0570;0580;0590;0600;0610; 0620;0630;0640	□
BPA0220	450 mm	0380;0390;0400;0580;0590;0600;0610;0620;0630; 0640	□
BPA0230	460 mm	0390;0400;0590;0600;0610;0620;0630;0640	□
BPA0240	470 mm	0300;0400;0600;0610;0620;0630;0640	□
BPA0250	480 mm	0300;0310;0610;0620;0630;0640	◆
BPA0260	490 mm	0300;0310;0320;0620;0630;0640	□
BPA0270	500 mm	0300;0310;0320;0330;0340;0630;0640	□
BPA0280	510 mm	0300;0310;0320;0330;0340;0640	□
BPA0290	520 mm	0300;0310;0320;0330;0340;0350	□
REAR SEAT TO FLOOR HEIGHT = SHr, SHr ≤ SHf			
BPA0300	380 mm	1380;1390;1400;1410;1415;1420;1430;1700;1710	□
BPA0310	390 mm	1380;1390;1400;1410;1415;1420;1430;1700;1710	□
BPA0320	400 mm	1430	□
BPA0330	410 mm		□
BPA0340	420 mm		□
BPA0350	430 mm		□
BPA0360	440 mm		□
BPA0370	450 mm	1370	◆
BPA0380	460 mm	1370	□
BPA0390	470 mm	1370	□
BPA0400	480 mm	1370;1380;1390;1400;1410;1415;1420	□
KNEE TO HEEL LENGTH, without seat cushion, UL ≤ SHf-40 mm			
BPA0430	270 mm (Only with Footrest high mounted)	0970;0980;0990	□
BPA0440	280 mm (Only with Footrest high mounted)	0970;0980;0990	□
BPA0450	290 mm (Only with Footrest high mounted)	0970;0980;0990	□
BPA0460	300 mm (Only with Footrest high mounted)	0970;0980;0990	□
BPA0470	310 mm (Only with Footrest high mounted)	0970;0980;0990	□
BPA0480	320 mm (Only with Footrest high mounted)	0970;0980;0990	□
BPA0490	330 mm (Only with Footrest high mounted)	0970;0980;0990	□
BPA0500	340 mm (Only with Footrest high mounted)	0970;0980;0990	□
BPA0510	350 mm (Only with Footrest high mounted)	0970;0980;0990	□
BPA0520	360 mm (Active frame only with Footrest high mounted)		□
BPA0530	370 mm (Active frame only with Footrest high mounted)		□
BPA0540	380 mm	1000	□
BPA0550	390 mm	1000	□
BPA0560	400 mm	1000	□
BPA0570	410 mm	1000	□
BPA0580	420 mm	1000	◆
BPA0590	430 mm	1000	□
BPA0600	440 mm	1000	□
BPA0610	450 mm	1000	□
BPA0620	460 mm (Only for Active frame)	1000	□
BPA0630	470 mm (Only for Active frame)	1000	□
BPA0640	480 mm (Only for Active frame)	1000	□

**SHf - Active frame/
Castor Size**3" : 400-480 mm
4" : 410-490 mm
5" : 430-510 mm
5,5" : 440-510 mm
6" : 450-520 mm
7" : 480-520 mm**SHf - Dynamic frame/
Castor Size**3" : 400-480 mm
4" : 410-490 mm
5" : 430-500 mm**SHr / Rear wheel size**22" : 380-450 mm
24" : 400-470 mm
25" : 400-470 mm
26" : 410-480 mm

BACKREST			
BACKREST TYPE			
BPA0650	Velcro tension adjustable, fixed height, infinitely variable angle adjustable, foldable and click in backrest (bent backrest tubes)	0720;0730;0740;0750;0760;0770;0780;0790;0800;0810;0820;0822;0824;1171;1172;1173;1174;1175;1176;1177;1180;1190;1195;1196	□
BPA0655	Velcro tension adjustable, fixed height, infinitely variable angle adjustable, foldable and click in backrest (straight backrest tubes)	0720;0730;0740;0750;0760;0770;0780;0790;0800;0810;0820;0822;0824;1171;1172;1173;1174;1175;1176;1177;1180;1190;1195;1196	□
BPA0660	Velcro tension adjustable, height adjustable, infinitely variable angle adjustable, foldable and click in backrest (bent backrest tubes)	0680;0690;0700;0710;0822;0824;1171;1172;1173;1174;1175;1176;1177	◆
BPA0663	Velcro tension adjustable, height adjustable, infinitely variable angle adjustable, foldable and click in backrest (straight backrest tubes)	0680;0690;0700;0710;0822;0824;1171;1172;1173;1174;1175;1176;1177	□
BPA0665	Narrow back, velcro tension adjustable, fixed height, infinitely variable angle adjustable, foldable and click in backrest	0680;0690;0700;0720;0730;0740;0750;0760;0770;0780;0790;0800;0810;0820;0822;0824;0850;0860;0870;0880;0890;1171;1172;1173;1174;1175;1176;1177;1180;1190;1195;1196	□
BPA0666	MatrX PB Elite std, incl hardware	0680;0690;0700;0710;0720;0730;0740;0750;0760;0770;0780;0790;0800;0810;0820;0840;0870;0880;0890;1170;1180;1190;1195;1196	□
BPA0667	MatrX PB Elite deep, incl hardware	0680;0690;0700;0710;0720;0730;0740;0750;0760;0770;0780;0790;0800;0810;0820;0840;0870;0880;0890;1170;1172;1173;1174;1175;1176;1177;1180;1190;1195;1196	□
BPA0669	Without Backrest, prepared for MatrX	0680;0690;0700;0710;0720;0730;0740;0750;0760;0770;0780;0790;0800;0810;0820;0822;0824;0840;0850;0860;0870;0880;0890;1170;1171;1172;1173;1174;1175;1176;1177;1180;1190;1195;1196	□
BACKREST HEIGHT (PLEASE NOTE THAT THE XLT CANNOT BE USED AS A SEAT DURING TRANSPORTATION IF THE BACK HEIGHT IS LOWER THAN 400 MM)			
FIXED BACKREST			
BPA0680	200 mm	0850;0860;0870;0880;0890	□
BPA0690	250 mm		□
BPA0700	300 mm		□
BPA0710	350 mm		□
HEIGHT ADJUSTABLE BACKREST			
BPA0720	300 mm (Height adjustable from 300 to 320 mm by 20 mm steps)		□
BPA0730	320 mm (Height adjustable from 300 to 320 mm by 20 mm steps)		□
BPA0740	350 mm (Height adjustable from 350 to 390 mm by 20 mm steps)		□
BPA0750	370 mm (Height adjustable from 350 to 390 mm by 20 mm steps)		◆
BPA0760	390 mm (Height adjustable from 350 to 390 mm by 20 mm steps)		□
BPA0770	400 mm (Height adjustable from 400 to 440 mm by 20 mm steps)		□
BPA0780	420 mm (Height adjustable from 400 to 440 mm by 20 mm steps)		□
BPA0790	440 mm (Height adjustable from 400 to 440 mm by 20 mm steps)		□
BPA0800	450 mm (Height adjustable from 450 to 490 mm by 20 mm steps)		□
BPA0810	470 mm (Height adjustable from 450 to 490 mm by 20 mm steps)		□
BPA0820	490 mm (Height adjustable from 450 to 490 mm by 20 mm steps)		□
MATR X BACKREST			
BPA0822	Height 310 mm		□
BPA0824	Height 410 mm		□
PUSH HANDLES			
BPA0830	Fixed push handles		◆
BPA0840	Foldable push handles		□
BPA0850	Height adjustable push handles, external, short grip		□
BPA0860	Height adjustable push handles, external, long grip for drum brakes		□
BPA0870	Height adjustable push handles, external, including neckrest, long grip for drum brakes		□
BPA0880	Height adjustable push handles, external, including headrest, long grip for drum brakes		□
BPA0890	Height adjustable push handles, external, including headrest with cheek support, long grip for drum brakes		□
BPA0900	Without push handles		□
SEAT			
SEAT TYPE			
BPA0910	Velcro tension- and depth adjustable, infinitely variable on 40 mm.		◆
SEAT CUSHIONS			
BPA0930	Standard seat cushion		□
BPA0960	Without seat cushion		◆
LEGRESTS / FOOTRESTS			
FOOTREST INCL. CALF STRAP			
BPA0970	Foot bow rigid		◆
BPA0980	Foot bow rigid with cover		□
BPA0990	Footrest one piece, flip-up, angle adjustable		□
BPA1000	Footrest, high mounted, for knee to heel length 270 to 370 mm		□
ARMRESTS/ SIDE GUARDS			
ARMRESTS/ SIDE GUARD LH TYPE			
BPA1010	Side guard, removable, LH	1730;1740;1750	◆
BPA1020	Mudguard, removable, LH	1180;1195;1570;1730;1740;1750	□
BPA1030	Armrest "Rio" height adjustable, short pad (280 mm), LH	1730;1740;1750	□
BPA1040	Armrest "Rio" height adjustable, long pad (380 mm), LH		□
BPA1050	Armrest Flip-up, height adjustable, short pad (280 mm), LH	1172;1174;1176;1180;1195;1690;1730;1740;1750	□
BPA1060	Armrest Flip-up, height adjustable, long pad (380 mm), LH	1172;1174;1176;1180;1195;1690	□
BPA1070	Armrest Hemi, LH	1730;1740;1750	□
BPA1080	Without armrest, LH	1690;1730;1740;1750	□
ARMRESTS/ SIDE GUARD RH TYPE			
BPA1090	Side guard, removable, RH	1730;1740;1750	◆
BPA1100	Mudguard, removable, RH	1190;1196;1570;1730;1740;1750	□
BPA1110	Armrest "Rio" height adjustable, short pad (280 mm), RH	1730;1740;1750	□
BPA1120	Armrest "Rio" height adjustable, long pad (380 mm), RH		□
BPA1130	Armrest Flip-up, height adjustable, short pad (280 mm), RH	1173;1175;1177;1190;1196;1690;1730;1740;1750	□
BPA1140	Armrest Flip-up, height adjustable, long pad (380 mm), RH	1173;1175;1177;1190;1196;1690	□
BPA1150	Armrest Hemi, RH	1730;1740;1750	□
BPA1160	Without armrest, RH	1690;1730;1740;1750	□
SEATING OPTIONS			
BACKREST OPTIONS			
BPA1170	Wedges to tension adjustable backrest		□
BPA1171	Headrest MatrX, 250 mm soft pad		□
BPA1172	Trunk support MatrX, fixed, offset, std, LH		□
BPA1173	Trunk support MatrX, fixed, offset, std, RH		□
BPA1174	Trunk support MatrX, swing away, large, LH		□
BPA1175	Trunk support MatrX, swing away, large, RH		□
BPA1176	Trunk support MatrX, swing away, offset, large, LH		□
BPA1177	Trunk support MatrX, swing away, offset, large, RH		□
BPA1180	Trunk support, LH		□
BPA1190	Trunk support, RH		□
BPA1195	Trunk support multi adjustable swingaway (LH)		□
BPA1196	Trunk support multi adjustable swingaway (RH)		□
FRAME			
FRAME COLOURS			
BPA1200	Anthracite Black		◆
BPA1210	Grey		□
BPA1220	Blue		□
BPA1230	Red		□
BPA1240	Orange		□
BPA1250	Blasted titanium		□

Backrest height	
Fixed	Height adjustable
200 mm	300-320 mm (20 mm step)
250 mm	350-390 mm (20 mm steps)
300 mm	400-440 mm (20 mm steps)
350 mm	450-490 mm (20 mm steps)

CASTOR ATTACHMENT TYPES			
BPA1260	Housing with angle adjustment		◆
FRONT CASTORS, 3" to 5" possible on Dynamic frame. 3" to 7" possible on Active frame.			
BPA1270	75x24 hard (3")		□
BPA1280	100x27 medium hard (4")		□
BPA1290	125x27 medium hard (5")		◆
BPA1300	140x40 soft (5,5")		□
BPA1315	150x30 medium hard (6")		□
BPA1320	150x30 soft (6")		□
BPA1330	150x32 pneumatic (6")		□
BPA1340	180x40 soft (7")		□
BPA1350	180x40 pneumatic (7")		□
FORKS			
BPA1360	Fixed steps		◆
REAR WHEELS SIZE/TYPE			
BPA1370	22" Standard	1490;1510;1530;1540;1820;1830	□
BPA1380	24" Standard		◆
BPA1390	24" Sporty, low profile (24x1"), aluminium handrims (<i>Complete wheel</i>)	1440;1450;1460;1470;1490;1510;1520;1530;1540;1600;1820;1830	□
BPA1400	24" Sporty, low profile (24x1"), black marathon evolution, aluminium handrims	1440;1450;1460;1470;1510;1520;1530;1540;1600;1820;1830	□
BPA1410	24" Spider wheel, low profile (24x1"), aluminium handrims (<i>Complete wheel.</i>)	1440;1450;1460;1470;1490;1500;1510;1520;1530;1540;1550;1600;1820;1830	□
BPA1415	24" Spinerger, low profile (24x1"), black marathon evolution, titanium handrims (<i>Complete wheel</i>)	1440;1450;1460;1470;1500;1510;1520;1530;1540;1600;1820;1830	□
BPA1420	25" Low profile (25x1"), aluminium handrims (<i>Complete wheel.</i>)	1450;1460;1470;1490;1510;1520;1530;1540;1600;1820;1830	□
BPA1430	26" Low profile (26x1"), aluminium handrims (<i>Complete wheel.</i>)	1450;1460;1470;1490;1510;1520;1530;1540;1600;1820;1830	□
TYRE TYPE, only for 22" and 24" standard rear wheels.			
BPA1440	Tyres pneumatic non marking		◆
BPA1450	Tyres pneumatic profiled		□
BPA1460	Tyres solid		□
BPA1470	Tyres punctureproof		□
BPA1490	Tyres black marathon evolution (<i>Only for 24" standard rear wheel and 24" Sporty</i>)		□
HANDRIM TYPE, only for 22" and 24" standard rear wheels.			
BPA1500	Handrim aluminium		◆
BPA1510	Handrim stainless		□
BPA1520	Handrim plastic coated		□
BPA1530	Handrim foam coated		□
BPA1540	Handrim black grip		□
REAR WHEEL FIXATION			
BPA1560	Standard rear wheel fixation		◆
BPA1570	Rear wheel extension (behind backrest frame)	1590;1620;1630;1700;1710;1820;1830	□
CAMBER ANGLE			
BPA1580	0° camber		□
BPA1590	3° camber	1820;1830	◆
WHEEL OPTIONS			
BPA1600	Spokeguard Frost (<i>Only for 22" and 24" standard rear wheels</i>)		□
SAFETY			
BRAKE TYPE			
BPA1610	Standard brake		□
BPA1620	Küschall performance brake	1670;1680;1820;1830	□
BPA1630	Küschall active brake	1670;1680;1820;1830	◆
BPA1640	One arm brake including brake lever extension, LH	1670;1680	□
BPA1650	One arm brake including brake lever extension, RH	1670;1680	□
SAFETY OPTIONS			
BPA1660	Pelvic belt		□
BPA1670	Brake lever extension, fold away, LH		□
BPA1680	Brake lever extension, fold away, RH		□
BPA1690	Auto lock for armrest		□
BPA1700	Anti-tip devices, swing away, LH (<i>For safety reasons we recommend the use of two anti-tip devices.</i>)	1711;1770	◆
BPA1701	Anti-tip devices, extended, LH (<i>For safety reasons we recommend the use of two anti-tip devices.</i>)	1710;1790	□
BPA1710	Anti-tip devices, swing away, RH (<i>For safety reasons we recommend the use of two anti-tip devices.</i>)	1780	◆
BPA1711	Anti-tip devices, extended, RH (<i>For safety reasons we recommend the use of two anti-tip devices.</i>)	1800	□
BPA1720	Reflectors kit (wheels and chassis)		□
ACCESSORIES			
BPA1730	Table tray including mounting kit (<i>Only with long armrest pad</i>)	1750	□
BPA1731	Table tray cushion		□
BPA1740	Locking device for table tray	1750	□
BPA1750	Mounting kit for table (if ordered separately)		□
BPA1760	Pump		□
BPA1770	Cane holder, LH		□
BPA1780	Cane holder, RH		□
BPA1790	Step tube / Tipper aid, LH		□
BPA1800	Step tube / Tipper aid, RH		□
BPA1810	Tool kit		□
BPA1815	Frame cover		□
BPA1820	Mounting bracket E-Motion		□
BPA1830	Mounting bracket E-Fix		□
SPECIAL			
FREIGHT & PRICE			
PACKAGING			
FREIGHT PRICE			
TOTAL PRICE			

SHf - Active frame/ Castor Size	
3"	: 400-480 mm
4"	: 410-490 mm
5"	: 430-510 mm
5,5"	: 440-510 mm
6"	: 450-520 mm
7"	: 480-520 mm

SHf - Dynamic frame/ Castor Size	
3"	: 400-480 mm
4"	: 410-490 mm
5"	: 430-500 mm

SHr / Rear wheel size	
22"	: 380-450 cm
24"	: 400-470 cm
25"	: 400-470 cm
26"	: 410-480 cm