

## **Certificate in relation to air transportation by passengers with walking disabilities**

We hereby certify that our product

**via GO V24**

including the lithium-ion batteries meets the requirements of

**UN Manual of Tests and Criteria" Part III subsection 38.3**

and therefore are suitable for transportation by passengers according

**IATA DGR Chapter 2.3.2.4**

The rated energy of the built-in drive unit battery (which is not intended for removal by the user) is **162 Wh**.

Since drive unit are portable and contain the battery, we recommend that you take them with you to the passenger cabin. The drive unit must be switched off and the control unit must be inserted into the holder provided on the drive unit to prevent it from unintentional activation.

The pilot in command must be informed about the location of the mobility aid.

**We strongly recommend contacting the airline (called operator) sufficiently in advance to request approval for transport.**

This certificate does not apply to damaged or defective devices or batteries.

Due to the frequent amendments of the respective regulations, which partly occur on an annual basis this certificate is valid until **31<sup>st</sup> December 2024**.

Albstadt, 13<sup>th</sup> of November 2023

Alber GmbH



Johannes Beltman  
Managing Director