

Cura 170

Please fill in completely!

The lift is individually made to fit your bathtub!
Please enter the dimensions of your bathtub to the picture that fits your bathtub most.

Request

Order

Purchaser _____

Client no. _____

Stamp _____

Date, signature _____

Reference# _____

Depth of bathtub* 385 cm

Length of supply hose 2 m

Drain _____ Headpart _____

right _____ right _____

left _____ left _____

centre (of the bathtub) _____

Free-standing bathtub _____

yes _____ no _____

Colour covermat _____

white _____ blue _____

User's weight _____

20 - 70 kg _____


50 - 130 kg _____

100 - 150 kg _____

* Please place a wood plate or the like across the bathtub and measure the depth in the centre of the bathtub.


The purchaser is responsible for complete and correct measurements!

1 Cura 170



A1	Bathtub bottom to drain	108 cm
A2	Drain to curve	5 cm
B	Bathtub width at the bottom	44 cm
C	Bathtub length at the top (inside)	153 cm
E	Bathtub width outside	69,5 cm
F	Bathtub width at the top (inside)	61 cm
X	Bathtub length outside	170 cm

2



A1	Bathtub bottom to drain	cm
A2	Drain to curve	cm
B1	Bathtub width at the bottom	cm
B2	Bathtub width at the bottom	cm
C	Bathtub length at the top (inside)	cm
E	Bathtub width outside	cm
F1	Bathtub width at the top (inside)	cm
F2	Bathtub width at the top (inside)	cm
X	Bathtub length outside	cm