



BATHLIFT TO LIE ON
CUSTOMER INFORMATION CUSTOMISATION

Customer information for placing a custom-made product on the market. Please fax the filled out form together with the dimension sheet to our fax no. +49 75 62 700-69!

The requested custom-made product is manufactured according to attached dimension sheet by the company:

AQUATEC GmbH
Alemannenstrasse 10
88316 Isny
Germany

To:
AQUATEC GmbH
(Customer service)
Alemannenstrasse 10
88316 Isny
Germany

Fax: +49 (0) 75 62 / 700-69

The custom-made product is for the patient:

Mr / Mrs _____

Street _____

Post code / City _____

By prescription of physician:

Name _____

Reachable at _____

This information is required for the production of a custom-made product according to article I, section 2d MDD.

Place _____

Stamp

Date, signature _____



Please fill in completely!

- Request
- Order

Purchaser _____

Client no. _____

Stamp _____

Date, signature _____

Reference _____

Depth of bathtub* _____ **cm**

Length of supply hose _____ **m**

Drain _____ **Headpart** _____

right _____ right _____

left _____ left _____

centre (of the bathtub) _____

Free-standing bathtub _____

yes _____ no _____

Colour covermat _____

white _____ blue _____

User's weight _____

20 - 70 kg _____

50 - 130 kg _____


100 - 150 kg _____

* Please place a wood plate or the like across the bathtub and measure the depth in the centre of the bathtub.

The purchaser is responsible for complete and correct measurements!


The lift is individually made to fit your bathtub!
Please enter the dimensions of your bathtub to the picture that fits your bathtub most.

1



A1	Bathtub bottom to drain	cm
A2	Drain to curve	cm
B	Bathtub width at the bottom	cm
C	Bathtub length at the top (inside)	cm
E	Bathtub width outside	cm
F	Bathtub width at the top (inside)	cm
X	Bathtub length outside	cm

2



A1	Bathtub bottom to drain	cm
A2	Drain to curve	cm
B1	Bathtub width at the bottom	cm
B2	Bathtub width at the bottom	cm
C	Bathtub length at the top (inside)	cm
E	Bathtub width outside	cm
F1	Bathtub width at the top (inside)	cm
F2	Bathtub width at the top (inside)	cm
X	Bathtub length outside	cm



Please fill in completely!

- Request
- Order

Purchaser _____

Client no. _____

Stamp _____

Date, signature _____

Reference _____

Depth of bathtub* _____ **cm**

Length of supply hose _____ **m**

Drain _____ **Headpart** _____

right right

left left

centre (of the bathtub)

Free-standing bathtub _____

yes no

Colour covermat _____

white blue

User's weight _____

20 - 70 kg

50 - 130 kg


100 - 150 kg

* Please place a wood plate or the like across the bathtub and measure the depth in the centre of the bathtub.

The purchaser is responsible for complete and correct measurements!


The lift is individually made to fit your bathtub!
Please enter the dimensions of your bathtub to the picture that fits your bathtub most.

3



A1 Bathtub bottom to drain	cm
A2 Drain to curve	cm
B1 Bathtub width at the bottom	cm
B2 Bathtub width at the bottom	cm
C Bathtub length at the top (inside)	cm
E Bathtub width outside	cm
F1 Bathtub width at the top (inside)	cm
F2 Bathtub width at the top (inside)	cm
X Bathtub length outside	cm

4



A1 Bathtub bottom to drain	cm
A2 Drain to curve	cm
B1 Bathtub width at the bottom	cm
B2 Bathtub width at the bottom	cm
C Bathtub length at the top (inside)	cm
E Bathtub width outside	cm
F1 Bathtub width at the top (inside)	cm
F2 Bathtub width at the top (inside)	cm
X Bathtub length outside	cm



Please fill in completely!

- Request
- Order

Purchaser

Client no. _____

Stamp _____

Date, signature _____

Reference

Depth of bathtub* _____ **cm**

Length of supply hose _____ **m**

Drain Headpart

right right

left left

centre (of the bathtub)

Free-standing bathtub

yes no

Colour covermat

white blue

User's weight


20 - 70 kg

50 - 130 kg

100 - 150 kg

The lift is individually made to fit your bathtub!
Please enter the dimensions of your bathtub to the picture that fits your bathtub most.

5



A1	Bathtub bottom to drain	cm
A2	Drain to curve	cm
B	Bathtub width at the bottom	cm
C	Bathtub length at the top (inside)	cm
E	Bathtub width outside	cm
F	Bathtub width at the top (inside)	cm
X	Bathtub length outside	cm

* Please place a wood plate or the like across the bathtub and measure the depth in the centre of the bathtub.

The purchaser is responsible for complete and correct measurements!



Please fill in completely!

- Request
 Order

Purchaser _____
Client no. _____
Stamp _____

Date, signature _____

Reference _____

Depth of bathtub* _____ **cm**
Length of supply hose _____ **m**

Drain **Headpart**
 right right
 left left
 centre (of the bathtub) _____

Free-standing bathtub
 yes no

Colour covermat
 white blue

User's weight
 20 - 70 kg
 50 - 130 kg
 100 - 150 kg


* Please place a wood plate or the like across the bathtub and measure the depth in the centre of the bathtub.

The purchaser is responsible for complete and correct measurements!

The lift is individually made to fit your bathtub!
Please enter the dimensions of your bathtub to the picture that fits your bathtub most.

6 Additional dimensioning of the handles, when they protrude into the bathtub


In case there is only one handle, viewer's line of sight:
 right left



A	Length of handle	cm
B	Outer edge bathtub to beginning of handle	cm
D	Distance between handles	cm

7 Additional dimensioning when the drain/overflow is in the centre of the bathtub

overflow at the wall side overflow at the access side



A1	Bathtub bottom to drain	cm
A2	Drain to curve	cm
E	Outer edge to inner edge (at the top)	cm
F	Outer edge to inner edge (at the bottom)	cm
G	Outer edge to inner edge (at the bottom)	cm
H	Outer edge to centre of drain/overflow	cm
I	Outer edge to outer edge of overflow	cm
K	Outer edge to centre of drain	cm



Please fill in completely!

- Request
- Order

Purchaser _____
 Client no. _____
 Stamp _____

Date, signature _____

Reference _____

Depth of bathtub* _____ **cm**
Length of supply hose _____ **m**

Drain **Headpart**
 right right
 left left
 centre (of the bathtub)

Free-standing bathtub
 yes no

Colour covermat
 white blue


User's weight
 20 - 70 kg
 50 - 130 kg
 100 - 150 kg

The lift is individually made to fit your bathtub!
 Please enter the dimensions of your bathtub to the picture that fits your bathtub most.

8 Special bathtub. Please add a photo for unusual bathtubs.

The sideflaps are not capable of bearing weight if there are openings here.

Sides are deepened:
 yes no



A	Leg length	cm
B	Leg length	cm
C	Bathtub length at the top (inside)	cm
D	Bathtub width at the top (inside)	cm
E	Bathtub bottom (curve/curve)	cm
F	Bathtub width at the bottom	cm
G	Outer edge to curve	cm
H	Outside dimension of bathtub	cm
I	Outer edge to outer edge of overflow	cm
K1	Outer edge to inner edge of bathtub (at the top)	cm
K2	Outer edge to inner edge of bathtub (at the top)	cm
L	Outer edge to inner edge of bathtub (at the top)	cm
M1	Outer edge to inner edge of bathtub (at the top)	cm
M2	Outer edge to inner edge of bathtub (at the top)	cm
P	Width of seat cut-out	cm

* Please place a wood plate or the like across the bathtub and measure the depth in the centre of the bathtub.

The purchaser is responsible for complete and correct measurements!