



Yes, you can.®

INVACARE France OPERATIONS SAS
Route de Saint Roch, 37230 Fondettes, France
Tel: 33 (0)247 626 466 Fax: 33 (0)247 626 488

Gjeldende fra juli 2018

Invacare® **Dahlia 30°**

POST 7
Allround rullestol med fastramme, tilt og
avtakbare benstøtter

Kunde/HMS:
Kontaktperson:
Adresse:

Telefon:
Mailadresse:
Bestillingnr./SF nr:

please make your selection with a cross in the box Standard feature Free of Charge (NC9999) non conformance with option 9999

→ Customer weight = 135 kg	
Code	Description
MODEL	
BVA0010	<input checked="" type="checkbox"/> DAHLIA 30°
SPECIAL	
BVA8000	<input type="checkbox"/> SPECIAL
PROPULSION TYPE	
BVA0020	<input type="checkbox"/> TRANSIT: 16" WHEELS WITH HUB BRAKES (NC : 0140, 0150, 0160, 0190, 1520, 1550, 1560, 1570, 1580, 1590, 1595, 1610, 1630, 1710, 1725, 1727)
BVA0030	<input checked="" type="checkbox"/> SELF PROPELLED: 24"/ 22" /20"WHEELS WITH USER BRAKES (NC : 0090, 0180, 1530, 1540)
BVA0040	<input type="checkbox"/> SELF- AND ASSISTANT PROPELLED: 24" / 22" WHEELS WITH HUB BRAKES (NC : 0090, 0120, 0140, 0180, 0190, 1530, 1540, 1550, 1660, 1680, 1690)
SEAT DIMENSIONS	
SEAT WIDTH	
BVA0050	<input type="checkbox"/> 390 MM (SMALL CHASSIS)
BVA0060	<input type="checkbox"/> 440 MM (MEDIUM CHASSIS)
BVA0070	<input type="checkbox"/> 490 MM (LARGE CHASSIS)
SEAT WIDTH EXTENSION (PRE ADJUSTMENT)	
BVA0080	<input type="checkbox"/> + 50 MM (NC : 0630)
BVA0090	<input type="checkbox"/> + 100 MM (NC : 0630)
BVA0100	<input checked="" type="checkbox"/> WITHOUT SEAT WIDTH EXTENSION (NC : 0460, 0470)
SEAT DEPTH	
BVA0110	<input checked="" type="checkbox"/> 430-500 MM
SEAT HEIGHT	
BVA0120	<input type="checkbox"/> 330 MM/ 13" (16" AND 20") (NC : 0637, 0670, 1440, 1450, 1460, 1470, 1480, 1490, 1500, 1510, 1560, 1570, 1580)
BVA0130	<input type="checkbox"/> 350 MM/ 14" (22", 20" AND 16") (NC : 0220, 0637, 0670, 1430, 1440, 1480, 1490, 1500, 1510, 1570, 1580)
BVA0140	<input type="checkbox"/> 380 MM/ 15" (20") (NC : 0220, 0637, 0670, 1430, 1440, 1480, 1490, 1500, 1510, 1540, 1560, 1570, 1580)
BVA0150	<input type="checkbox"/> 400 MM/ 16" (20" AND 22") (NC : 0670, 1430, 1440, 1480, 1490, 1500, 1510, 1540, 1570, 1580)
BVA0160	<input type="checkbox"/> 430 MM/ 17" (24" AND 22") (NC : 0670, 1430, 1440, 1480, 1490, 1500, 1510, 1540, 1550)
BVA0170	<input checked="" type="checkbox"/> 450 MM/ 18" (24" AND 16") (NC : 1430, 1440, 1550, 1560)
BVA0180	<input type="checkbox"/> 480 MM/ 19" (16") (NC : 1430, 1440, 1450, 1460, 1470, 1550, 1560, 1570, 1580, 1725, 1727)
BVA0190	<input type="checkbox"/> 500 MM/ 20" (20") (NC : 1420, 1430, 1440, 1450, 1460, 1470, 1540, 1560, 1570, 1580)
SEAT ANGLE LIMITATIONS	
BVA0200	<input type="checkbox"/> REDUCED TILT: POSITIVE SEAT ANGLE 10 DEGREE
BVA0210	<input type="checkbox"/> REDUCED TILT: POSITIVE SEAT ANGLE 20 DEGREE (NC : 0630, 0632)
BVA0220	<input checked="" type="checkbox"/> FULL TILT: POSITIVE SEAT ANGLE 30 DEGREE (NC : 0630, 0632)
RECLINE ANGLE LIMITATIONS	
BVA0225	<input checked="" type="checkbox"/> REDUCED BACKREST RECLINE: ANGLE 10 DEGREES
BVA0227	<input type="checkbox"/> REDUCED BACKREST RECLINE: ANGLE 20 DEGREES
BVA0230	<input type="checkbox"/> FULL BACKREST RECLINE: ANGLE 30 DEGREES
SEAT ANGLE - Negative	
BVA0240	<input type="checkbox"/> NEGATIVE SEAT ANGLE -1 (STANDARD IS + 0 DEGREE)
SEAT TYPE	
BVA0250	<input checked="" type="checkbox"/> DEPTH ADJUSTABLE 430-500 MM
SEAT CUSHION	
BVA0260	<input type="checkbox"/> CUSHION MISTRAL 2 LAYER SEAT UPHOLSTERY MANDATORY (NC : 0330, 0345)
BVA0270	<input type="checkbox"/> CUSHION FLOSHAPE SEAT UPHOLSTERY MANDATORY
BVA0310	<input type="checkbox"/> CUSHION VICAIR MULTIFUNCTIONAL (BLACK) (NC : 0330, 0340, 0345)
BVA0320	<input checked="" type="checkbox"/> WITHOUT SEAT CUSHION, (SEAT PLATE ONLY) (NC : 0330, 0340, 0345)
SEAT UPHOLSTERY	
BVA0330	<input type="checkbox"/> BLACK PLUSH TR35
BVA0340	<input type="checkbox"/> BLACK DARTEX TR26
BVA0345	<input type="checkbox"/> BLACK SPACER TR36
BACKREST	
BACKREST TUBE	
BVA0350	<input type="checkbox"/> BENT (20 DEGREES FROM 380 MM ABOVE SEATPLATE - ONLY MATRX) (NC: 0410, 0420, 0430, 0440, 0590, 1300)
BVA0360	<input checked="" type="checkbox"/> STRAIGHT (FLEX 3, MATRX AND LAGUNA)
BACKREST JOINT POSITION	
BVA0370	<input type="checkbox"/> + 50 MM (NC : 0470)
BVA0380	<input checked="" type="checkbox"/> WITHOUT BACKREST JOINT WIDTH EXTENSION
BACKREST TYPE	
BVA0390	<input type="checkbox"/> MATRX ELITE (INCL. HARDWARE MOUNTED TO FRAME, BACKREST SUPPLIED IN BOX) BACKREST HEIGHT MATRX MANDATORY (NC : 0460, 0470, 0490, 0495, 0500, 0510, 0520, 0530, 0540, 0555, 0560, 0570, 1140, 1150, 1160, 1170, 1180, 1190, 1200, 1210, 1220, 1230, 1240, 1250, 1260, 1265, 1270, 1280, 1290)
BVA0400	<input type="checkbox"/> MATRX ELITE DEEP (INCL. HARDWARE MOUNTED TO FRAME, BACKREST SUPPLIED IN BOX) BACKREST HEIGHT MATRX MANDATORY (NC : 0460, 0470, 0490, 0495, 0500, 0510, 0520, 0530, 0540, 0555, 0560, 0570, 1140, 1150, 1160, 1170, 1180, 1190, 1200, 1210, 1220, 1230, 1240, 1250, 1260, 1265, 1270, 1280, 1290)
BVA0405	<input type="checkbox"/> MATRX ELITE TR (INCL. HARDWARE MOUNTED TO FRAME, BACKREST SUPPLIED IN BOX) BACKREST HEIGHT MATRX MANDATORY (NC : 0460, 0470, 0490, 0495, 0500, 0510, 0520, 0530, 0540, 0555, 0560, 0570, 1140, 1150, 1160, 1170, 1180, 1190, 1200, 1210, 1220, 1230, 1240, 1250, 1260, 1265, 1270, 1280, 1290)
BVA0410	<input type="checkbox"/> LAGUNA 2 BACKREST PLATE (NC : 0455, 0457, 0460, 0470, 0490, 0495, 0510, 0520, 0530, 0540, 0555, 0610, 1112, 1265)
BVA0420	<input checked="" type="checkbox"/> LAGUNA BACKREST TENSION ADJUSTABLE, WITH BACKREST BOW AND NECKREST/HEADREST ATTACHMENT (NC : 0455, 0457, 0460, 0470, 0520, 0530, 0540, 0580, 0590, 0610, 1112, 1140, 1150, 1160, 1170, 1180, 1190, 1200, 1210, 1220, 1230, 1240, 1250, 1260, 1270)
BVA0430	<input type="checkbox"/> FLEX 3 BACKREST, FIXED PLATE (WIDTH ADJUSTABLE) (NC : 0455, 0457, 0490, 0495, 0555, 0610, 1112, 1265)
BVA0440	<input type="checkbox"/> FLEX 3 BACKREST, TENSION ADJUSTABLE(WIDTH ADJUSTABLE) (NC : 0455, 0457, 0470, 0495, 0610, 1112, 1290)
BVA0450	<input type="checkbox"/> BACKREST TUBES WITHOUT BACKREST (NC : 0455, 0457, 0460, 0470, 0490, 0500, 0510, 0520, 0530, 0540, 0555, 1140, 1150, 1160, 1170, 1180, 1190, 1200, 1210, 1220, 1230, 1240, 1250, 1260, 1265, 1270, 1280, 1290)
BACKREST HEIGHT MATRX	
BVA0455	<input type="checkbox"/> HEIGHT 41 CM
BVA0457	<input type="checkbox"/> HEIGHT 51 CM
BACKREST FLEX 3 ADJUSTMENT	
BVA0460	<input type="checkbox"/> FLEX 3 + 50 MM
BVA0470	<input type="checkbox"/> FLEX 3 + 100 MM (NC : 0630, 1725, 1727)
BVA0480	<input checked="" type="checkbox"/> WITHOUT BACKREST WIDTH EXTENSION
BACKREST COVER/ CUSHION	
BVA0490	<input type="checkbox"/> COVER, LATERAL BACKREST UPHOLSTERY MANDATORY
BVA0495	<input checked="" type="checkbox"/> COVER, LATERAL LOW (43CM) BACKREST UPHOLSTERY MANDATORY (NC : 0560)
BVA0500	<input type="checkbox"/> CUSHION, LAGUNA (LATERAL SUPPORT) BACKREST UPHOLSTERY MANDATORY (NC : 0555, 1265)
BVA0510	<input type="checkbox"/> CUSHION, MISTRAL 2 (WAIST SUPPORT) BACKREST UPHOLSTERY MANDATORY (NC : 0555, 0575, 1265)
BVA0520	<input type="checkbox"/> CUSHION, PASSAD 2 (SHOULDER SUPPORT) BACKREST UPHOLSTERY MANDATORY (NC : 0555, 1265)
BVA0530	<input type="checkbox"/> CUSHION, SHOULDER HIGH 05 (DARTEX) BACKREST UPHOLSTERY MANDATORY (NC : 0555, 0560, 0575, 1265)
BVA0540	<input type="checkbox"/> CUSHION, VICAIR MULTIFUNCTIONAL (BLACK) (NC : 0555, 0560, 0570, 0575, 1265)
BVA0550	<input type="checkbox"/> WITHOUT BACKREST COVER/ CUSHION (NC : 0555, 0560, 0570, 0575, 1265)
EXTRA LARGE BACKREST COVER / BACKREST CUSHION	
BVA0555	<input type="checkbox"/> XL BACKREST COVER/BACKREST CUSHION
BACKREST UPHOLSTERY	
BVA0560	<input type="checkbox"/> BLACK PLUSH TR35
BVA0570	<input checked="" type="checkbox"/> BLACK DARTEX TR26
BVA0575	<input type="checkbox"/> BLACK SPACER TR36

REA Dahlia 30°

Seat Height	Rear wheel size	Seat tilt as delivered	Backrest recline as delivered
330mm	16"	30	30
	20"	10	30
350mm	16"	30	30
	20"	20	20
	22"	10	30
380mm	20"	20	10
	20"	30	20
400mm	22"	20	10
	22"	30	20
430mm	24"	20	10
	16"	30	30
450mm	24"	30	20
	16"	30	30
480mm	16"	30	30
	20"	30	30
500mm	16"	30	30
	20"	30	30

PUSHHANDLES		
BVA0580	<input type="checkbox"/>	ANGLE ADJUSTABLE PUSHBAR
BVA0590	<input type="checkbox"/>	PUSH HANDLES DIVIDED
BVA0600	<input checked="" type="checkbox"/>	PUSH HANDLES CONNECTED,NO HEAD/NECKREST ATTACHMENT
BVA0610	<input type="checkbox"/>	PUSH HANDLES CONNECTED, WITH REA HEAD/NECKREST ATTACHMENT (NC : 1112)
SEATING		
SEAT TILT TYPE		
BVA0620	<input checked="" type="checkbox"/>	ASSISTANT OPERATED GAS PISTON FOR SEAT TILT BACKREST RECLINE MANDATORY (NC : 0915, 0920, 0930, 0950, 0975, 0980, 0990, 1010, 1750, 1755, 1756, 1757, 1760, 1765, 1770, 1775, 1780, 1790, 1850, 1890)
BVA0630	<input type="checkbox"/>	USER OPERATED TILT (SELF-TILT) WITH FIXATED ARMREST POLE BACKREST RECLINE MANDATORY (NC : 0915, 0920, 0930, 0950, 0975, 0980, 0990, 1010, 1750, 1755, 1756, 1757, 1760, 1765, 1770, 1775, 1780, 1790, 1850)
BVA0632	<input type="checkbox"/>	USER OPERATED TILT (SELF-TILT) WITH STANDARD ARMREST POLE BACKREST RECLINE MANDATORY (NC : 0915, 0920, 0930, 0950, 0975, 0980, 0990, 1010, 1750, 1755, 1756, 1757, 1760, 1765, 1770, 1775, 1780, 1790, 1850)
BVA0637	<input type="checkbox"/>	ELECTRIC OPERATED SEAT AND BACKREST TILT (NC : 0640, 0660, 0670, 1810)
BACKREST RECLINE OPERATION TYPE		
BVA0640	<input type="checkbox"/>	ASSISTANT OPERATED GAS PISTON FOR BACKREST RECLINE
BVA0660	<input checked="" type="checkbox"/>	FIXED BACKREST ANGLE (BARREL SCREW)
LEGRESTS / FOOTRESTS		
CENTRAL LEGRESTS TYPE		
BVA0670	<input type="checkbox"/>	CENTRAL LEGREST, INCL CALFPAD AND FOOTREST (NC : 0690, 0700, 0710, 0720, 0730, 0750, 0760, 0770, 0780, 0800, 0810, 1040, 1310, 1320, 1330, 1340, 1350)
LEGRESTS / HANGERS ADJUSTMENT		
BVA0680	<input checked="" type="checkbox"/>	LOW POSITION, STANDARD
BVA0690	<input type="checkbox"/>	HIGH POSITION FOR LONGER LEGS, 50 MM HIGHER
LEGRESTS / HANGERS LH TYPE		
BVA0700	<input type="checkbox"/>	HANGER ANGLE ADJUSTABLE - ALU 2008 LH. WITH KNEE PAD AND CALF SUPPORT (NC : 1330)
BVA0710	<input checked="" type="checkbox"/>	HANGER FIXED 80 DEG, SWING IN SWING OUT - ALU 2003, LH
BVA0720	<input type="checkbox"/>	HANGER FIXED 90 DEG, SWING IN SWING OUT - ALU 2003, LH (NC : 0820, 0840, 0850, 1310, 1480, 1490, 1500, 1510)
BVA0730	<input type="checkbox"/>	LEGREST AMPUTEE - ALU 2008 LH (NC : 0800, 0810, 0820, 0830, 0840, 0850, 1310, 1330, 1340, 1350)
BVA0740	<input type="checkbox"/>	WITHOUT LEGREST, LH (NC : 0800, 0810, 0820, 0830, 0840, 0850, 1310, 1330, 1340, 1350)
LEGRESTS / HANGERS RH TYPE		
BVA0750	<input type="checkbox"/>	HANGER ANGLE ADJUSTABLE - ALU 2008 RH. WITH KNEE PAD AND CALF SUPPORT (NC : 1330)
BVA0760	<input checked="" type="checkbox"/>	HANGER FIXED 80 DEG, SWING IN SWING OUT - ALU 2003, RH
BVA0770	<input type="checkbox"/>	HANGER FIXED 90 DEG, SWING IN SWING OUT - ALU 2003, RH (NC : 0860, 0880, 0890, 1320, 1480, 1490, 1500, 1510)
BVA0780	<input type="checkbox"/>	LEGREST AMPUTEE - ALU 2008 RH (NC : 0800, 0810, 0860, 0870, 0880, 0890, 1320, 1330, 1340, 1350)
BVA0790	<input type="checkbox"/>	WITHOUT LEGREST, RH (NC : 0800, 0810, 0860, 0870, 0880, 0890, 1320, 1330, 1340, 1350)
FOOTRESTS ONE PIECE		
BVA0800	<input type="checkbox"/>	FOOTREST ONE PIECE DEPTH AND ANGLE ADJUSTABLE SHORT TUBE. (NC : 0820, 0830, 0840, 0850, 0860, 0870, 0880, 0890)
BVA0810	<input type="checkbox"/>	FOOTREST ONE PIECE DEPTH AND ANGLE ADJUSTABLE LONG TUBE. (NC : 0820, 0830, 0840, 0850, 0860, 0870, 0880, 0890)
FOOTRESTS LH TYPE		
BVA0820	<input type="checkbox"/>	FOOTREST LOWER PART, FIXED WITH SHORT TUBE AND HEEL STRAP, LH. NOT TO COMBINE WITH FOOTPLATE EXTENDER. (NC : 1310, 1330, 1340)
BVA0830	<input checked="" type="checkbox"/>	FOOTREST LOWER PART, DEPTH AND ANGLE ADJUSTABLE WITH SHORT TUBE AND HEEL STRAP, LH (NC : 1330)
BVA0840	<input type="checkbox"/>	FOOTREST LOWER PART, FIXED WITH LONG TUBE AND HEEL STRAP, LH. NOT TO COMBINE WITH FOOTPLATE EXTENDER. (NC : 1310, 1330, 1340)
BVA0850	<input type="checkbox"/>	FOOTREST LOWER PART, DEPTH AND ANGLE ADJUSTABLE WITH LONG TUBE AND HEEL STRAP, LH (NC : 1330)
FOOTRESTS RH TYPE		
BVA0860	<input type="checkbox"/>	FOOTREST LOWER PART, FIXED WITH SHORT TUBE AND HEEL STRAP, RH. NOT TO COMBINE WITH FOOTPLATE EXTENDER. (NC : 1320, 1330, 1340)
BVA0870	<input checked="" type="checkbox"/>	FOOTREST LOWER PART, DEPTH AND ANGLE ADJUSTABLE WITH SHORT TUBE AND HEEL STRAP, RH (NC : 1330)
BVA0880	<input type="checkbox"/>	FOOTREST LOWER PART, FIXED WITH LONG TUBE AND HEEL STRAP, RH. NOT TO COMBINE WITH FOOTPLATE EXTENDER. (NC : 1320, 1330, 1340)
BVA0890	<input type="checkbox"/>	FOOTREST LOWER PART, DEPTH AND ANGLE ADJUSTABLE WITH LONG TUBE AND HEEL STRAP, RH (NC : 1330)
ARMRESTS (Heigh adjustable, Autolock attachment included)		
ARMRESTS LH TYPE		
BVA0900	<input type="checkbox"/>	STRAIGHT NARROW PAD, LH (NC : 1756)
BVA0910	<input checked="" type="checkbox"/>	STANDARD PAD, LH (NC : 1025, 1756)
BVA0915	<input type="checkbox"/>	HIGH ARMREST (+55MM), STANDARD PAD, LH (NC : 1025, 1756)
BVA0920	<input type="checkbox"/>	WIDE PAD, LH (NC : 1025, 1756)
BVA0930	<input type="checkbox"/>	HEMIPLEGIC ARMRESTS, LH. NOT TO COMBINE WITH TABLE ACCESSORIES. (NC : 1025, 1750, 1755, 1756, 1760, 1780, 1790)
BVA0950	<input type="checkbox"/>	WITHOUT ARMREST, LH (NC : 1020, 1025, 1750, 1755, 1760, 1780, 1790)
ARMRESTS RH TYPE		
BVA0960	<input type="checkbox"/>	STRAIGHT NARROW PAD, RH (NC : 1757)
BVA0970	<input checked="" type="checkbox"/>	STANDARD PAD, RH (NC : 1025, 1757)
BVA0975	<input type="checkbox"/>	HIGH ARMREST (+55MM), STANDARD PAD, RH (NC : 1025, 1757)
BVA0980	<input type="checkbox"/>	WIDE PAD RH (NC : 1025, 1757)
BVA0990	<input type="checkbox"/>	HEMIPLEGIC ARMRESTS, RH. NOT TO COMBINE WITH TABLE ACCESSORIES. (NC : 1025, 1750, 1755, 1756, 1760, 1770, 1775, 1780, 1790)
BVA1010	<input type="checkbox"/>	WITHOUT ARMREST, RH (NC : 1020, 1025, 1750, 1755, 1760, 1770, 1775, 1780, 1790)
SEATING OPTIONS		
ARMREST OPTIONS		
BVA1020	<input type="checkbox"/>	PAD, SIDE REST
BVA1025	<input type="checkbox"/>	PADDED ARMPAD COVER, BLACK PLUSH TR35
SEAT OPTIONS		
BVA1030	<input type="checkbox"/>	INCONTINENCE COVER
BVA1040	<input type="checkbox"/>	ABDUCTION CUSHION (NC : 1060)
BVA1050	<input type="checkbox"/>	SEAT PLATE EXTENSION IN WIDTH, FORWARD, REAR OR ASSYMETRIC (AS KIT, TWO UNITS OF PLATES) (NC : 1060)
BVA1060	<input type="checkbox"/>	2 CM SEAT HEIGHT INCREASED SEAT PLATE
BVA1070	<input type="checkbox"/>	TILT SCALE, SEATING
BACKREST OPTIONS		
BVA1080	<input type="checkbox"/>	HEADREST, SIDEWISE ADJUSTABLE (NC : 1110, 1112, 1120)
BVA1100	<input type="checkbox"/>	HEADREST WITH CHEEK SUPPORT, SIDEWISE ADJUSTABLE (NC : 1090, 1110, 1112, 1120)
BVA1090	<input type="checkbox"/>	COVER FOR HEADREST (NC : 1110, 1112, 1120)
BVA1110	<input checked="" type="checkbox"/>	NECKREST (NC : 1112)
BVA1112	<input type="checkbox"/>	NECKREST MATRX (NC : 1120, 1130)
BVA1120	<input type="checkbox"/>	COVER FOR NECKREST
BVA1125	<input type="checkbox"/>	UNIVERSAL HEAD/NECKREST HOLDER
BVA1130	<input type="checkbox"/>	HANDLES CLINIC KIT FOR HEAD-/NECKREST
TRUNK SUPPORT LH		
BVA1140	<input type="checkbox"/>	TRUNK SUPPORT MULTI ADJUSTABLE WITH STANDARD ARM (LH) (NC : 1260)
BVA1160	<input type="checkbox"/>	TRUNK SUPPORT MULTI ADJUSTABLE WITH EXTENDED ARM IN HEIGHT (LH) (NC : 1260)
BVA1180	<input type="checkbox"/>	TRUNK SUPPORT WITH STANDARD ARM (LH) (NC : 1260)
BVA1200	<input type="checkbox"/>	TRUNK SUPPORT WITH EXTENDED ARM IN HEIGHT (LH) (NC : 1260)
BVA1220	<input type="checkbox"/>	TRUNK SUPPORT SWINGAWAY (LH) (NC : 1260)
BVA1240	<input type="checkbox"/>	TRUNK SUPPORT MULTI ADJUSTABLE SWINGAWAY (LH) (NC : 1260)
TRUNK SUPPORT RH		
BVA1150	<input type="checkbox"/>	TRUNK SUPPORT MULTI ADJUSTABLE WITH STANDARD ARM (RH) (NC : 1260)
BVA1170	<input type="checkbox"/>	TRUNK SUPPORT MULTI ADJUSTABLE WITH EXTENDED ARM IN HEIGHT (RH) (NC : 1260)
BVA1190	<input type="checkbox"/>	TRUNK SUPPORT WITH STANDARD ARM (RH) (NC : 1260)
BVA1210	<input type="checkbox"/>	TRUNK SUPPORT WITH EXTENDED ARM IN HEIGHT (RH) (NC : 1260)
BVA1230	<input type="checkbox"/>	TRUNK SUPPORT SWINGAWAY (RH) (NC : 1260)
BVA1250	<input type="checkbox"/>	TRUNK SUPPORT MULTI ADJUSTABLE SWINGAWAY (RH) (NC : 1260)
OTHERS BACKREST OPTIONS		
BVA1260	<input type="checkbox"/>	TRUNK SUPPORT ATTACHMENT
BVA1265	<input type="checkbox"/>	LATERAL POSITIONING PADS (WITH VELCRO BETWEEN COVER AND BACKREST UPHOLSTERY), PAIR
BVA1270	<input type="checkbox"/>	LUMBAR SUPPORT CUSHION
BVA1280	<input type="checkbox"/>	WEDGES (ONE PAIR)
BVA1290	<input type="checkbox"/>	PRIVACY ATTACHMENT
BVA1300	<input type="checkbox"/>	TILT SCALE, BACKREST RECLINE
LEGREST OPTIONS		
BVA1310	<input type="checkbox"/>	FOOTPLATE EXTENDER LH (NC : 1340)
BVA1320	<input type="checkbox"/>	FOOTPLATE EXTENDER RH (NC : 1340)
BVA1330	<input type="checkbox"/>	CALF STRAP (FOR FOOTREST ONE PIECE)
BVA1340	<input type="checkbox"/>	FOOT PLATE PADDING
BVA1350	<input type="checkbox"/>	CALF STRAP PADDED

CHASSIS		
COLOUR STRAPS CHASSIS		
BVA1360	<input type="checkbox"/>	BLUE
BVA1370	<input type="checkbox"/>	SILVER
BVA1380	<input type="checkbox"/>	RED
BVA1390	<input type="checkbox"/>	CHAMPAGNE
BVA1400	<input type="checkbox"/>	PETROLIUM GREEN
BVA1405	<input checked="" type="checkbox"/>	WITHOUT COLOUR STRAPS
CASTOR ATTACHMENT TYPE		
BVA1410	<input checked="" type="checkbox"/>	FIXED
BVA1420	<input type="checkbox"/>	QR
FRONT CASTORS		
BVA1430	<input type="checkbox"/>	100X27 MEDIUM HARD (4")
BVA1440	<input type="checkbox"/>	125X27 MEDIUM HARD (5")
BVA1450	<input type="checkbox"/>	140X40 SOFT
BVA1460	<input checked="" type="checkbox"/>	150X27 MEDIUM HARD
BVA1470	<input type="checkbox"/>	150X30 SOFT
BVA1480	<input type="checkbox"/>	200X27 MEDIUM HARD
BVA1490	<input type="checkbox"/>	200X30 SOFT
BVA1500	<input type="checkbox"/>	200X32 PNEUMATIC
BVA1510	<input type="checkbox"/>	200X45 MEDIUM HARD
REAR WHEEL ATTACHMENT TYPE		
BVA1520	<input type="checkbox"/>	QR
BVA1530	<input type="checkbox"/>	FIXED (NC : 1725, 1727)
REAR WHEELS SIZE/TYRE		
BVA1540	<input type="checkbox"/>	16" (NC : 1590, 1595, 1600, 1610, 1630, 1650, 1660, 1670, 1680, 1690, 1710, 1725, 1727)
BVA1550	<input type="checkbox"/>	20" (NC : 1590, 1610, 1620, 1630, 1670, 1680, 1690, 1700)
BVA1560	<input type="checkbox"/>	22" (NC : 1590, 1670, 1700)
BVA1570	<input checked="" type="checkbox"/>	24" STANDARD REAR WHEELS
BVA1580	<input type="checkbox"/>	24" BASIC REAR WHEELS (NC : 1590, 1595, 1610, 1660, 1670, 1680, 1690, 1700)
TYRE TYPE		
BVA1590	<input type="checkbox"/>	LOW PROFILE PKT (24X1) (NC : 1650, 1660, 1680, 1690, 1700)
BVA1595	<input type="checkbox"/>	PNEUMATIC NON MARKING (NC : 1680, 1690, 1880)
BVA1600	<input type="checkbox"/>	PNEUMATIC BLOCK PATTERN (NC : 1880)
BVA1610	<input type="checkbox"/>	PNEUMATIC BLOCK PATTERN - PUNCTURE PROOF (NC : 1680, 1690, 1880)
BVA1620	<input checked="" type="checkbox"/>	TYRES SOLID
BVA1630	<input type="checkbox"/>	PNEUMATIC LOW PROFILE (NC : 1680, 1690, 1880)
HANDRIM TYPE		
BVA1650	<input checked="" type="checkbox"/>	ALUMINIUM 4-P
BVA1660	<input type="checkbox"/>	PLASTIC COATED 4-P
BVA1670	<input type="checkbox"/>	BRUSHED STAINLESS STEEL
BVA1680	<input type="checkbox"/>	ONE ARM DRIVE LH (NC : 1710, 1810)
BVA1690	<input type="checkbox"/>	ONE ARM DRIVE RH (NC : 1710, 1810)
BVA1700	<input type="checkbox"/>	NO HANDRIM
WHEEL OPTIONS		
SPOKEGUARD		
BVA1710	<input type="checkbox"/>	FROST
SAFETY		
BRAKE TYPE		
BVA1720	<input checked="" type="checkbox"/>	USER OPERATED TYRE BRAKES
BVA1725	<input type="checkbox"/>	ONE ARM BRAKE LH
BVA1727	<input type="checkbox"/>	ONE ARM BRAKE RH
SAFETY OPTIONS		
BVA1730	<input type="checkbox"/>	PELVIC BELT
BVA1740	<input type="checkbox"/>	REFLECTORS KIT (WHEELS AND CHASSIS)
ACCESSORIES		
TABLE TRAY		
BVA1750	<input type="checkbox"/>	TABLE TRAY INCL MOUNTING KIT (NC : 1765, 1790)
BVA1755	<input type="checkbox"/>	TABLE TRAY INCL MOUNTING KIT, WITH ELBOW PADS (NC : 1765, 1790)
BVA1756	<input type="checkbox"/>	HALF TABLE TRAY, CLEAR, FLIP UP, LH (ARMREST WITH STRAIGHT NARROW PAD INCL.) (NC : 1760, 1790)
BVA1757	<input type="checkbox"/>	HALF TABLE TRAY, CLEAR, FLIP UP, RH (ARMREST WITH STRAIGHT NARROW PAD INCL.) (NC : 1760, 1790)
BVA1770	<input type="checkbox"/>	SWING AWAY TABLE (NC : 1760, 1780, 1790)
BVA1775	<input type="checkbox"/>	SWING AWAY TABLE WITH ELBOW PADS (NC : 1760, 1780, 1790)
TABLE TRAY CUSHION		
BVA1760	<input type="checkbox"/>	TABLE CUSHION (NC : 1790)
BVA1765	<input type="checkbox"/>	PADDING (NC : 1790)
OTHERS ACCESSORIES		
BVA1780	<input type="checkbox"/>	LOCKING DEVICE FOR TABLE TRAY (NC : 1790)
BVA1790	<input type="checkbox"/>	MOUNTING KIT FOR TABLE (IF ORDERED SEPARATELY)
BVA1810	<input type="checkbox"/>	MOUNTING BRACKET VIAMOBIL (NC : 1870, 1900)
BVA1820	<input type="checkbox"/>	PUMP
BVA1830	<input type="checkbox"/>	CANE HOLDER
BVA1840	<input type="checkbox"/>	CLINIC KIT
BVA1850	<input type="checkbox"/>	HARNES ATTACHMENT BRACKET SEAT (PAIR), FRONT
BVA1860	<input type="checkbox"/>	HARNES ATTACHMENT BACKREST, UPPER
BVA1870	<input type="checkbox"/>	KERB LEVER (NC : 1900)
BVA1880	<input type="checkbox"/>	WASH SEALING KIT
BVA1890	<input type="checkbox"/>	TRANSPORTATION BRACKETS FOR OCCUPIED AND NON OCCUPIED TRANSPORT
BVA1900	<input type="checkbox"/>	VENTILATOR TRAY
BVA1920	<input type="checkbox"/>	DRIP STAND
BVA1930	<input type="checkbox"/>	RESTRICTOR KIT
SPECIAL		
FREIGHT & PRICE		
PACKAGING		
FREIGHT PRICE		
TOTAL PRICE		