

	POST 4		Kunde/HMS: Kontaktperson: Adresse: Telefon: Mailadresse: Bestillingnr./SF nr:	
--	---------------	--	--	---

☐ Please make your selection with a cross in the box : ◆ Standard feature ○ Free of Charge (NC9999) Non Conformance with ACTION5/MyOn HC Option 9999.

Gjeldende fra juli 2018

LPF-NO-MyOn HC-718.

LEGENDA B Both North America and Europe N North American option ONLY E European option ONLY	GPF DATAS
---	------------------

CODE	MODEL	EU DESCRIPTION	COMMENTS		(FOR EMEA)	(NC: 5600)
SJC0008	E	◆ MyOn HC				
SPECIAL						
N	E	◆ NO				(NC: 6099)
SJCSPEC	E	☐ SPECIAL			(MANDATORY FOR SPECIAL COLOUR)	
BASE MODEL FOR PRICE LIST						
SJC0100	E	◆ BASE MODEL				
FRAME TYPE						
SJC0200	B	◆ SWING AWAY FRAME				(NC: 3185, 3230, 3440, 3450)
SJC0210	B	☐ FIXED FRONT END 80 DEGREE ANGLE FRAME			(NC: 2800, 2810, 2820, 2830, 2840, 2850, 3000, 3010, 3020, 3030, 3040, 3050, 3100, 3110, 3120, 3130, 3140, 3150, 3400, 3410, 3420, 3430)	
PROPULSION TYPE						
SJC0300	B	◆ SELF PROPULSION				
SEAT						
SEAT WIDTH (SW)						
SJC1000	B	☐ 305 MM / 12"				Max user weight: 100kg / 220LBS (NC: 1150, 1160, 1190, 2460, 2465, 2470, 2475, 2480, 2840, 2850, 3400, 5005, 5006, 5007, 5008, 5009, 5010, 5020, 5030, 5040, 5042, 5044, 5090, 5600, 8170, 8180)
SJC1010	B	☐ 330 MM / 13"				Max user weight: 100kg / 220LBS (NC: 1150, 1160, 1190, 2460, 2465, 2470, 2475, 2480, 2840, 2850, 3400, 5005, 5006, 5007, 5008, 5009, 5010, 5020, 5030, 5040, 5042, 5044, 5090, 5600, 8170, 8180)
SJC1020	B	☐ 355 MM / 14"				Max user weight: 130kg / 290LBS (NC: 2460, 2465, 2470, 2475, 2480, 2840, 2850, 5005, 5006, 5007, 5008, 5009, 5010, 5020, 5030, 5040, 5042, 5044, 5090, 5600, 8170, 8180)
SJC1030	B	☐ 380 MM / 15"				Max user weight: 130kg / 290LBS (NC: 2850, 5005, 5008, 5090)
SJC1040	B	☐ 405 MM / 16"				Max user weight: 130kg / 290LBS (NC: 2850)
SJC1050	B	☐ 430 MM / 17"				Max user weight: 130kg / 290LBS
SJC1060	B	☐ 455 MM / 18"				Max user weight: 130kg / 290LBS
SJC1070	B	☐ 480 MM / 19"				Max user weight: 130kg / 290LBS (NC: 5005, 5006, 5007, 5090)
SJC1075	B	☐ 505 MM / 20"				Max user weight: 130kg / 290LBS - NON HEAVY DUTY (NC: 1100, 1110, 1170, 5005, 5006, 5007, 5008, 5009, 5090)
SJC1078	B	☐ 530 MM / 21"				Max user weight: 130kg / 290LBS - NON HEAVY DUTY (NC: 1100, 1110, 1170, 2460, 2465, 2470, 2475, 2480, 5005, 5006, 5007, 5008, 5009, 5010, 5020, 5030, 5040, 5042, 5044, 5090)
SEAT DEPTH (SD)						
SJC1140	B	☐ 450 MM / 18"				(NC: 5054)
SJC1170	B	☐ 350-400 MM / 14"-16" ADJUSTABLE			(Delivered at 375 mm / 15")	(NC: 1510, 3950, 4050, 4150, 5005, 5006, 5007, 5008, 5009, 5010, 5020, 5030, 5040, 5042, 5044, 5090, 7600)
SJC1180	B	☐ 400-450 MM / 16"-18" ADJUSTABLE			(Delivered at 425 mm / 17")	(NC: 1510, 5054)
SJC1190	B	☐ 450-500 MM / 18"-20" ADJUSTABLE			(Delivered at 475 mm / 19")	(NC: 1510, 5006, 5007)
SEAT HEIGHT IN FRONT (FSTF) OF THE CHAIR, MEASURED FROM FLOOR TO CENTER OF SEAT						
SJC1200	B	☐ 360 MM / 14"				(NC: 1310, 1320, 1330, 1340, 1350, 1360, 2840, 3040, 3140, 3420, 3430, 3450, 7040, 7050, 7055, 7060, 7070, 7080, 7090, 7100, 7110)
SJC1210	B	☐ 385 MM / 15"				(NC: 1320, 1330, 1340, 1350, 1360, 2840, 3040, 3140, 3430, 7000, 7020, 7030, 7040, 7070, 7080, 7090, 7100, 7110)
SJC1220	B	☐ 410 MM / 16"				(NC: 1330, 1340, 1350, 1360, 7080, 7090, 7100, 7110)
SJC1230	B	☐ 435 MM / 17"				(NC: 1300, 1340, 1350, 1360, 7110)
SJC1240	B	☐ 460 MM / 18"				(NC: 1300, 1310, 1350, 1360)
SJC1250	B	◆ 485 MM / 19"				(NC: 1300, 1310, 1320, 1360)
SJC1260	B	☐ 510 MM / 20"				(NC: 1300, 1310, 1320, 1330, 7000, 7020, 7030, 7040)
SJC1270	B	☐ 535 MM / 21"				(NC: 1300, 1310, 1320, 1330, 1340, 2830, 3030, 3130, 7000, 7020, 7030, 7040, 7050, 7055, 7060, 7070, 7080, 7090, 7100)
SEAT HEIGHT IN REAR (RSTF) OF THE CHAIR, MEASURED FROM FLOOR TO CENTER OF SEAT REAR						
SJC1300	B	☐ 360 MM / 14"				(NC: 7250, 7260, 8900)
SJC1310	B	☐ 385 MM / 15"				
SJC1320	B	☐ 410 MM / 16"				
SJC1330	B	☐ 435 MM / 17"				
SJC1340	B	◆ 460 MM / 18"				
SJC1350	B	◆ 485 MM / 19"				(NC: 7200, 7210, 7220)
SJC1360	B	☐ 510 MM / 20"				(NC: 7200, 7210, 7220, 7230, 7240, 7242, 8900)
SEAT UPHOLSTERY TYPE						
SJC1490	B	☐ WITHOUT SEAT UPHOLSTERY				
SJC1500	B	◆ NYLON SLING STEPLESS DEPTH ADJUSTABLE (-50 MM / -2") SEAT UPHOLSTERY			WITH PHONE HOLDER UNDER SEAT	
SJC1510	B	☐ NYLON BASIC FIXED SEAT UPHOLSTERY				
SJC1520	B	☐ TENSION ADJUSTABLE SEAT UPHOLSTERY & DEPTH ADJUSTABLE -50MM				

BACKREST			
BACKREST ATTACHEMENT TYPE			
SIC2000	E	<input type="checkbox"/>	FIXED ANGLE, 90° VS SEATRAIL
SIC2010	E	<input type="checkbox"/>	ANGLE ADJUSTABLE FROM -15° TO + 15°, SET AT 90° VS SEATRAIL
SIC2020	B	<input checked="" type="checkbox"/>	ANGLE ADJUSTABLE AND FOLD DOWN (-15° to + 15°)
FOLDS FORWARD ON SEAT & LOCKS IN POSITION			
BACKREST FIXED HEIGHT (UPHOLSTERY MOVES)			
MEASURED FROM THE TOP OF THE SEAT UPHOLSTERY TO THE TOP OF THE BACKREST UPHOLSTERY			
SIC2130	B	<input type="checkbox"/>	FIXED 7", HEIGHT 300MM / 12" WITH FIXED SHORT PUSHHANDLES
SIC2140	B	<input type="checkbox"/>	FIXED 7", HEIGHT 325MM / 13" WITH FIXED SHORT PUSHHANDLES
SIC2143	B	<input type="checkbox"/>	FIXED 7", HEIGHT 350MM / 14" WITH FIXED SHORT PUSHHANDLES
SIC2150	B	<input type="checkbox"/>	FIXED 7", HEIGHT 375MM / 15" WITH FIXED LONG PUSHHANDLES
SIC2153	B	<input type="checkbox"/>	FIXED 7", HEIGHT 400MM / 16" WITH FIXED LONG PUSHHANDLES
SIC2160	B	<input type="checkbox"/>	FIXED 7", HEIGHT 425MM / 17" WITH FIXED LONG PUSHHANDLES
SIC2163	B	<input type="checkbox"/>	FIXED 7", HEIGHT 450 MM / 18" WITH FIXED LONG PUSHHANDLES
SIC2170	B	<input type="checkbox"/>	FIXED 7", HEIGHT 475 MM / 19" WITH FIXED LONG PUSHHANDLES
SIC2173	B	<input type="checkbox"/>	FIXED 7", HEIGHT 500 MM / 20" WITH FIXED LONG PUSHHANDLES
HEIGHT 300 MM / 12" (NC: 2300, 2310, 2400, 2410, 2420, 2505, 8160, 8200) HEIGHT 325-350 MM / 13"-14" SET AT 325MM / 13" (2 positions, 25mm / 1" steps through uph) (NC: 2300, 2310, 2400, 2410, 2420, 2505, 8160, 8200) HEIGHT 375-400 MM / 15"-16" SET AT 375MM / 15" (2 positions, 25mm / 1" steps through uph) (NC: 2300, 2310, 2400, 2410, 2420, 2505, 8160, 8200) HEIGHT 375-400 MM / 15"-16" SET AT 400MM / 16" (2 positions, 25mm / 1" steps through uph) (NC: 2300, 2310, 2400, 2410, 2420, 8200) HEIGHT 425-450 MM / 17"-18" SET AT 425MM / 17" (2 positions, 25mm / 1" steps through uph) (NC: 2300, 2310, 2400, 2410, 2420) HEIGHT 425-450 MM / 17"-18" SET AT 450MM / 18" (2 positions, 25mm / 1" steps through uph) (NC: 2300, 2310, 2400, 2410, 2420) HEIGHT 475-500 MM / 19"-20" SET AT 475MM / 19" (2 positions, 25mm / 1" steps through uph) (NC: 2300, 2310, 2400, 2410, 2420) HEIGHT 475-500 MM / 19"-20" SET AT 500MM / 20" (2 positions, 25mm / 1" steps through uph) (NC: 2300, 2310, 2400, 2410, 2420)			
BACKREST HEIGHT ADJUSTABLE (BACKREST TUBE MOVES)			
MEASURED FROM THE TOP OF THE SEAT UPHOLSTERY TO THE TOP OF THE BACKREST UPHOLSTERY			
SIC2300	B	<input type="checkbox"/>	HEIGHT ADJUSTABLE FROM 360-400MM / 14.5"-16"
SIC2310	B	<input checked="" type="checkbox"/>	HEIGHT ADJUSTABLE FROM 410-490MM / 16"-19.5"
SET AT 400MM / 16" (NC: 2460, 2465, 2470, 2475, 2480, 2505) SET AT 450MM / 18" (NC: 2460, 2465, 2470, 2475, 2480)			
PUSHHANDLES TYPE (TO SELECT WHEN CHOSEN HEIGHT ADJUSTABLE BACKREST)			
SIC2400	B	<input type="checkbox"/>	WITHOUT PUSHHANDLES (7° degree bend)
SIC2410	B	<input checked="" type="checkbox"/>	REARSET (REA TYPE) (HEIGHT ADJUSTABLE + 150 MM / 6") (7° degree bend)
SIC2420	B	<input type="checkbox"/>	INTEGRATED (HEIGHT ADJUSTABLE + 135 MM / 5") (0° degree bend)
(NC: 8160, 9100, 9200) (NC: 2505)			
BACKREST MATRIX OPTIONS TYPE			
SIC2460	E	<input type="checkbox"/>	MATRIX ELITE BACK WITH HARDWARE 405 MM BACK HEIGHT
SIC2465	E	<input type="checkbox"/>	MATRIX ELITE BACK WITH HARDWARE 455 MM BACK HEIGHT
SIC2470	E	<input type="checkbox"/>	MATRIX ELITE DEEP BACK WITH HARDWARE 405 MM BACK HEIGHT
SIC2475	E	<input type="checkbox"/>	MATRIX ELITE DEEP BACK WITH HARDWARE 455 MM BACK HEIGHT
SIC2480	E	<input type="checkbox"/>	MATRIX HARDWARE EASY SET (WITHOUT MATRIX BACK)
ADDITIONAL DELIVERY TIME + 5 DAYS (NC: 2500, 2505, 2510, 5500, 5600) ADDITIONAL DELIVERY TIME + 5 DAYS (NC: 2500, 2505, 2510, 5500, 5600) ADDITIONAL DELIVERY TIME + 5 DAYS (NC: 2500, 2505, 2510, 5500, 5600) ADDITIONAL DELIVERY TIME + 5 DAYS (NC: 2500, 2505, 2510, 5500, 5600) ADDITIONAL DELIVERY TIME + 5 DAYS (NC: 2500, 2505, 2510, 5500, 5600)			
UPHOLSTERY TYPE			
SIC2500	B	<input type="checkbox"/>	WITHOUT BACKREST UPHOLSTERY / OMIT UPHOLSTERY
SIC2505	B	<input type="checkbox"/>	FIXED BACKREST, NYLON & MESH PADDED COVER
SIC2510	B	<input checked="" type="checkbox"/>	CONTOUR (TENSION) ADJUSTABLE STRAPS, NYLON & MESH COVER
WITH REFLECTIVE STRIPE & BACK POCKET (NC: 5500, 5600) WITH REFLECTIVE STRIPE & BACK POCKET			
HANGERS/FOOTPLATES			
HANGERS RH TYPE			
SIC3000	B	<input type="checkbox"/>	WITHOUT RH HANGER
SIC3010	B	<input type="checkbox"/>	70° STRAIGHT RH
SIC3020	B	<input checked="" type="checkbox"/>	80° STRAIGHT RH
SIC3030	B	<input type="checkbox"/>	90° STRAIGHT RH
SIC3040	B	<input type="checkbox"/>	ELEVATING LEGREST (0° TO 80°) RH
SIC3050	B	<input type="checkbox"/>	AMPUTEE LEGREST RH
(NC: 3170, 3175, 3185, 3200, 3210, 3230, 5400, 9400) (NC: 3170, 3175, 3185, 3440, 3450) (NC: 3170, 3175, 3185, 3440, 3450) (NC: 3170, 3175, 3185, 3440, 3450) (NC: 3170, 3175, 3185, 3200, 3440, 3450, 7080, 7090, 7100, 7110) (NC: 3170, 3175, 3185, 3230, 3400, 3440, 3450, 5400) (NC: 3170, 3175, 3185, 3200, 3210, 3230, 5400, 8100, 8120, 9400)			
HANGERS LH TYPE			
SIC3100	B	<input type="checkbox"/>	WITHOUT LH HANGER
SIC3110	B	<input type="checkbox"/>	70° STRAIGHT LH
SIC3120	B	<input checked="" type="checkbox"/>	80° STRAIGHT LH
SIC3130	B	<input type="checkbox"/>	90° STRAIGHT LH
SIC3140	B	<input type="checkbox"/>	ELEVATING LEGREST (0° TO 80°) LH
SIC3150	B	<input type="checkbox"/>	AMPUTEE LEGREST LH
(NC: 3170, 3175, 3185, 3300, 3310, 5400, 9400) (NC: 3170, 3175, 3185, 3440, 3450) (NC: 3170, 3175, 3185, 3440, 3450) (NC: 3170, 3175, 3185, 3200, 3440, 3450, 7080, 7090, 7100, 7110) (NC: 3170, 3175, 3185, 3400, 3440, 3450, 5400) (NC: 3170, 3175, 3185, 3300, 3310, 5400, 8100, 8120, 9400)			
FOOTPLATES RH TYPE			
SIC3200	B	<input checked="" type="checkbox"/>	FLIP UP RH FOOTPLATE
SIC3210	B	<input type="checkbox"/>	FLIP UP DEPTH & ANGLE ADJUSTABLE RH FOOTPLATE
SIC3230	E	<input type="checkbox"/>	ONE PIECE FOOTPLATE DEPTH & ANGLE ADJUSTABLE
WITH ADJUSTABLE HEEL STRAP (NC: 3410) WITH ADJUSTABLE HEEL STRAP (NC: 3400, 3410, 3420, 3430)			
FOOTPLATES LH TYPE			
SIC3300	B	<input checked="" type="checkbox"/>	FLIP-UP LH FOOTPLATE
SIC3310	B	<input type="checkbox"/>	FLIP-UP DEPTH & ANGLE ADJUSTABLE LH FOOTPLATE
WITH ADJUSTABLE HEEL STRAP (NC: 3410) WITH ADJUSTABLE HEEL STRAP			
FOOTREST TUBE (Given UL for adjustable footplate only. Add +15mm with fixed footplate)			
SIC3400	B	<input type="checkbox"/>	SWING AWAY FRAME, HIGH MOUNTED TUBE, FOR UL 160 to 375 MM
SIC3410	B	<input type="checkbox"/>	SWING AWAY FRAME, SHORT TUBE, FOR UL 285 TO 375 MM
SIC3420	B	<input checked="" type="checkbox"/>	SWING AWAY FRAME, STANDARD TUBE, FOR UL 295 to 420 MM
SIC3430	B	<input type="checkbox"/>	SWING AWAY FRAME, LONG TUBE, FOR UL 400 TO 525 MM
SIC3440	B	<input type="checkbox"/>	FIXED FRAME, STANDARD TUBE, FOR UL 320 TO 460 MM
SIC3450	B	<input type="checkbox"/>	FIXED FRAME, LONG TUBE, FOR UL 360 TO 500 MM
(NC: 9400)			
ARMRESTS			
MUDGUARDS / SIDE GUARDS (PAIR)			
SIC3700	B	<input checked="" type="checkbox"/>	WITHOUT SIDE GUARD/MUDGUARD
SIC3750	B	<input type="checkbox"/>	SIDE GUARDS
SIC3800	B	<input type="checkbox"/>	REMOVABLE MUDGUARD
(NC: 3910, 3920, 3930, 3940, 3960, 4010, 4020, 4030, 4040, 4060, 4110, 4120, 4130, 4140, 4160) (NC: 3910, 3920, 3930, 3940, 3960, 4010, 4020, 4030, 4040, 4060, 4110, 4120, 4130, 4140, 4160, 8160)			
ARMRESTS RH FORM			
SIC4000	B	<input type="checkbox"/>	WITHOUT ARMREST
SIC4010	B	<input checked="" type="checkbox"/>	REMOVABLE RH, HEIGHT ADJUSTABLE (N°.1) T-ARM
SIC4020	B	<input type="checkbox"/>	REMOVABLE RH, HEIGHT ADJUSTABLE (N°.1) T-ARM
SIC4030	B	<input type="checkbox"/>	SWING AWAY DETACHABLE HEIGHT ADJ RH (N°3)
SIC4040	B	<input type="checkbox"/>	SWING AWAY DETACHABLE HEIGHT ADJ RH (N°3)
SIC4050	B	<input type="checkbox"/>	TUBULAR ARMREST RH, SWING AWAY PADDED TUBE
SIC4060	B	<input type="checkbox"/>	REMOVABLE HEMI RH (N°4)
DESK ARMPAD (280 MM), ADJUSTABLE ON SEATRAIL BY 150 MM FULL LENGTH ARMPAD (380 MM), ADJUSTABLE ON SEATRAIL BY 50 MM DESK ARMPAD (280 MM) FULL LENGTH ARMPAD (380 MM) HEIGHT ADJUSTABLE 40 MM BY STEP OF 10 MM, PADDING 230 MM ARMPAD HEMI DARTEX 440 MM			
ARMRESTS LH FORM			
SIC4100	B	<input type="checkbox"/>	WITHOUT ARMREST
SIC4110	B	<input checked="" type="checkbox"/>	REMOVABLE LH, HEIGHT ADJUSTABLE (N°.1) T-ARM
SIC4120	B	<input type="checkbox"/>	REMOVABLE LH, HEIGHT ADJUSTABLE (N°.1) T-ARM
SIC4130	B	<input type="checkbox"/>	SWING AWAY DETACHABLE HEIGHT ADJ LH (N°3)
SIC4140	B	<input type="checkbox"/>	SWING AWAY DETACHABLE HEIGHT ADJ LH (N°3)
SIC4150	B	<input type="checkbox"/>	TUBULAR ARMREST LH, SWING AWAY PADDED TUBE
SIC4160	B	<input type="checkbox"/>	REMOVABLE HEMI LH (N°4)
DESK ARMPAD (280 MM), ADJUSTABLE ON SEATRAIL BY 150 MM FULL LENGTH ARMPAD (380 MM), ADJUSTABLE ON SEATRAIL BY 50 MM DESK ARMPAD (280 MM) FULL LENGTH ARMPAD (380 MM) HEIGHT ADJUSTABLE 40 MM BY STEP OF 10 MM, PADDING 230 MM ARMPAD HEMI DARTEX 440 MM			
SEATING OPTIONS TYPE			
CUSHION TYPE			
SIC5000	E	<input type="checkbox"/>	BASIC 50 MM / 2" FOAM CUSHION NYLON & ANTI-SLIP COVER
SIC5010	E	<input type="checkbox"/>	MATRIX FLO-TECH CONTOUR VISCO CUSHION
SIC5020	E	<input type="checkbox"/>	MATRIX FLO-TECH IMAGE CUSHION
SIC5030	E	<input type="checkbox"/>	MATRIX FLO-TECH LITE CUSHION
SIC5040	E	<input type="checkbox"/>	MATRIX FLO-TECH SOLUTION XTRA
SIC5042	E	<input type="checkbox"/>	MATRIX VI CUSHION
SIC5044	E	<input type="checkbox"/>	MATRIX LIBRA CUSHION
(NC: 5050, 5052, 5054) (NC: 5054) (NC: 5054) (NC: 5054) (NC: 5054) (NC: 5054) (NC: 5054)			
MATRIX CUSHION DEPTH			
(MANDATORY WITH CUSHION TYPE EXCEPT THE BASIC & CLASSIC ONE, SJA5000 & SJA5007)			
SIC5050	E	<input type="checkbox"/>	FOR 400MM SEAT DEPTH
SIC5052	E	<input type="checkbox"/>	FOR 450MM SEAT DEPTH
SIC5054	E	<input type="checkbox"/>	FOR 500MM SEAT DEPTH
ONLY WITH MATRIX VI & MATRIX LIBRA (NC: 5090)			
OTHER SEATING OPTIONS TYPE			
SIC5100	B	<input type="checkbox"/>	CLIP TO KEEP THE CHAIR FOLDED
SIC5400	B	<input type="checkbox"/>	ADJUSTABLE CALF STRAP
SIC5500	B	<input type="checkbox"/>	FOLDING BACKBRACE

COLOURS										
FRAME										
SJC6000	B	◆	PEARL WHITE							
SJC6010	B	□	ANTHRACITE BLACK							
SJC6020	B	□	POP BLUE							
SJC6030	B	□	DARK RED							
SJC6040	B	□	ELECTRIC GREEN							
SJC6050	B	□	HYPNOTIC PURPLE							
SJC6060	B	□	PEARL GREY							
SJC6070	B	□	MATT BLACK							
SJC6080	B	□	INTENSE YELLOW							
SJC6090	B	□	VIBRATING ORANGE							
SJC6099	E	□	SPECIAL COLOUR							ON DEMAND ONLY AND SPECIAL (SJSPEC) MANDATORY
ACCENT COLOURS (rear wheel plate, fork, disks on castorhousing, QR axle cap)										
SJC6100	B	◆	WITHOUT ACCENT							SILVER FORK & REARWHEEL SUPPORT
SJC6110	B	□	BLUE							OPTIONAL / UPCHARGE
SJC6200	B	□	PURPLE							OPTIONAL / UPCHARGE
FRONT WHEELS										
SJC7000	B	□	4" X 1.1/4" (100 X 30 MM) SOFTROLL							
SJC7020	B	□	4" X 3/4" (100 X 20 MM) SOLID							
SJC7030	B	□	4" X 1" (100 X 25 MM) SOLID, INTEGRATED "LIGHTS"							
SJC7040	B	□	4" X 1.3/8" (100 X 35 MM) ALU SOFTROLL							
SJC7055	B	◆	5" X 1" (125 X 25 MM) SOLID							
SJC7060	B	□	5" X 1" (125 X 25 MM) SOLID, INTEGRATED "LIGHTS"							
SJC7070	B	□	5" X 1.3/8" (125 X 35 MM) ALU SOFTROLL							
SJC7080	B	□	6" X 1" (150 X 25 MM) SOLID							
SJC7090	B	□	6" X 1.1/4" (150 X 30 MM) PNEUMATIC							
SJC7100	B	□	6" X 1.3/8" (150 X 35 MM) ALU SOFTROLL							
SJC7110	B	□	7" X 1.1/4" (180 X 30 MM) SOLID							
REAR WHEELS										
REAR WHEEL RIM TYPE (12mm QR)										
SJC7200	B	□	22" (560 MM) STANDARD, CROSS SPOKED WHEEL							(NC: 7310, 7320, 7330, 7410, 7425, 7430)
SJC7210	B	□	22" (560 MM) LIGHTWEIGHT, RADIAL SPOKED, BLACK RIM & HUB - 6 POINTS HANDRIM							(NC: 7310, 7320, 7330, 7410, 7425, 7430, 8000, 8160)
SJC7220	B	□	22" (560 MM) COMPOSITE WHEELS							(NC: 7310, 7320, 7410, 7425, 7430, 8000, 8160, 8600)
SJC7230	B	□	24" (610 MM) STANDARD, CROSS SPOKED WHEEL							(NC: 7320, 7330, 7425, 7430)
SJC7240	B	◆	24" (610 MM) LIGHTWEIGHT, RADIAL SPOKED, BLACK RIM & HUB - 6 POINTS HANDRIM							(NC: 7330, 7410, 8000, 8160)
SJC7242	B	□	24" (610 MM) COMPOSITE WHEELS							(NC: 7310, 7320, 7410, 7425, 7430, 8000, 8160)
SJC7244	B	□	24" (610 MM) SPINERGY WHEELS- BLACK							(NC: 7300, 7305, 7330, 7410, 7425, 7430, 8000, 8160)
SJC7250	B	□	25" (635 MM) LIGHTWEIGHT, RADIAL SPOKED, BLACK RIM & HUB - 6 POINTS HANDRIM							(NC: 7300, 7310, 7305, 7330, 7410, 8000, 8160)
SJC7260	B	□	25" (635 MM) SPINERGY WHEELS - BLACK							NC: 7300, 7310, 7305, 7330, 7410, 7420, 7425, 7430, 8000, 8160)
TIRE TYPE										
SJC7300	B	□	SOLID GREY 1.3/8" (35 MM)							
SJC7305	B	□	PNEUMATIC TYRE GREY 1.3/8" (35 MM) BLOCK PATTERN							
SJC7310	B	◆	HP PNEUMATIC TYRE GREY AND BLACK 1" (25 MM)							
SJC7320	B	□	SCHWALBE PNEUMATIC MARATHON PLUS EVO 1" (25 MM)							
SJC7330	B	□	FLAT FREE INSERTS 1.3/8" (35 MM)							
HANDRIM TYPE										
SJC7400	B	◆	ALUMINIUM HANDRIM							
SJC7410	E	□	STAINLESS STEEL							
SJC7420	B	□	HANDRIM COATED IN BLACK ANTI-SLIP							Only with Standard cross spoked wheel 24"
REAR WHEELS AXLE TYPE										
SJC7500	B	◆	QUICK RELEASE AXLE							
REAR WHEELS POSITION										
SJC7600	B	□	VERY ACTIVE WHEEL POSITION							From front of backrest to axle position + 4" (100 mm) (NC: 8900)
SJC7610	B	◆	ACTIVE REAR WHEEL POSITION							From front of backrest to axle position + 3" (75 mm)
SJC7620	B	□	MEDIUM ACTIVE WHEEL POSITION							From front of backrest to axle position + 2" (50 mm) (NC: 8900)
SJC7630	B	□	STANDARD REAR WHEELS POSITION							From front of backrest to axle position + 1" (25 mm)
CAMBER ANGLE										
SJC7700	B	□	0° CAMBER							STANDARD
SJC7710	B	□	1° CAMBER							STANDARD
SJC7720	B	◆	2° CAMBER							STANDARD
SJC7730	B	□	3° CAMBER							OPTIONAL / UPCHARGE
SJC7740	B	□	4° CAMBER							OPTIONAL / UPCHARGE
CHASSIS OPTIONS TYPE										
SJC8000	B	□	SPOKE PROTECTORS PAIR							Clear (except in black for the Standard cross spoked wheel 24" SJA7230)
SJC8020	B	□	PUMP - LOW PRESSURE							
SJC8030	B	□	PUMP - HIGH PRESSURE							
BRAKES TYPE										
SJC8100	B	◆	PUSH TO LOCK BRAKES							
SJC8110	B	□	PULL TO LOCK BRAKES							
SJC8120	B	□	PUSH TO LOCK BRAKES WITH EXTENTION							
SJC8130	B	□	PULL TO LOCK BRAKES WITH EXTENTION							
SJC8160	B	□	ATTENDANT HUB BRAKES (INCLUDES STANDARD USER PUSH/PULL TO LOCK BRAKE DEPENDING ON CONFIGURATION)							22" & 24" STANDARD CROSS SPOKED WHEELS ONLY
SAFETY OPTIONS TYPE										
SJC8200	E	□	TRANSPORTATION KIT INCLUDING BUCKLE BELT							(NC: 8300)
SJC8300	B	□	BUCKLE BELT							
SJC8500	B	◆	REAR ANTI-TIPPER ACTIVE RH							(NC: 9200)
SJC8550	B	◆	REAR ANTI-TIPPER ACTIVE LH							(NC: 9100)
SJC8600	B	□	REFLECTORS KIT (WHEELS)							
ATTENDANT OPTIONS TYPE										
SJC9000	B	□	CANE HOLDER							
SJC9100	B	□	TIPPING AID (LH)							
SJC9200	B	□	TIPPING AID (RH)							
ACCESSORIES TYPE										
SJC9300	B	□	USER TOOLKIT							
SJC9400	B	□	CASE HOLDER							