

FOR COMFORT AND PRESSURE RELIEF



The closed top and the open bottom make the cushion machine washable

Opening in the middle follows the contour of your shower and toilet commode chair

Optimal pressure redistribution by the 4 individual compartments filled with Vicair SmartCells.



Easily attached to a seat by use of the straps at each corner



- ▲ Perfect combination of comfort and skin protection in the shower and on the toilet
- ▲ Easy to clean because of the closed top side and the open structure at the bottom
- ▲ Equipped with individually adjustable attachment straps for optimal placement

The Vicair Hygiene consists of 4 individual compartments filled with Vicair SmartCells, making it a pressure-distributing and stable shower and toilet cushion. Because of the with SmartCells filled compartments, the cushion provides optimal positioning, redistributes pressure and minimizes shear forces. Because of the mesh material on the bottom, the cushion is, with the SmartCells inside, fully machine washable at 60°C. The cushion is easy to install because of the 4 adjustable straps at each corner. The Vicair Hygiene is suitable to use on the Invacare Ocean Ergo, but can also be used on other hygiene or shower commode chairs. For specific needs the filling of each of the 4 compartments can be easily adjusted.

PRODUCT PERFORMANCE

Skin protection / Anti decubitus

Pressure redistribution	▲▲▲▲△
Shear force reduction	▲▲▲△△
Microclimate control	▲▲▲▲△
Impact & vibration damping	▲▲▲▲▲

Positioning

Lateral stability	▲▲▲▲△
Frontal stability	▲▲▲▲△
Sliding reduction	▲▲▲▲△
Allows manual adjustments	▲▲▲▲▲
Increasing endurance	▲▲▲▲△

Maintenance & Safety

Durability	▲▲▲▲▲
Reliability	▲▲▲▲▲
Functionality when misused	▲▲▲△△
Easy to maintain	▲▲▲▲▲
Level of hygiene	▲▲▲▲▲

Functionality & Application

Comfort	▲▲▲▲▲
Ultra light weight	▲▲▲▲▲
Allows easy transfers	▲▲▲△△
Available standard sizes	▲▲△△△
Maximum user weight	▲▲▲▲▲

TECHNICAL SPECIFICATIONS

Thickness	8 cm
Width	45 cm
Depth	48
Average weight	600 grams
Max. user weight	250 kg
Machine washable	60°C max.