

Spare Parts Catalogue

Version: 2019/10/01a NO





Yes, you can.®

General Information

For any order of spare parts the following information is required:

- Product type*
- Product name & model*
- Serial Number*
- Part-Number*
- Part-Description*

Exploded view drawings and code numbers are not contractual and subject to modification, in relation to Invacare® continuous product improvement policy to meet end users needs.

The reproduction or use of this work, or of any substantial part of this work, in any form or by any electronic, mechanical or other means, including photocopying, scanning and recording, and in any information storage and retrieval system, including electronic transmission, is forbidden without the prior written permission of Invacare® International Sarl.



Yes, you can.®

TABLE OF CONTENTS

Informasjon

(User Manuals Storm4)	17
(User Manuals Ultra Low Maxx Seat)	18
(User Manuals LiNX Rehab Lite)	19
(User Manuals LiNX REM400)	20
(User Manuals LiNX REM500)	21
(User Manuals ACU/Compact Remote (LiNX))	22
(User Manuals IDC (LiNX))	23
(User Manuals ACS REM24SD)	24
(User Manuals ACS REM A & B)	25
(User Manuals ACS REM550)	26
(User Manuals ACS G90A)	27
(User Manuals ACS G91S)	28
(User Manuals REM 34 & REM 41 Tray/Centre Control)	29
(User Manuals IDC (ACS))	30
(User Manuals Compact Dual Control (ACS))	31
(User Manuals ACS Handprogrammer)	32
(User Manual Proximity Head Remote DSL-105)	33
(User Manuals Chin Control (2015-10-01 =>))	34
(User Manuals OptiCharge Charger)	35
(Change list Recaro 24V)	37

Chassis

(Storm4)	38
(Chassis Fittings)	39
(Chassis Fittings (with Ultra Low Maxx Seat))	40
(Front swingarm)	41
(Front Swing Arm Storm4 X-Plore)	42
(Motor Swingarm => 04.01.11)	43
(Motor Swingarm SSD)	44
(Motor Swingarm SSD (with Ultra Low Max Seat))	45
(Motor Swingarm SSD, Anti-Tippers flip up)	46
(Rear suspension => 03.08.12)	47
(Rear suspension 04.08.12 =>)	48
(Upper suspension, Storm4 X-Plore)	49
(Transport fixation => 03.08.12)	50
(Transport fixation 04.08.12 =>)	51
(Shrouds)	52
(Shrouds (with Ultra Low Maxx Seat))	54
(Shroud spring bracket/rear light)	55
(Storm4 TTPlus)	56
(Chassis Fittings)	39
(Front swingarm)	41
(Front Swing Arm Storm4 X-Plore)	42
(Motor Swingarm (GGB))	57
(Suspension)	58
(Upper suspension, Storm4 X-Plore)	49
(Transport Fixation)	59
(Shrouds)	60
(Shrouds Spring Bracket)	61



Yes, you can.®

TABLE OF CONTENTS

Sete

Sete	62
(Storm4)	63
(EasyAdapt Seat => 2017-08-10)	64
(Seat column / -raiser Storm4)	65
(Seat column / -raiser Storm4 X-Plore)	66
(Seat Height Adapter Storm4)	67
(Seat Height Adapter Storm4 X-Plore)	68
(Seat Tilt)	69
(Seat Frame)	70
(Seat Plate)	71
(Adapter Kit for Universal Seat)	72
(Seat Cushion)	73
(Modulite Seat)	75
(Seat column / -raiser Storm4)	65
(Seat column / -raiser Storm4 X-Plore)	66
(Seat Height Adapter Storm4)	67
(Seat Height Adapter Storm4 X-Plore)	68
(Seat Tilt)	76
(One Piece Seat Plate)	77
(Telescopic Seat Frame and Plate)	78
(Telescopic Seat Frame)	79
(Armrest bracket)	80
(Legrest Support)	81
(Seat Plate)	82
(Telescopic Seat Frame and Sling Seat)	83
(Telescopic Seat Frame)	84
(Armrest bracket)	80
(Legrest Support)	85
(Sling Seat)	86
(Telescopic Seat Frame and Universal Adapter)	87
(Telescopic Seat Frame)	79
(Armrest bracket)	80
(Legrest Support)	81
(Adapter Kit for Universal Seat)	88
(Seat Cushions "Standard" & "Matrx Flo-Tech")	89
(Seat Cushions "Matrx Libra")	92
(Seat Cushions "Matrx Vi")	93
(Ultra Low Maxx Seat)	94
(Seat Column)	95
(Seat Tilt)	96
(Seat Lifter with Tilt)	97
(Posterior-Tilt)	98
(Seat Plates & Side Rails)	99
(Seat Cushions "Matrx Vi")	100
(Seat Cushions "Matrx Flo-Tech")	101
(Seat Cushions "Matrx Libra")	102
(Recaro Seat)	103
(Seat column / -raiser Storm4)	65
(Seat column / -raiser Storm4 X-Plore)	66
(Seat Height Adapter Storm4)	67



Yes, you can.®

TABLE OF CONTENTS

(Seat Height Adapter Storm4 X-Plore)	68
(Seat Tilt)	69
(Seat Frame)	104
(Seat Plate)	105
(Recaro Seat)	106
(Change list Recaro 24V)	37
(Optimist Seat => 2018-12-31)	107
(Seat column / -raiser Storm4)	65
(Seat column / -raiser Storm4 X-Plore)	66
(Seat Height Adapter Storm4)	67
(Seat Height Adapter Storm4 X-Plore)	68
(Seat Tilt)	69
(Seat Plate)	108
(Seat Frame)	109
(Seat Cushions Optimist Seat)	110
(Storm4 TTPlus)	111
(EasyAdapt Seat => 2017-08-10)	112
(Seat column / Seatriser)	113
(Seat Height Adapter)	114
(Seat Tilt)	115
(Seat Frame)	70
(Seat Plate)	71
(Adapter Kit for Universal Seat)	72
(Seat Cushion)	73
(Recaro Seat)	116
(Seat column / Seatriser)	113
(Seat Height Adapter)	114
(Seat Tilt)	115
(Seat Frame)	104
(Seat Plate)	105
(Recaro Seat)	106
(Change list Recaro 24V)	37
(Optimist Seat => 2018-12-31)	117
(Seat column / Seatriser)	113
(Seat Height Adapter)	114
(Seat Tilt)	115
(Seat Plate)	108
(Seat Frame)	109
(Seat Cushions Optimist Seat)	110
(Modulite Seat)	118
(Seat column / Seatriser)	113
(Seat Height Adapter)	114
(Seat Tilt)	119
(One Piece Seat Plate)	77
(Telescopic Seat Frame and Plate)	78
(Telescopic Seat Frame)	79
(Armrest bracket)	80
(Legrest Support)	81
(Seat Plate)	82
(Telescopic Seat Frame and Sling Seat)	83
(Telescopic Seat Frame)	84
(Armrest bracket)	80



Yes, you can.®

TABLE OF CONTENTS

(Legrest Support)	85
(Sling Seat)	86
(Telescopic Seat Frame and Universal Adapter)	87
(Telescopic Seat Frame)	79
(Armrest bracket)	80
(Legrest Support)	81
(Adapter Kit for Universal Seat)	88
(Seat Cushions "Standard" & "Matrx Flo-Tech")	89
(Seat Cushions "Matrx Libra")	92
(Seat Cushions "Matrx Vi")	93

Rygg

Rygg	120
(Storm4)	121
(Backrest for EasyAdapt Seat => 2017-08-10)	122
(Backrestadjusting)	123
(Backrestcover for EasyAdapt)	124
(Backplate EasyAdapt)	125
(Backplate EasyAdapt for DAHL Docking Station)	126
(Backrest Cushions)	127
(Backrest Cushions Laguna & Onyx)	128
(Matrx backrest)	129
(Matrx Posture Back)	130
(Matrx Elite (Full) Back)	131
(Mounting Hardware for Matrx backrests)	132
(Push Handle EasyAdapt)	133
(Tension adjustable Sling Back)	134
("Universal"-Backrest frame)	135
(Backrest for Modulite Seat)	136
(Support for Backrest Frame)	137
(Backrest Frame)	138
(Sling Back, tension adjustable)	139
(Backrest support)	140
(Backposts)	141
(Connection tube)	142
(Back Recliner / Stabilization bracket)	143
(Backrest Cushions)	144
(Flex 3 / Laguna)	145
(Matrx Posture Back)	130
(Matrx Elite T/R Back)	147
(Matrx Elite Deep Back)	148
(Matrx Elite (Full) Back)	131
(Mounting Hardware for Matrx backrests)	132
(Flex3 back, tension adjustable)	149
(Backrest for Ultra Low Maxx Seat)	150
(Recliner)	151
(High Back)	152
(Backrest Hardware)	153
(Backrest Cushions)	154
(Matrx Backrests)	156
(Backrest Hardware)	157
(PB Back)	158



Yes, you can.®

TABLE OF CONTENTS

(Elite TR/TR HD Back)	159
(Elite Deep Back)	160
(Privacy Flaps)	161
(Backrest for Optimist Seat => 2018-12-31)	162
(Optimist Seats)	163
(Push Handle Optimist Seat)	164
(Storm4 TTPlus)	165
(Backrest for EasyAdapt Seat => 2017-08-10)	122
(Backrestadjusting)	123
(Backrestcover for EasyAdapt)	124
(Backplate EasyAdapt)	125
(Backplate EasyAdapt for DAHL Docking Station)	126
(Backrest Cushions)	127
(Backrest Cushions Laguna & Onyx)	128
(Matrx backrest)	129
(Matrx Posture Back)	130
(Matrx Elite (Full) Back)	131
(Mounting Hardware for Matrx backrests)	132
(Push Handle EasyAdapt)	133
(Tension adjustable Sling Back)	134
("Universal"-Backrest frame)	135
(Backrest for Modulite Seat)	166
(Support for Backrest Frame)	137
(Backrest Frame)	138
(Sling Back, tension adjustable)	139
(Backrest support)	140
(Backposts)	141
(Connection tube)	142
(Backrest Recliner / Stabilisation bracket)	167
(Backrest Cushions)	144
(Flex 3 / Laguna)	145
(Matrx Posture Back)	130
(Matrx Elite T/R Back)	147
(Matrx Elite Deep Back)	148
(Matrx Elite (Full) Back)	131
(Mounting Hardware for Matrx backrests)	132
(Flex3 back, tension adjustable)	149
(Backrest for Optimist Seat => 2018-12-31)	162
(Optimist Seats)	163
(Push Handle Optimist Seat)	164

Benstøtter og holdere

Benstøtter og holdere	168
(Storm4)	169
(Legrests for Ultra Low Maxx Seat)	170
(Central Legrests)	171
(Supports, Central Legrests)	172
(Central Legrest - powered LNX)	173
(Vari Legrests)	174
(Supports, Vari)	175
(Hanger Vari A (04/2009 =>))	176
(Hanger Vari F)	177



Yes, you can.®

TABLE OF CONTENTS

(Hanger Vari ADM man.)	178
(Hanger Vari ADE/ADP pow. (07/2005 =>))	179
(Footrests Vari, Flip-Up (01.08.10 =>))	180
(Footrests Vari, Flip-Up & angle adjustable (01.08.10 =>))	181
(Footrest Vari, Angle & Depth adj.)	182
(Pivot Plus Legrests)	183
(Legrest Support)	184
(Legrests)	185
(Legrest Support, Vari series 07/2005 ==>)	186
(Legrest Support, Vari series 07/2005 =>)	187
(Hanger Vari A (04/2009 =>))	176
(Hanger Vari F)	177
(Hanger Vari ADM man.)	178
(Hanger Vari ADE/ADP pow. (07/2005 =>))	179
(Central Legrest, manuell)	188
(Central Legrest, electric)	189
(Footrests Vari, Flip-Up (01.08.10 =>))	180
(Footrests Vari, Flip-Up & angle adjustable (01.08.10 =>))	181
(Support Tube with Footplate for central Legrest)	190
(Footboards)	191
(Lower Footrests)	192
(Storm4 TTPlus)	193
(Legrest Supports Vari (07/2005 =>))	194
(Hanger Vari A (04/2009 =>))	176
(Hanger Vari F)	177
(Hanger Vari ADM man.)	178
(Hanger Vari ADE/ADP pow. (07/2005 =>))	179
(Central Legrest, manuell)	188
(Central Legrest, powered)	195
(Footrests Vari, Flip-Up (01.08.10 =>))	180
(Footrests Vari, Flip-Up & angle adjustable (01.08.10 =>))	181
(Support Tube with Footplate for central Legrest)	190
(Footboards)	191
(Lower Footrests)	192

Armlen

Armlen	196
(Armrests for EasyAdapt Seat => 2017-08-10)	197
(Fixed Armrest, removable)	198
(Following armrests)	199
(Skirtguard for following armrest)	200
(Armrests & Hip Supports for Modulite Seat)	201
(T-Armrest)	202
(Padded Armrest)	203
(Inrigged Armrest)	204
(Flip-Up Armrest)	205
(Following Armrest)	206
(Hip Supports)	207
(Armrests for Ultra Low Maxx Seat)	208
(Two-post Reclining Arms, Flipup)	209
(Cantilever Arms, Flipup)	210
(Armpads, Desk & Full)	211



Yes, you can.®

TABLE OF CONTENTS

(Armpad, Ergonomic)	212
(Elbow Block (=> 2017-08-11))	213
(MACES Multi Axis Armpad Holder)	214
(Hip Supports)	207
(Armrest for Recaro Seat)	215
(Fixed Armrest, removable)	198
(Armrests for Optimist Seat => 2018-12-31)	216
(Armrests Arne Jensen for Optimist Seat)	217

Sving-og Drivhjul

(Front Fork)	218
(Front Wheels)	219

Drivhjul

(Storm4)	220
(Drive Wheel Trelleborg (8x3.00))	221
(Drive Wheel 14" (3.00 - 8))	222
(Storm4 TTPlus)	223
(Drive Wheels GGB)	224

Drivenhet

(Storm4)	225
(Propulsion => 04.01.11)	226
(Overview Propulsion)	227
(Storm4 TTPlus)	230
(GGB2 Motors)	231

Elektronikk

Elektronikk	232
(Storm4)	233
(ACS Electronic)	234
(ACS Electronic Modules)	235
(Power modules ACS 2)	236
(Holder Powermodule)	237
(Power module holder / Gyro)	238
(ACTs)	239
(Bracket ACT - Ultra Low Maxx Seat)	240
(G-Trac Sensor/Gyro)	241
(Holder Gyro)	242
(Light board)	243
(ACS Remotes & brackets)	244
(REM 24 SD with brackets)	245
(Remote REM 24SD)	246
(Remote Bracket Standard)	247
(Remote Bracket Standard - ULM and Soft Armpad on Modulite)	248
(Remote Bracket Swing away & height adjust.)	249
(Remote Bracket Swing away & height adjust. - ULM and Soft Armpad on Modulite)	250
(Remote Bracket Attendant)	251
(Remote Bracket Attendant Ultra Low Maxx Seat)	252



Yes, you can.®

TABLE OF CONTENTS

(Remote Bracket Attendant for Recaro-Seat)	253
(REM A & B with brackets)	254
(Remote Bracket Standard)	247
(Remote Bracket Swing away & height adjust.)	255
(Remote Bracket Attendant for Standard-/EasyAdapt-Seat)	256
(Remote Bracket Attendant for Recaro-Seat)	253
(REM 550 with brackets)	257
(Remote Bracket Standard)	247
(Remote Bracket Standard - ULM and Soft Armpad on Modulite)	248
(Remote Bracket Swing away & height adjust.)	249
(Remote Bracket Swing away & height adjust. - ULM and Soft Armpad on Modulite)	250
(Remote Bracket Attendant)	251
(Remote Bracket Attendant Ultra Low Maxx Seat)	252
(Remote bracket swingaway Quad Link - Ultra Low Maxx Seat)	258
(Remote Bracket Attendant for Recaro-Seat)	253
(REM 34 with brackets)	259
(Remote Bracket Standard)	260
(Remote Bracket Swing away & height adjust.)	261
(G90A with brackets)	262
(Remote Bracket Standard)	263
(Remote Bracket Standard - ULM and Soft Armpad on Modulite)	264
(Remote Bracket Swing away & height adjust.)	265
(Remote Bracket Swing away & height adjust. - ULM and soft Armpad on Modulite)	266
(Remote Bracket Attendant)	267
(Remote Bracket Attendant Ultra Low Maxx Seat)	268
(Remote Bracket Attendant for Recaro-Seat)	269
(Remote G91S & brackets)	270
(Remote Bracket Standard)	271
(Tray remote REM 41)	272
(Centre control REM 41)	273
(Chin Control, small (=> 16/12/31))	274
(Small Chincontrol)	275
(Framework for chincontrol)	276
(Chin Control - Modulite Seat (15/10/01 =>))	277
(User Manuals Chin Control (2015-10-01 =>))	34
(Chin Control DX-RJM)	278
(Chin Control powered - Details)	279
(Chin Control manual - Details)	280
(Mounting Hardware - Chin Control)	281
(Dual Control with brackets)	282
(Remote Bracket Dual Control for Standard-/EasyAdapt-Seat)	283
(Remote Bracket Dual Control for Modulite Seat)	284
(Remote Bracket Dual Control for Recaro Seat)	285
(Compact Dual Control with brackets)	286
(Compact Dual Control)	287
(Bracket, Compact Dual Control - Ultra Low Maxx Seat)	288
(Remote Bracket, Compact Dual Control/ACU - short)	289
(Remote Bracket, Compact Dual Control/ACU - long)	290
(Intuitive Dual Control (IDC) with brackets)	291
(IDC with bracket)	292
(User Manuals IDC (ACS))	30



Yes, you can.®

TABLE OF CONTENTS

(IDC with bracket (Ultra Low Maxx Seat))	293
(User Manuals IDC (ACS))	30
(Height adjustable Bracket, IDC)	294
(ACS Cables & Accessories)	295
(Bus-Cables ACS)	296
(T-Adapter)	297
(T-Adapter (Ultra Low Maxx Seat))	298
(Backrest Sensor - Ultra Low Maxx Seat)	299
(Tilt Sensor - Ultra Low Maxx Seat)	300
(Programming Kits)	301
(Magnetkey for Remotes)	302
(ACS Components)	303
(iPortal Modules and Accessories)	304
(Remote, Joystick, Control)	305
(Joystick Knobs)	306
(Buddy Button)	307
(Piko Button)	308
(Interface, Adapter, Other)	309
(LiNX Electronic)	310
(Modules)	311
(Power Modules)	312
(Bracket, Power Module)	313
(Actuator Modules (ACTs))	314
(G-Trac Sensor / Gyro)	315
(Remotes & Brackets)	316
(Remotes & Brackets)	317
(Remote REM110)	318
(Remote REM211)	319
(Remote REM216)	320
(Remote REM400)	321
(Toggle Switches for REM400)	322
(Remote REM500)	323
(Remote Bracket, Standard)	324
(Remote Bracket, Swing-away)	325
(Remote Bracket, Attendant)	326
(Remote Bracket, Nucleus)	327
(Remote Bracket, Swing-away REM500)	328
(Remote Bracket, Nucleus REM500)	329
(Dual Controls & Brackets)	330
(Dual Control REM050)	331
(Dual Control ACU200)	332
(Intuitive Dual Control (IDC))	333
(User Manuals IDC (LiNX))	23
(Remote Bracket, Attendant)	326
(Remote Bracket, Compact Dual Control/ACU - short)	289
(Remote Bracket, Compact Dual Control/ACU - long)	290
(Mounting Bracket, IDC)	334
(Height adjustable Bracket, IDC)	294
(Additional Remotes & Brackets)	335
(Compact Remote LiNX & Bracket)	336
(Compact Remotes LiNX)	337
(Bracket, Standard)	338



Yes, you can.®

TABLE OF CONTENTS

(Bracket, Nucleus)	339
(Chin Control, man.)	340
(Egg Switch ASL300)	341
(Chin Control, pow. - Modulite Seat)	342
(Chin Control, powered)	343
(Mounting Hardware, Chin Control (Modulite Seat))	344
(Head Controls (ASL104))	345
(ATOM Mini ASL104MLX Head Array LiNX)	346
(ATOM Proton ASL104PTLX Head Array LiNX)	347
(Egg Switch ASL300)	348
(Egg Switch ASL300)	341
(Wireless Accessories (ASL557/558))	349
(Micro Extremity Control (ASL130))	350
(Bracket, MEC - Ring Clamp Mount, Nucleus)	351
(Bracket, MEC - Tray Mount, Nucleus)	352
(Bracket, MEC - Tray Mount, Standard)	353
(Cable Mounting Sheets ASL - Modulite Back)	354
(Cables & Accessories)	355
(Bus Cables)	356
(Actuator Cable)	357
(Lifter/Tilt-Module Cables)	358
(Light Cable)	359
(Light Board)	360
(4-way Expansion Block)	361
(Slow-Stop-Sensor - Modulite Seat)	362
(Slow-Stop-Sensor - Recaro Seat)	363
(Accessory cables)	364
(Stereo-Mono-Splitter Cable)	365
(Programming Access Key)	366
(Components)	367
(Accessory cables)	364

(10-way switch & brackets)

(10-way switch & brackets)	368
(10-way switch)	369
(Bracket, front)	370
(Bracket, opposite side)	371
(Bracket, rear)	372
(USB Charger Port)	373
(Buddy Button)	307
(Egg Switch ASL300)	341
(Piko Button)	308
(Bodypoint Joystick Knobs)	374
(Storm4 TTPlus)	375
(ACS Modules)	376
(Powermodule GGB & Bracket)	377
(CLAM & Bracket)	378
(Relais & bracket)	379
(Lightboard GGB)	380
(ACS Remotes & Brackets)	381
(REM 24 SD with brackets)	382
(Remote REM 24SD)	246



Yes, you can.®

TABLE OF CONTENTS

(Remote Bracket Standard)	247
(Remote Bracket Swing away & height adjust.)	249
(Remote Bracket Attendant for Standard-/EasyAdapt-Seat)	256
(Remote Bracket Attendant for Recaro-Seat)	253
(REM A & B with brackets)	254
(Remote Bracket Standard)	247
(Remote Bracket Swing away & height adjust.)	255
(Remote Bracket Attendant for Standard-/EasyAdapt-Seat)	256
(Remote Bracket Attendant for Recaro-Seat)	253
(REM 550 with brackets)	383
(Remote Bracket Standard)	247
(Remote Bracket Swing away & height adjust.)	249
(Remote Bracket Attendant for Standard-/EasyAdapt-Seat)	256
(Remote Bracket Attendant for Recaro-Seat)	253
(REM 34 with brackets)	259
(Remote Bracket Standard)	260
(Remote Bracket Swing away & height adjust.)	261
(G90A with brackets)	384
(Remote Bracket Standard)	263
(Remote Bracket Swing away & height adjust.)	265
(Remote Bracket Attendant for Standard-/EasyAdapt-Seat)	385
(Remote Bracket Attendant for Recaro-Seat)	269
(Remote G91S & brackets)	270
(Remote Bracket Standard)	271
(Centre control REM 41)	273
(Tray remote REM 41)	272
(Chin Control, small (=> 16/12/31))	274
(Small Chincontrol)	275
(Framework for chincontrol)	276
(Chin Control - Modulite Seat (15/10/01 =>))	277
(User Manuals Chin Control (2015-10-01 =>))	34
(Chin Control DX-RJM)	278
(Chin Control powered - Details)	279
(Chin Control manual - Details)	280
(Mounting Hardware - Chin Control)	281
(Dual Control with brackets)	282
(Remote Bracket Dual Control for Standard-/EasyAdapt-Seat)	283
(Remote Bracket Dual Control for Modulite Seat)	284
(Remote Bracket Dual Control for Recaro Seat)	285
(Compact Dual Control with brackets)	386
(Compact Dual Control)	287
(Remote Bracket, Compact Dual Control/ACU - short)	289
(Remote Bracket, Compact Dual Control/ACU - long)	290
(Intuitive Dual Control (IDC) with bracket)	387
(IDC with bracket)	292
(User Manuals IDC (ACS))	30
(Height adjustable Bracket, IDC)	294
(ACS Cables & Accessories)	388
(Main cable harness)	389
(Plug Sheets el. legrest GGB)	390
(T-Adapter)	391
(T-Adapter Modulite Seat)	392



Yes, you can.®

TABLE OF CONTENTS

(Programming Kits)	301
(Magnetkey for Remotes)	302
(ACS Components)	303
(iPortal Modules and Accessories)	304
(Remote, Joystick, Control)	305
(Joystick Knobs)	306
(Buddy Button)	307
(Piko Button)	308
(Interface, Adapter, Other)	309

kabel og lader

(Storm4)	393
(Batteries & Accessories)	394
(Chargers)	395
(Chargers)	396
(Storm4 TTPlus)	397
(Batteries & Accessories)	398
(Chargers)	395
(Chargers)	396

Sikkerhedsudstyr

(Frontlight (German homologation))	399
(Frontlight - Seat Rail mounted (German homologation))	400
(Rearlight (German homologation))	401
(Rear Light Protector Kit)	402
(LED-Frontlight)	403
(LED-Frontlight - Seat Rail mounted)	404
(LED-Rearlight)	405
(Attachment Kit - Optimist Seat => 2018-12-31)	406
(Reflective Sign)	407
(Reflective Sign - ULM Seat)	408
(Kerb Rider)	409
(Lap Belts / H-Belt)	410
(Lap Belts - Modulite Seat)	411
(H-Belts - Modulite Seat)	412
(Bodypoint Pelvic Belts)	413
(Lap Belt - ULM Seat)	414
(Lap Belt - Max Seat)	415
(Rearview Mirror)	416
(Anti-Tippers)	417
(Storm4)	418
(Anti-Tipper Wheels => 04.01.11)	419
(Anti-Tipper for SSD Motors)	420
(Anti-Tippers for SSD Motors, flip up)	421
(Storm4 TTP)	422
(Anti-Tippers for GGB2 Motors)	423

Sete utstyr

(Headrests)	424
(Headrest with cheek support)	425
(Headrest & Bracket)	426



Yes, you can.®

TABLE OF CONTENTS

(Headrestholder EasyAdapt Seat)	427
(Headrest, sidewise adjustable)	428
(Elan Headrests - Modulite Seat)	429
(Headrests - Ultra Low Maxx Seat)	430
(Lateral Support)	431
(Lateral Supports - Modulite Seat)	432
(Lateral supports, Flex3 backrest)	433
(Lateral supports, Flex3 backrest - multi adjustable)	434
(Lateral supports, Flex3 backrest (tension adjustable) - multi adjustable)	435
(Lateral Support - Ultra Low Maxx Seat)	436
(Lateral Support - EasyAdapt Seat => 2017-08-10)	437
(Lateral Support (Dartex))	438
(Lateral Support, multi adjustable with flexible pad)	439
(Bag => 2017-08-10)	440

Annet Utstyr & Tilbehør

(Bag holder - Modulite Seat)	441
(Luggage Rack)	442
(Caneholder)	443
(Cane holder (Easy Adapt Seat & Std-Seat) =>2017-08-10)	444
(Cane holder (Recaro Seat))	445
(Cane holder (Modulite: backrest frame))	446
(Cane holder (Modulite: back posts man recliner))	447
(Cane holder (Modulite: back posts pow recliner))	448
(Fender for Castors)	449
(Fenders for Drive Wheels => 04.01.11)	450
(Fenders for Drive Wheels (SSD-Motors))	451
(Fenders (GGB Motors))	452
(Trays)	453
(Trays)	454
(Trays (Ultra Low Maxx Seat))	455
(km-Counter & Bracket)	456
(Mounting-kit km-Counter)	457
(Hour Counter)	458
(Bag)	459
(Toolkit)	460
(T-Nuts for Accessories Modulite & ULM Seat)	461

BG	Български	HE	עברית	PL	Polski
CS	Čeština	HR	Hrvatski	PT	Português
DA	Dansk	HU	Magyar	RU	Русский
DE	Deutsch	IT	Italiano	RO	Română
EN	English	JA	日本語	SA	العربية
ES	Español	LV	Latviešu	SK	Slovenčina
ET	Eesti	LT	Lietuvių	SL	Slovenščina
FI	Suomi	NL	Nederlands	SV	Svenska
FR	Français	NO	Norsk	UK	Українська

languages.iso

CS	1532570	(User Manual Storm4 Czech)	1
DA	1529704	(User Manual Storm4 Danish)	1
DE	1529687	(User Manual Storm4 German)	1
EN	1529686	(User Manual Storm4 English)	1
ES	1529690	(User Manual Storm4 Spanish)	1
ET	SP1543013	(User Manual Storm4 Estonian)	1
FI	1529702	(User Manual Storm4 Finsk)	1
FR	1529688	(User Manual Storm4 French)	1
HR	SP1650542	(User Manual Storm4 Croatian)	1
HU	SP1574045	(User Manual Storm4 Hungarian)	1
IT	1529689	(User Manual Storm4 Italian)	1
LT	1532635	(User Manual Storm4 Lithuanian)	1
NL	1529700	(User Manual Storm4 Dutch)	1
NO	1529703	BRUKSANVISNING STORM4 MANUAL NORSK	1
PL	1532571	(User Manual Storm4 Polish)	1
RU	SP1549066	(User Manual Storm4 Russian)	1
SA	SP1654738	(User Manual Storm4 Arabic)	1
SL	1532634	(User Manual Storm4 Slovenian)	1
SV	1529701	(User Manual Storm4 Swedish)	1

(User Manuals Ultra Low Maxx Seat)

BG	Български	HE	עברית	PL	Polski
CS	Čeština	HR	Hrvatski	PT	Português
DA	Dansk	HU	Magyar	RU	Русский
DE	Deutsch	IT	Italiano	RO	Română
EN	English	JA	日本語	SA	العربية
ES	Español	LV	Latviešu	SK	Slovenčina
ET	Eesti	LT	Lietuvių	SL	Slovenščina
FI	Suomi	NL	Nederlands	SV	Svenska
FR	Français	NO	Norsk	UK	Українська

languages.iso

DA	SP1585726	(User Manual Ultra Low Maxx Seat - Danish)	1
DE	SP1585727	(User Manual Ultra Low Maxx Seat - German)	1
EN	SP1585725	(User Manual Ultra Low Maxx Seat - English)	1
ES	SP1601887	(User Manual Ultra Low Maxx Seat - Spanish)	1
FI	SP1585728	(User Manual Ultra Low Maxx Seat - Finish)	1
FR	SP1585729	(User Manual Ultra Low Maxx Seat - French)	1
IT	SP1588266	(User Manual Ultra Low Maxx Seat - Italian)	1
NL	SP1585730	(User Manual Ultra Low Maxx Seat - Dutch)	1
NO	SP1585731	(User Manual Ultra Low Maxx Seat - Norwegian)	1
PT	SP1585732	(User Manual Ultra Low Maxx Seat - Portuguese)	1
SV	SP1585733	(User Manual Ultra Low Maxx Seat - Swedish)	1

(User Manuals LiNX Rehab Lite)

BG	Български	HE	עברית	PL	Polski
CS	Čeština	HR	Hrvatski	PT	Português
DA	Dansk	HU	Magyar	RU	Русский
DE	Deutsch	IT	Italiano	RO	Română
EN	English	JA	日本語	SA	العربية
ES	Español	LV	Latviešu	SK	Slovenčina
ET	Eesti	LT	Lietuvių	SL	Slovenščina
FI	Suomi	NL	Nederlands	SV	Svenska
FR	Français	NO	Norsk	UK	Українська

languages.iso

BG	SP1603635	(User Manual LiNX - Bulgarian)	1
CS	SP1603637	(User Manual LiNX - Czech)	1
DA	SP1603638	(User Manual LiNX - Danish)	1
DE	SP1603644	(User Manual LiNX - German)	1
EN	SP1602969	(User Manual LiNX - English)	1
ES	SP1603711	(User Manual LiNX - Spanish)	1
ET	SP1603641	(User Manual LiNX - Estonian)	1
FI	SP1603642	(User Manual LiNX - Finnish)	1
FR	SP1603643	(User Manual LiNX - French)	1
HR	SP1603636	(User Manual LiNX - Croatian)	1
HU	SP1603700	(User Manual LiNX - Hungarian)	1
IT	SP1603701	(User Manual LiNX - Italian)	1
LV	SP1603702	(User Manual LiNX - Latvian)	1
LT	SP1603703	(User Manual LiNX - Lithuanian)	1
NL	SP1603639	(User Manual LiNX - Dutch)	1
NO	SP1603704	(User Manual LiNX - Norwegian)	1
PL	SP1603705	(User Manual LiNX - Polish)	1
PT	SP1603706	(User Manual LiNX - Portuguese)	1
RU	SP1603707	(User Manual LiNX - Russian)	1
SK	SP1603799	(User Manual LiNX - Slovak)	1
SL	SP1603708	(User Manual LiNX - Slovenian)	1
SV	SP1603712	(User Manual LiNX - Swedisch)	1

(User Manuals LiNX REM400)

BG	Български	HE	עברית	PL	Polski
CS	Čeština	HR	Hrvatski	PT	Português
DA	Dansk	HU	Magyar	RU	Русский
DE	Deutsch	IT	Italiano	RO	Română
EN	English	JA	日本語	SA	العربية
ES	Español	LV	Latviešu	SK	Slovenčina
ET	Eesti	LT	Lietuvių	SL	Slovenščina
FI	Suomi	NL	Nederlands	SV	Svenska
FR	Français	NO	Norsk	UK	Українська

languages.iso

EN	SP1637423	(User Manual LiNX REM400 - English)	1
DA	SP1640738	(User Manual LiNX REM400 - Danish)	1
NL	SP1640739	(User Manual LiNX REM400 - Dutch)	1
FI	SP1640743	(User Manual LiNX REM400 - Finnish)	1
FR	SP1640744	(User Manual LiNX REM400 - French)	1
DE	SP1640745	(User Manual LiNX REM400 - German)	1
IT	SP1640746	(User Manual LiNX REM400 - Italian)	1
NO	SP1640747	(User Manual LiNX REM400 - Norwegian)	1
PT	SP1640748	(User Manual LiNX REM400 - Portuguese)	1
ES	SP1640749	(User Manual LiNX REM400 - Spanish)	1
SV	SP1640750	(User Manual LiNX REM400 - Swedish)	1

BG	Български	HE	עברית	PL	Polski
CS	Čeština	HR	Hrvatski	PT	Português
DA	Dansk	HU	Magyar	RU	Русский
DE	Deutsch	IT	Italiano	RO	Română
EN	English	JA	日本語	SA	العربية
ES	Español	LV	Latviešu	SK	Slovenčina
ET	Eesti	LT	Lietuvių	SL	Slovenščina
FI	Suomi	NL	Nederlands	SV	Svenska
FR	Français	NO	Norsk	UK	Українська

languages.iso

EN	SP1637425	(User Manual LiNX REM500 - EN)	1
DA	SP1640751	(User Manual LiNX REM500 - Danish)	1
NL	SP1640752	(User Manual LiNX REM500 - Dutch)	1
FI	SP1640753	(User Manual LiNX REM500 - Finnish)	1
FR	SP1640754	(User Manual LiNX REM500 - French)	1
DE	SP1640755	(User Manual LiNX REM500 - German)	1
IT	SP1640764	(User Manual LiNX REM500 - Italian)	1
NO	SP1640765	(User Manual LiNX REM500 - Norwegian)	1
PT	SP1640766	(User Manual LiNX REM500 - Portuguese)	1
ES	SP1640767	(User Manual LiNX REM500 - Spanish)	1
SV	SP1640768	(User Manual LiNX REM500 - Swedish)	1

(User Manuals ACU/Compact Remote (LiNX))

BG	Български	HE	עברית	PL	Polski
CS	Čeština	HR	Hrvatski	PT	Português
DA	Dansk	HU	Magyar	RU	Русский
DE	Deutsch	IT	Italiano	RO	Română
EN	English	JA	日本語	SA	العربية
ES	Español	LV	Latviešu	SK	Slovenčina
ET	Eesti	LT	Lietuvių	SL	Slovenščina
FI	Suomi	NL	Nederlands	SV	Svenska
FR	Français	NO	Norsk	UK	Українська

languages.iso

EN	SP1608055	(User Manual LiNX ACU/Compact Remote - English)	1
PL	SP1638981	(User Manual LiNX ACU/Compact Remote - Polish)	1
ET	SP1638982	(User Manual LiNX ACU/Compact Remote - Estonian)	1
FR	SP1638983	(User Manual LiNX ACU/Compact Remote - French)	1
RU	SP1638985	(User Manual LiNX ACU/Compact Remote - Russian)	1
SV	SP1638986	(User Manual LiNX ACU/Compact Remote - Swedish)	1
SK	SP1638987	(User Manual LiNX ACU/Compact Remote - Slovak)	1
SL	SP1638988	(User Manual LiNX ACU/Compact Remote - Slovenian)	1
NL	SP1638989	(User Manual LiNX ACU/Compact Remote - Dutch)	1
FI	SP1638990	(User Manual LiNX ACU/Compact Remote - Finnish)	1
IT	SP1638991	(User Manual LiNX ACU/Compact Remote - Italian)	1
DE	SP1638992	(User Manual LiNX ACU/Compact Remote - German)	1
HU	SP1639001	(User Manual LiNX ACU/Compact Remote - Hungarian)	1
BG	SP1639002	(User Manual LiNX ACU/Compact Remote - Bulgarian)	1
LV	SP1639003	(User Manual LiNX ACU/Compact Remote - Latvian)	1
LT	SP1639004	(User Manual LiNX ACU/Compact Remote - Lithuanian)	1
NO	SP1639005	(User Manual LiNX ACU/Compact Remote - Norwegian)	1
CS	SP1639006	(User Manual LiNX ACU/Compact Remote - Czech)	1
ES	SP1639007	(User Manual LiNX ACU/Compact Remote - Spanish)	1
DA	SP1639008	(User Manual LiNX ACU/Compact Remote - Danish)	1
HR	SP1639009	(User Manual LiNX ACU/Compact Remote - Croatian)	1
PT	SP1639019	(User Manual LiNX ACU/Compact Remote - Portuguese)	1

BG	Български	HE	עברית	PL	Polski
CS	Čeština	HR	Hrvatski	PT	Português
DA	Dansk	HU	Magyar	RU	Русский
DE	Deutsch	IT	Italiano	RO	Română
EN	English	JA	日本語	SA	العربية
ES	Español	LV	Latviešu	SK	Slovenčina
ET	Eesti	LT	Lietuvių	SL	Slovenščina
FI	Suomi	NL	Nederlands	SV	Svenska
FR	Français	NO	Norsk	UK	Українська

languages.iso

DA	SP1644921	(User Manual IDC LiNX - Danish)	1
DE	SP1644924	(User Manual IDC LiNX - German)	1
EN	SP1643411	(User Manual IDC LiNX - English)	1
ES	SP1644955	(User Manual IDC LiNX - Spanish)	1
FI	SP1644922	(User Manual IDC LiNX - Finnish)	1
FR	SP1644923	(User Manual IDC LiNX - French)	1
IT	SP1644925	(User Manual IDC LiNX - Italian)	1
NL	SP1644868	(User Manual IDC LiNX - Dutch)	1
NO	SP1644926	(User Manual IDC LiNX - Norwegian)	1
PT	SP1644927	(User Manual IDC LiNX - Portuguese)	1
SV	SP1644956	(User Manual IDC LiNX - Swedish)	1

(User Manuals ACS REM24SD)

BG	Български	HE	עברית	PL	Polski
CS	Čeština	HR	Hrvatski	PT	Português
DA	Dansk	HU	Magyar	RU	Русский
DE	Deutsch	IT	Italiano	RO	Română
EN	English	JA	日本語	SA	العربية
ES	Español	LV	Latviešu	SK	Slovenčina
ET	Eesti	LT	Lietuvių	SL	Slovenščina
FI	Suomi	NL	Nederlands	SV	Svenska
FR	Français	NO	Norsk	UK	Українська

languages.iso

CS	SP1533252	(User Manual REM24SD Czech)	1
DA	1436146	(User Manual REM24SD Danish)	1
DE	1436142	(User Manual REM24SD German)	1
EN	1436144	(User Manual REM24SD English)	1
ES	1436149	(User Manual REM24SD Spanish)	1
ET	SP1533253	(User Manual REM24SD Estonian)	1
FI	1436386	(User Manual REM24SD Finnish)	1
FR	1436143	(User Manual REM24SD French)	1
HR	SP1533251	(User Manual REM24SD Croatian)	1
HU	SP1533254	(User Manual REM24SD Hungarian)	1
IT	1436385	(User Manual REM24SD Italian)	1
LV	SP1543152	(User Manual REM24SD Latvian)	1
LT	SP1533255	(User Manual REM24SD Lithuanian)	1
NL	1436148	(User Manual REM24SD Dutch)	1
NO	1436147	(User Manual REM24SD Norwegian)	1
PL	SP1533256	(User Manual REM24SD Polish)	1
PT	1436150	(User Manual REM24SD Portuguese)	1
SL	SP1533257	(User Manual REM24SD Slovenian)	1
SV	1436145	(User Manual REM24SD Swedish)	1

(User Manuals ACS REM A & B)

BG	Български	HE	עברית	PL	Polski
CS	Čeština	HR	Hrvatski	PT	Português
DA	Dansk	HU	Magyar	RU	Русский
DE	Deutsch	IT	Italiano	RO	Română
EN	English	JA	日本語	SA	العربية
ES	Español	LV	Latviešu	SK	Slovenčina
ET	Eesti	LT	Lietuvių	SL	Slovenščina
FI	Suomi	NL	Nederlands	SV	Svenska
FR	Français	NO	Norsk	UK	Українська

languages.iso

IT	1529230	(User Manual REM A & B - Italian)	1
NO	1529232	(User Manual REM A & B - Norwegian)	1
ES	1529174	(User Manual REM A & B - Spanish)	1
FI	1529228	(User Manual REM A & B - Finnish)	1
EN	1529173	(User Manual REM A & B - English)	1
DA	1529171	(User Manual REM A & B - Danish)	1
NL	1529231	(User Manual REM A & B - Dutch)	1
DE	1529172	(User Manual REM A & B - German)	1
PT	1529233	(User Manual REM A & B - Portuguese)	1
SV	1529234	(User Manual REM A & B Swedish)	1
FR	1529229	(User Manual REM A & B - French)	1
CS	SP1533259	(User Manual REM A & B - Czech)	1
ET	SP1533265	(User Manual REM A & B - Estonian)	1
HR	SP1533258	(User Manual REM A & B - Croatian)	1
HU	SP1533266	(User Manual REM A & B - Hungarian)	1
LT	SP1533267	(User Manual REM A & B - Lithuanian)	1
PL	SP1533268	(User Manual REM A & B - Polish)	1
RO	SP1533270	(User Manual REM A & B - Romanian)	1
SL	SP1533269	(User Manual REM A & B - Slovenian)	1

(User Manuals ACS REM550)

BG	Български	HE	עברית	PL	Polski
CS	Čeština	HR	Hrvatski	PT	Português
DA	Dansk	HU	Magyar	RU	Русский
DE	Deutsch	IT	Italiano	RO	Română
EN	English	JA	日本語	SA	العربية
ES	Español	LV	Latviešu	SK	Slovenčina
ET	Eesti	LT	Lietuvių	SL	Slovenščina
FI	Suomi	NL	Nederlands	SV	Svenska
FR	Français	NO	Norsk	UK	Українська

languages.iso

IT	1528094	(User Manual REM550 Italiano)	1
NO	1528096	BRUKSANVISNING DX REM550 MANUAL NORSK	1
ES	1528091	(User Manual REM550 Spanish)	1
FI	1528092	(User Manual REM550 Finsk)	1
EN	1528090	BRUKSANVISNING DX-REM550 MANUAL	1
DA	1528088	(User Manual REM550 Danish)	1
NL	1528095	(User Manual REM550 Dutch)	1
DE	1528089	(User Manual REM550 German)	1
PT	1528097	(User Manual REM550 Portuguese)	1
SV	1528098	(User Manual REM550 Swedish)	1
FR	1528093	(User Manual REM550 French)	1
CS	SP1533280	(User Manual REM550 Czech)	1
ET	SP1533281	(User Manual REM550 Estonian)	1
HR	SP1533279	(User Manual REM550 Croatian)	1
HU	SP1533282	(User Manual REM550 Hungarian)	1
LT	SP1533283	(User Manual REM550 Lithuanian)	1
PL	SP1533284	(User Manual REM550 Polish)	1
SL	SP1533285	(User Manual REM550 Slovenian)	1

BG	Български	HE	עברית	PL	Polski
CS	Čeština	HR	Hrvatski	PT	Português
DA	Dansk	HU	Magyar	RU	Русский
DE	Deutsch	IT	Italiano	RO	Română
EN	English	JA	日本語	SA	العربية
ES	Español	LV	Latviešu	SK	Slovenčina
ET	Eesti	LT	Lietuvių	SL	Slovenščina
FI	Suomi	NL	Nederlands	SV	Svenska
FR	Français	NO	Norsk	UK	Українська

languages.iso

EN	SP1524441	(User Manual G90A - English)	1
DE	SP1536742	(User Manual G90A - German)	1
ES	SP1638870	(User Manual G90A - Spanish)	1
FR	SP1638871	(User Manual G90A - French)	1
NO	SP1638872	(User Manual G90A - Norwegian)	1
PT	SP1638873	(User Manual G90A - Portuguese)	1
SV	SP1638874	(User Manual G90A - Swedish)	1
DA	SP1638875	(User Manual G90A - Danish)	1
NL	SP1638876	(User Manual G90A - Dutch)	1

(User Manuals ACS G91S)

BG	Български	HE	עברית	PL	Polski
CS	Čeština	HR	Hrvatski	PT	Português
DA	Dansk	HU	Magyar	RU	Русский
DE	Deutsch	IT	Italiano	RO	Română
EN	English	JA	日本語	SA	العربية
ES	Español	LV	Latviešu	SK	Slovenčina
ET	Eesti	LT	Lietuvių	SL	Slovenščina
FI	Suomi	NL	Nederlands	SV	Svenska
FR	Français	NO	Norsk	UK	Українська

languages.iso

(User Manuals REM 34 & REM 41 Tray/Centre Control)

BG	Български	HE	עברית	PL	Polski
CS	Čeština	HR	Hrvatski	PT	Português
DA	Dansk	HU	Magyar	RU	Русский
DE	Deutsch	IT	Italiano	RO	Română
EN	English	JA	日本語	SA	العربية
ES	Español	LV	Latviešu	SK	Slovenčina
ET	Eesti	LT	Lietuvių	SL	Slovenščina
FI	Suomi	NL	Nederlands	SV	Svenska
FR	Français	NO	Norsk	UK	Українська

languages.iso

CS	SP1538031	(User Manual REM41 Tray Control Czech)	1
DA	SP1538022	(User Manual REM41 Tray Control Danish)	1
DE	SP1541402	(User Manual REM41 Tray Control German)	1
EN	SP1538023	(User Manual REM41 Tray Control English)	1
ES	SP1538024	(User Manual REM41 Tray Control Spanish)	1
FR	SP1538025	(User Manual REM41 Tray Control French)	1
IT	SP1538026	(User Manual REM41 Tray Control Italian)	1
NL	SP1538027	(User Manual REM41 Tray Control Dutch)	1
NO	SP1538028	(User Manual REM41 Tray Control Norwegian)	1
PT	SP1538029	(User Manual REM41 Tray Control Portuguese)	1
SL	SP1538032	(User Manual REM41 Tray Control Slovenian)	1
SV	SP1538030	(User Manual REM41 Tray Control Swedish)	1

(User Manuals IDC (ACS))

BG	Български	HE	עברית	PL	Polski
CS	Čeština	HR	Hrvatski	PT	Português
DA	Dansk	HU	Magyar	RU	Русский
DE	Deutsch	IT	Italiano	RO	Română
EN	English	JA	日本語	SA	العربية
ES	Español	LV	Latviešu	SK	Slovenčina
ET	Eesti	LT	Lietuvių	SL	Slovenščina
FI	Suomi	NL	Nederlands	SV	Svenska
FR	Français	NO	Norsk	UK	Українська

languages.iso

DE	1521480	(User Manual IDC - German)	1
EN	1521481	(User Manual IDC - English)	1
FR	1521482	(User Manual IDC - French)	1
NL	1521483	(User Manual IDC - Dutch)	1
NO	1521484	BRUKSANVISING INTUITIV LEDSAGERSTYRING	1
SV	1521485	(User Manual IDC - Swedish)	1
DA	1521486	(User Manual IDC - Danish)	1
IT	1522252	(User Manual IDC - Italian)	1
PT	1522254	(User Manual IDC - Portuguese)	1
ES	1522255	(User Manual IDC - Spanish)	1

(User Manuals Compact Dual Control (ACS))

BG	Български	HE	עברית	PL	Polski
CS	Čeština	HR	Hrvatski	PT	Português
DA	Dansk	HU	Magyar	RU	Русский
DE	Deutsch	IT	Italiano	RO	Română
EN	English	JA	日本語	SA	العربية
ES	Español	LV	Latviešu	SK	Slovenčina
ET	Eesti	LT	Lietuvių	SL	Slovenščina
FI	Suomi	NL	Nederlands	SV	Svenska
FR	Français	NO	Norsk	UK	Українська

languages.iso

BG	SP1589186	(User Manual Compact Dual Control - Bulgarian)	1
CS	SP1589313	(User Manual Compact Dual Control - Czech)	1
DA	SP1588488	(User Manual Compact Dual Control - Danish)	1
DE	SP1588489	(User Manual Compact Dual Control - German)	1
EN	SP1588462	(User Manual Compact Dual Control - English)	1
ES	SP1588496	(User Manual Compact Dual Control - Spanish)	1
ET	SP1589311	(User Manual Compact Dual Control - Estonian)	1
FI	SP1588490	(User Manual Compact Dual Control - Finnish)	1
FR	SP1588491	(User Manual Compact Dual Control - French)	1
HR	SP1589314	(User Manual Compact Dual Control - Croatian)	1
HU	SP1589315	(User Manual Compact Dual Control - Hungarian)	1
IT	SP1588492	(User Manual Compact Dual Control - Italian)	1
LV	SP1589317	(User Manual Compact Dual Control - Latvian)	1
LT	SP1589316	(User Manual Compact Dual Control - Lithuanian)	1
NL	SP1588493	(User Manual Compact Dual Control - Dutch)	1
NO	SP1588494	(User Manual Compact Dual Control - Norwegian)	1
PL	SP1589318	(User Manual Compact Dual Control - Polish)	1
PT	SP1588495	(User Manual Compact Dual Control - Portuguese)	1
RU	SP1589320	(User Manual Compact Dual Control - Russian)	1
RO	SP1589319	(User Manual Compact Dual Control - Romanian)	1
SK	SP1589321	(User Manual Compact Dual Control - Slovak)	1
SL	SP1589322	(User Manual Compact Dual Control - Slovenian)	1
SV	SP1588497	(User Manual Compact Dual Control - Swedish)	1
UK	SP1589323	(User Manual Compact Dual Control - Ukrainian)	1

(User Manuals ACS Handprogrammer)

BG	Български	HE	עברית	PL	Polski
CS	Čeština	HR	Hrvatski	PT	Português
DA	Dansk	HU	Magyar	RU	Русский
DE	Deutsch	IT	Italiano	RO	Română
EN	English	JA	日本語	SA	العربية
ES	Español	LV	Latviešu	SK	Slovenčina
ET	Eesti	LT	Lietuvių	SL	Slovenščina
FI	Suomi	NL	Nederlands	SV	Svenska
FR	Français	NO	Norsk	UK	Українська

languages.iso

DE	1419316	(User Manual Handheld Programmer German)	1
EN	1419317	(User Manual Handheld Programmer English)	1

(User Manual Proximity Head Remote DSL-105)

BG	Български	HE	עברית	PL	Polski
CS	Čeština	HR	Hrvatski	PT	Português
DA	Dansk	HU	Magyar	RU	Русский
DE	Deutsch	IT	Italiano	RO	Română
EN	English	JA	日本語	SA	العربية
ES	Español	LV	Latviešu	SK	Slovenčina
ET	Eesti	LT	Lietuvių	SL	Slovenščina
FI	Suomi	NL	Nederlands	SV	Svenska
FR	Français	NO	Norsk	UK	Українська

languages.iso

EN	SP1555731	(User Manual Proximity Head Remote DSL-105 English)	1
DE	SP1555732	(User Manual Proximity Head Remote DSL-105 German)	1
FR	SP1555733	(User Manual Proximity Head Remote DSL-105 French)	1
DA	SP1555839	(User Manual Proximity Head Remote DSL-105 Danish)	1
ES	SP1555840	(User Manual Proximity Head Remote DSL-105 Spanish)	1
FI	SP1555841	(User Manual Proximity Head Remote DSL-105 Finnish)	1
IT	SP1555842	(User Manual Proximity Head Remote DSL-105 Italian)	1
NL	SP1555843	(User Manual Proximity Head Remote DSL-105 Dutch)	1
NO	SP1555844	(User Manual Proximity Head Remote DSL-105 Norwegian)	1
PT	SP1555845	(User Manual Proximity Head Remote DSL-105 Portuguese)	1
SV	SP1555846	(User Manual Proximity Head Remote DSL-105 Swedish)	1

(User Manuals Chin Control (2015-10-01 =>))

BG	Български	HE	עברית	PL	Polski
CS	Čeština	HR	Hrvatski	PT	Português
DA	Dansk	HU	Magyar	RU	Русский
DE	Deutsch	IT	Italiano	RO	Română
EN	English	JA	日本語	SA	العربية
ES	Español	LV	Latviešu	SK	Slovenčina
ET	Eesti	LT	Lietuvių	SL	Slovenščina
FI	Suomi	NL	Nederlands	SV	Svenska
FR	Français	NO	Norsk	UK	Українська

languages.iso

CS	SP1589326	(User Manual Chin Control - Czech)	1
DA	SP1589394	(User Manual Chin Control - Danish)	1
DE	SP1589395	(User Manual Chin Control - German)	1
EN	SP1589325	(User Manual Chin Control - English)	1
ES	SP1589535	(User Manual Chin Control - Spanish)	1
FI	SP1589436	(User Manual Chin Control - Finnish)	1
FR	SP1589437	(User Manual Chin Control - French)	1
HR	SP1589438	(User Manual Chin Control - Croatian)	1
IT	SP1589439	(User Manual Chin Control - Italian)	1
LT	SP1589440	(User Manual Chin Control - Lithuanian)	1
NL	SP1589441	(User Manual Chin Control - Dutch)	1
NO	SP1589442	(User Manual Chin Control - Norwegian)	1
PL	SP1589443	(User Manual Chin Control - Polish)	1
PT	SP1589444	(User Manual Chin Control - Portuguese)	1
SL	SP1589445	(User Manual Chin Control - Slovenian)	1
SV	SP1589465	(User Manual Chin Control - Swedish)	1

(User Manuals OptiCharge Charger)





BG	Български	HE	עברית	PL	Polski
CS	Čeština	HR	Hrvatski	PT	Português
DA	Dansk	HU	Magyar	RU	Русский
DE	Deutsch	IT	Italiano	RO	Română
EN	English	JA	日本語	SA	العربية
ES	Español	LV	Latviešu	SK	Slovenčina
ET	Eesti	LT	Lietuvių	SL	Slovenščina
FI	Suomi	NL	Nederlands	SV	Svenska
FR	Français	NO	Norsk	UK	Українська

languages.iso

EN	SP1636951	(User Manual OptiCharge Charger - English)	1
BG	SP1639049	(User Manual OptiCharge Charger - EDO (BG/CS/EN/ET/HR/HU/LT/LV/PL/RO/RU/SK/SL))	1
CS	SP1639049	(User Manual OptiCharge Charger - EDO (BG/CS/EN/ET/HR/HU/LT/LV/PL/RO/RU/SK/SL))	1
ET	SP1639049	(User Manual OptiCharge Charger - EDO (BG/CS/EN/ET/HR/HU/LT/LV/PL/RO/RU/SK/SL))	1
HR	SP1639049	(User Manual OptiCharge Charger - EDO (BG/CS/EN/ET/HR/HU/LT/LV/PL/RO/RU/SK/SL))	1
HU	SP1639049	(User Manual OptiCharge Charger - EDO (BG/CS/EN/ET/HR/HU/LT/LV/PL/RO/RU/SK/SL))	1
LT	SP1639049	(User Manual OptiCharge Charger - EDO (BG/CS/EN/ET/HR/HU/LT/LV/PL/RO/RU/SK/SL))	1
LV	SP1639049	(User Manual OptiCharge Charger - EDO (BG/CS/EN/ET/HR/HU/LT/LV/PL/RO/RU/SK/SL))	1
PL	SP1639049	(User Manual OptiCharge Charger - EDO (BG/CS/EN/ET/HR/HU/LT/LV/PL/RO/RU/SK/SL))	1
RO	SP1639049	(User Manual OptiCharge Charger - EDO (BG/CS/EN/ET/HR/HU/LT/LV/PL/RO/RU/SK/SL))	1
RU	SP1639049	(User Manual OptiCharge Charger - EDO (BG/CS/EN/ET/HR/HU/LT/LV/PL/RO/RU/SK/SL))	1
SK	SP1639049	(User Manual OptiCharge Charger - EDO (BG/CS/EN/ET/HR/HU/LT/LV/PL/RO/RU/SK/SL))	1
SL	SP1639049	(User Manual OptiCharge Charger - EDO (BG/CS/EN/ET/HR/HU/LT/LV/PL/RO/RU/SK/SL))	1
DA	SP1640325	(User Manual OptiCharge Charger - Danish)	1
DE	SP1640326	(User Manual OptiCharge Charger - German)	1
ES	SP1640327	(User Manual OptiCharge Charger - Spanish)	1
FI	SP1640328	(User Manual OptiCharge Charger - Finnish)	1
FR	SP1640329	(User Manual OptiCharge Charger - French)	1
IT	SP1640330	(User Manual OptiCharge Charger - Italian)	1
NL	SP1640331	(User Manual OptiCharge Charger - Dutch)	1
NO	SP1640332	(User Manual OptiCharge Charger - Norwegian)	1
PT	SP1640333	(User Manual OptiCharge Charger - Portuguese)	1

(User Manuals OptiCharge Charger)

SV	SP1640334	(User Manual OptiCharge Charger - Swedish)	1
SA	SP1654740	(User Manual OptiCharge - Arabic)	1

				
Storm4	O	O	P	P
G50	O	O	P	special
FDX	O	O	P	special
Typhoon	O	P	O	O
Typhoon II	O	P	O	O
Zephyr	O	P	O	O
Tornado	O	P	O	O
Storm3	O	P	O	O


- -

(For further information please contact your customer service.)

1

(Storm4)


-	(Chassis Fittings)	1
-	(Chassis Fittings (with Ultra Low Maxx Seat))	1
-	(Front swingarm)	1
-	(Front Swing Arm Storm4 X-Plore)	1
-	(Motor Swingarm => 04.01.11)	1
-	(Motor Swingarm SSD)	1
-	(Motor Swingarm SSD (with Ultra Low Maxx Seat))	1
-	(Motor Swingarm SSD, Anti-Tippers flip up)	1
-	(Rear suspension => 03.08.12)	1
-	(Rear suspension 04.08.12 =>)	1
-	(Upper suspension, Storm4 X-Plore)	1
-	(Transport fixation => 03.08.12)	1
-	(Transport fixation 04.08.12 =>)	1
-	(Shrouds)	1
-	(Shrouds (with Ultra Low Maxx Seat))	1
-	(Shroud spring bracket/rear light)	1





1	SP1530317	(Main Bolt)	1
2	SP1540337	(Main Bolt (Storm4 X-Plore))	1
3	SP1526347	BATTERILUKE STORM4 STORM4	1
4	SP1524458	(Cable support, cpl.)	1
5	SP1561820	(Plug for M8)	3

(Chassis Fittings (with Ultra Low Maxx Seat))

1602106.kim




1	SP1530317	(Main Bolt)	1
2	SP1540337	(Main Bolt (Storm4 X-Plore))	1
3	SP1543040	(Battery door)	1
4	SP1533336	(Cable support, cpl.)	1
5	SP1561820	(Plug for M8)	3




1	SP1527905	SVINGARM HØ STORM4	1
2	SP1527904	SVINGARM VE STORM4	1

(Front Swing Arm Storm4 X-Plore)

1537881X.dwg




1	SP1540359	(Swingarm, Storm4 X-Plore, RH)	1
2	SP1540360	(Swingarm, Storm4 X-Plore, LH)	1
3	SP1540361	(Suspension bracket, Storm4 X-Plore, RH)	1
4	SP1541357	(Suspension bracket, Storm4 X-Plore, LH)	1
5	SP1541358	(Cover suspension bracket, Storm4 X-Plore, RH)	1
6	SP1541359	(Cover suspension bracket, Storm4 X-Plore, LH)	1
7	SPL01043	GLIDER SVINGARM STORM4	1



1	-	(No longer available.)	1
2	SP1541799	(Swingarm RH)	1
3	SP1530581	FORINGSSETT STORM4	1

(Motor Swingarm SSD)


1533674X.dwg



1	SP1541363	MONTERINGSKIT SSD MOTOR STORM4 KPL HØ OG VE	1
2	SP1541367	LAGERKIT SVINGARM STORM4	1
3	SPL01043	GLIDER SVINGARM STORM4	1
4	SP1541761	SVINGARM SSD MOTOR HØ STORM4	1
4	SP1541762	SVINGARM SSD MOTOR VE STORM4	1
5	SP1533974	(Bridge)	1

(Motor Swingarm SSD (with Ultra Low Max Seat))


1602356.kim



1	SP1533974	(Bridge)	1
2	SP1541363	MONTERINGSKIT SSD MOTOR STORM4 KPL HØ OG VE	1
3	SP1541367	LAGERKIT SVINGARM STORM4	1
4	SP1587251	(Swingarm LH)	1
5	SP1587252	(Swingarm RH)	1
6	SPL01043	GLIDER SVINGARM STORM4	1

(Motor Swingarm SSD, Anti-Tippers flip up)


1567106.no



1	SP1541363	MONTERINGSKIT SSD MOTOR STORM4 KPL HØ OG VE	1
2	SP1552980	(Swingarm LH)	1
3	SP1552979	(Swingarm RH)	1
4	SPL01043	GLIDER SVINGARM STORM4	1


(Rear suspension => 03.08.12)

152896X.no



1	SP1530318	DEMPER SOFT KPL STORM4	1
1	SP1530337	DEMPER MEDIUM KPL STORM4	1
1	SP1530338	DEMPER HARD KPL STORM4	1
1	SP1540343	(Shock absorber, Storm4 X-Plore, cpl.)	1
2	SP1530444	FESTEKIT STØTDEMPERE STORM4	1


(Rear suspension 04.08.12 =>)



1	SP1564530	(Shock absorber, soft, cpl.)	1
2	SP1564531	(Shock absorber, middle, cpl.)	1
3	SP1564532	(Shock absorber, hard, cpl.)	1
-	SP1564636	(Shock absorber, Storm4 X-Plore, cpl.)	1
4	SP1528950	(Shock absorber support, cpl.)	1
5	SP1564637	(Mounting-Kit)	1

Invacare® Storm4 | Chassis
(Upper suspension, Storm4 X-Plore)


1540327X.iso




1	SP1540348	(Shock absorber front swing arm, soft, cpl.)	1
1	SP1540349	(Shock absorber front swing arm, hard, cpl.)	1
2	SP1540350	(Mountingkit)	1

(Transport fixation => 03.08.12)

156419.dwg




1	SP1563104	TRANSPORTFESTER BAK PAR STORM4 <i>not for passenger transportation</i>	1
2	SP1528974	TRANSPORTFESTE FORAN PAR STORM4 <i>not for passenger transportation</i>	1



1	SP1528951	TRANSPORTFESTE BAK STORM4 <i>not for passenger transportation</i>	1
2	SP1528974	TRANSPORTFESTE FORAN PAR STORM4 <i>not for passenger transportation</i>	1

(Shrouds)

SH4Shrouds.no




	SP1553166	LAKKSTIFT BLÅ STORM4 STORM4		1
	SP1553148	LAKKSTIFT SVART STORM4 STORM4		1
	SP1553167	(Touch-up paint - Champagne)		1
	SP1553169	(Touch-up paint - Red)		1
	SP1553168	(Touch-up paint - White (RAL 9016))		1
	SP1553210	(Touch-up paint - Green (RAL 6003))		1
	SP1553211	LAKKSTIFT SØLV STORM4 STORM4		1
1	SP1530339	FRONTKÅPE BLÅ STORM4	→ 12.04.2012	1
1	SP1559801	DEKSEL FORAN BLÅ STORM4 FRA APRIL 2012	13.04.2012 →	1
2	SP1530340	(Front shroud, cpl., Black (RAL 9005))	→ 09.04.2012	1
2	SP1559802	DEKSEL FORAN SORT STORM4 FRA APRIL 2012	10.04.2012 →	1
3	SP1530341	(Front shroud, cpl., Champagne)	→ 12.04.2012	1
3	SP1559803	(Front shroud, cpl. - champagne)	13.04.2012 →	1
4	SP1530342	(Front shroud, cpl., Red)	→ 12.04.2012	1
4	SP1559804	(Front shroud, cpl. - red)	13.04.2012 →	1
5	SP1530343	(Front shroud, cpl., White (RAL 9016))	→ 12.04.2012	1
5	SP1559805	(Front shroud, cpl. - white (RAL 9016))	13.04.2012 →	1
6	SP1530344	(Front shroud, cpl., Green (RAL 6003))	→ 12.04.2012	1
6	SP1559806	(Front shroud, cpl. - green (RAL 6003))	13.04.2012 →	1
7	SP1530345	DEKSEL STORT SØLV V/HEIS STORM4		1
8	SP1530346	DEKSEL BAK BLÅ STORM4		1
9	SP1530347	DEKSEL BAK SORT STORM4		1
10	SP1530348	(Rear shroud, cpl., Champagne)		1
11	SP1530427	(Rear shroud, cpl., Red)		1
12	SP1530428	(Rear shroud, cpl., White (RAL 9016))		1
13	SP1530429	(Rear shroud, cpl., Green (RAL 6003))		1
14	1419274	REFLEKTOR RØD TYPHOON		1
15	0572-06266	INVACARE LOGO STORM4		1
16	1440769	(Hand screw M6x15)		1
17	SP1522428	DECAL KIT STORM4 STORM4		1
18	SP1523651	DECAL KIT STORM4 SORT STORM4		1

19	SP1529145	(Decal-Kit Champagne)	1
20	SP1529393	(Decal-Kit Red (RAL 3007))	1
21	SP1529395	(Decal-Kit White (RAL 9016))	1
22	SP1529397	(Decal-Kit Green (RAL 6003))	1

(Shrouds (with Ultra Low Maxx Seat))


SH4Shrouds iso



	SP1553166	LAKKSTIFT BLÅ STORM4 STORM4	1
	SP1553148	LAKKSTIFT SVART STORM4 STORM4	1
	SP1553167	(Touch-up paint - Champagne)	1
	SP1553169	(Touch-up paint - Red)	1
	SP1553168	(Touch-up paint - White (RAL 9016))	1
	SP1553210	(Touch-up paint - Green (RAL 6003))	1
	SP1553211	LAKKSTIFT SØLV STORM4 STORM4	1
1	SP1559801	DEKSEL FORAN BLÅ STORM4 FRA APRIL 2012	1
2	SP1559802	DEKSEL FORAN SORT STORM4 FRA APRIL 2012	1
3	SP1559803	(Front shroud, cpl. - champagne)	1
4	SP1559804	(Front shroud, cpl. - red)	1
5	SP1559805	(Front shroud, cpl. - white (RAL 9016))	1
6	SP1559806	(Front shroud, cpl. - green (RAL 6003))	1
7	SP1546754	(Middle shroud - White (RAL 9006))	1
8	SP1530346	DEKSEL BAK BLÅ STORM4	1
9	SP1530347	DEKSEL BAK SORT STORM4	1
10	SP1530348	(Rear shroud, cpl., Champagne)	1
11	SP1530427	(Rear shroud, cpl., Red)	1
12	SP1530428	(Rear shroud, cpl., White (RAL 9016))	1
13	SP1530429	(Rear shroud, cpl., Green (RAL 6003))	1
14	1419274	REFLEKTOR RØD TYPHOON	1
15	0572-06266	INVACARE LOGO STORM4	1
16	1440769	(Hand screw M6x15)	1
17	SP1522428	DECAL KIT STORM4 STORM4	1
18	SP1523651	DECAL KIT STORM4 SORT STORM4	1
19	SP1529145	(Decal-Kit Champagne)	1
20	SP1529393	(Decal-Kit Red (RAL 3007))	1
21	SP1529395	(Decal-Kit White (RAL 9016))	1
22	SP1529397	(Decal-Kit Green (RAL 6003))	1

(Shroud spring bracket/rear light)

1528810X.iso




1	SP1528810	(Shroud for spring bracket without light LH)	1
2	SP1528811	(Shroud for spring bracket without light RH)	1
3	SP1552927	(Shroud for spring bracket, with baggage carrier, without light, LH)	1
4	SP1552928	(Shroud for spring bracket, with baggage carrier, without light, RH)	1
5	SP1527836	DEKSEL BAK LED LYKT VE STORM4	1
6	SP1527835	DEKSEL BAK LED LYKT HØ STORM4	1
7	SP1552929	(Shroud for spring bracket, with baggage carrier, for LED-rearlight, LH)	1
8	SP1552930	(Shroud for spring bracket, with baggage carrier, for LED-rearlight, RH)	1
9	SP1528627	(Shroud for spring bracket for rearlight Germany, LH)	1
10	SP1528626	(Shroud for spring bracket for rearlight Germany, RH)	1

Invacare® Storm4 | Chassis
(Storm4 TTPlus)

-	(Chassis Fittings)	1
-	(Front swingarm)	1
-	(Front Swing Arm Storm4 X-Plore)	1
-	(Motor Swingarm (GGB))	0
-	(Suspension)	0
-	(Upper suspension, Storm4 X-Plore)	1
-	(Transport Fixation)	0
-	(Shrouds)	0
-	(Shrouds Spring Bracket)	0


1581904.iso




1	SP1582232	(Motor swingarm (GGB2), cpl. LH)	1
2	SP1582233	(Motor swingarm (GGB2), cpl. RH)	1
3	SP1574130	(Bridge (GGB2), cpl.)	1
4	SP1576525	(Shockabsorber Support X-Plore (GGB2), LH)	1
5	SP1576526	(Shockabsorber Support X-Plore (GGB2), RH)	1
6	SP1582218	(Brake Lever Kit (GGB2), cpl.)	1
7	SP1541367	LAGERKIT SVINGARM STORM4	1
8	5007791	HÄNDTK GUMMI T HANDBRO	1
9	SP1574283	(Rubber Stop)	1

(Suspension)

1581876 iso




1	SP1582237	(Spring Element <= 80kg)	1
2	SP1582363	(Spring Element 70-100kg)	1
3	SP1582364	(Spring Element 90-150kg)	1
4	SP1577234	(Suspension Support, LH)	1
5	SP1574142	(Suspension Support, RH)	1




1	SP1582234	(Transport fixation, rear)	1
2	SP1528974	TRANSPORTFESTE FORAN PAR STORM4	1

(Shrouds)

1581874.000



	SP1553166	LAKKSTIFT BLÅ STORM4 STORM4	1
	SP1553148	LAKKSTIFT SVART STORM4 STORM4	1
	SP1553167	(Touch-up paint - Champagne)	1
	SP1553169	(Touch-up paint - Red)	1
	SP1553168	(Touch-up paint - White (RAL 9016))	1
	SP1553210	(Touch-up paint - Green (RAL 6003))	1
	SP1553211	LAKKSTIFT SØLV STORM4 STORM4	1
1	SP1559801	DEKSEL FORAN BLÅ STORM4 FRA APRIL 2012	1
2	SP1559802	DEKSEL FORAN SORT STORM4 FRA APRIL 2012	1
3	SP1559803	(Front shroud, cpl. - champagne)	1
4	SP1559804	(Front shroud, cpl. - red)	1
5	SP1559805	(Front shroud, cpl. - white (RAL 9016))	1
6	SP1559806	(Front shroud, cpl. - green (RAL 6003))	1
7	SP1530345	DEKSEL STORT SØLV V/HEIS STORM4	1
8	SP1581950	(Rear Shroud, cpl. - blue (RAL 5026))	1
9	SP1581951	DEKSEL BAK SORT STORM4 MODULITE	1
10	SP1581952	(Rear Shroud, cpl. - champagne)	1
11	SP1582000	(Rear Shroud, cpl. - transparent cherry red)	1
12	SP1582001	(Rear Shroud, cpl. - white (RAL 9016))	1
13	SP1582002	(Rear Shroud cpl. - olive green (RAL 6003))	1
14	1419274	REFLEKTOR RØD TYPHOON	1
15	0572-06266	INVACARE LOGO STORM4	1
16	1440769	(Hand screw M6x15)	1
17	SP1522428	DECAL KIT STORM4 STORM4	1
18	SP1523651	DECAL KIT STORM4 SORT STORM4	1
19	SP1529145	(Decal-Kit Champagne)	1
20	SP1529393	(Decal-Kit Red (RAL 3007))	1
21	SP1529395	(Decal-Kit White (RAL 9016))	1
22	SP1529397	(Decal-Kit Green (RAL 6003))	1



1	SP1577034	(Shroud for spring bracket without light LH GGB2)	1
2	SP1577035	(Shroud for spring bracket without light RH GGB2)	1
3	SP1577056	(Shroud for spring bracket with luggage-rack, w/o light LH GGB2)	1
4	SP1577057	(Shroud for spring bracket with luggage-rack, w/o light RH GGB2)	1
5	SP1574138	(Shroud for spring bracket with LED light LH GGB2)	1
6	SP1574139	(Shroud for spring bracket with LED light RH GGB2)	1
7	SP1577054	(Shroud for spring bracket with luggage-rack, with LED light LH GGB2)	1
8	SP1577055	(Shroud for spring bracket with luggage-rack, with LED light RH GGB2)	1
9	SP1577036	(Shroud for spring bracket with light LH GGB2)	1
10	SP1577037	(Shroud for spring bracket with light RH GGB2)	1

Sete

ST4setel.no




1	-	(Storm4)	1
2	-	(Storm4 TTPlus)	1

2	-	(EasyAdapt Seat => 2017-08-10)	1
3	-	(Modulite Seat)	1
4	-	(Ultra Low Maxx Seat)	1
5	-	(Recaro Seat)	1
6	-	(Optimist Seat => 2018-12-31)	1


(EasyAdapt Seat => 2017-08-10)

1573129.iss




1	-	(Seat column / -raiser Storm4) [1025/1026/1030/1040]	1
2	-	(Seat column / -raiser Storm4 X-Plore) [1025/1026/1030/1040]	1
3	-	(Seat Height Adapter Storm4) [1040]	1
4	-	(Seat Height Adapter Storm4 X-Plore) [1040]	1
5	-	(Seat Tilt) [1002][1014]	1
6	-	(Seat Frame)	1
7	-	(Seat Plate)	1
8	-	(Adapter Kit for Universal Seat) [1064]	1
9	-	(Seat Cushion)	1

(Seat column / -raiser Storm4)



1	SP1530459	SETELØFT TIL SEPT 2013 STORM4 KPL	→ 14.05.2012	1
-	SP1568362	(Mounting kit for seatbracket)	→ 14.05.2012	1
2	SP1556263	SETELØFT FRA SEPT 2013 STORM 4 KPL	15.05.2012 →	1
-	SP1568361	(Mounting Kit lifter for seat bracket)	15.05.2012 →	1
3	SP1548536	(Seat column, fixed, height 45cm and 48cm, adjustable)		1
4	SP1529390	KABEL SETELØFT STORM4	→ 14.05.2012	1
5	SP1555994	KABEL FARTSREDUKSJON STORM4	15.05.2012 →	1
6	SP1559850	KABEL SETEHEIS LINAK STORM4	15.05.2012 →	1

(Seat column / -raiser Storm4 X-Plore)




1510800X.kso

1	SP1530459	SETELØFT TIL SEPT 2013 STORM4 KPL	→ 15.09.2013	1
-	SP1568362	(Mounting kit for seatbracket)	→ 15.09.2013	1
2	SP1556263	SETELØFT FRA SEPT 2013 STORM 4 KPL	16.09.2013 →	1
-	SP1568361	(Mounting Kit lifter for seat bracket)	16.09.2013 →	1
3	SP1548536	(Seat column, fixed, height 45cm and 48cm, adjustable)		1
4	SP1529390	KABEL SETELØFT STORM4	→ 15.09.2013	1
5	SP1555994	KABEL FARTSREDUKSJON STORM4	16.09.2013 →	1
6	SP1559850	KABEL SETEHEIS LINAK STORM4	16.09.2013 →	1

(Seat Height Adapter Storm4)


1560782.iso




1	SP1540334	SETEHØYDE ADAPTER STORM4	→ 19.04.2012	1
2	SP1556447	(Seatheightadapter cpl.)	20.04.2012 → 14.05.2012	1
3	SP1555071	SETEHØYDE ADAPTER STORM4 FRA 15/5-2012	15.05.2012 →	1

(Seat Height Adapter Storm4 X-Plore)


1560782.iso




1	SP1540334	SETEHØYDE ADAPTER STORM4	→ 19.04.2012	1
2	SP1556447	(Seatheightadapter cpl.)	20.04.2012 → 15.09.2013	1
3	SP1555071	SETEHØYDE ADAPTER STORM4 FRA 15/5-2012	16.09.2013 →	1



1	SP1551991	TILT BÆRER STORM4 [1002][1014]	20.04.2012 →	1
2	SP1559711	TILT KPL U BÆRER STORM4 [1002][1014]	20.04.2012 →	1
3	-	(No longer available - please refer to Replacement Info for which part is to be used instead.) [1002][1014] - SP1568271	→ 19.04.2012	0
4	SP1537803	TILT HOLDER STORM4 [1002][1014]	→ 19.04.2012	1
5	SP1551010	(Tilt spindle, cpl.) [1002]	18.03.2011 →	1
6	SP1551009	(Mounting-Kit)	18.03.2011 →	1
7	SP1532482	(Manual tilt, cpl.) [1002]	→ 17.03.2011	1
8	SP1532483	(Mounting-Kit for manual tilt)	→ 17.03.2011	1
9	SP1532480	MOTOR TILT STORM4 [1014] ACS		1
9	SP1607900	(REAC Motor RE35 stroke=50 (w/o switchbox) - LiNX) [1014] LiNX	01.01.2017 →	1
10	SP1532481	MONTERINGSKIT TILT STORM4		1
11	SP1588880	(Module Cover, cpl.)	01.01.2016 →	1
12	SP1603273	(Cable Support Sheets, cpl.) [1030][1040]		1



1	SP1536856	FRONTRAMME SETE RECARO STORM4	1
2	SP1542967	BAKRAMME EASYADAPT SETE STORM4	1




1573205.iso

1	SP1530486	SETEPLATE KORT VE STORM4	1
2	SP1530487	SETEPLATE KORT HØ STORM4	1
1	SP1530488	SETEPLATE MID VE STORM4	1
2	SP1530489	SETEPLATE MID HØ STORM4	1
1	SP1530490	SETEPLATE LANG VE STORM4	1
2	SP1530491	SETEPLATE LANG HØ STORM4	1
3	SP1530492	DEKKPLATE SETE STORM4	1

(Adapter Kit for Universal Seat)

156173.no




1 SP1563052


(Adapter Kit for Universal Seat)

1

Matrix Flo-tech



Onyx Cushion




Matrix Flo-tech			Onyx Cushion		
Wheelchair	Wheelchair	Color	Wheelchair	Wheelchair	Color
1	380mm	410mm	25	380mm	410mm
2	380mm	460mm	26	380mm	460mm
3	380mm	510mm	27	380mm	510mm
4	430mm	410mm	28	430mm	410mm
5	430mm	460mm	29	430mm	460mm
6	430mm	510mm	30	430mm	510mm
7	480mm	410mm	31	480mm	410mm
8	480mm	460mm	32	480mm	460mm
9	480mm	510mm	33	480mm	510mm
10	530mm	430mm	34	530mm	430mm
11	530mm	460mm	35	530mm	460mm
12	530mm	510mm	36	530mm	510mm
13	380mm	410mm	37	380mm	410mm
14	380mm	460mm	38	380mm	460mm
15	380mm	510mm	39	380mm	510mm
16	430mm	410mm	40	430mm	410mm
17	430mm	460mm	41	430mm	460mm
18	430mm	510mm	42	430mm	510mm
19	480mm	410mm	43	480mm	410mm
20	480mm	460mm	44	480mm	460mm
21	480mm	510mm	45	480mm	510mm
22	530mm	430mm	46	530mm	430mm
23	530mm	460mm	47	530mm	460mm
24	530mm	510mm	48	530mm	510mm
		⇒ 30.06.14			⇒ 30.06.14
		Black Smooth			Grey Plush
		⇒ 30.11.14			Black orange seam
		Dartex Grey			Black
		01.12.14 ⇒			
		Dartex Black			

1573212.ISO

1	SPFLB15/16	SETEPUTE SB38 SD41 SORT STORM4	1	
2	SPFLB15/18	SETEPUTE SB38 SD46 SORT STORM4	1	
3	SPFLB15/20	SETEPUTE SB38 SD51 SORT STORM4	1	
4	SPFLB17/16	SETEPUTE SB43 SD41 SORT STORM4	1	
5	SPFLB17/18	SETEPUTE SB43 SD46 STD STORM4	1	
6	SPFLB17/20	SETEPUTE SB43 SD51 SORT STORM4	1	
7	SPFLB19/16	SETEPUTE SB48 SD41 SORT STORM4	1	
8	SPFLB19/18	SETEPUTE SB48 SD46 SORT STORM4	1	
9	SPFLB19/20	SETEPUTE SB48 SD51 SORT STORM4	1	
10	SPFCB21/17	SETEPUTE SB53 SD43 SORT STORM4	1	
11	SPFCB21/18	SETEPUTE SB53 SD46 SORT STORM4	1	
12	SPFCB21/20	SETEPUTE SB53 SD51 SORT STORM4	1	
13	SP1525250-0232	\$ BRUK SP1529550-TR26 STORM4 DARTEX	→ 30.11.2014	1
13	SP1529546-TR26	(Seat Cushion Flo-tech Lite 15x16" Dartex black)	01.12.2014 →	1
14	-	(No longer available - please choose Dartex black.)	→ 30.11.2014	1
		Replacement: SP1529547-TR26		
14	SP1529547-TR26	(Seat Cushion Flo-tech Lite 15x18" Dartex black)	01.12.2014 →	1
15	-	(No longer available - please choose Dartex black.)	→ 30.11.2014	1
		Replacement: SP1529548-TR26		
15	SP1529548-TR26	(Seat Cushion Flo-tech Lite 15x20" Dartex black)	01.12.2014 →	1
16	SP1525254-0232	SETEPUTE SB43SD41 EASYAD STORM4 DARTEX	→ 30.11.2014	1
16	SP1529549-TR26	(Seat Cushion Flo-tech Lite 17x16" Dartex black)	01.12.2014 →	1
17	-	(No longer available - please choose Dartex black.)	→ 30.11.2014	1
		Replacement: SP1529550-TR26		
17	SP1529550-TR26	SETEPUTE SB43SD46 DARTEX STORM4 TDX	01.12.2014 →	1
18	-	(No longer available - please choose Dartex black.)	→ 30.11.2014	1
		Replacement: SP1529551-TR26		
18	SP1529551-TR26	SETEPUTE SB43SD51 DARTEX STORM4	01.12.2014 →	1
19	SP1525257-0232	SETEPUTE SB48SD41 EASYAD STORM4 DARTEX	→ 30.11.2014	1
19	SP1529552-TR26	(Seat Cushion Flo-tech Lite 19x16" Dartex black)	01.12.2014 →	1
20	-	(No longer available - please choose Dartex black.)	→ 30.11.2014	1
		Replacement: SP1529553-TR26		
20	SP1529553-TR26	(Seat Cushion Flo-tech Lite 19x18" Dartex black)	01.12.2014 →	1

(Seat Cushion)


21	-	(No longer available - please choose Dartex black.) <i>Replacement: SP1529554-TR26</i>	→ 30.11.2014	1
21	SP1529554-TR26	SETEPUTE SB48SD51 DARTEX STORM4	01.12.2014 →	1
22	SP1525260-0232	SETEPUTE SB53SD43 EASYAD STORM4 DARTEX	→ 30.11.2014	1
22	SP1529555-TR26	(Seat Cushion Flo-tech Contour 21x17" Dartex black)	01.12.2014 →	1
23	SP1525261-0232	SETEPUTE SB53SD46 EASYAD STORM4 DARTEX	→ 30.11.2014	1
23	SP1529556-TR26	(Seat Cushion Flo-tech Contour 21x18" Dartex black)	01.12.2014 →	1
24	SP1525262-0232	SETEPUTE SB53SD51 EASYAD STORM4 DARTEX	→ 30.11.2014	1
24	SP1529570-TR26	(Seat Cushion Flo-tech Contour 21x20" Dartex black)	01.12.2014 →	1
25	-	(No longer available - please choose Black Chenille.) <i>Replacement: SP1529546-TR35</i>		1
25	SP1529546-TR35	(Seat Cushion Flo-tech Lite 15x16" Black Chenille)		1
26	-	(No longer available - please choose Black Chenille.) <i>Replacement: SP1529547-TR35</i>		1
26	SP1529547-TR35	(Seat Cushion Flo-tech Lite 15x18" Black Chenille)		1
27	SP1525253-0237	(Seat Cushion Flo-tech Lite 15x20" Grey Plush)		1
28	SP1525254-0237	(Seat Cushion Flo-tech Lite 17x16" Grey Plush)		1
29	SP1525255-0237	(Seat Cushion Flo-tech Lite 17x18" Grey Plush)		1
30	-	(No longer available - please choose Black Chenille.) <i>Replacement: SP1529551-TR35</i>		1
30	SP1529551-TR35	(Seat Cushion Flo-tech Lite 17x20" Black Chenille)		1
31	SP1525257-0237	(Seat Cushion Flo-tech Lite 19x16" Grey Plush)		1
32	SP1525258-0237	(Seat Cushion Flo-tech Lite 19x18" Grey Plush)		1
33	SP1525259-0237	(Seat Cushion Flo-tech Lite 19x20" Grey Plush)		1
34	SP1525260-0237	(Seat Cushion Flo-tech Contour 21x17" Grey Plush)		1
35	SP1525261-0237	(Seat Cushion Flo-tech Contour 21x18" Grey Plush)		1
36	SP1525262-0237	SETEPUTE SB53SD51 EASYAD STORM4		1
37	SP1529546	SETEPUTE EASY-A SB38SD41 SORT M ORANG SØM		1
38	SP1529547	SETEPUTE EASY-A SB38SD46 SORT M ORANG SØM		1
39	SP1529548	SETEPUTE EASY-A SB38SD51 SORT M ORANG SØM		1
40	SP1529549	SETEPUTE EASY-A SB43SD41 SORT M ORANG SØM		1
41	SP1529550	SETEPUTE EASY-A SB43SD46 SORT M ORANG SØM		1
42	SP1529551	SETEPUTE EASY-A SB43SD51 SORT M ORANG SØM		1
43	SP1529552	SETEPUTE EASY-A SB48SD41 SORT M ORANG SØM		1
44	SP1529553	SETEPUTE EASY-A SB48SD46 SORT M ORANG SØM		1
45	SP1529554	SETEPUTE EASY-A SB48SD51 SORT M ORANG SØM		1
46	SP1529555	SETEPUTE EASY-A SB53SD43 SORT M ORANG SØM		1
47	SP1529556	SETEPUTE EASY-A SB53SD46 SORT M ORANG SØM		1
48	SP1529570	SETEPUTE EASY-A SB53 SD5 SORT M ORANG SØM		1
49	SP1532058	(Seat Cushion Onyx 16x16" Black)		1
50	SP1534080	(Seat Cushion Onyx 16x18" Black)		1
51	SP1532059	(Seat Cushion Onyx 18x16" Black)		1
52	SP1532060	(Seat Cushion Onyx 18x18" Black)		1
53	SP1534136	(Seat Cushion Onyx 20x16" Black)		1
54	SP1534199	(Seat Cushion Onyx 20x18" Black)		1




1	-	(Seat column / -raiser Storm4) [1023/1025/1030/1040]	1
2	-	(Seat column / -raiser Storm4 X-Plore) [1023/1025/1030/1040]	1
3	-	(Seat Height Adapter Storm4) [1040]	1
4	-	(Seat Height Adapter Storm4 X-Plore) [1040]	1
5	-	(Seat Tilt) [1002][1014]	1
6	-	(One Piece Seat Plate) [2110]	1
7	-	(Telescopic Seat Frame and Plate) [2111]	1
8	-	(Telescopic Seat Frame and Sling Seat) [2112]	1
9	-	(Telescopic Seat Frame and Universal Adapter) [2109]	1
10	-	(Seat Cushions "Standard" & "Matrx Flo-Tech") [2141/2142/2143]	1
10	-	(Seat Cushions "Matrx Libra") [2147]	1
10	-	(Seat Cushions "Matrx Vi") [2144]	1

(Seat Tilt)

1583167.indd



1	SP1551991	TILT BÆRER STORM4 [1002][1014]		1
2	SP1559711	TILT KPL U BÆRER STORM4 [1002][1014]		1
3	SP1551010	(Tilt spindle, cpl.) [1002]		1
4	SP1551009	(Mounting-Kit)		1
5	SP1532480	MOTOR TILT STORM4 [1014] ACS		1
5	SP1607900	(REAC Motor RE35 stroke=50 (w/o switchbox) - LiNX) [1014] LiNX	01.01.2017 →	1
6	SP1532481	MONTERINGSKIT TILT STORM4		1
7	SP1573217	SETE SUPPORT MODULITE [2001]		1
8	SP1588880	(Module Cover, cpl.)	01.01.2016 →	1
9	SP1603273	(Cable Support Sheets, cpl.) [1030][1040]		1




1571864.iso


1	SP1571297	(One Piece Seat Plate)	1
1	SP1574030	(One Piece Seat Plate, pre-drilled for ACT)	1
2	SPU21087	(Adhesive hook fastener strip 38x400mm)	2

(Telescopic Seat Frame and Plate)

1574269 iso



1	-	(Telescopic Seat Frame)	1
2	-	(Armrest bracket)	1
3	-	(Legrest Support)	1
4	-	(Seat Plate)	1




1572009.iso

1	SP1572977	(Seat Frame 510mm, pre set up to SW 430mm, cpl.)	1
2	SP1582216	SETERAMME KPL MODULITE MODULITE	1
3	SP1571810	(Sealing Strap with plugs)	1
4	SP1571811	(Groove Stone with mounting parts)	1
5	SP1554646	RØRBRAKETT SETEBUNN MODULITE	1
6	SP1588443	C-PROFIL SETE 51 CM LANG MODULITE	1
7	SP1588444	C-PROFIL SETE 44 CM LANG MODULITE	1

(Armrest bracket)


1571834 Rev



1 SP1570189

ARMLENSBRAKETT MODULITE


1



1572021.iso

1	SP1558935	BENSTØTTE FESTE SB38+43 MODULITE	1
2	SP1559294	BENSTØTTE FESTE SB48+53 MODULITE	1
3	SP1558134	MELLOMSTYKKE FESTE BENST MODULITE	1
4	SP1558135	BRAKETT SETEPLATE BENST MODULITE SB48-53	1

Invacare® Storm4 | Sete
(Seat Plate)




1572310.iso

1	SP1572199	SETEPLATE MODULITE	1
-	SP1656888	(Screw Kit Modulite Seat Plate) <i>Including 6 screws and 6 washers</i>	1
2	SPU21087	(Adhesive hook fastener strip 38x400mm)	2


(Telescopic Seat Frame and Sling Seat)

1574270.00




1	-	(Telescopic Seat Frame)	1
2	-	(Armrest bracket)	1
3	-	(Legrest Support)	1
4	-	(Sling Seat)	1

(Telescopic Seat Frame)




1571950.ISO

1	SP1572978	(Seat Frame 510mm, pre set up to SW 430mm, cpl.)	1
2	SP1582217	(Seat Frame 545mm, pre set up to SW 430mm, cpl.)	1
3	SP1571810	(Sealing Strap with plugs)	1
4	SP1571811	(Groove Stone with mounting parts)	1
5	SP1555400	BRACKET C-PROFIL MODULIT TDX SP2	1
6	SP1588443	C-PROFIL SETE 51 CM LANG MODULITE	1
7	SP1588444	C-PROFIL SETE 44 CM LANG MODULITE	1



1	SP1559305	FIRKANTRØR FESTE BENST MODULITE	1
2	SP1559307	FIRKANTRØR FESTE BENSTØT SB48+SB53 KPL MODULITE	1
3	SP1558134	MELLOMSTYKKE FESTE BENST MODULITE	1
4	SP1558135	BRAKETT SETEPLATE BENST MODULITE SB48-53	1

Invacare® Storm4 | Sete
(Sling Seat)




1571762.iso

1	SP1557877	VELCRO SETE SB38-44 STORM4 TDX SP2	1
2	SP1555911	VELCRO SETEBUNN SB46 53 STORM4 TDX SP2	1

(Telescopic Seat Frame and Universal Adapter)


1581866.iso



1	-	(Telescopic Seat Frame)	1
2	-	(Armrest bracket)	1
3	-	(Legrest Support)	1
4	-	(Adapter Kit for Universal Seat)	1

(Adapter Kit for Universal Seat)

1581063.iso



1 SP1580750

(Adapter Kit for Universal Seat 510mm)

1


(Seat Cushions "Standard" & "Matrix Flo-Tech")

Standard



1	410mm	Black
2	460mm	
3	510mm	

Matrix Flo-tech Contur Visco



64	430mm	430mm	Black
----	-------	-------	-------

Matrix Flo-tech Lite & Contour



4	380mm	410mm	Black smooth	
5	380mm	460mm		
6	380mm	510mm		
7	430mm	410mm		
8	430mm	460mm		
9	430mm	510mm		
10	480mm	410mm		
11	480mm	460mm		
12	480mm	510mm		
13	530mm	430mm		
14	530mm	460mm		
15	530mm	510mm		
16	380mm	410mm		Dartex Black
17	380mm	460mm		
18	380mm	510mm		
19	430mm	410mm		
20	430mm	460mm		
21	430mm	510mm		
22	480mm	410mm		
23	480mm	460mm		
24	480mm	510mm		
25	530mm	430mm		
26	530mm	460mm		
27	530mm	510mm		

28	380mm	410mm	Black Plush	
29	380mm	460mm		
30	380mm	510mm		
31	430mm	410mm		
32	430mm	460mm		
33	430mm	510mm		
34	480mm	410mm		
35	480mm	460mm		
36	480mm	510mm		
37	530mm	430mm		
38	530mm	460mm		
39	530mm	510mm		
40	380mm	410mm		Black orange seam
41	380mm	460mm		
42	380mm	510mm		
43	430mm	410mm		
44	430mm	460mm		
45	430mm	510mm		
46	480mm	410mm		
47	480mm	460mm		
48	480mm	510mm		
49	530mm	430mm		
50	530mm	460mm		
51	530mm	510mm		

52	380mm	410mm	Black Spacer fabric
53	380mm	460mm	
54	380mm	510mm	
55	430mm	410mm	
56	430mm	460mm	
57	430mm	510mm	
58	480mm	410mm	
59	480mm	460mm	
60	480mm	510mm	
61	530mm	430mm	
62	530mm	460mm	
63	530mm	510mm	

1606764 ISO

1	1420460	SETEPUTE SD41 STD DRAGON VERTIC <i>[2141]</i>	1
2	1420461	SETEPUTE SD46 STD DRAGON VERTIC <i>[2141]</i>	1
3	1420462	SETEPUTE SD51 STD DRAGON VERTIC <i>[2141]</i>	1
4	SPFLB15/16	SETEPUTE SB38 SD41 SORT STORM4 <i>[2142+2151]</i>	1
5	SPFLB15/18	SETEPUTE SB38 SD46 SORT STORM4 <i>[2142+2151]</i>	1
6	SPFLB15/20	SETEPUTE SB38 SD51 SORT STORM4 <i>[2142+2151]</i>	1
7	SPFLB17/16	SETEPUTE SB43 SD41 SORT STORM4 <i>[2142+2151]</i>	1
8	SPFLB17/18	SETEPUTE SB43 SD46 STD STORM4 <i>[2142+2151]</i>	1
9	SPFLB17/20	SETEPUTE SB43 SD51 SORT STORM4 <i>[2142+2151]</i>	1
10	SPFLB19/16	SETEPUTE SB48 SD41 SORT STORM4 <i>[2142+2151]</i>	1
11	SPFLB19/18	SETEPUTE SB48 SD46 SORT STORM4 <i>[2142+2151]</i>	1
12	SPFLB19/20	SETEPUTE SB48 SD51 SORT STORM4 <i>[2142+2151]</i>	1
13	SPFCB21/17	SETEPUTE SB53 SD43 SORT STORM4 <i>[2142+2151]</i>	1
14	SPFCB21/18	SETEPUTE SB53 SD46 SORT STORM4 <i>[2142+2151]</i>	1
15	SPFCB21/20	SETEPUTE SB53 SD51 SORT STORM4 <i>[2142+2151]</i>	1
16	SP1529546-TR26	(Seat Cushion Flo-tech Lite 15x16" Dartex black) <i>[2142+2152]</i>	1
17	SP1529547-TR26	(Seat Cushion Flo-tech Lite 15x18" Dartex black) <i>[2142+2152]</i>	1
18	SP1529548-TR26	(Seat Cushion Flo-tech Lite 15x20" Dartex black) <i>[2142+2152]</i>	1

(Seat Cushions "Standard" & "Matrix Flo-Tech")


19	SP1529549-TR26	(Seat Cushion Flo-tech Lite 17x16" Dartex black) <i>[2142+2152]</i>	1
20	SP1529550-TR26	SETEPUTE SB43SD46 DARTEX STORM4 TDX <i>[2142+2152]</i>	1
21	SP1529551-TR26	SETEPUTE SB43SD51 DARTEX STORM4 <i>[2142+2152]</i>	1
22	SP1529552-TR26	(Seat Cushion Flo-tech Lite 19x16" Dartex black) <i>[2142+2152]</i>	1
23	SP1529553-TR26	(Seat Cushion Flo-tech Lite 19x18" Dartex black) <i>[2142+2152]</i>	1
24	SP1529554-TR26	SETEPUTE SB48SD51 DARTEX STORM4 <i>[2142+2152]</i>	1
25	SP1529555-TR26	(Seat Cushion Flo-tech Contour 21x17" Dartex black) <i>[2142+2152]</i>	1
26	SP1529556-TR26	(Seat Cushion Flo-tech Contour 21x18" Dartex black) <i>[2142+2152]</i>	1
27	SP1529570-TR26	(Seat Cushion Flo-tech Contour 21x20" Dartex black) <i>[2142+2152]</i>	1
28	SP1529546-TR35	(Seat Cushion Flo-tech Lite 15x16" Black Chenille) <i>[2142+2154]</i>	1
29	SP1529547-TR35	(Seat Cushion Flo-tech Lite 15x18" Black Chenille) <i>[2142+2154]</i>	1
30	SP1529548-TR35	(Seat Cushion Flo-tech Lite 15x20" Black Chenille) <i>[2142+2154]</i>	1
31	SP1529549-TR35	(Seat Cushion Flo-tech Lite 17x16" Black Chenille) <i>[2142+2154]</i>	1
32	SP1529550-TR35	(Seat Cushion Flo-tech Lite 17x18" Black Chenille) <i>[2142+2154]</i>	1
33	SP1529551-TR35	(Seat Cushion Flo-tech Lite 17x20" Black Chenille) <i>[2142+2154]</i>	1
34	SP1529552-TR35	(Seat Cushion Flo-tech Lite 19x16" Black Chenille) <i>[2142+2154]</i>	1
35	SP1529553-TR35	SETEPUTE SB48 SD46 SORT STORM4 <i>[2142+2154]</i>	1
36	SP1529554-TR35	SETEPUTE SB48SD51 PLYSJ STORM4 <i>[2142+2154]</i>	1
37	SP1529555-TR35	(Seat Cushion Flo-tech Contour 21x17" Black Chenille) <i>[2142+2154]</i>	1
38	SP1529556-TR35	(Seat Cushion Flo-tech Contour 21x18" Black Chenille) <i>[2142+2154]</i>	1
39	SP1529570-TR35	(Seat Cushion Flo-tech Contour 21x20" Black Chenille) <i>[2142+2154]</i>	1
40	SP1529546	SETEPUTE EASY-A SB38SD41 SORT M ORANG SØM <i>[2142+2153]</i>	1
41	SP1529547	SETEPUTE EASY-A SB38SD46 SORT M ORANG SØM <i>[2142+2153]</i>	1
42	SP1529548	SETEPUTE EASY-A SB38SD51 SORT M ORANG SØM <i>[2142+2153]</i>	1
43	SP1529549	SETEPUTE EASY-A SB43SD41 SORT M ORANG SØM <i>[2142+2153]</i>	1
44	SP1529550	SETEPUTE EASY-A SB43SD46 SORT M ORANG SØM <i>[2142+2153]</i>	1
45	SP1529551	SETEPUTE EASY-A SB43SD51 SORT M ORANG SØM <i>[2142+2153]</i>	1
46	SP1529552	SETEPUTE EASY-A SB48SD41 SORT M ORANG SØM <i>[2142+2153]</i>	1
47	SP1529553	SETEPUTE EASY-A SB48SD46 SORT M ORANG SØM <i>[2142+2153]</i>	1
48	SP1529554	SETEPUTE EASY-A SB48SD51 SORT M ORANG SØM <i>[2142+2153]</i>	1
49	SP1529555	SETEPUTE EASY-A SB53SD43 SORT M ORANG SØM <i>[2142+2153]</i>	1
50	SP1529556	SETEPUTE EASY-A SB53SD46 SORT M ORANG SØM <i>[2142+2153]</i>	1

(Seat Cushions "Standard" & "Matrx Flo-Tech")

51	SP1529570	SETEPUTE EASY-A SB53 SD5 SORT M ORANG SØM <i>[2142+2153]</i>	1
52	1629216	(Seat Cushion Flo-tech Lite 15x16" Black Spacer fabric)	1
53	1629217	(Seat Cushion Flo-tech Lite 15x18" Black Spacer fabric)	1
54	1629218	(Seat Cushion Flo-tech Lite 15x20" Black Spacer fabric)	1
55	1629219	(Seat Cushion Flo-tech Lite 17x16" Black Spacer fabric)	1
56	1629220	(Seat Cushion Flo-tech Lite 17x18" Black Spacer fabric)	1
57	1629221	(Seat Cushion Flo-tech Lite 17x20" Black Spacer fabric)	1
58	1629222	(Seat Cushion Flo-tech Lite 19x16" Black Spacer fabric)	1
59	1629223	(Seat Cushion Flo-tech Lite 19x18" Black Spacer fabric)	1
60	1629224	(Seat Cushion Flo-tech Lite 19x20" Black Spacer fabric)	1
61	1629225	(Seat Cushion Flo-tech Contour 21x17" Black Spacer fabric)	1
62	1629226	(Seat Cushion Flo-tech Contour 21x18" Black Spacer fabric)	1
63	1629227	(Seat Cushion Flo-tech Contour 21x20" Black Spacer fabric)	1
64	SP1501113	SETEPUTE 43*43 TDX FLOTEC CONTUR VISCO <i>[2143]</i>	1

(Seat Cushions "Matrx Libra")


1603356.iss



Matrx Libra

		mm			
		380	430	480	530
mm	410	1	4	7	-
	450	2	5	8	10
	510	3	6	9	11
	560	-	-	-	12

1	LC1516	(Matrx Libra cushions 380x410 mm)	1
2	LC1518	(Matrx Libra cushions 380x460 mm)	1
3	LC1520	(Matrx Libra cushions 380x510 mm)	1
4	LC1716	(Matrx Libra cushions 430x410 mm)	1
5	LC1718	(Matrx Libra cushions 430x460 mm)	1
6	LC1720	(Matrx Libra cushions 430x510 mm)	1
7	LC1916	(Matrx Libra cushions 480x410 mm)	1
8	LC1918	(Matrx Libra cushions 480x460 mm)	1
9	LC1920	(Matrx Libra cushions 480x510 mm)	1
10	LC2118	(Seat Cushion Matrx Libra 21x18")	1
11	LC2120	(Seat Cushion Matrx Libra 21x20")	1
12	LC2122	(Seat Cushion Matrx Libra 21x22")	1




Matrx Vi

		mm			
		380	430	480	530
410		1	4	7	-
460	mm	2	5	8	10
510		3	6	9	11

1	1579786	(Matrx Flo-Tech VI cushions 410x410 mm)	1
2	1579787	(Matrx Flo-Tech VI cushions 410x460 mm)	1
3	1579788	(Matrx Flo-Tech VI cushions 410x510 mm)	1
4	1579920	(Matrx Flo-Tech VI cushions 460x410 mm)	1
5	1579921	(Matrx Flo-Tech VI cushions 460x460 mm)	1
6	1579922	(Matrx Flo-Tech VI cushions 460x510 mm)	1
7	1579923	(Matrx Flo-Tech VI cushions 510x410 mm)	1
8	1579924	(Matrx Flo-Tech VI cushions 510x460 mm)	1
9	1579925	(Matrx Flo-Tech VI cushions 510x510 mm)	1
10	1579926	(Matrx Flo-Tech VI cushions 560x460 mm)	1
11	1579927	(Matrx Flo-Tech VI cushions 560x510 mm)	1

(Ultra Low Maxx Seat)

-	(Seat Column)	1
-	(Seat Tilt)	1
-	(Seat Lifter with Tilt)	1
-	(Posterior-Tilt)	1
-	(Seat Plates & Side Rails)	1
-	(Seat Cushions "Matrx Vi")	1
-	(Seat Cushions "Matrx Flo-Tech")	1
-	(Seat Cushions "Matrx Libra")	1




1 SP1589681

(Seat column, cpl.)


1

(Seat Tilt)


1602021.km



1	SP1589710	(Seat Mounting Adapter)	1
2	SP1589711	(Seat Adapter Rails)	1
3	SP1585523	(Tilt Module, cpl.)	1
4	SP1588440	(Tilt Actuator, cpl. with hardware)	1
-	SPTRE0755	(Cable, Tilt Actuator)	1



1	SP1589710	(Seat Mounting Adapter)	1
2	SP1589711	(Seat Adapter Rails)	1
3	SP1585524	(Lifter/Tilt Module, cpl.)	1
4	SP1588441	(Tilt actuator, cpl. with hardware)	1
5	SP1588442	(Lifter Actuator, cpl. with mounting hardware)	1
6	SPTRC2071	(Cover, Rear Wiring Lifter, cpl.)	1
7	SP1604199	ENERGIKJEDE TDX-SP2	1




1 SPTRC1923

(Pre-Tilt Mounting Kit 5°)
Posterior-Tilt

1


1588152 no



1	SPTRC2187	(Seat Pans - SW 405-505mm (16-20"))	1
2	SPTRC2189	(Seat Pans - SW 505-610mm (19-24"))	1
3	SPTRC2141	(Seat Rail Kit SD 380-480mm, cpl.)	1
4	SPTRC2142	(Seat Rail Kit SD 480-580mm, cpl.)	1
5	SPTRC2090	(End cap kit (2 pcs))	1

(Seat Cushions "Matrx Vi")

1588153.iss



Matrx Vi



		mm									
		405	430	455	480	505	530	555	580	610	
mm	380	1	4	4	8	8	-	-	-	-	
	405	1	4	4	8	8	-	-	-	-	
	430	1	4	4	8	8	-	-	-	-	
	455	2	5	5	9	9	12	12	15	15	
	480	2	5	5	9	9	12	12	15	15	
	505	3	6	6	10	10	13	13	16	16	
	530	3	6	6	10	10	13	13	16	16	
	555	-	7	7	11	11	14	14	17	17	
	580	-	7	7	11	11	14	14	17	17	


1	1579786	(Matrx Flo-Tech VI cushions 410x410 mm)	1
2	1579787	(Matrx Flo-Tech VI cushions 410x460 mm)	1
3	1579788	(Matrx Flo-Tech VI cushions 410x510 mm)	1
4	1579920	(Matrx Flo-Tech VI cushions 460x410 mm)	1
5	1579921	(Matrx Flo-Tech VI cushions 460x460 mm)	1
6	1579922	(Matrx Flo-Tech VI cushions 460x510 mm)	1
7	1584183	(Matrx Vi 18x22")	1
8	1579923	(Matrx Flo-Tech VI cushions 510x410 mm)	1
9	1579924	(Matrx Flo-Tech VI cushions 510x460 mm)	1
10	1579925	(Matrx Flo-Tech VI cushions 510x510 mm)	1
11	1584184	(Matrx Vi 20x22")	1
12	1579926	(Matrx Flo-Tech VI cushions 560x460 mm)	1
13	1579927	(Matrx Flo-Tech VI cushions 560x510 mm)	1
14	1584247	(Matrx Vi 22x22")	1
15	1584384	(Matrx Vi 24x18")	1
16	1584417	(Matrx Vi 24x20")	1
17	1584418	(Matrx Vi 24x22")	1

(Seat Cushions "Matrix Flo-Tech")



1588508.jpg

Matrix Flo-Tech Lite

							
		405	430	455	480	505	
Matrix Flo-Tech Lite		380	1	1	5	5	9
		405	2	2	6	6	9
		430	2	2	6	6	10
		455	3	3	7	7	10
		480	3	3	7	7	-
		505	4	4	8	8	-
		530	4	4	8	8	-




Matrix Flo-Tech Contour

											
		405	430	455	480	505	530	555	580	610	
Matrix Flo-Tech Contour		380	11	11	15	15	19	-	-	-	-
		405	12	12	16	16	19	-	-	-	-
		430	12	12	16	16	20	20	24	24	27
		455	13	13	17	17	21	21	24	24	27
		480	13	13	17	17	21	21	25	25	28
		505	14	14	18	18	22	22	25	25	28
		530	14	14	18	18	22	22	26	26	29
		555	-	-	-	-	23	23	26	26	29

1	FLB17/15	(Flo-Tech Lite 17X15")	1
2	FLB17/16	(Flo-Tech Lite 17X16")	1
3	FLB17/18	(Flo-Tech Lite 17X18")	1
4	FLB17/20	(Flo-Tech Lite 17X20")	1
5	FLB19/15	(Flo-Tech Lite 19x15")	1
6	FLB19/16	(Flo-Tech Lite 19x16")	1
7	FLB19/18	(Flo-Tech Lite 19x18")	1
8	SPFLB19/20	SETEPUTE SB48 SD51 SORT STORM4	1
9	SP1584714	(Flo-Tech Lite 20x15")	1
10	SP1584713	(Flo-Tech Lite 20x17")	1
11	FCB17/15	(Flo-Tech Contour 17x15")	1
12	SPFCB17_16	(Flo-Tech Contour 17x16")	1
13	SPFCB17_18	(Flo-Tech Contour 17x18")	1
14	SPFCB17_20	(Flo-Tech Contour 17x20")	1
15	SPFCB19_15	(Flo-Tech Contour 19x15")	1
16	SPFCB19_16	(Flo-Tech Contour 19x16")	1
17	SPFCB19_18	(Flo-Tech Contour 19x18")	1
18	SPFCB19_20	(Flo-Tech Contour 19x20")	1
19	SPFCB20_15	(Flo-Tech Contour 20x15")	1
20	SPFCB21/17	SETEPUTE SB53 SD43 SORT STORM4	1
21	SPFCB21/18	SETEPUTE SB53 SD46 SORT STORM4	1
22	SPFCB21/20	SETEPUTE SB53 SD51 SORT STORM4	1
23	SP1584739	(Flo-Tech Contour 21x22")	1
24	SP1584740	(Flo-Tech Contour 23x17")	1
25	SP1584741	(Flo-Tech Contour 23x19")	1
26	SP1584742	(Flo-Tech Contour 23x21")	1
27	SP1584743	(Flo-Tech Contour 24x17")	1
28	SP1584770	(Flo-Tech Contour 24x19")	1
29	SP1584771	(Flo-Tech Contour 24x21")	1

(Seat Cushions "Matrx Libra")


1603355.iss



Matrx Libra

		mm						
		405	430	455	480	505	530	555
mm	405	1	4	7	10	13	-	-
	430	1	4	7	10	13	-	-
	455	2	5	8	11	14	16	19
	480	2	5	8	11	14	16	19
	505	3	6	9	12	15	17	20
	530	3	6	9	12	15	17	20
	555	-	-	-	-	-	18	21
	580	-	-	-	-	-	18	21


1	LC1616	(Matrx Libra cushions 410x410 mm)	1
2	LC1618	(Matrx Libra cushions 410x460 mm)	1
3	LC1620	(Matrx Libra cushions 410x510 mm)	1
4	LC1716	(Matrx Libra cushions 430x410 mm)	1
5	LC1718	(Matrx Libra cushions 430x460 mm)	1
6	LC1720	(Matrx Libra cushions 430x510 mm)	1
7	LC1816	(Matrx Libra cushions 460x410 mm)	1
8	LC1818	(Matrx Libra cushions 460x460 mm)	1
9	LC1820	(Matrx Libra cushions 380x510 mm)	1
10	LC1916	(Matrx Libra cushions 480x410 mm)	1
11	LC1918	(Matrx Libra cushions 480x460 mm)	1
12	LC1920	(Matrx Libra cushions 480x510 mm)	1
13	LC2016	(Matrx Libra cushions 510x410 mm)	1
14	LC2018	(Matrx Libra cushions 510x460 mm)	1
15	LC2020	(Matrx Libra cushions 510x510 mm)	1
16	LC2118	(Seat Cushion Matrx Libra 21x18")	1
17	LC2120	(Seat Cushion Matrx Libra 21x20")	1
18	LC2122	(Seat Cushion Matrx Libra 21x22")	1
19	LC2218	(Matrx Libra cushions 560x460 mm)	1
20	LC2220	(Matrx Libra cushions 560x510 mm)	1
21	LC2222	(Seat Cushion Matrx Libra 22x22")	1



1	-	(Seat column / -raiser Storm4) [1025/1026/1030/1040]	1
2	-	(Seat column / -raiser Storm4 X-Plore) [1025/1026/1030/1040]	1
3	-	(Seat Height Adapter Storm4) [1040]	1
4	-	(Seat Height Adapter Storm4 X-Plore) [1040]	1
5	-	(Seat Tilt) [1002][1014]	1
6	-	(Seat Frame)	1
7	-	(Seat Plate)	1
8	-	(Recaro Seat)	1


Invacare® Storm4 | Sete
(Seat Frame)

1524555X.dwg



1	SP1536856	FRONTRAMME SETE RECARO STORM4	1
2	SP1542966	SETERAMME BAKRE AJ/RECAR STORM4	1

1573206.iso




1 SP1530519

SETEPLATE RECARO STORM4


1

(Recaro Seat)



1422546.ISO


i	-	(Recaro Mobility 24V is no longer available. Please pay attention to the change list!)	1
i	-	(Change list Recaro 24V)	1
	F367002	SÄKERHETSÅLÅTE FÖRLÅNG	1
	Z163001A	HOFTEBELTE MED HURTIÅLÅS RECARO SETE STORM3	1
	SP1423173	(RECARO SEAT MOUNTING KIT)	1
1	SP1552597	(Recaroseat Mobilty, 12V, cpl.)	1
1	SP1422546	SITTEENHET RECARO N-JOY SB49-53 RH77-83 EL	1
2	N89999	HODESTÅTTE RECARO, STOR	1
3	SP1423177	SKYVEHÅNDTAK RECARO	1
4	1435421	HÅNDTAK Å22	1
5	SP1552596	RYGG RECARO ELEKTRISK RECARO <i>Recaro Mobility 24V is no longer available. Please pay attention to the compatibility chart!</i>	1
5	SP1552781	(Backrest Recaro Mobility, mechanic)	1
6	SP1419207	SETEDEL RECARO STORM4	1



1	-	(Seat column / -raiser Storm4) [1025/1026/1030/1040]	1
2	-	(Seat column / -raiser Storm4 X-Plore) [1025/1026/1030/1040]	1
3	-	(Seat Height Adapter Storm4) [1040]	1
4	-	(Seat Height Adapter Storm4 X-Plore) [1040]	1
5	-	(Seat Tilt) [1002][1014]	1
6	-	(Seat Plate)	1
7	-	(Seat Frame)	1
8	-	(Seat Cushions Optimist Seat)	1


(Seat Plate)

1573207.iso




1	SP1524647	(Seatplate Arne Jensen 38cm)	1
2	SP1524643	(Seatplate Arne Jensen 45cm)	1
3	SP1524646	SETEPLATE AJ SB50 OPTIMIST SETE	1

1573478.iso



1	SP1428490	(Seat support Optimist Seat, SW 38cm)	1
2	SP1428110	KIT OPTIMIST SETE INKL AKTUATOR RYGG STORM3	1
3	SP1428491	SETE BUNN OPTIMIST SB50 OPTIMIST	1
4	1419417	AKTUATOR RYGG OPTIMIST STORM/TWISTER	1

(Seat Cushions Optimist Seat)




1	1419632	(Seat Cushion positive, SW 38cm, black)	1
2	5313953	OPTIMIST SETEPUTE SB 45 VOKSEN POSITIV	1
3	5314038	OPTIMIST SETEPUTE SB 50 VOKSEN POSITIV	1
4	1419633	(Seat Cushion negative, SW 38cm, black)	1
5	5313958	OPTIMIST SETEPUTE SB 45 VOKSEN NEGATIV	1
6	5314041	(Seat Cushion negative, SW 50cm, black)	1
7	1419634	(Seat Cushion paraplegic, SW 38cm, black)	1
8	5313962	OPTIMIST SETEPUTE SB 45 VOKSEN PARAPLEGI	1
9	5314044	OPTIMIST SETEPUTE SB 50 VOKSEN PARAPLEGI	1


1	-	(EasyAdapt Seat => 2017-08-10) [1180]	1
2	-	(Recaro Seat) [1152]	1
3	-	(Optimist Seat => 2018-12-31) [1171]	1
4	-	(Modulite Seat) [2100]	1

(EasyAdapt Seat => 2017-08-10)

1503147.kso




1	-	(Seat column / Seatriser) <i>[1025/1026/1030/1040]</i>	1
2	-	(Seat Height Adapter) <i>[1040]</i>	1
3	-	(Seat Tilt) <i>[1002][1014]</i>	1
4	-	(Seat Frame)	1
5	-	(Seat Plate)	1
6	-	(Adapter Kit for Universal Seat) <i>[1064]</i>	1
7	-	(Seat Cushion)	1



1	SP1556263	SETELØFT FRA SEPT 2013 STORM 4 KPL	1
-	SP1568361	(Mounting Kit lifter for seat bracket)	1
2	SP1548536	(Seat column, fixed, height 45cm and 48cm, adjustable)	1
3	SP1579553	(Slow Stop cable lifter GGB2)	1
4	SP1579637	(Lifter cable LINAK)	1

(Seat Height Adapter)


1583151.no




1 SP1555071

SETEHØYDE ADAPTER STORM4 FRA 15/5-2012


1




1	SP1551991	TILT BÆRER STORM4 [1002][1014]	1
2	SP1559711	TILT KPL U BÆRER STORM4 [1002][1014]	1
3	SP1551010	(Tilt spindle, cpl.) [1002]	1
4	SP1551009	(Mounting-Kit)	1
5	SP1588045	(Tilt Motor Reac 23.3 (incl. Hardware)) [1014]	1
6	SP1532481	MONTERINGSKIT TILT STORM4	1
7	SP1580197	(Relais cover, cpl.)	1
8	SP1603273	(Cable Support Sheets, cpl.) [1030][1040]	1




1	-	(Seat column / Seatriser) [1025/1026/1030/1040]	1
2	-	(Seat Height Adapter) [1040]	1
3	-	(Seat Tilt) [1002][1014]	1
4	-	(Seat Frame)	1
5	-	(Seat Plate)	1
6	-	(Recaro Seat)	1





1	-	(Seat column / Seatriser) [1025/1026/1030/1040]	1
2	-	(Seat Height Adapter) [1040]	1
3	-	(Seat Tilt) [1002][1014]	1
4	-	(Seat Plate)	1
5	-	(Seat Frame)	1
6	-	(Seat Cushions Optimist Seat)	1



1	-	(Seat column / Seatriser) [1025/1026/1030/1040]	1
2	-	(Seat Height Adapter) [1040]	1
3	-	(Seat Tilt) [1002][1014]	1
4	-	(One Piece Seat Plate) [2110]	1
5	-	(Telescopic Seat Frame and Plate) [2111]	1
6	-	(Telescopic Seat Frame and Sling Seat) [2112]	1
7	-	(Telescopic Seat Frame and Universal Adapter) [2109]	1
8	-	(Seat Cushions "Standard" & "Matrx Flo-Tech") [2141/2142/2143]	1
8	-	(Seat Cushions "Matrx Libra") [2147]	1
8	-	(Seat Cushions "Matrx Vi") [2144]	1



1	SP1551991	TILT BÆRER STORM4 [1002][1014]	1
2	SP1559711	TILT KPL U BÆRER STORM4 [1002][1014]	1
3	SP1551010	(Tilt spindle, cpl.) [1002]	1
4	SP1551009	(Mounting-Kit)	1
5	SP1588045	(Tilt Motor Reac 23.3 (incl. Hardware)) [1014]	1
6	SP1532481	MONTERINGSKIT TILT STORM4	1
7	SP1573217	SETE SUPPORT MODULITE [2001]	1
8	SP1580196	(Seat Height Adapter GGB2) [2001]+[1290]/[1275]/[1483]	1
9	SP1580197	(Relais cover, cpl.)	1
10	SP1603273	(Cable Support Sheets, cpl.) [1030][1040]	1




1	-	(Storm4)	0
2	-	(Storm4 TTPlus)	0


-	(Backrest for EasyAdapt Seat => 2017-08-10)	1
-	(Backrest for Modulite Seat)	1
-	(Backrest for Ultra Low Maxx Seat)	1
-	(Backrest for Optimist Seat => 2018-12-31)	1

(Backrest for EasyAdapt Seat => 2017-08-10)

1573229.no




1	-	(Backrestadjusting)	1
2	-	(Backrestcover for EasyAdapt)	1
3	-	(Backplate EasyAdapt)	1
3	-	(Backplate EasyAdapt for DAHL Docking Station)	1
4	-	(Backrest Cushions) [1985/1986/1987/1988]	1
5	-	(Push Handle EasyAdapt) [1307]	1
6	-	(Tension adjustable Sling Back) [1185/1186]	1
7	-	("Universal"-Backrest frame) [0904]	1




-	SP1529390	KABEL SETELØFT STORM4	1
-	SP1529529	ROMSENSOR MED KONTAKT STORM4	1
-	SP1562940	(Adaptercable Shark to reverse actuator) <i>used for TT-Plus</i>	1
-	SP1440726	(ACTUATORCABLE RECLINER ACS2 / SHARK) <i>used for TT-Plus</i>	1
1	-	(No longer available.)	1
2	-	(No longer available.)	1
3	-	(No longer available.)	1
4	-	(No longer available.)	1
5	1548506	(Microswitch, Short lever arm)	1
5	SP1526445	KABEL RYGGJUSTERINGBRYTE STORM4	1
6	SP1530520	MOTOR RYGG STORM4	1
7	SP1420483	PLASTRATT STORM4	1
8	SP1421696	REFLEX STORM4	1
9	-	(No longer available.)	1
10	-	(No longer available.)	1
11	SP1535255	MONTERINGSSATS ARML OPPF STORM4 EASY ADAPT	1
12	SP1537863	BRO FOR LENGDEKOMP STORM4	1
13	SP1532897	(Backrest sliding plate)	1
14	SP1516962	MOTOR AKSLING RYGG EA STORM4	1

(Backrestcover for EasyAdapt)




1	SP1528487	RYGGDEKSEL EASYADAPT STORM4	1
2	SP1537706	DEKSEL RYGG M ARMLEN STORM4	1



1	SP1530525	RYGGPL EA SENT STORM4	1
2	SP1530526	RYGGPL EA HØ STORM4	1
3	SP1530527	RYGGPL EA VE STORM4	1

(Backplate EasyAdapt for DAHL Docking Station)

1579892 iso






1	SP1541384	RYGGPLATE MODULITE DAHL STORM4 EASYADAPT KPL	1
2	SP1579843	RYGGPLATE MODULITE DAHL STORM4 EASYADAPT MIDT	1
3	SP1579844	RYGGPLATE MODULITE DAHL STORM4 EASYADAPT VE	1
4	SP1579845	RYGGPLATE MODULITE DAHL STORM4 EASYADAPT HØ	1

-	(Backrest Cushions Laguna & Onyx)	1
-	(Matrix backrest)	1

(Backrest Cushions Laguna & Onyx)




Laguna



			
1	380mm		
2	430mm		grey plush
3	480mm		
4	530mm		
5	380mm	⇒ 30.11.14	01.12.14 ⇒ dartex grey dartex black
6	430mm		
7	480mm		
8	530mm		
9	380mm		black
10	430mm		
11	480mm		
12	530mm		
13	380mm		black orange seam
14	430mm		
15	480mm		
16	530mm		


Onyx



				
17	400mm	500mm		
18	400mm	630mm		
19	430mm	500mm		black
20	430mm	630mm		
21	450mm	500mm		
22	450mm	630mm		
23	480mm	500mm		
24	480mm	630mm		

1526905X/ISO

1	1444339	RYGGPUTE LAGUNA TR32 372X550 BELLIS		1
2	-	(No longer available - please choose Dartex black.)		1
3	1444340	RYGGPUTE LAGUNA 46 TR32		1
4	1444341	RYGGPUTE LAGUNA 50 TR32		1
5	-	(No longer available - please choose Dartex black.)	→ 30.11.2014	1
5	SP1427118-TR26	RYGGPUTE SB39 LAGUNA 372 X550 TR26	01.12.2014 →	1
6	-	(No longer available - please choose Dartex black.)	→ 30.11.2014	1
6	SP1427119-TR26	RYGGPUTE LAGUNA 402X550 TR26	01.12.2014 →	1
7	-	(No longer available - please choose Dartex black.)	→ 30.11.2014	1
7	SP1427120-TR26	RYGGPUTE LAGUNA 452X550 SORT	01.12.2014 →	1
8	-	(No longer available - please choose Dartex black.)	→ 30.11.2014	1
8	SP1427121-TR26	RYGGPUTE S LAGUNA 492X550 TR26	01.12.2014 →	1
9	SP1529585	RYGGPUTE EASYADAPT SB38 STORM4 M SORT		1
10	SP1529586	RYGGPUTE EASYADAPT SB43 STORM4 STD SORT		1
11	SP1529587	RYGGPUTE EASYADAPT SB48 STORM4 SORT		1
12	SP1529596	RYGGPUTE EASYADAPT SB53 STORM4 M SORT		1
13	SP1529538	RYGGPUTE LAGUNA 37CM SORT M ORANG SØM		1
14	SP1529539	RYGGPUTE LAGUNA 43CM SORT M ORANG SØM		1
15	SP1529540	RYGGPUTE LAGUNA 48CM SORT M ORANG SØM		1
16	SP1529545	RYGGPUTE LAGUNA 50CM SORT M ORANG SØM		1
17	SP1532061	(ONYX Back 40x50cm (16" x 20"))		1
18	SP1532062	(ONYX Back 40x63cm (16" x 25"))		1
19	SP1532064	(ONYX Back 43x50cm (17" x 20"))		1
20	SP1532065	(ONYX Back 43x63cm (17" x 25"))		1
21	SP1532118	(ONYX Back 45x50cm (18" x 20"))		1
22	SP1532194	(ONYX Back 45x63cm (18" x 25"))		1
23	SP1532211	(ONYX Back 48x50cm (19" x 20"))		1
24	SP1532212	(ONYX Back 48x63cm (19" x 25"))		1




1	-	(Matrix Posture Back)	1
2	-	(Matrix Elite (Full) Back)	1
3	-	(Mounting Hardware for Matrix backrests)	1

(Matrx Posture Back)


158885 no

Posture Back



1	1544896	Ryggtrekk Posture Back s b40,5 rh40,5	1
2	1544899	Ryggtrekk Posture Back t rekk sb43 rh40,5	1
3	1544995	Ryggtrekk Posture Back s b45,5 rh40,5	1
4	1544897	Ryggtrekk Posture Back s b40,5 rh50,5	1
5	1544982	Ryggtrekk Posture Back t rekk sb43 rh50,5	1
6	1544996	Ryggtrekk Posture Back s b45,5 rh50,5	1
7	1515727	RYGG Posture Back 40.5w x 40.5h No H/ware <i>[2174]</i>	1
8	1515729	Posture Back 43w x 40.5 No H/ware <i>[2174]</i>	1
9	1515731	RYGG Posture Back 45.5w x 40.5h No H/ware <i>[2174]</i>	1
10	1515728	RYGG Posture Back 40.5w x 50.5h No H/ware <i>[2174]</i>	1
11	1515730	Posture Back 43w x 50.5h No H/ware <i>[2175]</i>	1
12	1515732	Rygg Posture Back b45,5 h50,5 ufeste <i>[2175]</i>	1


Elite Full Back




1	SP1532094	(Matrx Elite Back 15x20") [2176]	1
2	SP1517755	(Matrx Elite Back 17x20", cpl.) [2176]	1
3	SP1517756	(Matrx Elite Back 18x20", cpl.) [2176]	1

(Mounting Hardware for Matrix backrests)

163/807/00






1 SP1529001

KJØREHÅNDTAK STORM4

1


(Tension adjustable Sling Back)

1533409X.dwg



1	SP1531825	RYGG VELCRO SB43 STORM4	1
1	SP1533409	RYGG VELCRO SB53 STORM4	1
2	SP1531830	RYGGPOLSTER SB43 STORM4	1
2	SP1533405	RYGGPOLSTER SB53 STORM4	1
3	SP1537258	VELCROBÅND SB43CM STORM4	1
3	SP1537296	VELCROBÅND SB53CM STORM4	1
4	SP1537297	RYGGBÆRER SB43 STORM4	1
4	SP1537298	RYGGBÆRER SB53 STORM4	1

1531052X.dwg




1 SP1531052

RAMME RYGG UNIVERSAL EL-STOLER PORTA


1

(Backrest for Modulite Seat)

1566095.no




1	-	(Support for Backrest Frame) [0975+2161]	1
2	-	(Backrest Frame) [2161]	1
3	-	(Sling Back, tension adjustable) [2171]	1
4	-	(Backrest support)	1
5	-	(Backposts) [2162/2163]	1
6	-	(Connection tube) [1307]	1
7	-	(Back Recliner / Stabilization bracket) [0977]	1
8	-	(Backrest Cushions) [2172/2173/2174/2175/2176/2177/2178]	1
9	-	(Flex3 back, tension adjustable) [2179]	1



1571833.iso

1	SP1571542	(Backrest Support RH)	1
2	SP1571553	(Backrest Support LH)	1
3	5824117	STJERNEGREP FESTE NY MANØVERBOKS MOOVER <i>(only swedish homologation)</i>	2


Invacare® Storm4 | Rygg
(Backrest Frame)



1571831.iso

1	SP1567810	RYGGRØR TDX MODULITE	1
2	0021-70002	FINGERSKRUE M8 ITALIA RYGGRØR	1

(Sling Back, tension adjustable)




1571829.iso


1	SP1567529	RYGG VELCRO SB38 RH48 MODULITE	1
2	SP1567465	VELCRO POLSTER HØ48 SB43 MODULITE	1
3	SP1567530	RYGG VELCRO SB48 RH48 MODULITE	1
4	SP1567531	(Sling Back, tension adjustable & upholstery, BW 530mm / BH 480mm, cpl.)	1
5	SP1567627	(Sling Back, tension adjustable & upholstery, BW 380mm / BH 540mm. cpl.)	1
6	SP1567746	RYGG VELCRO SB43 RH43 MODULITE	1
7	SP1567747	RYGG VELCRO SB48 RH54 MODULITE	1
8	SP1567748	VELCRORYGG KPL RH54 RB53 MODULITE	1
9	SP1563563	(Backrest Upholstery BW 380)	1
10	SP1563653	(Backrest Upholstery BW 430)	1
11	SP1563654	(Backrest Upholstery BW 480)	1
12	SP1563655	(Backrest Upholstery BW 530)	1
13	SP1571492	(Sling back)	1
14	SP1571493	SLING BACK TDX SP2	1
15	SP1571494	(Sling back)	1
16	SP1571510	(Sling back)	1

(Backrest support)

1572928 Rev




-	-	(Design change on 2016-05-02 - We suggest to change older supports on both sides.)		1
1	SP1582709	FESTE RYGGRØR HØ MODULITE <i>Compatible with Flip-Up & Following Armrest [2193/2192/2194] - replaces SP1571560</i>		1
2	SP1582708	FESTE RYGGRØR VE MODULITE <i>Compatible with Flip-Up & Following Armrest [2193/2192/2194] - replaces SP1571561</i>		1
3	SP1582952	RYGGFESTE 53 HØ MODULITE <i>Compatible with Flip-Up & Following Armrest [2193/2192/2194] - replaces SP1571568</i>		1
4	SP1582951	RYGGFESTE 53 VE MODULITE <i>Compatible with Flip-Up & Following Armrest [2193/2192/2194] - replaces SP1571569</i>		1
5	SP1571676	(Covers)		1
6	SP1579877	(Cross Handle M6) <i>(Swedish Homologation only!)</i>	→ 31.07.2015	2
6	SP1589324	(Cross handle M8) <i>(Swedish Homologation only!)</i>	01.08.2015 →	2
7	SP1590656	(Hole plug (for powered backrest only))	01.07.2016 →	1



1	SP1587177	RYGGRØR HØYDE 48 HØ MODULITE <i>Compatible with Flip-Up & Following Armrest [2193/2192/2194] - replaces SP1563102</i>	1
2	SP1587184	RYGGRØR HØ 54 MODULITE <i>Compatible with Flip-Up & Following Armrest [2193/2192/2194] - replaces SP1563103</i>	1
3	SP1587176	RYGGRØR HØYDE 48 VE MODULITE <i>Compatible with Flip-Up & Following Armrest [2193/2192/2194] - replaces SP1560576</i>	1
4	SP1587183	RYGGRØR VE 54 MODULITE <i>Compatible with Flip-Up & Following Armrest [2193/2192/2194] - replaces SP1559578</i>	1
5	SP1587163	(Backpost man. BH 480mm - RH) <i>Compatible with Flip-Up & Following Armrest [2193/2192/2194] - replaces SP1572489</i>	1
6	SP1587182	(Backpost man. BH 540mm - RH) <i>Compatible with Flip-Up & Following Armrest [2193/2192/2194] - replaces SP1572490</i>	1
7	SP1587162	(Backpost man. BH 480mm - LH) <i>Compatible with Flip-Up & Following Armrest [2193/2192/2194] - replaces SP1572491</i>	1
8	SP1587178	(Backpost man. BH 540mm - LH) <i>Compatible with Flip-Up & Following Armrest [2193/2192/2194] - replaces SP1572492</i>	1
9	SP1586122	HÅNDTAK HOLK MODULITE <i>replaces item SP1563651</i>	1
10	0021-70002	FINGERSKRUE M8 ITALIA RYGGRØR	1
11	1567999	(Bolt M8x16 8.8 black DIN912)	1
12	SP5816477	SRØRBEN TIL RØR SORT MOOVER 85 (GL25)	2
13	SP1606879	(Backpost el. BH 480mm - RH w/o Pushhandle)	1
14	SP1606893	(Backpost el. BH 540mm - RH w/o Pushhandle)	1
15	SP1606880	(Backpost el. BH 480mm - LH w/o Pushhandle)	1
16	SP1606894	(Backpost el. BH 540mm - LH w/o Pushhandle)	1
17	SP1606877	(Backpost man. BH 480mm - RH w/o Pushhandle)	1
18	SP1606881	(Backpost man. BH 540mm - RH w/o Pushhandle)	1
19	SP1606878	(Backpost man. BH 480mm - LH w/o Pushhandle)	1
20	SP1606882	(Backpost man. BH 540mm - LH w/o Pushhandle)	1

(Connection tube)




1571835.ISO

1	SP1570370	KJØREBØYLE MODULITE SB38 MODULITE	1
2	SP1570371	KJØREBØYLE SB43 MODULITE STORM4	1
3	SP1570372	KJØREBØYLE SB48 LS1 MODULITE	1
4	SP1570373	KJØREBØYLE SB53 MODULITE STORM4 TDXSP2	1

(Back Recliner / Stabilization bracket)


1571097.kim




1	SP1567169	ACTUATOR BRAKETT MODULITE	1
2	SP1567168	BRAKETT STABILISATOR RYG STORM4 MODULITE	1
6	SP1607105	(LiNX Actuator cable LA23 - 950mm)	1
3	SP1571702	(Cover)	1
4	SP1571704	(Locking bolt)	1
5	SP1557820	ACTUATOR RYGG MODULITE	1
6	SP1572166	KABEL ACTUATOR MODULITE	1
6	SP1607124	(LiNX Actuator cable LA23 - 1500mm)	1

(Backrest Cushions)

1588915 no




1	-	(Flex 3 / Laguna)		1
2	-	(Matrx Posture Back)		1
3	-	(Matrx Elite T/R Back)	01.07.2015 →	1
4	-	(Matrx Elite Deep Back)	01.07.2015 →	1
5	-	(Matrx Elite (Full) Back)		1
6	-	(Mounting Hardware for Matrx backrests)		1



Flex 3

Laguna




1	SP1568489	FLEX3 RYGG TIL MODULITE STORM4 TDX [2172]/[2173]	1
2	SP1572321	(Clamps (cpl. kit))	1
3	1497411	RATT M12 LANG KPL NAKKE/HODESTØTTE	1
4	SP1529585	RYGGPUTE EASYADAPT SB38 STORM4 M SORT [2173+2151+2121]	1
5	SP1529586	RYGGPUTE EASYADAPT SB43 STORM4 STD SORT [2173+2151+2122]	1
6	SP1529587	RYGGPUTE EASYADAPT SB48 STORM4 SORT [2173+2151+2123]	1
7	SP1529596	RYGGPUTE EASYADAPT SB53 STORM4 M SORT [2173+2151+2124]	1
8	SP1427118-TR26	RYGGPUTE SB39 LAGUNA 372 X550 TR26 [2173+2152+2121]	1
9	SP1427119-TR26	RYGGPUTE LAGUNA 402X550 TR26 [2173+2152+2122]	1
10	SP1427120-TR26	RYGGPUTE LAGUNA 452X550 SORT [2173+2152+2123]	1
11	SP1427121-TR26	RYGGPUTE S LAGUNA 492X550 TR26 [2173+2152+2124]	1
12	SP1529538-TR35	(Backrest Cushion Laguna, Plush black - SW 380mm) [2173+2154+2121]	1
13	SP1529539-TR35	(Backrest Cushion Laguna, Plush black - SW 430mm) [2173+2154+2122]	1
14	SP1529540-TR35	(Backrest Cushion Laguna, Plush black - SW 480mm) [2173+2154+2123]	1
15	SP1529545-TR35	(Backrest Cushion Laguna, Plush black - SW 530mm) [2173+2154+2124]	1
16	SP1529538	RYGGPUTE LAGUNA 37CM SORT M ORANG SØM [2173+2153+2121]	1
17	SP1529539	RYGGPUTE LAGUNA 43CM SORT M ORANG SØM [2173+2153+2122]	1
18	SP1529540	RYGGPUTE LAGUNA 48CM SORT M ORANG SØM [2173+2153+2123]	1
19	SP1529545	RYGGPUTE LAGUNA 50CM SORT M ORANG SØM [2173+2153+2124]	1

(Flex 3 / Laguna)

20	SP1427118-TR36	(CUSHION, BACKREST LAGUNA 372X550, SPACER BLACK) <i>[2173+2155+2121]</i>	1
21	SP1427119-TR36	Ryggpute Laguna b40,2 h5 5 Spacer sort <i>[2173+2155+2122]</i>	1
22	SP1427120-TR36	Ryggpute Laguna b45,2 h5 5 Spacer sort <i>[2173+2155+2123]</i>	1
23	SP1427121-TR36	Ryggpute Laguna b49,2 h5 5 Spacer sort <i>[2173+2155+2124]</i>	1

1580006 no

Elite T/R Back




1	1580110	Ryggtrekk PB Elite TiltR ecline sb48 rh46	1
2	1580297	Ryggtrekk PB Elite TiltR ecline sb48 rh51	1
3	1580298	Ryggtrekk PB Elite TiltR ecline sb51 rh41	1
4	SP1580029	(Elite TR Back 16x20", cpl.) [2177+2121]	1
5	SP1580030	(Elite TR Back 17x20", cpl.) [2177+2122]	1
6	SP1580031	(Elite TR Back 18x20", cpl.) [2177+2123]	1

(Matrx Elite Deep Back)

158913.no


Elite Deep Back



2	1544077	Ryggtrekk Elite dyp Back sb43 rh50,5	1
3	1544081	Ryggtrekk Elite dyp Back sb45,5 rh50,5	1
4	1543511	RYGG ELITE DEEP SB38 X50,5 UTEN FESTE [2178+2121]	1
5	1543519	RYGG ELITE DEEP SB43RH50 MATRIX [2178+2122]	1
6	1543523	RYGG ELITE DEEP SB46 H50 MATRIX [2178+2123]	1


(Flex3 back, tension adjustable)

1602719.iss




1	SP1579675	RAMME STILLBAR RYGG SB39 AZALEA	1
2	SP1579692	RAMME STILLBAR RYGG SB44 AZALEA	1
3	SP1579693	RAMME STILLBAR RYGG SB49 AZALEA	1
4	SP1543667	Ryggtrekk Lateral sb44 D artex sort	1
5	SP1543669	Ryggtrekk Lateral sb49 D artex sort	1
6	SP1543671	Ryggtrekk LATERAL sb55 D artex sort	1
7	SP1580557	Festesett feste av rygg mot ryggamme (sett à 4)	1

(Backrest for Ultra Low Maxx Seat)



1588115.iso



-	(Recliner)	1
-	(High Back)	1
-	(Matrix Backrests)	1




-	SPTRE0322	(Actuator Cable, Backrest Recliner)	1
1	SPTRC2191	(Backrest Recliner with Hardware)	1
2	SPTRC2138	(Backrest shroud with mounting hardware)	1
3	SPTRC1926	(Extended Sheer Reduction (ESR))	1
		[2060/2063/2066]	

(High Back)

High Back




1	-	(Backrest Hardware)	1
2	-	(Backrest Cushions)	1



1	SPTRC2136	(High Back Interface Kit SW 405-505mm (16-20"))	1
2	SPTRC2137	(High Back Interface Kit SW 505-610mm (19-24"))	1
3	SPTRC2188	(Push Bars SW 16-20")	1
4	SPTRC2190	(Push Bars SW 20-24")	1
5	SPTRX8497	(High Back Center Plate SW 405-505mm (16-20"))	1
7	SPTRX7883	(Linear Bearing)	4

(Backrest Cushions)

158201 no




1	SP1583769	(High Back 405x405mm (16x16"), cpl.)	1
2	SP1583772	(High Back 430x405mm (17x16"), cpl.)	1
3	SP1583774	(High Back 455x405mm (18x16"), cpl.)	1
4	SP1583857	(High Back 480x405mm (19x16"), cpl.)	1
5	SP1583930	(High Back 505x405mm (20x16"), cpl.)	1
6	SP1583932	(High Back 530x405mm (21x16"), cpl.)	1
7	SP1584053	(High Back 555x405mm (22x16"), cpl.)	1
8	SP1584055	(High Back 580x405mm (23x16"), cpl.)	1
9	SP1584081	(High Back 610x405mm (24x16"), cpl.)	1
10	SP1583771	(High Back 405x455mm (16x18"), cpl.)	1
11	SP1583773	(High Back 430x455mm (17x18"), cpl.)	1
12	SP1583775	(High Back 455x455mm (18x18"), cpl.)	1
13	SP1583858	(High Back 480x455mm (19x18"), cpl.)	1
14	SP1583931	(High Back 505x455mm (20x18"), cpl.)	1
15	SP1583934	(High Back 530x455mm (21x18"), cpl.)	1
16	SP1584054	(High Back 555x455mm (22x18"), cpl.)	1
17	SP1584056	(High Back 580x455mm (23x18"), cpl.)	1
18	SP1584118	(High Back 610x455mm (24x18"), cpl.)	1
19	SPTRS0425-1616	(Backrest Cushion - High Back 16x16")	1
20	SPTRS0425-1716	(Backrest Cushion - High Back 17x16")	1
21	SPTRS0425-1816	(Backrest Cushion - High Back 18x16")	1
22	SPTRS0425-1916	(Backrest Cushion - High Back 19x16")	1
23	SPTRS0425-2016	(Backrest Cushion - High Back 20x16")	1
24	SPTRS0425-2116	(Backrest Cushion - High Back 21x16")	1
25	SPTRS0425-2216	(Backrest Cushion - High Back 22x16")	1
26	SPTRS0425-2316	(Backrest Cushion - High Back 23x16")	1
27	SPTRS0425-2416	(Backrest Cushion - High Back 24x16")	1
28	SPTRS0425-1618	(Backrest Cushion - High Back 16x18")	1
29	SPTRS0425-1718	(Backrest Cushion - High Back 17x18")	1
30	SPTRS0425-1818	(Backrest Cushion - High Back 18x18")	1
31	SPTRS0425-1918	(Backrest Cushion - High Back 19x18")	1

32	SPTRS0425-2018	(Backrest Cushion - High Back 20x18")	1
33	SPTRS0425-2118	(Backrest Cushion - High Back 21x18")	1
34	SPTRS0425-2218	(Backrest Cushion - High Back 22x18")	1
35	SPTRS0425-2318	(Backrest Cushion - High Back 23x18")	1
36	SPTRS0425-2418	(Backrest Cushion - High Back 24x18")	1


(Matrx Backrests)

1588242 Rev

Matrx Back



1	-	(Backrest Hardware)	1
2	-	(PB Back)	1
3	-	(Elite TR/TR HD Back)	1
4	-	(Elite Deep Back)	1
-	-	(Privacy Flaps)	1




1	SPTRC1951	(Matrx Back Interface Kit SW 405-505mm (16-20"))	1
2	SPTRC2095	(Matrx Back Interface Kit SW 505-610mm (19-24"))	1
3	SPTRX8299	(Matrx Center Plate SW 405-505mm (16-20"))	1
4	SPTRX8779	(Matrx Center Plate SW 505-610mm (19-24"))	1
5	SPTRC2188	(Push Bars SW 16-20")	1
6	SPTRC2190	(Push Bars SW 20-24")	1
7	SPTRX7883	(Linear Bearing)	1

(PB Back)


158204 no

Posture Back




1	1515727	RYGG Posture Back 40.5w x 40.5h No H/ware	1
2	1515729	Posture Back 43w x 40.5 No H/ware	1
3	1515731	RYGG Posture Back 45.5w x 40.5h No H/ware	1
4	1515728	RYGG Posture Back 40.5w x 50.5h No H/ware	1
5	1515730	Posture Back 43w x 50.5h No H/ware	1
6	1515732	Rygg Posture Back b45,5 h50,5 ufeste	1
7	1544896	Ryggtrekk Posture Back s b40,5 rh40,5	1
8	1544899	Ryggtrekk Posture Back t rekk sb43 rh40,5	1
9	1544995	Ryggtrekk Posture Back s b45,5 rh40,5	1
10	1544897	Ryggtrekk Posture Back s b40,5 rh50,5	1
11	1544982	Ryggtrekk Posture Back t rekk sb43 rh50,5	1
12	1544996	Ryggtrekk Posture Back s b45,5 rh50,5	1

Elite T/R Back




1	SP1580029	(Elite TR Back 16x20", cpl.)	1
2	SP1580030	(Elite TR Back 17x20", cpl.)	1
3	SP1580031	(Elite TR Back 18x20", cpl.)	1
4	SP1580032	(Elite TR Back 19x20", cpl.)	1
5	SP1580033	ELITE TR RYGG 505x500mm 20x20"	1
6	SP1584160	(Elite TR Back HD 21x20", cpl.)	1
7	SP1584178	(Elite TR Back HD 22x20", cpl.)	1
8	SP1584180	(Elite TR Back HD 23x20", cpl.)	1
9	SP1584182	(Elite TR Back HD 24x20", cpl.)	1
10	1580110	Ryggtrekk PB Elite TiltR ecline sb48 rh46	1
11	1580297	Ryggtrekk PB Elite TiltR ecline sb48 rh51	1
12	1580298	Ryggtrekk PB Elite TiltR ecline sb51 rh41	1
13	1580299	Ryggtrekk PB Elite TiltR ecline sb51 rh46	1
14	1580300	Ryggtrekk PB Elite TiltR ecline sb51 rh51	1
15	SPTRS0375-2120	(Backrest Cushion - Elite TR HD Back 21x20")	1
16	SPTRS0375-2220	(Backrest Cushion - Elite TR HD Back 22x20")	1
17	SPTRS0375-2320	(Backrest Cushion - Elite TR HD Back 23x20")	1
18	SPTRS0375-2420	(Backrest Cushion - Elite TR HD Back 24x20")	1

Elite Deep Back



1	1543515	RYGG ELITE DEEP SB41 H51 MATRIX	1
2	1543519	RYGG ELITE DEEP SB43RH50 MATRIX	1
3	1543523	RYGG ELITE DEEP SB46 H50 MATRIX	1
4	1543987	RYGG Elite Deep Back 48w x 50.5hNo H/ware	1
5	1543991	RYGG Elite Deep SB50.5 H50.5h UTEN FESTE	1
6	1584177	RYGG Posture Back 530 x 500 Deep Back No H/Ware	1
7	1584179	(Elite Deep HD Back 22x20", cpl.)	1
8	1584181	(Elite Deep HD Back 23x20", cpl.)	1
9	1583933	(Elite Deep HD Back 24x20", cpl.)	1
10	1544073	Ryggtrekk Elite dyp Back sb40,5 rh50,5	1
11	1544077	Ryggtrekk Elite dyp Back sb43 rh50,5	1
12	1544081	Ryggtrekk Elite dyp Back sb45,5 rh50,5	1
13	1544085	Ryggtrekk Elite dyp Back sb48 rh50,5	1
14	1544124	Ryggtrekk Elite dyp Back sb50,5 rh50,5	1
15	SPTRS0344-2120	(Backrest Cushion - Elite Deep HD Back 21x20")	1
16	SPTRS0344-2220	(Backrest Cushion - Elite Deep HD Back 22x20")	1
17	SPTRS0344-2320	(Backrest Cushion - Elite Deep HD Back 23x20")	1
18	SPTRS0344-2420	(Backrest Cushion - Elite Deep HD Back 24x20")	1




405/430mm	①
455/480mm	②
505/530mm	③
555/580mm	④
605/630mm	⑤

1	1551800	Trekk intim strliten b33 I41		1
2	1551903	Trekk intim strmedium b3 8 I41		1
3	1551904	Trekk intim strstor b43 I41		1
4	SPPFLAP-XL	(Privacy Flap SW 555/580mm)	11.08.2017 →	1
5	SPPFLAP-XXL	(Privacy Flap SW 605/630mm)	11.08.2017 →	1

(Backrest for Optimist Seat => 2018-12-31)

1573230.no




1	-	(Optimist Seats)	1
2	-	(Push Handle Optimist Seat)	1

OPTIMIST 1			OPTIMIST 2			OPTIMIST 3			OPTIMIST 5		
1	38cm	46cm	5	38cm	46cm	11	45cm	60cm	13	45cm	70cm
2	45cm	46cm	6	38cm	46cm	12	50cm	70cm	14	50cm	70cm
3	45cm	60cm	7	45cm	46cm						
4	50cm	60cm	8	45cm	60cm						
			9	45cm	60cm						
			10	50cm	60cm						

1	1419630	OPTIMISTSETE 1 SB38 H46 UTEN AKTUATOR	1
2	5313888	(Backrest cushion Optimist Typ1, SW45cm)	1
3	5313889	OPTIMISTSETE 1 SB45 H60 UTEN AKTUATOR	1
4	1419625	OPTIMISTSETE 1 SB50 H46 UTEN AKTUATOR	1
5	1419631	(Backrest cushion Optimist Typ2, SW38cm)	1
6	5313867	OPTIMISTSETE 2 EL SELE SB38 H60 UTEN AKTUATOR	1
7	5313882	OPTIMISTSETE 2 SB45 H47 UTEN AKTUATOR	1
8	5313883	OPTIMISTSETE 2 SB45 H60 UTEN AKTUATOR	1
9	5313895	OPTIMISTSETE VOKSEN M EL SELE UTEN AKTUATOR	1
10	1419628	(Backrest cushion Optimist Typ2, SW50cm)	1
11	5313878	OPTIMISTSETE 3 SB45 H47 UTEN AKTUATOR STORM4	1
12	5313877	OPTIMISTSETE 3 SB50 H60 UTEN AKTUATOR	1
13	5313872	OPTIMISTSETE 5 SB45 H60 UTEN AKTUATOR	1
14	5313871	OPTIMISTSETE 5 SB50 H60 KPL	1

(Push Handle Optimist Seat)

1441556X.dwg



1 1441553


(Push Handle Arne Jensen)

1

-	(Backrest for EasyAdapt Seat => 2017-08-10)	1
-	(Backrest for Modulite Seat)	0
-	(Backrest for Optimist Seat => 2018-12-31)	1

(Backrest for Modulite Seat)


1566095.no



1	-	(Support for Backrest Frame) [0975+2161]	1
2	-	(Backrest Frame) [2161]	1
3	-	(Sling Back, tension adjustable) [2171]	1
4	-	(Backrest support)	1
5	-	(Backposts) [2162/2163]	1
6	-	(Connection tube) [1307]	1
7	-	(Backrest Recliner / Stabilisation bracket) [0977]	1
8	-	(Backrest Cushions) [2172/2173/2174/2175/2176/2177/2178]	1
9	-	(Flex3 back, tension adjustable) [2179]	1

(Backrest Recliner / Stabilisation bracket)



1583079.iss



1	SP1567169	ACTUATOR BRAKETT MODULITE	1
2	SP1567168	BRAKETT STABILISATOR RYG STORM4 MODULITE	1
3	SP1571702	(Cover)	1
4	SP1571704	(Locking bolt)	1
5	SP1557820	ACTUATOR RYGG MODULITE	1
6	SP1572166	KABEL ACTUATOR MODULITE	1
7	SP1579551	(Actuator Cable LA23 GGB2 (~850mm))	1

Benstøtter og holdere

ST46621-01



1	-	(Storm4)	0
2	-	(Storm4 TTPlus)	0

-	(Legrests for Ultra Low Maxx Seat)	1
-	(Legrest Support, Vari series 07/2005 ==>)	1
-	(Legrest Support, Vari series 07/2005 =>)	1
-	(Hanger Vari A (04/2009 =>))	1
-	(Hanger Vari F)	1
-	(Hanger Vari ADM man.)	1
-	(Hanger Vari ADE/ADP pow. (07/2005 =>))	1
-	(Central Legrest, manuell)	1
-	(Central Legrest, electric)	1
-	(Footrests Vari, Flip-Up (01.08.10 =>))	1
-	(Footrests Vari, Flip-Up & angle adjustable (01.08.10 =>))	1
-	(Support Tube with Footplate for central Legrest)	1
-	(Footboards)	1
-	(Lower Footrests)	1


(Legrests for Ultra Low Maxx Seat)

-	(Central Legrests)	1
-	(Vari Legrests)	1
-	(Pivot Plus Legrests)	1


-	(Supports, Central Legrests)	1
-	(Central Legrest - powered LNX)	1

Invacare® Storm4 | Benstøtter og holdere
(Supports, Central Legrests)

1588-405.no




1	SPTRC1945	(Legrest Support - SW 405-505mm (16-20"))	1
2	SPTRC2096	(Legrest Support - SW 505-610mm (19-24"))	1



-	SPTRC1946	SENERBENSTØTTE LNX TDX SP2		1
-	SPTRE0756	(Actuator cable)		1
1	SPTRC2124	(Legrest actuator, cpl. with mounting hardware)		1
2	SPTRC2125_C01	(Foot plates Medium (240mm depth), Pair, cpl.)		1
2	SPTRC2125_C03	FOTPLATER STORE KPL PAR TDX SP2	01.01.2016 →	1
3	SPTRC2126	(Cover, cpl. with mounting hardware)		1
4	SPTRC2035	(Calf Pads (pair), cpl. with mounting hardware)		1
5	SPTRC2262	(Mounting Hardware, foot plates)		1
6	SPTRC2222	(Footboard 305x240mm, cpl. with mounting hardware)	01.01.2016 →	1

(Vari Legrests)


-	(Supports, Vari)	1
-	(Hanger Vari A (04/2009 =>))	1
-	(Hanger Vari F)	1
-	(Hanger Vari ADM man.)	1
-	(Hanger Vari ADE/ADP pow. (07/2005 =>))	1
-	(Footrests Vari, Flip-Up (01.08.10 =>))	1
-	(Footrests Vari, Flip-Up & angle adjustable (01.08.10 =>))	1
-	(Footrest Vari, Angle & Depth adj.)	1



1588410.ISO


1	SP1586651	BENSTØTTE OVERDEL MAN KPL VE	1
2	SP1587187	BENSTØTTE OVERDEL EL KPL VE	1
3	SP1586706	BENSTØTTE OVERDEL MAN KPL HØ	1
4	SP1587186	(Legrest Support Vari RH el. cpl.)	1
5	SPTRC2276	(Cross Brace 16-22" W)	1
6	SPTRC2521	(Cross Brace 20-24" W)	1
7	SPTRC2061	(End cap kit)	1
8	SP1584416	(Legrest Power Supply Cable 900mm - LH)	1
9	SP1584302	(Legrest Power Supply Cable 700mm - RH)	1

(Hanger Vari A (04/2009 =>))



1528842.iso


-	SP1527820	(Rework-Kit Vari A Legrest Version 1 => Version 2) <i>Vari A Version 1 => Version 2 Upgrade Kit</i>	1
1	SP1527498	BENSTØTTE OVERDEL VE MAN VREG LKOMP FRA 07/2005	1
1	SP1527492	BENST OVERDEL MAN HØ STORM4	1
2	SP1434955	(Cover Legrest cpl.)	1
3	L34011	(Handknob M8 Tristarfemale)	1
4	SP1450447	LEGGPLATEHOLDER KPL VE NY MAI 05	1
4	SP1450446	LEGGPLATEHOLDER KPL HØ NY MAI 05	1
5	6301333	LEGGSTØTTE EL MAN VREG BENST STORM/TWISTER	1



1449347X.iso

1	SP1449347	BENSTØTTE OVERDEL VARI F STORM TDX BORA VE	1
1	SP1449346	BENSTØTTE OVERDEL VARI F STORM TDX BORA HØ	1
2	SP1434955	(Cover Legrest cpl.)	1
3	1434956	(Mechanic)	1
4	SPZ173178	FJÆR STORM4 BORA	1


(Hanger Vari ADM man.)



1449409X.iso

1	6301333	LEGGSTØTTE EL MAN VREG BENST STORM/TWISTER	1
2	SP1450446	LEGGPLATEHOLDER KPL HØ NY MAI 05	1
2	SP1450447	LEGGPLATEHOLDER KPL VE NY MAI 05	1
3	SPZ173178	FJÆR STORM4 BORA	1
4	SPZ174134	(SLIDING BUSHING)	1
5	SP1450410	KNEPOLSTER BENST EL VE 07.2005->	1
6	SP1450448	KNEPOLSTER BENSTG EL HØ 07.2005->	1
7	SP1449293	BENST HOLDER OVERDEL PAR EL VREG 07.2005->	1
8	SPZ174130	BENST ERGO REGULERINGSMEK KPL	1
9	SPN12196	FJÆR STORM4 BORA	1
10	SP1006490	(Brake Lever Rubber -LT2000)	1
11	SP1449409	BENSTØTTE OVERDEL ADM VE STORM4 MAN VREG	1
11	SP1449410	BENSTØTTE OVERDEL ADM HØ STORM4 MAN VREG	1

(Hanger Vari ADE/ADP pow. (07/2005 =>))




1449407X.iso

1	6301333	LEGGSTØTTE EL MAN VREG BENST STORM/TWISTER	1
2	SP1450447	LEGGPLATEHOLDER KPL VE NY MAI 05	1
2	SP1450446	LEGGPLATEHOLDER KPL HØ NY MAI 05	1
3	SPZ173178	FJÆR STORM4 BORA	1
4	SPZ174134	(SLIDING BUSHING)	1
5	SP1450448	KNEPOLSTER BENSTG EL HØ 07.2005->	1
6	SP1450410	KNEPOLSTER BENST EL VE 07.2005->	1
7	SP1449293	BENST HOLDER OVERDEL PAR EL VREG 07.2005->	1
8	SPZ151006	MOTOR EL BENST STORM4	1
9	SPN12196	FJÆR STORM4 BORA	1
10	SP1449408	BENST OVERDEL ADP HØ STORM4	1
10	SP1449407	BENST OVERDEL ADP VE STORM4	1


(Footrests Vari, Flip-Up (01.08.10 =>))

152964X.lin



1	SP1541575	FORPLATE OPPFELLBAR KPL VENSTRE	1
2	SP1541574	FOTPLATE OPPFELLBAR KPL HØYRE	1
3	SP1529645	FOTPLATE KLP SB43-45 VE STORM4	1
4	SP1529644	FOTPLATE KLP SB43-45 HØ STORM4	1
5	SP1543045	(Plug-in tube LH/RH)	1
6	SP1442467	LEGGBÅND HELKAPPE PAR MIRAGE, STORM4, TDX	1


(Footrests Vari, Flip-Up & angle adjustable (01.08.10 =>))



1	SP1541580	FOTPLATE M/INNERRØR KPL VE VREG	1
2	SP1541579	FOTPLATE M/INNERRØR KPL	1
3	SP1529645	FOTPLATE KLP SB43-45 VE STORM4	1
4	SP1529644	FOTPLATE KLP SB43-45 HØ STORM4	1
5	SP1423204	(FOOT PLATE TUBE, ANGLE ADJUSTMENT LH CPL.)	1
6	SP1423203	(FOOT PLATE TUBE, ANGLE ADJUSTMENT RH CPL.)	1
7	SP1442467	LEGGBÅND HELKAPPE PAR MIRAGE, STORM4, TDX	1

(Footrest Vari, Angle & Depth adj.)


1583898.no



1	SP1419344	FOTPLATE VE VREG OG DREG KPL MED INNERRØR	1
2	SP1419345	FOTPLATE HØ VREG OG DREG KPL MED INNERRØR	1

-	(Legrest Support)	1
-	(Legrests)	1


(Legrest Support)



1645204.iso

1	SPTRC2276	(Cross Brace 16-22" W)	1
2	SPTRC2521	(Cross Brace 20-24" W)	1
3	SPTRC2277	(Legrest Support (pair))	1
4	SPTRC2061	(End cap kit)	1
5	SPTRC2570	(Phono Connector Kit (TDX-SP2 only))	1


11.08.2017 →



1	SPTRC2070	(Pivot Plus Legrest man. LH, cpl. with support)	1
2	SPTRC2069	(Pivot Plus Legrest man. RH, cpl. with support)	1
3	SPTRC2068	(Pivot Plus Legrest pow. LH, cpl. with support)	1
4	SPTRC2067	(Pivot Plus Legrest pow. RH, cpl. with support)	1
5	SPTRC2127-C02	(Calf Pad RH/LH with mounting hardware)	1
6	SPTRC2127-C02	(Calf Pad RH/LH with mounting hardware)	1
7	SPTRC2129	(Wheel, Footplate - cpl.)	1
8	SPTRC2128	(Legrest Actuator, cpl.)	1
9	SPTRC2109-C01	(Pivot Mount Assy LH/RH)	1
10	SPTRC2108-C09	(Footplate Assy 5.25x8", Heavy Duty - LH)	1
11	SPTRC2108-C10	(Footplate Assy 5.25x8", Heavy Duty - RH)	1
-	SPTRE0756	(Actuator cable) LH / links	1
-	SPTRE0754	(Actuator cable) RH / rechts	1


(Legrest Support, Vari series 07/2005 ==>)

1449734 iso



1	SP1449732	BENSTØTTEHOLDER MAN HØ 5CM STORM	1
2	SP1449736	BENSTØTTEHOLDER MAN HØ 8CM STORM TDX BORA	1
3	SP1449733	BENSTØTTEHOLDER MAN VE 5CM STORM	1
4	SP1449737	BENSTØTTEHOLDER MAN VE 8CM STORM TDX BORA	1
5	SP1531071	BENSTØTTEHOLDER EL HØ STORM4 BORA 5CM	1
6	SP1531073	BENSTØTTEHOLDER EL 8CM STORM4 HØ	1
7	SP1531072	BENSTØTTEHOLDER EL VE STORM4 BORA 5CM	1
8	SP1531082	BENSTØTTEHOLDER EL 8CM STORM4 VE	1
9	SP1530795	KABEL I HOLDER F BENST STORM4 BORA TDX	1

(Legrest Support, Vari series 07/2005 =>)




1653784.ISO

1	SP60109510	(Legrest Support Vari series, 5cm - LH)	1
2	SP1449737	BENSTØTTEHOLDER MAN VE 8CM STORM TDX BORA	1
3	SP60109511	(Legrest Support Vari series, 5cm - RH)	1
4	SP1449736	BENSTØTTEHOLDER MAN HØ 8CM STORM TDX BORA	1
5	SP1530795	KABEL I HOLDER F BENST STORM4 BORA TDX	1

Invacare® Storm4 | Benstøtter og holdere
(Central Legrest, manuell)


1529387X.dwg



1	SP1529387	SENERBENST VREG MAN STORM4	1
2	1436070	LEGGPLATE HØ KPL M/HOLDER SENTERBENST	1
3	1436069	LEGGPLATE VE KPL M/HOLDER SENTERBENST	1

Invacare® Storm4 | Benstøtter og holdere
(Central Legrest, electric)


152938X.no




1	SP1529388	SENERBENST VREG EL STORM4 1 MOTOR	1
2	1436076	LEGGPLATE VE KPL EL SENTERBENST	1
3	1440594	LEGGPLATE HØ KPL EL SENTERBENST	1
4	1440593	GASSFJÆR STORM4	1
5	SP1440600	ACTUATOR IVC SENTERBENST TDX SP2	1
6	SP1574298	FESTEBRAKETT IVC SENTER BENSTØTTE STOR TYPE	1
7	SP1576334	(Mounting Hardware Arm 2-E)	1
8	SP1524744	(Adapter Cable AMP -> DX2 (L=50mm))	1

(Support Tube with Footplate for central Legrest)

1440290X.NO






1	SP1440290	FOTPLATE OPPFELL INNERØR 32-38 HØ/VE STORM4 T2	1
1	SP1440291	FOTPLATE OPPFELLB MEDIUM INNERRØR SENTER BENST	1
2	SP1440442	INNERRØR SENTERBENSTØTTE STORM 4 EXTRA LITEN	1
2	SP1440443	INNERØR SENTERBENSTØTTE STORM 4 EXTRA MEDIUM	1
3	SPF00002	(Hexagon Nut M8 bk - DIN 985)	2
4	SP1440517	(Darning)	1
5	SP1546720	(Footboardrubber)	1



1645695.iso

1	SP1423156	FOTPLATE HEL KPL STORM3	1
2	SP1645738	(Footboard bracket RH)	1
3	SP1645737	(Footboard bracket LH)	1
4	Z174167	(Calf strap L=930mm)	1
5	SP1533540	FOTPLATE SB53 STORM4	1
6	SP1552599	(Footboard bracket RH - SW530mm)	1
7	SP1552600	(Footboard bracket LH - SW530mm)	1
8	SP1552767	(Calf strap L=1130mm)	1
9	Z174159	FOTBRETT UNDERLAG HEL FOTPLA STORM/TWISTER	1

(Lower Footrests)




F3_BSTU.iso

1	SP1421662	BENSTØTTEUNDERDEL HØ STORM3	1
1	SP1421663	BENSTØTTEUNDERDEL VE STORM3	1
2	SP1421664	BENSTØTTEUNDERDEL HØ VINKELREG. STORM	1
2	SP1421665	BENSTØTTEUNDERDEL VE VINKELREG. STORM	1
3	SP1419344	FOTPLATE VE VREG OG DREG KPL MED INNERRØR	1
3	SP1419345	FOTPLATE HØ VREG OG DREG KPL MED INNERRØR	1


-	(Legrest Supports Vari (07/2005 =>))	0
-	(Hanger Vari A (04/2009 =>))	0
-	(Hanger Vari F)	0
-	(Hanger Vari ADM man.)	0
-	(Hanger Vari ADE/ADP pow. (07/2005 =>))	0
-	(Central Legrest, manuell)	0
-	(Central Legrest, powered)	0
-	(Footrests Vari, Flip-Up (01.08.10 =>))	1
-	(Footrests Vari, Flip-Up & angle adjustable (01.08.10 =>))	1
-	(Support Tube with Footplate for central Legrest)	1
-	(Footboards)	1
-	(Lower Footrests)	1

(Legrest Supports Vari (07/2005 =>))

1582374.00



1	SP1449732	BENSTØTTEHOLDER MAN HØ 5CM STORM	1
2	SP1449736	BENSTØTTEHOLDER MAN HØ 8CM STORM TDX BORA	1
3	SP1449733	BENSTØTTEHOLDER MAN VE 5CM STORM	1
4	SP1449737	BENSTØTTEHOLDER MAN VE 8CM STORM TDX BORA	1
5	SP1449314	BENSTØTTEHOLDER EL HØ 5CM STORM3 NY 010705	1
6	SP1449315	BENSTØTTEHOLDER EL HØ 8CM STORM3 NY 010705 <i>EasyAdapt</i>	1
7	SP1449735	BENSTØTTEHOLDER EL VE 5CM STORM3 TO TY 010705 <i>EasyAdapt</i>	1
7	SP1537398	(Footrest/Legrest receptacle, powered, LH, 5cm - R-Net w/ lifter) <i>Modulite</i>	1
8	SP1449739	BENSTØTTEHOLDER EL VE 8CM ST03 TO TY NY 010705 <i>EasyAdapt</i>	1
9	SP1450232	KABEL MED KONTAKT BENST TDX SP	1
10	SP1450228	BENST FESTEDEL TDX	1



1	SP1529388	SENERBENST VREG EL STORM4 1 MOTOR	1
2	1436076	LEGGPLATE VE KPL EL SENTERBENST	1
3	1440594	LEGGPLATE HØ KPL EL SENTERBENST	1
4	1440593	GASSFJÆR STORM4	1
5	SP1440600	ACTUATOR IVC SENTERBENST TDX SP2	1
6	SP1574298	FESTEBRAKETT IVC SENTER BENSTØTTE STOR TYPE	1
7	SP1576334	(Mounting Hardware Arm 2-E)	1
8	SP1579552	(Adapter Cable, Centr. Legrest (GGB2))	1

Armlen




SPCSTORM440.iso

1	-	(Armrests for EasyAdapt Seat => 2017-08-10)	1
2	-	(Armrests & Hip Supports for Modulite Seat)	1
3	-	(Armrests for Ultra Low Maxx Seat)	1
4	-	(Armrest for Recaro Seat)	1
5	-	(Armrests for Optimist Seat => 2018-12-31)	1


(Armrests for EasyAdapt Seat => 2017-08-10)

1573233.no




1	-	(Fixed Armrest, removable) <i>[0635/0636]</i>	1
2	-	(Following armrests) <i>[0620]</i>	1
3	-	(Skirtguard for following armrest) <i>[0623]</i>	1

(Fixed Armrest, removable)



1421692X.iso


1	SP1425846	(Armrest cpl. 38cm LH Recaro man. Back)	1
2	SP1425845	(Armrest cpl. 38cm RH Recaro man. Back)	1
3	SP1436534	ARMLN DHREG VE SB39 SD41 H34 FLEX TOR/TYPH	1
3	SP1436536	ARMLN DHREG VE SB39 SD41 H38 FLEX TOR TYPH	1
4	SP1436533	ARMLN DHREG HØ SB39 SD41 H34 FLEX TOR/TYPH	1
4	SP1436535	ARMLN DHREG HØ SB39 SD41 H38 FLEX TOR/TYPH	1
5	SP1421694	ARMLN VE KPL HREG AVTGB 340MM STORM TDX	1
5	SP1421695	ARMLN VE KPL DHREG 380MM STORM3	1
6	SP1421692	ARMLN HØ KPL HREG AVTGB 340MM STORM TDX	1
6	SP1421693	ARMLN HØ KPL DHREG 380MM STORM3	1
7	SP1443004	ARMLNSPOLSTER VE STORM STORM TDX BORA	1
7	SP1449043	ARMLNSPOLSTER HØ STORM STORM TDX BORA	1
7	SP1528846	ARMLNSPOLSTER VE LANG STORM TDX BORA	1
7	SP1528845	ARMLNSPOLSTER HØ LANG STORM4 TDX BORA	1
8	SP1510857	ARMLNSPOLSTER HOLDER VE STORM3/4 2009-> <i>For chairs bevor 01/01/09 also order a new armpad.</i>	1
8	SP1510856	ARMLNSPOLSTER HOLDER HØ STORM3/4 2009-> <i>For chairs bevor 01/01/09 also order a new armpad.</i>	1
9	SP1423225	HOLDER ARMLN KPL VENSTR FAST ARMLN STORM3	1
9	SP1423223	HOLDER ARMLN KPL HØYRE FAST ARMLN STORM3	1
10	SP1428078	KLESBESKYTTER STORM4 TDX	1
11	SP1420471	KLESBESKYTTER TDX	1



1	SP1536036	ARMLÉN OPPFELL VE STORM4	1
1	SP1536100	ARMLÉN LANG VE OPPFELLB STORM4	1
2	SP1536035	ARMLÉN OPPFELL HØ STORM4	1
2	SP1536099	ARMLÉN LANG HØ OPPFELLB STORM4	1
3	SP1443004	ARMLÉNSPOLSTER VE STORM STORM TDX BORA	1
3	SP1528846	ARMLÉNSPOLSTER VE LANG STORM TDX BORA	1
4	SP1449043	ARMLÉNSPOLSTER HØ STORM STORM TDX BORA	1
4	SP1528845	ARMLÉNSPOLSTER HØ LANG STORM4 TDX BORA	1
5	SP1535504	(Knurled screw M8x30)	1
6	SP1535244	FESTE ARMLÉN RYGGMONT STORM4	1

(Skirtguard for following armrest)


1425860X.dwg



1	SP1426303	KLESBESKYTTER OPPF ARML STORM4	1
2	SP1442549	FESTE ARMLEN VE TDX	1
3	SP1442550	FESTE ARMLEN HØ TDX	1

(Armrests & Hip Supports for Modulite Seat)


1571738.000




1	-	(T-Armrest) [2191+2195/2196/2197]	1
2	-	(Padded Armrest) [2198]	1
3	-	(Inrigged Armrest) [2183+2195/2196/2197/2199]	0
4	-	(Flip-Up Armrest) [2193+2195/2196/2197/2199]	0
5	-	(Following Armrest) [2192/2194+2195/2196/2197/2199]	0
6	-	(Hip Supports) [2184/2185/2186/2187/2188/2189]	0

(T-Armrest)

1572570 iso




1	SP1567750	FESTE POSTER ARMLEN TDX	1
2	SP1570319	FESTE ARMLEN ERGONOMISK MODULITE	1
3	SP1585088	(Armrest Bracket for Soft Armpad)	1
4	SP1443004	ARMLENSPOLSTER VE STORM STORM TDX BORA	1
5	SP1449043	ARMLENSPOLSTER HØ STORM STORM TDX BORA	1
6	SP1528846	ARMLENSPOLSTER VE LANG STORM TDX BORA	1
7	SP1528845	ARMLENSPOLSTER HØ LANG STORM4 TDX BORA	1
8	SP1571589	ARMLENSPOLSTER ERGONO VE STORM4 TDX	1
9	SP1571588	ARMLENSPOLSTER ERGONO HØ STORM4 TDX	1
10	SP1585068	ARMLENSPOLSTER VE MYKT TDX SP2	1
11	SP1585069	ARMLENSPOLSTER HØ MYKT TDX SP2	1
12	SP1588374	ARMLEN U POLSTER VE TDX SP2 MODULITE	1
13	SP1588375	ARMLEN U POLSTER HØ TDX SP2 MODULITE	1
14	SP1601882	SIDESTØTTE HØY VE UTEN FESTE	1
15	SP1601883	SIDESTØTTE HØY HØ UTEN FESTE	1
16	SP1572940	KLESBESKYTTER VE MODULITE	1
17	SP1572948	KLESBESKYTTER HØ MODULITE	1
18	SP1567749	BRAKETT KLESBESKYTTER MODULITE	1
19	SPP31298	(Reflector, yellow)	1
20	SPL10024	(Wing knobs M6x16)	1




1	SP1571987	ARMLEN PADDED KPL U/FEST MODULITE HØ/VE	1
2	SP1571962	ARMLENSPOLSTER PADDED MODULITE	1
3	SPL10024	(Wing knobs M6x16)	1
4	SP1567749	BRAKETT KLESBESKYTTER MODULITE	1
5	SPP31298	(Reflector, yellow)	1

(Inrigged Armrest)




1	SP1567750	FESTE POSTER ARMLÉN TDX	1
2	SP1570319	FESTE ARMLÉN ERGONOMISK MODULITE	1
3	SP1585088	(Armrest Bracket for Soft Ampad)	1
4	SP1443004	ARMLÉNSPOLSTER VE STORM STORM TDX BORA	1
5	SP1449043	ARMLÉNSPOLSTER HØ STORM STORM TDX BORA	1
6	SP1528846	ARMLÉNSPOLSTER VE LANG STORM TDX BORA	1
7	SP1528845	ARMLÉNSPOLSTER HØ LANG STORM4 TDX BORA	1
8	SP1571589	ARMLÉNSPOLSTER ERGONO VE STORM4 TDX	1
9	SP1571588	ARMLÉNSPOLSTER ERGONO HØ STORM4 TDX	1
10	SP1585068	ARMLÉNSPOLSTER VE MYKT TDX SP2	1
11	SP1585069	ARMLÉNSPOLSTER HØ MYKT TDX SP2	1
12	SP1585037	(Siderest inrigged LH)	1
13	SP1585038	(Siderest inrigged RH)	1
14	SP1605654	(Siderest inrigged LH (SD410mm))	1
15	SP1605653	(Siderest inrigged RH (SD410mm))	1
16	SP1578119	KLESBESKYTTER HØ/VE TDX	1
17	SP1603118	(Skirt Guard LH/RH (SD410mm))	1
18	SP1567749	BRACKETT KLESBESKYTTER MODULITE	1
19	SPL10024	(Wing knobs M6x16)	1
20	SPP31298	(Reflector, yellow)	1
21	SP1582989	(Reflector yellow D=43/38mm)	1




-	-	(Retrofitting the armrests on chairs built before 15-03-27 also requires new Pos. A & B and Sling Back with respective cut-outs!) !!!		0
1	SP1582668	ARMLLEN FLIP UP VE KPL MUDULITE	19.01.2015 →	1
2	SP1582644	ARMLLEN FLIP UP HØ KPL MODULITE	19.01.2015 →	1
3	SP1585412	HOLDER VE POLSTER OPPFEL MODULITE		1
4	SP1585413	HOLDR HØ POLSTER OPPFELL MODULITE		1
5	SP1585414	FESTE MYKT POLSTER HØ MODULITE		1
6	SP1585415	FESTE MYKT POLSTER VE MODULITE		1
7	SP1528846	ARMLLENSPOLSTER VE LANG STORM TDX BORA		1
8	SP1528845	ARMLLENSPOLSTER HØ LANG STORM4 TDX BORA		1
9	SP1443004	ARMLLENSPOLSTER VE STORM STORM TDX BORA		1
10	SP1449043	ARMLLENSPOLSTER HØ STORM STORM TDX BORA		1
11	SP1571589	ARMLLENSPOLSTER ERGONO VE STORM4 TDX		1
12	SP1571588	ARMLLENSPOLSTER ERGONO HØ STORM4 TDX		1
13	SP1585068	ARMLLENSPOLSTER VE MYKT TDX SP2		1
14	SP1585069	ARMLLENSPOLSTER HØ MYKT TDX SP2		1
15	SP1590599	(Cap)		1
16	SP1582989	(Reflector yellow D=43/38mm)		1

(Following Armrest)



1585182.iso

-	-	(Retrofitting the armrests on chairs built before 15-03-27 also requires new Pos. A & B and Sling Back with respective cut-outs!) !!!		0
1	SP1585263	ARMLEN OPPFELLBAR VE MODULITE KORT	19.01.2015 →	1
2	SP1585264	ARMLEN OPPFELLBAR HØ MODULITE KORT	19.01.2015 →	1
3	SP1585267	ARMLEN OPPFELL VE LANG MODULITE FOLLOWING	19.01.2015 →	1
4	SP1585268	ARMLEN OPPFELLB HØ LANG MODULITE FOLLOWING	19.01.2015 →	1
5	SP1585412	HOLDER VE POLSTER OPPFEL MODULITE		1
6	SP1585413	HOLDR HØ POLSTER OPPFELL MODULITE		1
7	SP1585414	FESTE MYKT POLSTER HØ MODULITE		1
8	SP1585415	FESTE MYKT POLSTER VE MODULITE		1
9	SP1443004	ARMLENSPOLSTER VE STORM STORM TDX BORA		1
10	SP1449043	ARMLENSPOLSTER HØ STORM STORM TDX BORA		1
11	SP1528846	ARMLENSPOLSTER VE LANG STORM TDX BORA		1
12	SP1528845	ARMLENSPOLSTER HØ LANG STORM4 TDX BORA		1
13	SP1571589	ARMLENSPOLSTER ERGONO VE STORM4 TDX		1
14	SP1571588	ARMLENSPOLSTER ERGONO HØ STORM4 TDX		1
15	SP1585068	ARMLENSPOLSTER VE MYKT TDX SP2		1
16	SP1585069	ARMLENSPOLSTER HØ MYKT TDX SP2		1
17	SP1583156	(Amrest Support, long, LH)		1
18	SP1583155	(Amrest Support, long, RH)		1
19	SP1583154	STØTTE ARMLLEN VE KORT		1
20	SP1582639	(Amrest Support, short, RH)		1
21	SP1590599	(Cap)		1
22	SP1582989	(Reflector yellow D=43/38mm)		1
23	SP1583166	(Bearing Plate, cpl.)		1




1585183.ISO

1	SP1584411	FESTEBRAKETT LÅRSTØTTE TDX SP2 MODULITE VE	1
1	SPTRC2235	(Hip Support, Quick Release bracket, ext. (+50mm height), LH)	1
2	SP1584410	FESTEBRAKETT LÅRSTØTTE T STORM4 TDX LOW-RIDER HØ	1
2	SPTRC2234	(Hip Support, Quick Release bracket, ext. (+50mm height), RH)	1
3	SP1584414	LÅRSTØTTE PUTE 33X10 MODULITE	1
4	SP1584413	LÅRSTØTTE PUTE 23*10 CM MODULITE	1
5	SPTRC1864-C02	(Hip Support Pad 180x150mm)	1
6	SP1584412	LÅRSTØTTE PUTE 14*10 CM MODULITE	1

(Armrests for Ultra Low Maxx Seat)

1	-	(Two-post Reclining Arms, Flipup)	0
2	-	(Cantilever Arms, Flipup)	0
3	-	(Armpads, Desk & Full)	0
4	-	(Armpad, Ergonomic)	0
5	-	(Elbow Block (=> 2017-08-11))	0
6	-	(MACES Multi Axis Armpad Holder)	0
7	-	(Hip Supports)	0


(Two-post Reclining Arms, Flipup)




1	SPTRC1933	(Two-post Reclining Arms, Flipback 240-330mm)	1
2	SPTRC1934	(Two-post Reclining Arms, Flipback 320-405mm)	1
3	SPTRC2135	(Flip Back Lever with connector LH)	1
4	SPTRC2134	(Flip Back Lever with connector RH)	1
5	SPTRM0130	(Armrest-pipe cap)	1
6	SPTRC2060	(Clamp with square tube adapter)	1
7	SP1602015	(Reflector Kit ULM - LH)	1
7	SP1602016	(Reflector Kit ULM - RH)	1

(Cantilever Arms, Flipup)

1602155.iso




1	SPTRC1936	(Cantilever Armrests Flipback 240-330mm)	1
2	SPTRC1937	(Cantilever Armrests Flipback 320-405mm)	1
3	SPTRM0130	(Armrest-pipe cap)	1
4	SPTRC2060	(Clamp with square tube adapter)	1
5	SP1602015	(Reflector Kit ULM - LH)	1
5	SP1602016	(Reflector Kit ULM - RH)	1



1	SPTRX8038	(Armpad, Desk Amrest, LH)		1
2	SPTRX8037	(Armpad, Desk Amrest, RH)		1
3	SPTRX8036	(Armpad, Full Amrest, LH)		1
4	SPTRX8035	(Armpad, Full Amrest, RH)		1
5	SPTRC2133	(Desk Armrest Body LH with mounting hardware)		1
6	SPTRC2132	(Desk Armrest Body RH with mounting hardware)		1
7	SPTRC2131	(Full Armrest Body LH with mounting hardware)		1
8	SPTRC2130	(Full Armrest Body RH with mounting hardware)		1
9	SPTRC1935	(Desk to Full Armrest Extension Kit) <i>Needed for Pos. 7 & 8!</i>	→ 10.08.2017	1
9	SPTRC2526	(Armrest Extension Kit (Desk -> Full)) <i>Needed for Pos. 7 & 8!</i>	11.08.2017 →	1
10	SP1602292	(Screw set for armpad fixation) <i>[1494]</i>		1

(Armpad, Ergonomic)


1588506 Rev



1	SPTRC2552-C01	(Ergonomic Arm Trough LH)	1
2	SPTRC2552-C02	(Ergonomic Arm Trough RH)	1
3	SPTRC2526	(Armrest Extension Kit (Desk -> Full)) <i>Needed for Pos. 1 & 2!</i>	1
4	SPTRC2102	(Mounting Hardware for Otto Bock)	1
5	SPTRC2553	(Flat hand pad)	1

(Elbow Block (=> 2017-08-11))


1603428.iss



1	SPTRM1059	(Elbow Block pad 3x10 LH)	1
2	SPTRM1060	(Elbow Block pad 3x10 RH)	1
3	SPTRM1058	(Elbow Block pad 4x4)	1
4	SPTRM1055	(Elbow Block Hardware, cpl. LH)	1
5	SPTRM1056	(Elbow Block Hardware, cpl. RH)	1

(MACES Multi Axis Armpad Holder)


1588515.00



1 SPTRC1938

(MACES Armpad Holder Kit, Multi Axis)

1




1 -


(Fixed Armrest, removable)
[0636]

1

(Armrests for Optimist Seat => 2018-12-31)




(Armrests Arne Jensen for Optimist Seat)




1	1419657	ARMLEN OPPFELL 30CM HØ SNAPLOKK OPTIMIST KPL	1
2	1419658	ARMLEN OPPFELL 30CM VE SNAPLOKK OPTIMIST KPL	1
3	1419659	ARMLEN OPPFELL 35CM HØ SNAPLOKK OPTIMIST KPL	1
4	1419660	ARMLEN OPPFELL 35CM VE SNAPLOKK OPTIMIST KPL	1
5	1419661	ARMLEN OPPFELL 40 CM HØ SNAPLOKK OPTIMIST KPL	1
6	1419662	ARMLEN OPPFELL 40CM VE SNAPLOKK OPTIMIST KPL	1
7	1419647	ARMLENSPOLSTER HØYRE VOKSEN ANATOMISK,KORT 30	1
8	1419648	ARMLENSPOLSTER VENSTRE VOKSEN ANATOMISK,KORT 30	1
9	5314011	ARMLENSPOLSTER HØYRE VOKSEN ANATOMISK,KORT 35	1
10	5314010S	ARMLENSPOLSTER VENSTRE VOKSEN ANATOMISK,KORT 35	1
11	5314011L	ARMLENSPOLSTER HØYRE VOKSEN ANATOMISK,LANG 40	1
12	5314010	ARMLENSPOLSTER VENSTRE VOKSEN ANATOMISK,LANG 40	1

(Front Fork)

152890X.lin



1	SP1588447	GAFFEL KPL VE STORM4 <i>We recommend to exchange the front forks by pairs, because of new design. - replaces SP1552920</i>	30.07.2015 →	1
2	SP1588448	GAFFEL KPL HØ STORM4 <i>We recommend to exchange the front forks by pairs, because of new design. - replaces SP1552926</i>	30.07.2015 →	1
3	SP1530553	HJULSKJERM SVINGHJUL VE STORM4		1
4	SP1530530	HJULSKJERM SVINGHJUL HØ STORM4		1
5	SP1530554	NAV DEKSEL KPL STORM3 STORM3 SVINGHJUL		1
-	-	(Please refer to Service Manual for correct bearing orientation!)		1
6	SP1530529	GAFFELHUSLAGER KPL STORM4		1




1	SP1418531	SVINGHJUL KPL 3.00-6 KPL GRÅ LUFT STORM 3	1
1	SP1422932	SVINGHJUL 3.00-6 KPL GRÅ PUNKTERINGSFRI	1
1	SP1430231	SVINGHJUL 3.00-6 GRÅ KOMPLETT STORM 3	1
2	SP1423182	KULELAGER KIT SVINGHJUL STORM 3	1
3	1415175	SLANGE 3.00-6 90GR TR87 STORM3	1
3	SP1557438	(TUBE 3.00 X 6 90 ° TR 87, PUNCTURE PROTECTED)	1
4	1418532	DEKK 3.00-6 GRÅ STORM3	1
4	1430229	DEKK 3.00-6 GRÅ PUNKTERINGSSIKKER STORM3	1
5	SP1423181	HJUL MONTERINGSKIT STORM4	1

(Storm4)

-	(Drive Wheel Trelleborg (8x3.00))	1
-	(Drive Wheel 14" (3.00 - 8))	1


(Drive Wheel Trelleborg (8x3.00))

1421832X.dwg



1	SP1421800	DRIVHJUL TRELLEBORG STORM4 KPL		1
1	SP1421832	DRIVHJUL TRELLEBORG KPL STORM4 PUNKTERINGSIKKER		1
2	SP1423184	FELG KPL.DRIVHJUL STORM 3		1
3	M12016	SLANGE 3.00 X 8 T. STORM/ G40		1
4	M11968	DEKK 3.00-8 TRELLEBORG GRÅ STORM STORM 3 STD.		1
5	1509054-0002	NAV DRIVHJUL,STORM 3 <i>Replaced by SP1541458</i>	→ 04.01.2011	1
5	SP1541458	NAV 4 HULL STORM4 <i>Replaces 1509054-0002</i>	05.01.2011 →	1

(Drive Wheel 14" (3.00 - 8))



1	SP1532414	DRIVHJUL KPL M LUFT STORM4		1
2	SP1532415	DRIV HJUL PUNKT SIKKERT STORM4 KPL		1
3	SP1532416	HJUL KOMPAKT 3.00-8 GRÅ STORM4		1
4	SP1541373	FELG KPL DRIVHJUL STORM4		1
-	SP1522975	Bolt 5/16-18 x 7/8" TDX		5
5	SP1522940	DEKK LUFT 14X 3 (3,00-8) STORM4		1
6	SP1141342	DEKK KOMPACT GB2 14X3 STORM4 TDX		1
7	SP1039651	SLANGE 14" STORM4 MAX		1
8	SP1541374	(Tube 14", puncture protected)		1
9	SP1163728	(Wheel hub 5-hole)	05.01.2011 →	1
		<i>05.01.2011 => - Replaces SP1536413</i>		
10	SP1536413	HJULNAV KPL STORM4	→ 04.01.2011	1
		<i>=>05.01.2011 - Replaced by SP1163728</i>		


-

(Drive Wheels GGB)

0


(Drive Wheels GGB)

1581908.iso




1	SP1582235	(Wheel mounting kit)	1
2	SP1582236	(Outer Wheel Rim 14" GGB2 incl. Sticker)	1
3	SP1522940	DEKK LUFT 14X 3 (3,00-8) STORM4	1
4	SP1141342	DEKK KOMPAKT GB2 14X3 STORM4 TDX	1
5	SP1039651	SLANGE 14" STORM4 MAX	1
6	SP1541374	(Tube 14", puncture protected)	1

-	(Propulsion => 04.01.11)	0
-	(Overview Propulsion)	0



1	SP1530578	(Drivingmotor, 6km/h, cpl., LH) => S/N 061122812867	1
1	SP1530577	(Drivingmotor, 6km/h, cpl., RH) => S/N 061122812867	1
1	SP1530580	\$MOTOR STORM4 VE UTGÅTT => S/N 061122610150	1
1	SP1530579	\$MOTOR STORM4 HØ UTGÅTT => S/N 061122610150	1
2	SP1101725	PROPP FOR FRIKPL HENDEL STORM4	1
3	SP1522967	BØRSTER 4 STK KIT TDX STORM 4	1
4	SP1551712	(Labelset "Push/Drive")	1



Bora	6km/h	c	ACS LINX Shark	1	2	33	34	35	37	---	---	51	45
			VR2	3	4	33	34	35	41	---	---	51	45
			R-Net	5	6	33	34	35	43	---	---	51	45
	10km/h	d	ACS LINX Shark	7	8	31	32	35	36	---	---	51	45
		VR2	9	10	31	32	35	40	---	---	51	45	
		R-Net	11	12	31	32	35	42	---	---	51	45	
BoraHD	6km/h		ACS LINX Shark	54	55	33	34	35	39	---	---	51	46
	10km/h		Shark	56	57	31	32	35	38	---	---	51	46
	6km/h		R-Net	58	59	33	34	35	62	---	---	51	46
	10km/h			60	61	31	32	35	63	---	---	51	46
FDX	6km/h	e		13	14	33	34	35	37	52	53	---	46
	6km/h	f	ACS Shark	54	55	33	34	35	39	52	53	---	46
	10km/h			15	16	31	32	35	38	52	53	---	46
Kite	6km/h		ACS LINX Shark	17	18	33	34	35	37	---	---	51	45
	10km/h			19	20	31	32	35	36	---	---	51	45
Kite HD	6km/h		ACS LINX Shark	21	22	33	34	35	39	---	---	51	46
	10km/h			23	24	31	32	35	38	---	---	51	46
	12km/h			72	73	33	34	35	78	---	---	51	46
Storm4	6km/h		ACS Shark	47	48	34	33	35	39	---	---	51	46
	10km/h			27	28	32	31	35	38	---	---	51	46
	12km/h			74	75	34	33	35	78	---	---	51	46
TDX-SP	6km/h	a		25 / 66	26 / 67	34	33	35	39 / 65	52	53	---	46
	8km/h	a	ACS LINX Shark	29 / 68	30 / 69	34	33	35	38 / 64	52	53	---	46
TDX-SP2	10km/h	b		49 / 70	50 / 71	32	31	35	38 / 64	52	53	---	46
	12km/h			76	77	34	33	35	78	52	53	---	46
TDX-SP NB	6km/h		ACS LINX Shark	25 / 66	26 / 67	34	33	35	39 / 65	52	53	---	46
TDX-SP2 NB	10km/h			49 / 70	50 / 71	32	31	35	38 / 64	52	53	---	46

a) 28.09.2011 →
 b) 07.09.2011 →
 c) 18.04.2011 →
 d) 01.08.2011 →
 e) → 31.05.2013
 f) 01.06.2013 →
 (S/N 13F06900078)

1163099X.iso


-	N23004	KIL KRYSS DRIVAXEL (MO	1
1	SP1524195	DRIVENHET 6KMT VE BORA	1
2	SP1524194	DRIVENHET 6KMT HØ BORA	1
3	SP1541408	(Propulsion cpl. LH, Gearbox i = 23:1, Motor, 2-pole, 250W, slow, 600mm cable, VR2)	1
4	SP1541410	(Propulsion cpl. RH, Gearbox i = 23:1, Motor, 2-pole, 250W, slow, 600mm cable, VR2)	1
5	SP1541409	(Propulsion cpl. LH, Gearbox i = 23:1, Motor, 2-pole, 250W, slow, 600mm cable, R-Net)	1
6	SP1541411	(Propulsion cpl. RH, Gearbox i = 23:1, Motor, 2-pole, 250W, slow, 600mm cable, R-Net)	1
7	SP1524197	(Propulsion cpl. LH, Gearbox i = 18:1, Motor, 2-pole, 250W, fast, 600mm cable)	1
8	SP1524196	(Propulsion cpl. RH, Gearbox i = 18:1, Motor, 2-pole, 250W, fast, 600mm cable)	1
9	SP1541404	(Propulsion cpl. LH, Gearbox i = 18:1, Motor, 2-pole, 250W, fast, 600mm cable, VR2)	1
10	SP1541406	(Propulsion cpl. RH, Gearbox i = 18:1, Motor, 2-pole, 250W, fast, 600mm cable, VR2)	1
11	SP1541405	(Propulsion cpl. LH, Gearbox i = 18:1, Motor, 2-pole, 250W, fast, 600mm cable, R-Net)	1
12	SP1541407	(Propulsion cpl. RH, Gearbox i = 18:1, Motor, 2-pole, 250W, fast, 600mm cable, R-Net)	1
13	SP1168983	(Propulsion cpl. LH, Gearbox i = 23:1, Motor, 2-pole, 250W, slow, 600mm cable)	1
14	SP1168984	(Propulsion cpl. RH, Gearbox i = 23:1, Motor, 2-pole, 250W, slow, 600mm cable)	1
15	SP1168989	(Propulsion cpl. LH, Gearbox i = 18:1, Motor, 4-pole, 340W, fast, 760mm cable)	1
16	SP1168988	(Propulsion cpl. RH, Gearbox i = 18:1, Motor, 4-pole, 340W, fast, 760mm cable)	1
17	SP1568147	(Propulsion cpl. LH, Gearbox i = 23:1, Motor, 2-pole, 250W, slow, 600mm cable)	1

Declutching knob included!



(Overview Propulsion)

18	SP1568148	(Propulsion cpl. RH, Gearbox i = 23:1, Motor, 2-pole, 250W, slow, 600mm cable) <i>Declutching knob included!</i>	1
19	SP1568149	(Propulsion cpl. LH, Gearbox i = 18:1, Motor, 2-pole, 250W, fast, 600mm cable) <i>Declutching knob included!</i>	1
20	SP1568150	(Propulsion cpl. RH, Gearbox i = 18:1, Motor, 2-pole, 250W, fast, 600mm cable) <i>Declutching knob included!</i>	1
21	SP1568151	(Propulsion cpl. LH, Gearbox i = 23:1, Motor, 4-pole, 340W, slow, 760mm cable)	1
22	SP1568152	(Propulsion cpl. RH, Gearbox i = 23:1, Motor, 4-pole, 340W, slow, 760mm cable)	1
23	SP1568153	(Propulsion cpl. LH, Gearbox i = 18:1, Motor, 4-pole, 340W, fast, 760mm cable)	1
24	SP1568208	(Propulsion cpl. RH, Gearbox i = 18:1, Motor, 4-pole, 340W, fast, 760mm cable)	1
25	SP1602974	(Propulsion cpl. LH, Gearbox i = 23:1, Motor 4-pole 340W, slow, 760mm cable (incl. declutching lever))	1
26	SP1602975	(Propulsion cpl. RH, Gearbox i = 23:1, Motor 4-pole 340W, slow, 760mm cable (incl. declutching lever))	1
27	SP1541366	DRIVMOTOR KPL VE STORM4 SSD <i>S/N 061122610151 =></i>	1
28	SP1541365	DRIVMOTOR KPL HØ STORM4 SSD <i>S/N 061122610151 =></i>	1
29	SP1170307	DRIVMOTOR KPL VE TDX TDX SSD TDX SP2	1
30	SP1170306	DRIVMOTOR KPL HØ TDX TDX SSD TDX SP2	1
31	SP1163413	(Gearbox, i = 18:1, LH (Storm4 & TDX-SP/SP2 RH))	1
32	SP1163412	(Gearbox, i = 18:1, RH (Storm4 & TDX-SP/SP2 LH))	1
33	SP1163415	(Gearbox, i = 23:1, LH (Storm4 & TDX-SP/SP2 RH))	1
34	SP1163414	(Gearbox, i = 23:1, RH (Storm4 & TDX-SP/SP2 LH))	1
35	SP1542953	(Clamp band with clutch)	1
36	SP1163099	(Motor, 2-pole, 250W, fast, 600mm cable, ACS)	1
37	SP1164046	(Motor, 2-pole, 250W, slow, 600mm cable, ACS)	1
38	SP1165863	(Motor, 4-pole, 340W, fast, 760mm cable, ACS)	1
39	SP1165865	(Motor, 4-pole, 340W, slow, 760mm cable, ACS)	1
40	SP1544840	(Motor, 2-pole, 250W, fast, 600mm cable, VR2)	1
41	SP1544795	(Motor, 2-pole, 250W, slow, 600mm cable, VR2)	1
42	SP1544841	(Motor, 2-pole, 250W, fast, 600mm cable, R-Net)	1
43	SP1544796	(Motor, 2-pole, 250W, slow, 600mm cable, R-Net)	1
45	SP1540358	(Kit, brushes with caps, 2-pole (Set of 2))	1
46	SP1522967	BØRSTER 4 STK KIT TDX STORM 4	1
47	SP1536022	(Propulsion cpl. LH, G: i = 23:1, M: 4-pole 340W slow, 760mm cable (w/ declutching knob)) <i>S/N 061122812868 =></i>	1
48	SP1536021	(Propulsion cpl. RH, G: i = 23:1, M: 4-pole 340W slow, 760mm cable (w/ declutching knob)) <i>S/N 061122812868 =></i>	1
49	SP1602976	DRIVENHET 10KM/T VE 4PL-18 TDX-SP/SP2	1
50	SP1602998	DRIVENHET 10KM/T HØ 4PL-18 TDX-SP/SP2	1
51	SP1579595	FRIKOPLINGS HÅNDTAK STORM4 OG ALLE SSD MOTOR <i>Single part. Consider to buy another one for the other side if required. - Replaces SP1168668, SP1168669 & SP1169005</i>	1
52	SP1168667	(Declutching lever, cpl. - LH)	1
53	SP1168666	(Declutching lever, cpl. - RH)	1
54	SP1552782	(Propulsion cpl. LH, Gearbox i = 23:1, Motor 4-pole, 340W, slow, 760mm cable, ACS)	1
55	SP1552783	(Propulsion cpl. RH, Gearbox i = 23:1, Motor 4-pole, 340W, slow, 760mm cable, ACS)	1
56	SP1552872	(Propulsion cpl. LH, Gearbox i = 18:1, Motor 4-pole, 340W, fast, 760mm cable, ACS)	1


57	SP1552873	(Propulsion cpl. RH, Gearbox i = 18:1, Motor 4-pole, 340W, fast, 760mm cable, ACS)		1
58	SP1552832	(Propulsion cpl. LH, Gearbox i = 23:1, Motor 4-pole, 340W, slow, 760mm cable, R-Net)		1
59	SP1552833	(Propulsion cpl. RH, Gearbox i = 23:1, Motor 4-pole, 340W, slow, 760mm cable, R-Net)		1
60	SP1552836	(Propulsion cpl. LH, Gearbox i = 18:1, Motor 4-pole, 340W, fast, 760mm cable, R-Net)		1
61	SP1552837	(Propulsion cpl. RH, Gearbox i = 18:1, Motor 4-pole, 340W, fast, 760mm cable, R-Net)		1
62	SP1562498	(Motor, 4-Pol, 340W, slow, 760mm cable, R-Net)		1
63	SP1562499	(Motor, 4-Pol, 340W, fast, 760mm cable, R-Net)		1
64	SP1188842	(DuraWatt Motor, 4-pole, 340W, fast, 760mm cable) <i>TDX-SP2 only!</i>	11.08.2017 →	1
65	SP1188843	(DuraWatt Motor, 4-pole, 340W, slow, 760mm cable) <i>TDX-SP2 only!</i>	11.08.2017 →	1
66	SP1639953	(DuraWatt Propulsion, cpl. - 6 km/h - LH) <i>TDX-SP2 only!</i>	11.08.2017 →	1
67	SP1639952	(DuraWatt Propulsion, cpl. - 6 km/h - RH) <i>TDX-SP2 only!</i>	11.08.2017 →	1
68	SP1639957	(DuraWatt Propulsion, cpl. - 8 km/h - LH) <i>TDX-SP2 only!</i>	11.08.2017 →	1
69	SP1639956	(DuraWatt Propulsion, cpl. - 8 km/h - RH) <i>TDX-SP2 only!</i>	11.08.2017 →	1
70	SP1639959	(DuraWatt Propulsion, cpl. - 10 km/h - LH) <i>TDX-SP2 only!</i>	11.08.2017 →	1
71	SP1639958	(DuraWatt Propulsion, cpl. - 10 km/h - RH) <i>TDX-SP2 only!</i>	11.08.2017 →	1
72	SP1647029	(DuraWatt Propulsion 5K, cpl. - 12km/h - LH (Kite HD))		1
73	SP1647030	(DuraWatt Propulsion 5K, cpl. - 12km/h - RH (Kite HD))		1
74	SP1647214	(DuraWatt Propulsion 5K, cpl. - 12 km/h - LH (Storm4))		1
75	SP1647215	(DuraWatt Propulsion 5K, cpl. - 12 km/h - RH (Storm4))		1
76	SP1647221	(DuraWatt Propulsion 5K, cpl. - 12 km/h - LH (TDX-SP2)) <i>TDX-SP2 only!</i>		1
77	SP1647220	(DuraWatt Propulsion 5K, cpl. - 12 km/h - RH (TDX-SP2)) <i>TDX-SP2 only!</i>		1
78	SP1192560	(DuraWatt Motor 5K, 4-pole, 340W, 760mm cable)		1



1	SP1559356	TT2 MOTOR VE STORM4	1
2	SP1559359	TT2 MOTOR HØ STORM4	1
3	SP1582231	(Mounting Hardware GGB2 Motor)	1







1	-	(Storm4)	0
2	-	(Storm4 TTPlus)	0



1	-	(ACS Electronic)	1
2	-	(LiNX Electronic)	1

(ACS Electronic)

FDX79.no




1	-	(ACS Electronic Modules)	1
2	-	(ACS Remotes & brackets)	1
3	-	(ACS Cables & Accessories)	1
4	-	(ACS Components)	1


-	(Power modules ACS 2)	1
-	(Holder Powermodule) <i>NOT for Ultra Low Maxx Seat!</i>	1
-	(Power module holder / Gyro) <i>only with Ultra Low Maxx Seat!</i>	1
-	(ACTs)	1
-	(Bracket ACT - Ultra Low Maxx Seat)	1
-	(G-Trac Sensor/Gyro)	1
-	(Holder Gyro)	1
-	(Light board)	1

(Power modules ACS 2)

1524244X.dwg



-	SP1555702	(Gasket kit for ACS2, Powermodule and ACT)	1
1	SP1524244	POWERMODUL DX2 PMA90L STORM4/CHASSWHEEL <i>Lock washer between module and bracket!</i>	1
1	SP1537975	POWERMODUL DX2 PMA90LG FOR GYRO <i>Lock washer between module and bracket!</i>	1



1 SP1543129


(Holder Powermodule)

Lock washer between module and bracket!

1

(Power module holder / Gyro)

154261X.kso




1 SP1546769

(Power module holder / Gyro)

Lock washer between module and bracket!


1



1	SP1528542	ACTUATORMODUL ACT2 DX2 STORM4 TDX BORA	1
2	SP1529382	ACTUATORMODUL ACT4 DX2 STORM4 TDX BORA	1

(Bracket ACT - Ultra Low Maxx Seat)


158505 iso



1 SP1584343

(Bracket ACTs (Ultra Low Maxx Seat))

1



1537848X.iso


1 SP1541535

(G-Trac (Gyro) Sensor, cpl.)

1

(Holder Gyro)



1537974X.dwg



1 SP1541625


GYRO MODUL HOLDER STORM4

1




-	SP1524237	KABEL LYS TDX	1
1	SP1528890	LYSFORDELER STORM4	1

(ACS Remotes & brackets)




1	-	(REM 24 SD with brackets) [1601/1602]	1
2	-	(REM A & B with brackets) [1638/1639/1643/1644]	1
3	-	(REM 550 with brackets) [1617/1618]	1
4	-	(REM 34 with brackets) [1611]	1
5	-	(G90A with brackets) [1614/1616]	1
6	-	(Remote G91S & brackets) [1615]	1
7	-	(Tray remote REM 41) [1607]	1
7	-	(Centre control REM 41) [1608]	1
8	-	(Chin Control, small (=> 16/12/31)) [1609]	1
9	-	(Chin Control - Modulite Seat (15/10/01 =>)) [1719]	1
10	-	(Dual Control with brackets) [1291/1292]	1
11	-	(Compact Dual Control with brackets)	1
12	-	(Intuitive Dual Control (IDC) with brackets) [1645]	1


1602274 Rev



1	-	(Remote REM 24SD) [1601/1602]	0
2	-	(Remote Bracket Standard) [1601+1207/1208]	1
3	-	(Remote Bracket Standard - ULM and Soft Armpad on Modulite) [1601+1207/1208+2020]	1
4	-	(Remote Bracket Swing away & height adjust.) [1601+1707/1708]	1
5	-	(Remote Bracket Swing away & height adjust. - ULM and Soft Armpad on Modulite) [1601+1707/1708+2020]	1
6	-	(Remote Bracket Attendant) [1602+1207/1208]	1
7	-	(Remote Bracket Attendant Ultra Low Maxx Seat) [1602+1207/1208+2020]	1
8	-	(Remote Bracket Attendant for Recaro-Seat) [1602+1207/1208]	1



1	1436541	(Remote REM24SD)	1
2	SPZ183229	JOYSTIKK KOMPLETT	1




1559851.iso

1	SPS183040	STYREBOKSHOLDER STD KPL STORM3 G40	1
2	SP1436542	HOLDER DX-REM 24 SD RETT STORM SPORT/3 OPPFL.ARML	1
3	L31001	VINGSKRUV M 6X10 GARANT G24 S	1


(Remote Bracket Standard - ULM and Soft Armpad on Modulite)

1588513.no



1	SP1587309	(Remote bracket Std.)	1
2	SP1587693	(Remote bracket MACES Std - REM 24SD/A/B/550)	1
3	SPS183040	STYREBOKSHOLDER STD KPL STORM3 G40	1
4	L31001	VINGSKRUV M 6X10 GARANT G24 S	1

(Remote Bracket Swing away & height adjust.)




1559852.iso


1	SP1498222	FESTE FOR STYREBOKST TIL SWING AWAY	1
2	SP1496249	SWING AWAY VE NY MOD2008 KPL	1
3	SP1496248	SWING AWAY HØ NY MOD2008 KPL	1

(Remote Bracket Swing away & height adjust. - ULM and Soft Armpad on Modulite)

1588514.no




1	SP1587329	(Remote bracket, swingaway)	1
2	SP1587695	(Remote bracket MACES swingaway - REM 24SD/A/B/550)	1
3	SP1496249	SWING AWAY VE NY MOD2008 KPL	1
4	SP1496248	SWING AWAY HØ NY MOD2008 KPL	1



1	1422806	(Remote Support LH)	1
2	1422810	(Remote Support RH)	1
3	SP1422937	FESTE STYREBOKS RYGGRØR STORM4 TDX SP2	1
4	SP1422936	(Bracket for attendant control RH)	1
5	L31001	VINGSKRUV M 6X10 GARANT G24 S	1


(Remote Bracket Attendant Ultra Low Maxx Seat)

1588437 no



1	SP1588427	BRAKETT LEDSAGERSTYRING STYREBOKS BAK	1
2	SP1587228	HOLDER STYREBOKS LEDS REM 24/A/B/550 ULM	1

(Remote Bracket Attendant for Recaro-Seat)




1563480.ISO

1	1422806	(Remote Support LH)	1
2	1422810	(Remote Support RH)	1
3	SP1422814	(Assembly kit for dual control cpl., Recaro, LH)	1
4	SP1422813	(Assembly kit for dual control cpl., Recaro, RH)	1
5	L31001	VINGSKRUV M 6X10 GARANT G24 S	1
6	1435421	HÅNDTAK Ø22	1

(REM A & B with brackets)


1503649 no



1	SP1523440	(Remote REM A) [1638/1643]	1
2	SP1523445	(Remote REM B) [1639/1644]	1
3	-	(Remote Bracket Standard) [1638/1639+1207/1208]	1
4	-	(Remote Bracket Swing away & height adjust.) [1638/1639+1707/1708]	1
5	-	(Remote Bracket Attendant for Standard-/EasyAdapt-Seat) [1643/1644+1207/1208]	1
6	-	(Remote Bracket Attendant for Recaro-Seat) [1643/1644+1207/1208]	1

(Remote Bracket Swing away & height adjust.)


1559857.iso



1	SP1536180	HOLDER STYREBOKS SHARK/REM A/REM B	1
2	SP1496249	SWING AWAY VE NY MOD2008 KPL	1
3	SP1496248	SWING AWAY HØ NY MOD2008 KPL	1


(Remote Bracket Attendant for Standard-/EasyAdapt-Seat)

1560837.no



1	1422806	(Remote Support LH)	1
2	1422810	(Remote Support RH)	1
3	SP1422937	FESTE STYREBOKS RYGGRØR STORM4 TDX SP2	1
4	SP1422936	(Bracket for attendant control RH)	1
5	L31001	VINGSKRUV M 6X10 GARANT G24 S	1


1602275 Rev






1	SP1523444	STYREBOKS REM550 STORM4 TDX BORA [1617/1618]	1
2	-	(Remote Bracket Standard) [1617+1207/1208]	1
3	-	(Remote Bracket Standard - ULM and Soft Armpad on Modulite) [1617+1207/1208+2020]	1
4	-	(Remote Bracket Swing away & height adjust.) [1617+1707/1708]	1
5	-	(Remote Bracket Swing away & height adjust. - ULM and Soft Armpad on Modulite) [1617+1707/1708+2020]	1
6	-	(Remote Bracket Attendant) [1618+1207/1208]	1
7	-	(Remote Bracket Attendant Ultra Low Maxx Seat) [1618+1207/1208+2020]	1
8	-	(Remote bracket swingaway Quad Link - Ultra Low Maxx Seat) [1617+2075/2076+2020]	1
9	-	(Remote Bracket Attendant for Recaro-Seat) [1618+1207/1208]	1

(Remote bracket swingaway Quad Link - Ultra Low Maxx Seat)

1638916.no




1	SP1606020	(Quad Link, metric - LH)	1
2	SP1591028	(Quad Link, metric - RH)	1
3	SPTRC2060	(Clamp with square tube adapter)	1
4	SP1640530	(Mounting support Quad Link)	1
5	SP1640529	(Screws and magnets kit Quad Link)	1
6	SP1603276	(Remote Mounting Quad Link)	1



1560846.iso

1	SP1440677	(Remote REM 34) [1611]	1
2	-	(Remote Bracket Standard) [1611+1207/1208]	1
3	-	(Remote Bracket Swing away & height adjust.) [1611+1707/1708]	1


(Remote Bracket Standard)



1560844.iso


1	SP1434921	HOLDER FOR REM 34 STORM4	1
2	SPS183040	STYREBOKSHOLDER STD KPL STORM3 G40	1
3	L31001	VINGSKRUV M 6X10 GARANT G24 S	1

(Remote Bracket Swing away & height adjust.)




1560845.iso

1	SP1527125	FESTE FOR REM34 STORM4	1
2	SP1496249	SWING AWAY VE NY MOD2008 KPL	1
3	SP1496248	SWING AWAY HØ NY MOD2008 KPL	1



1	SP1504720	STYREBOKS DX-REM G90A STORM4 [1614/1616]	1
2	-	(Remote Bracket Standard) [1614+1207/1208]	1
3	-	(Remote Bracket Standard - ULM and Soft Armpad on Modulite) [1614+1207/1208+2020]	1
4	-	(Remote Bracket Swing away & height adjust.) [1614+1707/1708]	1
5	-	(Remote Bracket Swing away & height adjust. - ULM and soft Armpad on Modulite) [1614+1707/1708+2020]	1
6	-	(Remote Bracket Attendant) [1616+1207/1208]	1
7	-	(Remote Bracket Attendant Ultra Low Maxx Seat) [1616+1207/1208+2020]	1
8	-	(Remote Bracket Attendant for Recaro-Seat) [1618+1207/1208]	1




1562939.ISO

1	SP1504822	(Remote Support Standard)	1
2	SPS183040	STYREBOKSHOLDER STD KPL STORM3 G40	1
3	L31001	VINGSKRUV M 6X10 GARANT G24 S	1


(Remote Bracket Standard - ULM and Soft Armpad on Modulite)

1588518.no



1	SP1587310	(Remote bracket Std G90)	1
2	SP1587694	(Remote bracket MACES Std - G90)	1
3	SPS183040	STYREBOKSHOLDER STD KPL STORM3 G40	1
4	L31001	VINGSKRUV M 6X10 GARANT G24 S	1

(Remote Bracket Swing away & height adjust.)




1560310.iso

1	SP1517352	FESTE FOR G90 SWINGAWAY STORM4 TDX BORA	1
2	SP1496249	SWING AWAY VE NY MOD2008 KPL	1
3	SP1496248	SWING AWAY HØ NY MOD2008 KPL	1

(Remote Bracket Swing away & height adjust. - ULM and soft Armpad on Modulite)


1588519.no



1	SP1587331	(Remote bracket G90)	1
2	SP1587696	(Remote bracket MACES, swingaway - G90)	1
3	SP1496249	SWING AWAY VE NY MOD2008 KPL	1
4	SP1496248	SWING AWAY HØ NY MOD2008 KPL	1

(Remote Bracket Attendant)


1560853.iso



1	SP1509435	(Remote Support attendant control G90A, LH)	1
2	SP1509439	(Remote Support attendant control G90A, RH)	1
3	SP1422937	FESTE STYREBOKS RYGGRØR STORM4 TDX SP2	1
4	SP1422936	(Bracket for attendant control RH)	1
5	L31001	VINGSKRUV M 6X10 GARANT G24 S	1


(Remote Bracket Attendant Ultra Low Maxx Seat)

1588438 No



1	SP1588427	BRAKETT LEDSAGERSTYRING STYREBOKS BAK	1
2	SP1587229	(Adapter G90 (with screws))	1

(Remote Bracket Attendant for Recaro-Seat)





1563503.iso

1	SP1509435	(Remote Support attendant control G90A, LH)	1
2	SP1509439	(Remote Support attendant control G90A, RH)	1
3	SP1422814	(Assembly kit for dual control cpl., Recaro, LH)	1
4	SP1422813	(Assembly kit for dual control cpl., Recaro, RH)	1
5	L31001	VINGSKRUV M 6X10 GARANT G24 S	1
6	1435421	HÅNDTAK Ø22	1


(Remote G91S & brackets)

1560680 iso




1	SP1505187	(Remote Control G91S) [1615]	1
2	-	(Remote Bracket Standard) [1615+1207/1208]	1

(Remote Bracket Standard)




1	SP1509690	HOLDER FOR DX-REMG90 HØ JUSTERBAR HØYDE/DYBDE		1
2	SP1527504	(Swan neck remote bracket with adapter plate)		1
3	SP1603801	(Mounting bracket, cpl.) <i>Old bracket no longer available!</i>	05.08.2013 →	1
4	SP1603847	(Mounting Adapter G91) <i>Old bracket no longer available!</i>	05.08.2013 →	1

(Tray remote REM 41)



1	L36999	JOYSTICK KULE DIA32mm DX-REM41	1
2	Z183001A	JOYSTICK GUMMI ACS	1
3	SP1527930	BORDSTYRING DX-REM41-D DOLPHIN MASTER	1
4	Z177026	(SUPPORT BRACKET)	1
5	N04109	KULELEDD HØ SERVO G40	1
6	SP4100136	(SUPPORT ROD, CPL..)	1
7	U89010A.99	(ROD. M. L=560 MM)	1
8	D20109	(Countersunk bolt DIN966-M05x12 bk)	1
9	SPZ177011	HOLDER FOR BORD SVINGBAR KPL STORM3/STORM4	1
10	Z177025	FESTE KIT BOR STORM4	1
11	SP1541444	(Sliding bar for tray 540mm)	1
12	SP1541413	(Sliding bar for tray 620mm)	1
13	SP1541414	(Sliding bar for tray 660mm)	1
14	5009880	(Table 540mm, without cut swing away)	1
15	5009682	(Table 620mm, without cut swing away)	1
16	Z17024A	BORD HEL PLATE STORM3	1





1	SP1530588	DX-REM41 MED SENTERFESTE TDX	1
1	SP1530589	(CENTRAL STEERING CPL., SWING-AWAY LH, REM41-D)	1
2	Z177024.00	FESTE F MITTSTYRNING,	1
3	F60208	(Cap nut DIN1587-M08 zc)	1
4	L01037	(Sks Nut Cap M12)	1
5	F20106	MUTTER SEKSKANT SELVLÅSENDE G40	1
6	G00104	SKIVE 13X24X2,5 G40	1
7	Z17032A	(Bush Rd. 19 X 3.5 X 25)	1
8	L01055	ENDESTYKKE FESTEBR.	1
9	L01008	(Tube End Plug SI 22)	1
10	Z183216.P40	FESTE F MANBOX MITTST,	1
11	1559862	KABEL BUS 1,7 M	1
12	Z183089.00	BREDD JUST ARM MITTST,	1
13	Z183086.00	ÄNDRÖR HÖ, STORM	1
13	Z183092.00	ÄNDRÖR VÅ, STORM	1
14	H281005A.00	(Support rail, chrome plated)	1
15	5006720	KLEMME 22MM TDX	1
16	G00106	SKIVE M8 DRIVHJUL MIRAGE	1
17	SP1527930	BORDSTYRING DX-REM41-D DOLPHIN MASTER	1
23	L36999	JOYSTICK KULE DIA32mm DX-REM41	1
25	Z183001A	JOYSTICK GUMMI ACS	1
26	Z183289	SENTERSTYRT STØTTE RØR KPL	1
27	L31001	VINGSKRUV M 6X10 GARANT G24 S	1
28	Z183289	SENTERSTYRT STØTTE RØR KPL	1

(Chin Control, small (=> 16/12/31))


-	(Small Chincontrol)	1
-	(Framework for chincontrol)	1

1552295X.iso



1	1448705	\$HAKESTYRING KPL MONT EKSTRA JOYSTICK	1
1	1421770	(Chin control G80/G90, cpl.)	1
1	SP1544255	(Chin control G91, cpl.)	1
2	Z183282	(Joystick, DX-RJM (with Binder plug, 4-pole))	1
3	Z183265	BRUK GJS26826 DYNAMICS	1
4	1442943	(Switchjoystick for REM 24 SD, cpl.)	1
4	Z183283	(Switchjoystick G80/G90, cpl.)	1
4	SP1544250	(Switchjoystick for G91S, cpl.)	1
5	SP1552294	(Switchjoystick bracket)	1
6	SP1552295	(Joystick clamp)	1
7	L32901	SPAK M8 STORM4	1
8	SP1442944	(Harness Chin Control for REM 24 SD)	1
8	1442963	\$KABELSTAMME VERTIKAL SWINGAWAY HAKESTYRING	1
9	Z183275	(Harness Chin Control for G80/G90)	1
10	SP1541715	(Harness Chin Control for G91S)	1

(Framework for chincontrol)





1422942X.iso



2	-	(No longer available.) <i>SP1422944 / SP1429945 / SP1422946 no longer available</i>	1
1	SP1422943	HAKESTYRING RAMMEVERK HORIZONTAL SWING AWAY VE	1
1	SP1422942	(Framework for chin control, manual, horizontal swing away, RH, cpl.)	1
-	SP1552335	(Framework for chin control, powered, horizontal swing away, LH (w/o actuator & switch)) <i>Not shown / Nicht abgebildet / Non illustré</i>	1
-	SP1552336	(Framework for chin control, powered, horizontal swing away, RH (w/o actuator & switch)) <i>Not shown / Nicht abgebildet / Non illustré</i>	1
3	SP1552337	(Actuator for horizontal swing away chin control, cpl.)	1
3	SP1552418	(Actuator for vertical swing away chin control, cpl.) <i>Not shown / Nicht abgebildet / Non illustré</i>	1

(Chin Control - Modulite Seat (15/10/01 =>))

1591328.lin




-	-	(If retrofitting, the user manual must be ordered separately!)	1
-	-	(User Manuals Chin Control (2015-10-01 =>))	1
1	-	(Chin Control DX-RJM)	1
2	-	(Mounting Hardware - Chin Control)	1



1	-	(Chin Control powered - Details)	1
2	-	(Chin Control manual - Details)	1
-	SP1601389	(Kit, Swing-away mounting arm, powered, LH, cpl. with stereo-switch (no joystick nor harness included)) <i>Please order respective harness & joystick separately!</i>	1
-	SP1601390	(Kit, Swing-away mounting arm, powered, RH, cpl. with stereo-switch (no joystick nor harness included)) <i>Please order respective harness & joystick separately!</i>	1
-	SP1601391	(Kit, Swing-away mounting arm, manual, LH, cpl. (no joystick nor harness included)) <i>Please order respective harness & joystick separately!</i>	1
-	SP1601408	(Kit, Swing-away mounting arm, manual, RH, cpl. (no joystick nor harness included)) <i>Please order respective harness & joystick separately!</i>	1
-	SP1601528	(Kit, Joysticks & Harness for REM 24 SD)	1
-	SP1601529	(Kit, Joysticks & Harness for G91S)	1


(Chin Control powered - Details)

1591297.jpg



-	SP1605661	(Fixation Kit (cable ties, spiral coiled tube & end caps))	01.04.2016 →	1
1	SP1580094	(Actuator cpl. LH)		1
2	SP1591292	(Actuator cpl. RH)		1
3	SP1605659	(Remote Rod Assy) <i>New designed bracket (2016-04-01 =>) for Modulite Seat needed!</i>	01.04.2016 →	1
4	SP1589144	(Rod cylinder connection chin, cpl.)		1
5	SP1605660	(Joystick Rod Assy)	01.04.2016 →	1
6	SP1589146	(Mounting arm for Switch Joystick)		1
7	SP1589142	(Mounting hardware for DX-RJM, cpl.)		1
8	Z183282	(Joystick, DX-RJM (with Binder plug, 4-pole))		1
9	1442943	(Switchjoystick for REM 24 SD, cpl.)		1
9	SP1544250	(Switchjoystick for G91S, cpl.)		1
10	Z183265	BRUK GJS26826 DYNAMICS		1
11	SP1442944	(Harness Chin Control for REM 24 SD)		1
11	SP1541715	(Harness Chin Control for G91S)		1
12	SP1580701	(Stereo-switch box, cpl. with mounting hardware)		1
13	SP1589145	(Adjustment Screw)		1
14	SP1601528	(Kit, Joysticks & Harness for REM 24 SD)		1
14	SP1601529	(Kit, Joysticks & Harness for G91S)		1
15	SP1602721	(End cap A12, black)	01.04.2016 →	3


(Chin Control manual - Details)



-	SP1605661	(Fixation Kit (cable ties, spiral coiled tube & end caps))	01.04.2016 →	1
1	SP1577456	(ASM Chin Control man. LH, cpl.)		1
2	SP1591293	(ASM Chin Control man. RH, cpl.)		1
3	SP1589144	(Rod cylinder connection chin, cpl.)		1
4	SP1605660	(Joystick Rod Assy)	01.04.2016 →	1
5	SP1589146	(Mounting arm for Switch Joystick)		1
6	SP1589142	(Mounting hardware for DX-RJM, cpl.)		1
7	Z183282	(Joystick, DX-RJM (with Binder plug, 4-pole))		1
8	1442943	(Switchjoystick for REM 24 SD, cpl.)		1
8	SP1544250	(Switchjoystick for G91S, cpl.)		1
9	Z183265	BRUK GJS26826 DYNAMICS		1
10	SP1442944	(Harness Chin Control for REM 24 SD)		1
10	SP1541715	(Harness Chin Control for G91S)		1
11	SP1580666	(Adjustment Handle)		1
12	SP1601528	(Kit, Joysticks & Harness for REM 24 SD)		1
12	SP1601529	(Kit, Joysticks & Harness for G91S)		1
13	SP1602721	(End cap A12, black)	01.04.2016 →	2

(Mounting Hardware - Chin Control)


1591309 iso



1	SP1577173	(Mounting Hardware, Chin Control - Modulite Seat) <i>Design change 2016-04-01 =></i>	1
2	SP1604355	(Clamp Lever M8x40)	1

(Dual Control with brackets)


1564180 no



1	-	(No longer available.) [1291/1292]	1
2	-	(Remote Bracket Dual Control for Standard-/EasyAdapt-Seat) [1291/1292]	1
3	-	(Remote Bracket Dual Control for Modulite Seat)	1
4	-	(Remote Bracket Dual Control for Recaro Seat) [1291/1292]	1

(Remote Bracket Dual Control for Standard-/EasyAdapt-Seat)

1560885.iso




1 SP1422818

FESTEBRAKETT LEDSAGERSTY RYGGMONT ENEKL

1


(Remote Bracket Dual Control for Modulite Seat)



1574201.iso

1	SP1571799	(Remote Bracket Dual Control LH)	1
2	SP1571800	(Remote Bracket Dual Control RH)	1

(Remote Bracket Dual Control for Recaro Seat)




1563504.iso


1	SP1422819	(Assembly kit for dual control Recaro)	1
2	1435421	HÅNDTAK Ø22	1

(Compact Dual Control with brackets)

1588511.000



1	-	(Compact Dual Control) [1720]		1
2	-	(Bracket, Compact Dual Control - Ultra Low Maxx Seat) [1207/1208 + 1720 for 2020]	03.06.2015 →	1
3	-	(Remote Bracket, Compact Dual Control/ACU - short)	01.07.2015 →	1
4	-	(Remote Bracket, Compact Dual Control/ACU - long)	01.07.2015 →	1




1 SP1586988

(Remote DX2-ACU)

1

(Bracket, Compact Dual Control - Ultra Low Maxx Seat)

1588512.no




1 SP1586795

HOLDER STYREBOKS DUAL CONTROL

1

(Remote Bracket, Compact Dual Control/ACU - short)

1588857 .ind




1 SP1587904

(Remote Attachment - short)

1

(Remote Bracket, Compact Dual Control/ACU - long)

158856 iso



1 SP1587119

(Remote Attachment - long)


1

(Intuitive Dual Control (IDC) with brackets)

-	(IDC with bracket)	1
-	(IDC with bracket (Ultra Low Maxx Seat))	1
-	(Height adjustable Bracket, IDC)	1

(IDC with bracket)


1602413.iss



-	(The IDC incl. mounting hardware has been revised, the old hardware is no longer available.)	01.01.2016 →	0
-	(If retrofitting, the user manual must be ordered separately!)		1
-	(User Manuals IDC (ACS))		1
1	SP1603425 INTUITIV LEDSAGERSTYRING ALLE STOLER MED ACS		1
2	SP1603426 HOLDER INTUITIV LEDSAGERSTYRING NY TYPE		1
3	SP1601919 GUMMIFORING FESTE IDC STORM4 MODULITE		1

(IDC with bracket (Ultra Low Maxx Seat))


1602412.lin



-	(The IDC incl. mounting hardware has been revised, the old hardware is no longer available.)	01.01.2016 →	0
- -	(If retrofitting, the user manual must be ordered separately!)		1
- -	(User Manuals IDC (ACS))		1
1	SP1603425 INTUITIV LEDSAGERSTYRING ALLE STOLER MED ACS		1
2	SP1603427 (IDC Bracket, cpl. with mounting hardware for Ultra Low Maxx Seat)		1
3	SP1603426 HOLDER INTUITIV LEDSAGERSTYRING NY TYPE		1
4	SP1601919 GUMMIFORING FESTE IDC STORM4 MODULITE		1

(Height adjustable Bracket, IDC)

1605484 iso



1 SP1604195


(Height adjustable bracket for IDC, cpl.)

1

(ACS Cables & Accessories)

-	(Bus-Cables ACS)	1
-	(T-Adapter)	1
-	(T-Adapter (Ultra Low Maxx Seat))	1
-	(Backrest Sensor - Ultra Low Maxx Seat)	1
-	(Tilt Sensor - Ultra Low Maxx Seat)	1
-	(Programming Kits)	1
-	(Magnetkey for Remotes)	1


(Bus-Cables ACS)



Position	L [mm]
1	300
2	500
3	700
4	1100
5	1700
6	2000
7	2500
8	2700
9	1060

1	GSM63003	KABEL BUS 0,3M ACS,	1
2	SPGSM63005	(Bus cable ACS - 500mm) <i>Replace GSM63004</i>	1
3	GSM63007	ACS-BUSKABEL L=700 MM	1
4	GSM63011	KABEL-BUS 1,1M ACS, ST utgått bruk 63010	1
5	1559862	KABEL BUS 1,7 M <i>Replace GSM63016</i>	1
6	GSM63020	DX BUS CABLE 2.0M STORM <i>Replace 1419386</i>	1
7	1440424	(Bus cable ACS, 2500mm)	1
8	SP1559873	(Bus cable ACS - 2700mm) <i>Replace GSM63026</i>	1
9	SP1550402	(Extension cable ACS - 1060mm)	1
-	GSM9887	DX-BUS GURLY CABLE	1


1534069X.dwg



1	SP1559812	HOLDER DX-SKT4 STORM4 FRA 010512	20.04.2012 →	1
1	SP1537966	\$HOLDER BRUK SP1559812 STORM4	→ 19.04.2012	1
2	SP1537967	FORGRENING BUSS SKT4 ACS		1

(T-Adapter (Ultra Low Maxx Seat))


1586599 No



1	SP1586798	(T-Adapter with mounting hardware)	1
2	P10072	(Bus 4-way socket (T-Adapter))	1

(Backrest Sensor - Ultra Low Maxx Seat)


1602018.kim




1	SP1589753	(Switch bracket)	1
2	SPTRE0753	(Tilt Switch)	2

(Tilt Sensor - Ultra Low Maxx Seat)

1602019.kim




1	SP1602017	(Switch bracket)	1
2	SPTRE0753	(Tilt Switch)	1



1	1449280	(USW wizard dongle kit, EDL-USW, cpl.)	1
2	6001450	(Hand Programmer, Instructions in French, German, Spanish)	1
	SP1441551	(Programming Adapter Shark DK-ADAPT (necessary to program Shark))	1

(Magnetkey for Remotes)



1564806.iso


1 SP1419695 (Magnet key with bag)

1


-	(iPortal Modules and Accessories)	1
-	(Remote, Joystick, Control)	1
-	(Joystick Knobs)	1
-	(Buddy Button)	1
-	(Piko Button)	1
-	(Interface, Adapter, Other)	1

(iPortal Modules and Accessories)

1603-435 no



1	-	(No longer available.) [1629] obsolete	0
2	-	(No longer available.) [1624] obsolete	1
3	-	(No longer available.) [1662] obsolete	1
4	-	(No longer available.) [1634] obsolete	1
5	SP1545531	(Gooseneck short for iPod touch & iPhone) [1635]	1
6	SP1545532	(Gooseneck long for iPod touch & iPhone) [1636]	1




1	-	(No longer available.) [1641/1642]	1
2	SP1586989	(Remote Joystick) [1621]	1
3	SP1586990	STYREBOKS LEDSAGER DX2 LOW-FORCE [1722]	1
4	-	(No longer available.) [1620]	1
5	P10175	(Mini-Joystick with Push function) [1622]	1
-	SP1442679	(Head-Remote) [1613]	1
-	SP1603442	(Cable for Head-Remote ultrasonic porximity sensor)	1

(Joystick Knobs)

1603437.no




1	L36999	JOYSTICK KULE DIA32mm DX-REM41 <i>[1650]</i>	1
2	P10397	PROPP HAKESTYRING <i>[1651]</i>	1
3	P39302	JOYSTICK KNOTT CARROT ACS <i>[1652]</i>	1
4	P39303	T-GREP ACS STORM4 MODULITE <i>[1653]</i>	1
5	P39305	(Joystick Knob Flat T-Plate) <i>[1654]</i>	1
6	P10306	\$PROGRAMMERINGSSENHET DM4 GARANT 63E PRO <i>[1655]</i>	1
7	Z183265	BRUK GJS26826 DYNAMICS <i>[1656]</i>	1
8	P39306	(Joystick Knob Straight Bar) <i>[1657]</i>	1
9	P39307	(Joystick Knob Finger Tube) <i>[1658]</i>	1
10	P39304	(Joystick Knob Rount T-Bar) <i>[1659]</i>	1




1	SP5813964	(Buddy Button - red) <i>[1685]</i>	1
2	SP5813966	(Buddy Button - yellow) <i>[1686]</i>	1
3	SP5813965	(Buddy Button - green) <i>[1687]</i>	1
4	SP5813981	BUDDY BUTTOM BLÅ EA STORM4 EASY ADAPT <i>[1690]</i>	1

(Piko Button)

1603436.iss






1	SP1535717	PIKO BUTTON RØD DYNAMIC <i>[1711]</i>	1
2	SP1535718	Piko Button Yellow <i>[1712]</i>	1
3	SP1535723	PIKO BUTTON GRØNN DYNAMIC <i>[1713]</i>	1
4	SP1535724	PICO BRYTER ORANGE STORM4 <i>[1714]</i>	1
5	SP1535725	PIKO BUTTON BLÅ EA STORM4 EASYADAPT <i>[1715]</i>	1
6	SP1535726	PIKO BUTTON SORT DYNAMIC <i>[1716]</i>	1



1	SP1544036	(4 Switch Interface) [1630]	1
2	P10163	(Sip & Puff Module(DX-SNP)) [1633]	1
-	SPZ183267	(Bracket for Sipp and Puff) [1668]	1
3	P10072	(Bus 4-way socket (T-Adapter)) [1680]	1
4	P30080	(5-1 Switch Adaptor) [1681]	1
5	P40107	(Charger Socket Modul) [1661]	1

(LiNX Electronic)


1502714.no



1	-	(Modules)	1
2	-	(Remotes & Brackets)	1
3	-	(Cables & Accessories)	1
4	-	(Components)	1

1	-	(Power Modules)	1
2	-	(Bracket, Power Module)	1
3	-	(Actuator Modules (ACTs))	1
4	-	(G-Trac Sensor / Gyro)	1

(Power Modules)




① 120A - Drive only

② 120A - Drive, Actuators & Light

1639232.iso

1	SP1602834	(Power Module LiNX DLX-PM120 120A - Drive only)	1
2	SP1602837	(Power Module LiNX DLX-PM120AL 120A - Drive, Actuators & Light)	1





1 SP1639173

(Bracket, Power Module LiNX)


1

(Actuator Modules (ACTs))

1604610 iso




1	SP1602883	(LiNX 2-Channel Actuator Module (DLX-ACT200))	1
2	SP1602884	(LiNX 4-Channel Actuator Module (DLX-ACT400))	1



1	SP1607817	(G-Trac Sensor (Gyro) LiNX, cpl. with bracket)	1
2	SP1602885	(G-Trac Sensor (Gyro) LiNX DLX-GYR100)	1

(Remotes & Brackets)


-	(Remotes & Brackets)	1
-	(Dual Controls & Brackets)	1
-	(Additional Remotes & Brackets)	1



1642452.iso


1	-	(Remote REM110)	1
2	-	(Remote REM211)	1
3	-	(Remote REM216)	1
4	-	(Remote REM400)	1
5	-	(Remote REM500)	1
6	-	(Remote Bracket, Standard)	1
7	-	(Remote Bracket, Swing-away)	1
8	-	(Remote Bracket, Attendant)	1
9	-	(Remote Bracket, Nucleus)	1
10	-	(Remote Bracket, Swing-away REM500)	1
11	-	(Remote Bracket, Nucleus REM500)	1

Invacare® Storm4 | Elektronikk
(Remote REM110)



1605135.iso


1	SP1602806	(LiNX Remote DLX-REM110 - Drive Only)	1
2	SP1607776	(LiNX Joystick Knob REM050/110/211/216)	1
3	SP1653457	(LiNX Gaiter REM050/110/211/216)	1



1605136.iso


1	SP1602825	(LiNX Remote DLX-REM211 - Drive and Actuators)	1
2	SP1607776	(LiNX Joystick Knob REM050/110/211/216)	1
3	SP1653457	(LiNX Gaiter REM050/110/211/216)	1

Invacare® Storm4 | Elektronikk
(Remote REM216)



1605137.iso


1	SP1602826	(LiNX Remote DLX-REM216 - Drive, Actuators and Lights)	1
2	SP1607776	(LiNX Joystick Knob REM050/110/211/216)	1
3	SP1653457	(LiNX Gaiter REM050/110/211/216)	1



1641559.iso

1	SP1647646	(LiNX Remote DLX-REM400-B (Rev. B)) <i>replaces: SP1602827 DLX-REM400</i>	1
2	SP60109657	(Joystick Knob REM400)	1
3	-	(Toggle Switches for REM400) <i>DLX-REM400-B is mandatory for Toggle Switches!</i>	1
4	SP1647187	(LiNX Remote DLX-REM400-B with Toggle Switches, cpl.)	1
5	SP1654025	(LiNX DLX-REM400 Gaiter)	1

(Toggle Switches for REM400)




1647650.iso

1 SP1647186

(DLX-SW400TGL Toggle Switches for DLX-REM400-B)
Not to be used with old DLX-REM400!

1




1 SP1602828

(LiNX Remote DLX-REM500)

1


(Remote Bracket, Standard)



1605138.iso

1	SP1603417	(LiNX Remote bracket, Standard)	1
2	SP1603419	(Remote support tube)	1

(Remote Bracket, Swing-away)




1605139.iso

1	SP1496249	SWING AWAY VE NY MOD2008 KPL	1
2	SP1496248	SWING AWAY HØ NY MOD2008 KPL	1
3	SP1603418	(LiNX Remote support for height adj. & swing-away bracket)	1
4	SP60110244	(LiNX Remote Support for height adj. & swing-away bracket (for soft armpad))	1

(Remote Bracket, Attendant)

1605185 iso




1 SP1603363

(LiNX Remote bracket, Attendant)

1


(Remote Bracket, Nucleus)



1654655.iso

1	SP1582573	(Nucleus Center Mount ASL660C)	1
2	SP1637863	(Adapter Kit REM110-400 on Nucleus)	1
3	SP1640858	(Adapter Ø22mm round to square)	1


(Remote Bracket, Swing-away REM500)



1641875.iso

1	SP1640294	(Remote Bracket REM500, multi-adjustable - RH)	1
2	SP1639960	(Remote Bracket REM500, multi-adjustable - LH)	1
3	SP1603419	(Remote support tube)	1


(Remote Bracket, Nucleus REM500)



1641876.iso


1	SP1582573	(Nucleus Center Mount ASL660C)	1
2	SP1640564	(Mounting Bracket for REM500)	1
3	SP1640858	(Adapter Ø22mm round to square)	1

(Dual Controls & Brackets)



1642453.iso


1	-	(Dual Control REM050)		1
2	-	(Dual Control ACU200)	11.08.2017 →	1
3	-	(Intuitive Dual Control (IDC))	01.01.2018 →	1
4	-	(Remote Bracket, Attendant)		1
5	-	(Remote Bracket, Compact Dual Control/ACU - short)	11.08.2017 →	1
6	-	(Remote Bracket, Compact Dual Control/ACU - long)	11.08.2017 →	1
7	-	(Mounting Bracket, IDC)	01.01.2018 →	1
8	-	(Height adjustable Bracket, IDC)	01.01.2018 →	1



1605186.iso

1	SP1602792	(Dual Control LiNX DLX-REM050)	1
2	SP1607776	(LiNX Joystick Knob REM050/110/211/216)	1
3	SP1653457	(LiNX Gaiter REM050/110/211/216)	1

Invacare® Storm4 | Elektronikk
(Dual Control ACU200)




1641652.iso

1 SP1607393

(Compact Dual Control LiNX DLX-ACU200)


1



1645696.iso

-	(If retrofitting, the user manual must be ordered separately!)	1
-	(User Manuals IDC (LiNX))	1
1	SP1645692 (Intuitive Dual Control LiNX, cpl.)	1

(Mounting Bracket, IDC)




1645697.ISO

1	SP1603426	HOLDER INTUITIV LEDSAGERSTYRING NY TYPE	1
2	SP1601919	GUMMIFORING FESTE IDC STORM4 MODULITE	1

(Additional Remotes & Brackets)



-	(Compact Remote LiNX & Bracket)	1
-	(Chin Control, man.)	1
-	(Chin Control, pow. - Modulite Seat)	1
-	(Head Controls (ASL104))	1
-	(Wireless Accessories (ASL557/558))	1
-	(Micro Extremity Control (ASL130))	1
-	(Cable Mounting Sheets ASL - Modulite Back)	1

(Compact Remote LiNX & Bracket)




1642454.iso

1	-	(Compact Remotes LiNX)	1
2	-	(Bracket, Standard)	1
3	-	(Bracket, Nucleus)	1




1	SP1602829	(Proportional Compact Remote DLX-CR400)	1
2	SP1607394	(Proportional Compact Remote Low Force DLX-CR400LF)	1

(Bracket, Standard)




1641879.iso

1	SP1637891	(Remote Bracket Standard)	1
2	SPS183040	STYREBOKSHOLDER STD KPL STORM3 G40	1




1641878.iso

1	SP1582573	(Nucleus Center Mount ASL660C)	1
2	SP1637862	(Remote Bracket for Nucleus)	1
3	SP1640858	(Adapter Ø22mm round to square)	1



1641774.iso



1	SP1640066	(Swing Away Mechanism, Atom Head Array)	1
2	-	(Egg Switch ASL300)	1
3	SP1640527	(Chin mount for ASL138)	1
4	SP1640702	(Switch Mount for Swing Away Mechanism)	1
5	SP1607766	(Extremity Control Joystick ASL 138 LX)	1
6	SP1640291	(Assy, ASL Chinmount for ASL138 Joystick (incl. ASL300))	1
7	SP1654755	(Foam Joystick knobs for ASL-138 (pack of 3))	1



1642049.iso


1	SP1640086	(Egg Switch ASL300, black)	1
2	SP1646840	(Extension Cable ASL807 (Mono) 72"/183cm lenght)	1

(Chin Control, pow. - Modulite Seat)



1645849.iso


1	-	(Chin Control, powered)	1
2	-	(Mounting Hardware, Chin Control (Modulite Seat))	1



1645360.ISO


-	SP1602891	(LiNX Bus cable 1000mm (GLM-BUS100))	1
1	SP1580094	(Actuator cpl. LH)	1
2	SP1591292	(Actuator cpl. RH)	1
3	SP1605659	(Remote Rod Assy)	1
4	SP1589144	(Rod cylinder connection chin, cpl.)	1
5	SP1605660	(Joystick Rod Assy)	1
6	SP1642589	(Flip Switch Box LiNX)	1
7	SP1602721	(End cap A12, black)	1
8	SP1642566	(Switch Module, cpl.)	1
9	SP1535717	PIKO BUTTON RØD DYNAMIC	1
9	SP1535725	PIKO BUTTON BLÅ EA STORM4 EASYADAPT	1
10	SP1607766	(Extremity Control Joystick ASL 138 LX)	1
11	SP1645740	(Bracket, ASL 138 LX)	1
12	SP1589145	(Adjustment Screw)	1
13	SP1654755	(Foam Joystick knobs for ASL-138 (pack of 3))	1

(Mounting Hardware, Chin Control (Modulite Seat))



1645361.iso


1	SP1577173	(Mounting Hardware, Chin Control - Modulite Seat)	1
2	SP1642569	(Bracket, ASL Interface Box - LH (Modulite))	1
3	SP1642045	(Bracket, ASL Interface Box - RH (Modulite))	1



1653436.iso

1	-	(ATOM Mini ASL104MLX Head Array LiNX)	1
2	-	(ATOM Proton ASL104PTLX Head Array LiNX)	1
3	-	(Egg Switch ASL300)	1

(ATOM Mini ASL104MLX Head Array LiNX)




1643205.iso

1 SP1640089

(Atom Head Array LiNX - Mini ASL104M LX)

1


(ATOM Proton ASL104PTLX Head Array LiNX)



1643206.iso

1	SP1607755	(Atom Head Array LiNX - Proton ASL104PT LX)	1
2	SP1640066	(Swing Away Mechanism, Atom Head Array)	1
3	SP1640702	(Switch Mount for Swing Away Mechanism)	2



Invacare® Storm4 | Elektronikk
(Egg Switch ASL300)



1653459.iso

1	SPTC1273	(Multi-Mount Hardware Kit)	1
2	-	(Egg Switch ASL300)	1


(Wireless Accessories (ASL557/558))



1642369.iso

1	SP1640061	(Wireless Triple Switch Receiver ASL557-3X)	1
2	SP1640062	(ATOM Wireless Mouse Emulator ASL558X)	1


(Micro Extremity Control (ASL130))



1653437.ISO

1	SP1607756	(MEC Proportional Joystick ASL130 (LiNX), cpl. w/ Interface Box)	1
2	SP1642569	(Bracket, ASL Interface Box - LH (Modulite))	1
3	SP1642045	(Bracket, ASL Interface Box - RH (Modulite))	1
4	-	(Bracket, MEC - Ring Clamp Mount, Nucleus)	1
5	-	(Bracket, MEC - Tray Mount, Nucleus)	1
6	-	(Bracket, MEC - Tray Mount, Standard)	1


(Bracket, MEC - Ring Clamp Mount, Nucleus)



1641725.iso

1	SP1582573	(Nucleus Center Mount ASL660C)	1
2	SP1582572	(Ring Clamp ASL644 for MEC ASL130)	1
3	SP1640858	(Adapter Ø22mm round to square) <i>not for ULM Seat!</i>	1


(Bracket, MEC - Tray Mount, Nucleus)



1641726.iso

1	SP1582573	(Nucleus Center Mount ASL660C)	1
2	SP1640454	(Hitch Mount for Tray)	1
3	SP1640290	(Tray for MEC ASL130)	1
4	SP1640858	(Adapter Ø22mm round to square) <i>not for ULM Seat!</i>	1


(Bracket, MEC - Tray Mount, Standard)



1653451.iso

1	SP1640290	(Tray for MEC ASL130)	1
2	SP1637858	(Adapter Kit, square tube)	1
3	SP1603419	(Remote support tube)	1
4	SP1642370	(Clamp Lever M6x20)	1

(Cable Mounting Sheets ASL - Modulite Back)



1647092.ISO

1 SP1647087

(Cable Sheets on Modulite Back for ASL components)



1

-	(Bus Cables)	1
-	(Actuator Cable)	1
-	(Lifter/Tilt-Module Cables)	1
-	(Light Cable)	1
-	(Light Board)	1
-	(4-way Expansion Block)	1
-	(Slow-Stop-Sensor - Modulite Seat)	1
-	(Slow-Stop-Sensor - Recaro Seat)	1
-	(Accessory cables)	1
-	(Stereo-Mono-Splitter Cable)	1
-	(Programming Access Key)	1

(Bus Cables)

LiNX

1605413.jpg



Bus cable

Lenght [mm]


①	300
②	500
③	700
④	1000
⑤	1200
⑥	1500
⑦	1700
⑧	2000
⑨	2500

Extension cable

Lenght [mm]

⑩	640
⑪	900


1	SP1602888	(LiNX Bus cable 300mm (GLM-BUS030))	1
2	SP1602889	(LiNX Bus cable 500mm (GLM-BUS050))	1
3	SP1602890	(LiNX Bus cable 700mm (GLM-BUS070))	1
4	SP1602891	(LiNX Bus cable 1000mm (GLM-BUS100))	1
5	SP1603489	(LiNX Bus cable 1200mm (GLM-BUS120))	1
6	SP1602892	(LiNX Bus cable 1500mm (GLM-BUS150))	1
7	SP1602893	(LiNX Bus cable 1700mm (GLM-BUS170))	1
8	SP1602894	(LiNX Bus cable 2000mm (GLM-BUS200))	1
9	SP1602904	(LiNX Bus cable 2500mm (GLM-BUS250))	1
10	SP1602905	(LiNX Extension cable 640mm (GLM-EXT064))	1
11	SP1602906	(LiNX Extension cable 900mm (GLM-EXT090))	1




!	-	(No longer available - please refer to Replacement Info for which part is to be used instead.)	0
1	SP1607105	(LiNX Actuator cable LA23 - 950mm) <i>Only Backrest Recliner - Replaces SP1605680</i>	1
1	SP1607124	(LiNX Actuator cable LA23 - 1500mm) <i>Replaces SP1605815</i>	1

(Lifter/Tilt-Module Cables)

150665 iso




1	SP1607899	(Sensor Cable LiNX - 300mm)	1
2	SP1573204	ADAPTER KABEL TYCO ACT MODULITE LNX	1
2	SP1576062	(Extension Cable (1100mm))	1



1 SP1603336 (LiNX Light cable)


1



1 1415004

LYSFORDELER PRINT G40


1



1	SP1607842	(4-way Expansion Block LiNX, cpl. with mounting hardware)	1
2	SP1602907	(4-way Expansion Block, LiNX, GLM-CONX4 (single part))	1

(Slow-Stop-Sensor - Modulite Seat)


1609178/00





1	SP1607887	(Slow-Stop-Sensor LiNX, cpl. with bracket for Modulite Seat)	1
2	SP1604482	(Slow-Stop-Sensor LiNX)	1

(Slow-Stop-Sensor - Recaro Seat)


1639179.00



1	SP1607888	(Slow-Stop-Sensor LiNX, cpl. with mounting hardware for Recaro Seat)	1
2	SP1603338	(LiNX Slow-Stop-Sensor, cpl.)	1



1	SP1637339	(Accessory cable LiNX (w/o lifter))	1
2	SP1637563	(Accessory cable LiNX (with lifter))	1




1653460.iso

1 SP1644958

(Stereo-Mono-Splitter Cable)

1

(Programming Access Key)




1641880.iso

1 SP1602886

LINX PROGRAMMERINGS KEY LINX

1







1642451.ISO


1	-	(Accessory cables)	1
2	-	(10-way switch & brackets)	1
3	-	(USB Charger Port)	1
4	-	(Buddy Button)	1
5	-	(Egg Switch ASL300)	1
6	-	(Piko Button)	1
7	-	(Bodypoint Joystick Knobs)	1

(10-way switch & brackets)

1608155.kso



1	-	(10-way switch)	1
2	-	(Bracket, front)	1
3	-	(Bracket, opposite side)	1
4	-	(Bracket, rear)	1




1 SP1602838

(LiNX 10-way switch (DLX-SW100))

1

(Bracket, front)




1600157.160

1 SP1606388

(Bracket 10-way switch - front)


1



1	SP1606388	(Bracket 10-way switch - front)	1
2	SP1603419	(Remote support tube)	1
3	1440437-9005	(Joystick support Standard)	1
4	1417938	(Countersunk bolt DIN7991-M05x20-8.8 zc)	2
5	1422890	(Washer DIN125-Ø5,3 bk)	2
6	F20111	(Hexagon nut DIN985-M05 bk)	2
7	SP1606371	(Bracket, opposite side, cpl.)	1

(Bracket, rear)


160165.no



1 SP1606389

(Bracket 10-way switch - rear)

1







1641896.iso

1 SP1607868

(USB Charger Port)


1

(Bodypoint Joystick Knobs)



1653500.iso

1	PC101A	(Bodypoint Joystick Handle - 3" U-Shaped, ergonomic)	1
2	PC102A	(Bodypoint Joystick Handle - 4" U-Shaped, ergonomic)	1
3	PC107A	(Bodypoint Joystick Handle - Domed Rubber)	1
4	PC110A	(Bodypoint Joystick Handle - Grooved Mushroom)	1




1	-	(ACS Modules)	0
2	-	(ACS Remotes & Brackets)	0
3	-	(ACS Cables & Accessories)	0
4	-	(ACS Components)	0

(ACS Modules)

-	(Powermodule GGB & Bracket)	0
-	(CLAM & Bracket)	0
-	(Relais & bracket)	0
-	(Lightboard GGB)	0

(Powermodule GGB & Bracket)


1581910 iso




1	SP1568490	POWERMODUL DX-GB2 STORM4	1
2	SP1582003	(Bracket Powermodule GGB2)	1

Invacare® Storm4 | Elektronikk
(CLAM & Bracket)


1581911.no



1	P10081	CLAM FOR LYS FUNKSJONER STORM4	1
2	SP1582004	(Bracket CLAM)	1




1	SP1580133	(Relais, bracket & cover, cpl.)	1
2	SP1580197	(Relais cover, cpl.)	1




1	SP1587248	LYSFORDELER STORM4 TT STORM4 TT <i>LED only! [1437]</i>	1
1	SP1528890	LYSFORDELER STORM4 <i>German homologation only! [1367]</i>	1

(ACS Remotes & Brackets)



1	-	(REM 24 SD with brackets) [1601/1602]		1
2	-	(REM A & B with brackets) [1638/1639/1643/1644]		1
3	-	(REM 550 with brackets) [1617/1618]		1
4	-	(REM 34 with brackets) [1611]		1
5	-	(G90A with brackets) [1614/1616]		1
6	-	(Remote G91S & brackets) [1615]		1
7	-	(Centre control REM 41) [1608]		1
7	-	(Tray remote REM 41) [1607]		1
8	-	(Chin Control, small (=> 16/12/31)) [1609]		1
9	-	(Chin Control - Modulite Seat (15/10/01 =>)) [1719]		1
10	-	(Dual Control with brackets) [1291/1292]	→ 30.06.2015	1
11	-	(Compact Dual Control with brackets) [1720]	01.07.2015 →	1
12	-	(Intuitive Dual Control (IDC) with bracket) [1645]		1


(REM 24 SD with brackets)



1563648.iso


1	-	(Remote REM 24SD) [1601/1602]	0
2	-	(Remote Bracket Standard) [1601+1207/1208]	1
3	-	(Remote Bracket Swing away & height adjust.) [1601+1707/1708]	1
4	-	(Remote Bracket Attendant for Standard-/EasyAdapt-Seat) [1602+1207/1208]	1
5	-	(Remote Bracket Attendant for Recaro-Seat) [1602+1207/1208]	1

1563647 no



1	SP1523444	STYREBOKS REM550 STORM4 TDX BORA [1617/1618]	1
2	-	(Remote Bracket Standard) [1617+1207/1208]	1
3	-	(Remote Bracket Swing away & height adjust.) [1617+1707/1708]	1
4	-	(Remote Bracket Attendant for Standard-/EasyAdapt-Seat) [1618+1207/1208]	1
5	-	(Remote Bracket Attendant for Recaro-Seat) [1618+1207/1208]	1

Invacare® Storm4 | Elektronikk
(G90A with brackets)




1563650.iso

1	SP1504720	STYREBOKS DX-REM G90A STORM4 [1614/1616]	1
2	-	(Remote Bracket Standard) [1614+1207/1208]	1
3	-	(Remote Bracket Swing away & height adjust.) [1614+1707/1708]	1
4	-	(Remote Bracket Attendant for Standard-/EasyAdapt-Seat) [1616+1207/1208]	1
5	-	(Remote Bracket Attendant for Recaro-Seat) [1616+1207/1208]	1

(Remote Bracket Attendant for Standard-/EasyAdapt-Seat)


1560853.iso



1	SP1509435	(Remote Support attendant control G90A, LH)	1
2	SP1509439	(Remote Support attendant control G90A, RH)	1
3	SP1422937	FESTE STYREBOKS RYGGRØR STORM4 TDX SP2	1
4	SP1422936	(Bracket for attendant control RH)	1
5	L31001	VINGSKRUV M 6X10 GARANT G24 S	1

(Compact Dual Control with brackets)

1588802 Rev




1	-	(Compact Dual Control)	1
2	-	(Remote Bracket, Compact Dual Control/ACU - short)	1
3	-	(Remote Bracket, Compact Dual Control/ACU - long)	1

(Intuitive Dual Control (IDC) with bracket)

-	-	(IDC with bracket)	1
-	-	(Height adjustable Bracket, IDC)	1

(ACS Cables & Accessories)

-	(Main cable harness)	0
-	(Plug Sheets el. legrest GGB)	0
-	(T-Adapter)	0
-	(T-Adapter Modulite Seat)	0
-	(Programming Kits)	0
-	(Magnetkey for Remotes)	0




1 SP1579550

(Main Cable Harness Storm4 TT)


1

(Plug Sheets el. legrest GGB)

1566930 Iso



1	SP1577086	(Plug Sheet GGB - Modulite Seat)	1
2	SP1581723	(Plug Sheet GGB - EasyAdapt Seat)	1




1 SP1579548

(T-Adapter GGB2, cpl.)

1

(T-Adapter Modulite Seat)

1581913 IRO



1 SP1579547


(T-Adapter Modulite GGB2, cpl.)

1


-	(Batteries & Accessories)	1
-	(Chargers)	1
-	(Chargers)	1

(Batteries & Accessories)

1581915.iso









-	(Important advice! Replace batteries in pairs. Never mix: manufactures, date codes or technologies (Gel with AGM batteries))		1
1	SP1531651 (Battery Strap)		1
2	SP1524214 STØTTEPLATE BATTERI STORM4		1
3	SP1518073 (Battery strap)		1
4	SP1525580 SIKRINGSHOLDER MED SIKRING KPL		1
-	SP1588075 (Cable Tie Set V0 (10 of each T18R/T30R/T50R))		1
5	1433734 SIKRING 80A TIL BATTERI		1
6	SP1524666 BATTERIKABEL ACS2 ACS		1
6	SP1607814 (Battery Cable) LiNX	01.01.2017 →	1
7	P52996 POLSKO DEKSEL		4
8	SPP50016 BATTERI 12V/73,6AH MK		1



164-1060.iso

1	-	(No longer available.)	1
2	P21079	(Power Cord - EU Plug (IEC 60210-C13 / CEE 7/7 Type F 'Schuko'))	1
3	P21080	(Power Cord - UK Plug (IEC 60320-C13 / BS 1363 Type G))	1
4	1443841	(Power Cord - Swiss Plug (IEC 60320-C13 / SEV 1011 Type J))	1
5	SP1553849	(Power Cord - Israel Plug (IEC 60320-C13 / SI 32 Type H))	1
6	SP1606541	(Charger OptiCharge 8A - EU Plug) [1343]	1
6	SP1606543	(Charger OptiCharge 8A - UK Plug) [1343]	1
6	SP1606544	LADER OPTICHARGE 10A OPTICHARGE [1313]	1
6	SP1606545	(Charger OptiCharge 10A - UK-Plug) [1313]	1
7	-	(No longer available.)	1

(Chargers)




1653735.iso

1	-	(No longer available.)	1
2	P21079	(Power Cord - EU Plug (IEC 60210-C13 / CEE 7/7 Type F 'Schuko'))	1
3	SP1606541	(Charger OptiCharge 8A - EU Plug) [1343]	1
3	SP1606544	LADER OPTICHARGE 10A OPTICHARGE [1313]	1
4	5813987	(Charger GI Popular 8A (IP54) - EU Plug) [1388]	1
5	-	(No longer available.)	1
6	SP1653465	(Charger WPC300 10A (IP54) - EU Plug) [1906]	1

-	(Batteries & Accessories)	1
-	(Chargers)	1
-	(Chargers)	1

(Batteries & Accessories)


1581915.iso



-	(Important advice! Replace batteries in pairs. Never mix: manufactures, date codes or technologies (Gel with AGM batteries))	1
1	SP1531651 (Battery Strap)	1
2	SP1524214 STØTTEPLATE BATTERI STORM4	1
3	SP1518073 (Battery strap)	1
4	SP1579554 (Fuse Cable (GGB2))	1
-	SP1588075 (Cable Tie Set V0 (10 of each T18R/T30R/T50R))	1
5	1433734 SIKRING 80A TIL BATTERI	1
6	SP1579555 (Battery Cable ACS2 (GGB2))	1
7	P52996 POLSKO DEKSEL	4
8	SPP50016 BATTERI 12V/73,6AH MK <i>Replaces P50016</i>	1

Invacare® Storm4 | Sikkerhetsutstyr
(Frontlight (German homologation))


1528320K.lin



1	SP1528319	(Frontlight, complete with bracket, LH (Germany))	1
1	SP1528547	BRUK SP1527958 STORM4 <i>not shown</i>	1
2	SP1528320	(Frontlight, complete with bracket, RH (Germany))	1
2	SP1528548	BRUK SP1527957 STORM4 <i>not shown</i>	1
3	SP1530571	HOLDER LYSAN KPL VE STORM4	1
3	SP1530570	HOLDER LYSAN KPL HØ STORM4	1
4	SP1419271	(Light Bulb 24V/21W, yellow (Turn Signal front/rear))	1
5	SPP31108	(Light Bulb 24V/10W (BA15s), white (Front))	1

(Frontlight - Seat Rail mounted (German homologation))


1573323.iss



-	(Bracket design changed on 2015-01-10. For Retrofitting please order SP1556591.)		1
1	SP1573237	(Frontlight, complete with bracket, LH (Germany))	1
2	SP1573252	(Frontlight, complete with bracket, RH (Germany))	1
3	SP1573253	(Frontlight bracket, LH, cpl.)	1
4	SP1573320	(Frontlight bracket, RH, cpl.)	1
5	SP1556591	(Lamp bracket) <i>New design!</i>	10.01.2015 → 2
-	SP1419271	(Light Bulb 24V/21W, yellow (Turn Signal front/rear))	1
-	SPP31108	(Light Bulb 24V/10W (BA15s), white (Front))	1

Invacare® Storm4 | Sikkerhetsutstyr
(Rearlight (German homologation))


1528769X.dwg



!!!	-	(Rear Light Protector Kit)	0
1	SP1530572	BAKLYS KPL STORM4	1
2	SP1530574	HOLDER LAMPE BAK VE STORM4	1
3	SP1530573	HOLDER LAMPE BAK VE STORM4	1
4	SP1419271	(Light Bulb 24V/21W, yellow (Turn Signal front/rear))	1
5	P31191	PÆRE 24V/5W BAKLYS STORM 3	1

Invacare® Storm4 | Sikkerhetsutstyr
(Rear Light Protector Kit)


1585420.000



1 SP1581827

(Rear Light Protector Exchange Kit)


1



-	SP1541713	(Adapter cable LED light, LH) <i>with ULM Seat only!</i>	1
-	SP1541654	ADAPTERKABEL LED-LYS HØ <i>with ULM Seat only!</i>	1
1	SP1527958	LYKT VE KPL LED STORM4	1
2	SP1527957	LYKT HØ KPL LED STORM4	1
3	SP1560611	FESTE LYKT FORAN VE STORM4	1
4	SP1562245	FESTE LYKT FORAN HØ STORM4	1


(LED-Frontlight - Seat Rail mounted)

1573222.no



-	SP1541713	(Adapter cable LED light, LH) <i>FDX only!</i>	1
-	SP1541654	ADAPTERKABEL LED-LYS HØ <i>FDX only!</i>	1
-		(Bracket design changed on 2015-01-10. For Retrofitting please order SP1556591.)	1
1	SP1573223	LYKT KPL VE MED FESTE STORM4 MODULITE	1
2	SP1573224	LYKT KPL MED FESTE HØ STORM4 MODULITE	1
3	SP1573226	FESTE LYS FORAN HØ STORM4 MODULITE	1
4	SP1573227	FESTE LYS FORAN HØ MODULITE	1
5	SP1556591	(Lamp bracket) <i>New design!</i>	2


10.01.2015 →



-	SP1529712	BAKLYS CABLE STORM4	1
1	SP1530287	BLINKLYS BAK STORM4	1
2	SP1530288	BAKLYS STORM4	1

(Attachment Kit - Optimist Seat => 2018-12-31)

1567678.iso




SP1563883

BRAKETT LYKTHOLDER STORM4 OPTIMIST

1

(Reflective Sign)


SP1430410.no



1	SP1430410	(Reflective Sign ECE-R 69, complete)	1
2	1430834	(Rear Marking Sign ECE-R 69)	1
3	Z12008A.P40	(Gripping clip with hexagon hole, RAL 9005 black)	2
4	D50004	(Cylinder Bolt M6x30-8.8 bk - DIN 912)	4
5	F20109	MUTTER 6 KT M6	6
6	D00102	SEKSKANT SKRUE M8X45	1
7	G00107	STOPPSKIVE 8,4X17X1,6	2
8	F00002	SEKSKANTMUTTER M8	1
9	D90187	(Fillister Head Screw M6x20-10.9 bk - DIN 7380)	2
10	G00115	SKIVE 6,4X12,5X1,6 MISTRAL	2
11	1430418-9005	(Mounting Bracket for Reflective Sign - RAL 9005 black)	1
12	SP1580671	(Bracket, Reflective Sign, cpl.)	1

(Reflective Sign - ULM Seat)


1508516.no



1 1430834

(Rear Marking Sign ECE-R 69)


1




1	SP1528952	KANTKLATRER STORM4	→ 12.04.2012	1
1	SP1551945	KANTKLATRER NY SØYLE 3HU STORM4	13.04.2012 →	1
2	SP1530562	(Kerb Climber)		1
3	B11037B	GUMMIFOT T KANTKLÄTTRA		1

(Lap Belts / H-Belt)

1529460X.dwg




1	SP1530575	H-SELE FESTE EASYADAPT STORM4	1
2	SP1530576	H-SELE EASYADAPT STORM4	1
3	SP1529460	(Lap Belt hook & loop, cpl.)	1
4	SP1529461	HOFTEBELTE M HURTIGLÅS STORM4	1



1571836.ISO


1	SP1570381	HOFTEBELTE MODULITE FESTET PÅ RYGG	1
2	SP1571183	HOFTEBELTE C TIL SKINNE MODULITE	1
3	1422845	HOFTEBELTE MED SNAP LÅS STORM 3/ DRAGON	1

(H-Belts - Modulite Seat)




1573208.iso

1	SP1570385	(H-Belt for Backrest Frame)	1
2	SP1570380	(H-Belt for Back Posts)	1



1647489.iso

1	HB205-M46-B2	Hoftebelte Ctr-Pull PB m edium l46 Pads C-Mt	1
2	HB205-L62-B2	Hoftebelte Ctr-Pull PB s tor l62 Pads C-Mt	1
3	HB216-M46-B2	Hoftebelte Dual-Pull SR medium l46 Pads C-Mt	1
4	HB216-L62-B2	Hoftebelte Dual-Pull SR stor l62 Pads C-Mt	1
5	EB225M	Hoftebelte Evoflex Pelvi c Stabilizer medium Latc	1
6	EB225L	Hoftebelte Evoflex Pelvi c Stabilizer stor Latch	1
7	HKA10-1	T-Slot feste Kit M6 14mm til 10mm T-Slot Pr	1




1588510.ISO

1 SP1571183

HOFTEBELTE C TIL SKINNE MODULITE

1




1582371.ISO

1 1551067

HOFTEBELTE, LANG MAX

1

(Rearview Mirror)




1422955X.iso

1	SP1422955	SPEIL KPL KUN SPEIL STORM 3	1
2	F276005.00	HOLDER SPEIL DRAGON STORM4 HUSK SP1527151	1
3	SP1423435	HOLDER SPEIL STORM3	1
4	L31001	VINGSKRUV M 6X10 GARANT G24 S	1
5	SP1422957	SPEIL FESTE BRAKETT STORM4 PÅ SWINGAWAY	1
6	SP1527151	FESTE FOR SPEILFESTE STORM4 TDX BORA	1

-	(Storm4)	0
-	(Storm4 TTP)	0

(Storm4)


-	(Anti-Tipper Wheels => 04.01.11)	0
-	(Anti-Tipper for SSD Motors)	0
-	(Anti-Tippers for SSD Motors, flip up)	0



1	SP1530555	ANTI-TIP-HJUL Ø50 KPL STORM4	1
2	SP1530556	ANTI-TIP-HJUL Ø80 KPL STORM4	1

(Anti-Tipper for SSD Motors)


1533673X.dwg



1	SP1555703	(Anti-Tip-Wheel cpl., Ø50)	1
2	SP1555704	HJUL TIPPESIKRING STORM4	1

(Anti-Tippers for SSD Motors, flip up)

1552985X.dwg




1	SP1552420	(Anti-Tipper, flip up, cpl.)	1
2	SP1550631	(Mini indexing plunger)	1
3	SP1556002	(Anti-Tip-Wheel, Ø80)	1

(Storm4 TTP)

-

(Anti-Tippers for GGB2 Motors)

0



1 SP1582215


(Anti-Tip-Wheel Kit, Ø80mm)

1

(Headrests)

-	(Headrest with cheek support) <i>[1427]</i>	1
-	(Headrest & Bracket) <i>[1350]</i>	1
-	(Headrestholder EasyAdapt Seat)	1
-	(Headrest, sidewise adjustable) <i>[1432]</i>	1
-	(Elan Headrests - Modulite Seat)	1
-	(Headrests - Ultra Low Maxx Seat) <i>[2020]+[2092/2093/2094]</i>	1

(Headrest with cheek support)




1446560X.iso

1	1446560	HODE/SKINN STØTTE KPL TR26 MED SIDEFORFLYTNING	1
1	SP1541658	HODESTØTTE KINN EA STORM4 EASY ADAPT	1
3	SP1444304	(Cushion, cheek support)	1
5	SP1529572	HOLDER HODESTØTTE EASYADAPT	1
6	5716497	HYLSE MESSINGSKRUE FOR LOMME REA	1
7	SP1448661	(stop neck and headrest)	1
8	1445838	SKINNE HODESTØTTE KPL BE LLIS AZAL FOC SPI3	1
9	1485538	FESTE HODESTØTTE STORM4	1

Invacare® Storm4 | Sete utstyr
(Headrest & Bracket)


1422925FFH.dwg



1	SP1422925	HODESTØTTE KPL CONTOUR SORT STORM3/DRAG U/FESTE	1
2	SP1422865	HODESTØTTEPUTE U FESTE STORM4	1
3	SP1422914	HOLDER HODESTØTTE RYGG STORM4 TDX	1

(Headrestholder EasyAdapt Seat)


1529463X.dwg



1	SP1529463	NAKKESTØTTE FESTEBRACKETT STORM4 EASYADAPT	1
2	1434623	LÅSERATT 75MM STORM4	1


(Headrest, sidewise adjustable)

153091X.no



1	SP1541657	NAKKESTØTTE SIDEV JUST BRUK SP1530119	1
2	SP1530567	(Holder, cpl.)	1


(Elan Headrests - Modulite Seat)



1649324.iso

1	SP1603255	(ELAN Cushion, 150mm Standard) <i>[2323]</i>	1
1	SP1603254	(ELAN Cushion, 255mm Standard) <i>[2324]</i>	1
1	SP1603253	(ELAN Cushion, 355mm Standard) <i>[2322]</i>	1
2	SP1603291	(ELAN Cushion, 230mm x 180mm Occipital) <i>[2325]</i>	1
2	SP1603290	(ELAN Cushion, 305mm x 200mm Occipital) <i>[2326]</i>	1
3	SP1626823	(ELAN Cushion, 280mm x 255mm 4-point) <i>[2327]</i>	1
4	SP1603250	(ELAN Headrest Hardware, adjustable, w/o cushion)	1
5	SPTRC2480	(Headrest Clamp Kit, Low Profile Knob)	1
6	1603252	(ELAN Headrest Adapter for Flex3, Backrest Plate) <i>Flex3 Backrest Plate only</i>	1
7	1629412	(ELAN Headrest Adapter for Flex3, Tension adjustable) <i>Flex3 Tension adjustable only</i>	1

(Headrests - Ultra Low Maxx Seat)




1643207.iso

1	SPTRC1929	(Onyx Auto Style Headrest)	→ 10.08.2017	1
2	-	(No longer available - please refer to Replacement Info for which part is to be used instead.) <i>Replacement: SP1603250 & SP1603254</i>	→ 11.02.2018	1
3	-	(No longer available - please refer to Replacement Info for which part is to be used instead.) <i>Replacement: SP1603250 & SP1603291</i>	→ 11.02.2018	1
4	SPTRX2685	(Headrest Shim Kit Posture Back)		1
4	1544188	FESTE ELITE FOR HODESTØTTEFESTE 36-43cm		1
4	1544187	FESTE Elite FOR HODESTØTTE FESTE 46-51CM		1
5	SPTRC2480	(Headrest Clamp Kit, Low Profile Knob)		1
6	1603252	(ELAN Headrest Adapter for Flex3, Backrest Plate) <i>NOT for tension adjustable Flex3 Back!</i>		1
7	SPTRC2308-C02	(Onyx Auto Style Headrest Cushion)	11.08.2017 →	1
8	SP1603255	(ELAN Cushion, 150mm Standard)		1
8	SP1603254	(ELAN Cushion, 255mm Standard)		1
8	SP1603253	(ELAN Cushion, 355mm Standard)		1
9	SP1603291	(ELAN Cushion, 230mm x 180mm Occipital)		1
9	SP1603290	(ELAN Cushion, 305mm x 200mm Occipital)		1
10	SP1626823	(ELAN Cushion, 280mm x 255mm 4-point)		1
11	SP1603249	(Multi-axis Headrest Hardware, horizontally adjustable, w/o cushion)		1
12	SP1603250	(ELAN Headrest Hardware, adjustable, w/o cushion)		1

-	(Lateral Supports - Modulite Seat)	1
-	(Lateral Support - Ultra Low Maxx Seat) <small>[2020]+[2043/2044/2045/2046]</small>	1
-	(Lateral Support - EasyAdapt Seat => 2017-08-10)	1

(Lateral Supports - Modulite Seat)


1627167.no



1	-	(Lateral supports, Flex3 backrest)	1
2	-	(Lateral supports, Flex3 backrest - multi adjustable)	1
3	-	(Lateral supports, Flex3 backrest (tension adjustable) - multi adjustable)	1

(Lateral supports, Flex3 backrest)


1007155.00



-	-	(Design change on 2016-07-01 - we suggest to change older supports on both sides.)	1
1	SP1552678	Sidestøtte forlenget ve	1
2	SP1552679	Sidestøtte forlenget hø	1
3	1543574	SIDESTØTTEFESTE KPL LAGUNA 2, CLEMATIS	1

(Lateral supports, Flex3 backrest - multi adjustable)


160205.00



-	-	(Design change on 2016-07-01 - we suggest to change older supports on both sides.)	1
1	SP1555692	Sidestøtte multi multi r eg forlenget ve	1
2	SP1555691	Sidestøtte multi multi r eg forlenget hø	1
3	1543574	SIDESTØTTEFESTE KPL LAGUNA 2, CLEMATIS	1

(Lateral supports, Flex3 backrest (tension adjustable) - multi adjustable)


1607219.no



1	SP1562920	SIDESTØTTE VELCRORYGG VE TDX SP2 MODULITE	1
2	SP1562919	(Trunk support 2012, ri, multi adj, extended , Flex3 tension adj.)	1
3	SP1554143	SIDESTØTTE FESTE FLEX3 TIL VELCRO RYGG	1

(Lateral Support - Ultra Low Maxx Seat)

1588596.iss




1	SPTRC1806	(Lateral Hardware Assy LH (w/o sleeve))	1
2	SPTRC1805	(Lateral Hardware Assy RH (w/o sleeve))	1
3	SPTRC2123-C03	(Lateral Pad Kit 145x95 mm & sleeve for hardware)	1
4	SPTRC2123-C04	(Lateral Pad Kit 165x110 mm & sleeve for hardware)	1
5	SPTRC2091	(Extension Kit, Lateral Support (Set of 2 pcs) incl. longer sleeves)	1

(Lateral Support - EasyAdapt Seat => 2017-08-10)

1	-	(Lateral Support (Dartex)) [1446]	1
2	-	(Lateral Support, multi adjustable with flexible pad) [1445]	1

(Lateral Support (Dartex))

1570269 iso




SP1561878

SIDESTØTTE FAST DARTEX STORM4

1


(Lateral Support, multi adjustable with flexible pad)

1486065 iso




1	SP1486065	SIDESTØTTE STOFF VE DRAGON VERTIC	1
2	SP1486066	SIDESTØTTE STOFF HØ DRAGON VERTIC	1
3	1451822	SIDESTØTTE VENSTRE U/FES REA BELLIS	1
4	1451821	SIDESTØTTE HØYRE U/FESTE REA BELLIS	1
5	1524903	HOLDER SIDESTØTTE AZALEA	1
6	5253170	SIDESTØTTE FESTE KPL REA	1

1536178X.dwg




1	SP1536178	VESKE STORM4 HUSK SP1537708	1
2	SP1537708	FESTE FOR RYGGSEKK STORM4	1



1645632.ISO

1	SP1644222	(Kit, Bag holder, cpl.)	1
2	SP1644704	(Kit, IDC Adapter)	1
3	SP1644702	(Kit, Headrest Adapter)	1
4	SP1644703	(Kit, Bag hooks)	1

(Luggage Rack)




1546314X.iso

1	SP1546314	FESTE FOR BAGASJEBÆRER STORM4	1
2	SP1551619	BAGASJEBÆRER STORM4 (BE- STILL FESTE, SP1546314)	1
3	SP1556240	(Brackets (pair), Luggage Rack)	1
4	SP1551497	(Luggage Rack)	1
5	F384020	MONT SATS BATTERI BOKS,	2


-	(Cane holder (Easy Adapt Seat & Std-Seat) =>2017-08-10)	1
-	(Cane holder (Recaro Seat))	1
-	(Cane holder (Modulite: backrest frame))	1
-	(Cane holder (Modulite: back posts man recliner))	1
-	(Cane holder (Modulite: back posts pow recliner))	1

(Cane holder (Easy Adapt Seat & Std-Seat) =>2017-08-10)

1567015.iso




1	SP1528444	KRYKKEHOLDER STORM4	1
2	SP1529004	KRYKKEHOLDER HØ EASYADA STORM4	1
3	SP1420493	BÅND TIL KRYKKEHOLDER STORM4	1



1	SP1529535	(Caneholder Recaro cpl)	1
2	SP1529256	(Caneholder Recaro cpl)	1
3	SP1420493	BÅND TIL KRYKKEHOLDER STORM4	1

(Cane holder (Modulite: backrest frame))


1571962.jpg



1	SP1570396	KRYKKEHOLDER RUNDT RØR STORM3	1
2	SP1420493	BÅND TIL KRYKKEHOLDER STORM4	1

(Cane holder (Modulite: back posts man recliner))


1571903.no



1	SP1570443	(Cane holder, manual backrest)	1
2	SP1420493	BÅND TIL KRYKKEHOLDER STORM4	1


(Cane holder (Modulite: back posts pow recliner))

1571837.no



1	SP1420493	BÅND TIL KRYKKEHOLDER STORM4	1
2	SP1570382	KRYKKEHOLDER STORM4 STORM4 MODULITE	1


1567879.dwg



1	SP1530530	HJULSKJERM SVINGHJUL HØ STORM4	1
1	SP1530553	HJULSKJERM SVINGHJUL VE STORM4	1

(Fenders for Drive Wheels => 04.01.11)

	RH		LH			RH		LH			RH		LH	
	1	2	3	4		1	2	3	4		1	2	3	4
Kurze Sitzplatte Short Seatplate					Mittlere Sitzplatte Middle Seatplate					Lange Sitzplatte Long Seatplate				
T=43 B=38	ü	ü	ü	ü	T=48 B=38	ü	ü	ü	ü	T=53 B=38	ü	ü	ü	ü
T=43 B=43	ü	ü	ü	ü	T=48 B=43	ü	ü	ü	ü	T=53 B=43	ü	ü	ü	ü
T=43 B=48	ü	ü	ü	ü	T=48 B=48	ü	ü	ü	ü	T=53 B=48	ü	ü	ü	ü
T=43 B=53	ü	ü	ü	ü	T=48 B=53	ü	ü	ü	ü	T=53 B=53	ü	ü	ü	ü
T=41 B=38	ü	ü	ü	ü	T=46 B=38	ü	ü	ü	ü	T=51 B=38	ü	ü	ü	ü
T=41 B=43	ü	ü	ü	ü	T=46 B=43	ü	ü	ü	ü	T=51 B=43	ü	ü	ü	ü
T=41 B=48	ü	ü	ü	ü	T=46 B=48	ü	ü	ü	ü	T=51 B=48	ü	ü	ü	ü
T=41 B=53	ü	ü	ü	ü	T=46 B=53	ü	ü	ü	ü	T=51 B=53	ü	ü	ü	ü
T=38 B=38	ü	ü	ü	ü	T=43 B=38	ü	ü	ü	ü	T=48 B=38	ü	ü	ü	ü
T=38 B=43	ü	ü	ü	ü	T=43 B=43	ü	ü	ü	ü	T=48 B=43	ü	ü	ü	ü
T=38 B=48	ü	ü	ü	ü	T=43 B=48	ü	ü	ü	ü	T=48 B=48	ü	ü	ü	ü
T=38 B=53	ü	ü	ü	ü	T=43 B=53	ü	ü	ü	ü	T=48 B=53	ü	ü	ü	ü
Recaro Seat	ü	ü	ü	ü	Optimist Seat	ü	ü	ü	ü					




1	SP1530557	SKJERM BAK HØ STORM4 BRUKSP1540339	1
2	SP1530558	SKJERM BAK VE STORM4 BRUK SP1540340	1
3	SP1530559	HJULSKJERM KORT HØ STORM4 brukSP1540339	1
4	SP1530560	HJULSKJERM KORT VE STORM4 brukSP1540340	1

Invacare® Storm4 | Annet Utstyr & Tilbehør
(Fenders for Drive Wheels (SSD-Motors))

1534858X.no


	RH		LH			RH		LH			RH		LH	
	1	2	3	4		1	2	3	4		1	2	3	4
Kurze Sitzplatte Short Seatplate					Mittlere Sitzplatte Middle Seatplate					Lange Sitzplatte Long Seatplate				
T=430 B=380	✓	✓	✗	✗	T=480 B=380	✓	✓	✗	✗	T=530 B=380	✓	✓	✗	✗
T=430 B=430	✓	✓	✗	✗	T=480 B=430	✓	✓	✗	✗	T=530 B=430	✓	✓	✗	✗
T=430 B=480	✓	✓	✗	✗	T=480 B=480	✓	✓	✗	✗	T=530 B=480	✓	✓	✗	✗
T=430 B=530	✓	✓	✗	✗	T=480 B=530	✗	✗	✓	✓	T=530 B=530	✗	✗	✓	✓
T=410 B=380	✓	✓	✗	✗	T=460 B=380	✓	✓	✗	✗	T=510 B=380	✓	✓	✗	✗
T=410 B=430	✓	✓	✗	✗	T=460 B=430	✓	✓	✗	✗	T=510 B=430	✓	✓	✗	✗
T=410 B=480	✓	✓	✗	✗	T=460 B=480	✓	✓	✗	✗	T=510 B=480	✓	✓	✗	✗
T=410 B=530	✓	✓	✗	✗	T=460 B=530	✗	✗	✓	✓	T=510 B=530	✗	✗	✓	✓
T=380 B=380	✓	✓	✗	✗	T=430 B=380	✓	✓	✗	✗	T=480 B=380	✓	✓	✗	✗
T=380 B=430	✓	✓	✗	✗	T=430 B=430	✓	✓	✗	✗	T=480 B=430	✓	✓	✗	✗
T=380 B=480	✓	✓	✗	✗	T=430 B=480	✓	✓	✗	✗	T=480 B=480	✓	✓	✗	✗
T=380 B=530	✓	✓	✗	✗	T=430 B=530	✗	✗	✓	✓	T=480 B=530	✗	✗	✓	✓
Recaro Seat	✓	✓	✗	✗	Optimist Seat	✗	✗	✓	✓					



- | | | | |
|---|-----------|------------------------|---|
| 1 | SP1540339 | SKJERM BAK HØ STORM4 | 1 |
| 2 | SP1540340 | SKJERM BAK VE STORM4 | 1 |
| 3 | SP1540341 | (Mud Guard, short, RH) | 1 |
| 4 | SP1540342 | (Mud Guard, short, LH) | 1 |



(Fenders (GGB Motors))

1583913.iso




1	SP1583917	SKJERM TT2 VE BAK STORM4 [1275/1290/1483]+[1150/1152/1158/1199/2100] / [1275/1290/1483]+[[1180]-[1978/1980/1981]]	1
2	SP1583918	SKJERM TT2 HØ BAK STORM4 [1275/1290/1483]+[1150/1152/1158/1199/2100] / [1275/1290/1483]+[[1180]-[1978/1980/1981]]	1
3	SP1583919	(Fender LH, short) [1275/1290/1483]+[[1180/1171]+[1978/1980/1981]]	1
4	SP1583920	(Fender RH, short) [1275/1290/1483]+[[1180/1171]+[1978/1980/1981]]	1

1588630 iso




1	-	(Trays)	1
2	-	(Trays (Ultra Low Maxx Seat))	1


(Trays)




1	SP1548198	BORD UTEN UTSPARR SB38 STORM4	1
2	SP1548201	BORD UTEN UTSPARRING STORM4	1
3	SP1422954	BORD INKL FESTE KPL STORM UTEN UTSPARRING	1
4	SP1584644	(Swing-Out Tray, 760mm, Seatwidth 480mm to 580mm, no Cutout, cpl. with holder)	1
5	SP1535050	BORD 540MM SB33-43 UTEN UTSKJÆRING STORM4	1
6	SP1535010	BORD 620MM SB38-48 UTEN UTSKJÆRING STORM4	1
7	SP1535011	BORD 660MM SB44-53 UTEN UTSKJÆRING STORM4	1
8	SP1586049	(Tray table, 760mm, no Cutout)	1
9	SP1548138	(Swing-Out Tray, 540mm, Seatwidth 330mm to 430mm, Cutout RH, cpl. with holder)	1
10	SP1548199	BORD STYREBOKS HØ SIDE STORM4	1
11	SP1422953	BORD NEDFELLBART VE STYREBOKS HØ STORM3	1
12	SP1584663	(Swing-Out Tray, 760mm, Seatwidth 480mm to 580mm, Cutout RH, cpl. with holder)	1
13	SP1535051	BORD 540MM SB33-43 STORM4 UTSKJÆRING HØ	1
14	SP1535012	BORD 620 MM SB38-48 UTSK HØ STORM 4 UTEN FESTE	1
15	SP1535049	BORD 660MM SB 44-53 UTSK HØ SIDE UTEN FESTE	1
16	SP1586048	(Tray table, 760mm, Cutout RH)	1
17	SP1548197	(Swing-Out Tray, 540mm, Seatwidth 330mm to 430mm, Cutout LH, cpl. with holder)	1
18	SP1548200	BORD STYREBOKS VE SIDE STORM 4	1
19	SP1422952	BORD NEDFELLBART HØ STYREBOKS VE STORM3	1
20	SP1584662	(Swing-Out Tray, 760mm, Seatwidth 480mm to 580mm, Cutout LH, cpl. with holder)	1
21	SP1535009	BORD 540MM SB33-43 STORM4 UTSKJÆRING VE	1
22	SP1534974	BORD 620 MM SB38-48 UTSK VE STORM 4 UTEN FESTE	1
23	SP1534975	BORD 660MM SB38-48 UT VE STORM4 UTSKJÆRING VE	1
24	SP1586047	(Tray table, 760mm, Cutout LH)	1
25	SPZ177011	HOLDER FOR BORD SVINGBAR KPL STORM3/STORM4	1
26	Z177025	FESTE KIT BOR STORM4	1
27	SP1590090	(Tray Support bumper for Soft Armpad (increased height))	1




1	SPTRC2065	(Flip and fold away Tray LH)	1
2	SPTRC2066	(Flip and fold away Tray RH)	1



1	SP1530563	(Speedometer holder REM A & REM B)	1
2	SP1530564	(Speedometer holder REM 550)	1
3	SP1546263	(Speedometer holder REM24)	1
4	Z183259	(Speedometer)	1
5	1421733	(Reed-sensor cable l=2500mm)	1




1	SP1548705	(Mounting-kit Speedometer sensor)	1
2	SP1540353	(Mounting-kit Speedometer sensor SSD-motors)	1




-	(Not for retrofitting - only spares!)	1
1	SP1530275 (Hour counter)	1
1	SP1608033 (Hour counter (GGB Motors))	1

1567256.no



1	SP1430417	ARMLENSVESKE VE KPL STORM 3	1
2	SP1430855	HOLDER FOR ARMLENSVESKE STORM4 TDX BORA	1
3	SP5314967	VESKE ARMLEN STORM4	1
4	F370032-9005	HOLDER DOBBEL U ARMLEN STORM4	1




1 SP1422868

VERKTØY SETT

1

(T-Nuts for Accessories Modulite & ULM Seat)

157-4067 no



1	SP1556401	MUTTER SKINNE 80MM 2*M8 MODULITE	1
2	SP1572444	MUTTER SKINNE 55MM 2*M8 MODULITE	1
3	SP1555168	MUTTER SKINNE 45MM 2*M8 MODULITE	1
4	SP1557341	MUTTER SKINNE M8 BRED MODULITE	1
5	SP1573158	MUTTER SKINNE VRIDBAR M8 MODULITE	1
6	SP1573721	MUTTER SKINNE VRIDBAR M6 MODULITE	1



Yes, you can.®



Belgium & Luxemburg:

Invacare nv, Autobaan 22, B-8210 Loppem
Tel: (32) (0)50 83 10 10, Fax: (32) (0)50 83 10 11
belgium@invacare.com / www.invacare.be

Denmark:

Invacare A/S, Sdr. Ringvej 37, DK-2605 Brøndby
Tel: (45) (0)36 90 00 00, Fax: (45) (0)36 90 00 01
denmark@invacare.com / www.invacare.dk

Deutschland:

Invacare Aquatec GmbH, Alemannenstraße 10, D-88316 Isny
Tel: (49) (0)75 62 7 00 0, Fax: (49) (0)75 62 7 00 66
info@invacare-aquatec.com / www.invacare-aquatec.de

European Distributor Organisation:

Invacare, Kleiststraße 49, D-32457 Porta Westfalica
Tel: (49) (0)57 31 754 540, Fax: (49) (0)57 31 754 541
edo@invacare.com / www.invacare.eu.com

España:

Invacare SA, c/Areny s/n, Polígon Industrial de Celrà, E-17460 Celrà (Girona)
Tel: (34) (0)972 49 32 00, Fax: (34) (0)972 49 32 20
contactsp@invacare.com / www.invacare.es

France:

Invacare Poirier SAS, Route de St Roch, F-37230 Fondettes
Tel: (33) (0)2 47 62 64 66, Fax: (33) (0)2 47 42 12 24
contactfr@invacare.com / www.invacare.fr

Ireland:

Invacare Ireland Ltd, Unit 5 Seatown Business Campus, Seatown Road, Swords, County Dublin - Ireland
Tel: (353) 1 810 7084, Fax: (353) 1 810 7085
ireland@invacare.com / www.invacare.ie

Italia:

Invacare Mecc San s.r.l., Via dei Pini 62, I-36016 Thiene (VI)
Tel: (39) 0445 38 00 59, Fax: (39) 0445 38 00 34
italia@invacare.com / www.invacare.it

Nederland:

Invacare BV, Celsiusstraat 46, NL-6716 BZ Ede
Tel: (31) (0)318 695 757, Fax: (31) (0)318 695 758
nederland@invacare.com, csede@invacare.com / www.invacare.nl

Norge:

Invacare AS, Grensesvingen 9, Postboks 6230, Etterstad, N-0603 Oslo
Tel: (47) (0)22 57 95 00, Fax: (47) (0)22 57 95 01
norway@invacare.com, island@invacare.com / www.invacare.no

Österreich:

Invacare Austria GmbH, Herzog Odilostrasse 101, A-5310 Mondsee
Tel: (43) 6232 5535 0, Fax: (43) 6232 5535 4
info-austria@invacare.com / www.invacare.at

Portugal:

Invacare Lda, Rua Estrada Velha, 949, P-4465-784 Leça do Balio
Tel: (351) (0)225 1059 46/47, Fax: (351) (0)225 1057 39
portugal@invacare.com / www.invacare.pt

Sverige & Suomi:

Invacare AB, Fagerstagatan 9, S-163 91 Spånga
Tel: (46) (0)8 761 70 90, Fax: (46) (0)8 761 81 08
sweden@invacare.com, finland@invacare.com / www.invacare.se

Switzerland:

Invacare AG, Benkenstrasse 260, CH-4108 Witterswil
Tel: (41) (0)61 487 70 80, Fax: (41) (0)61 487 70 81
switzerland@invacare.com / www.invacare.ch

United Kingdom:

Invacare Limited, Pencoed Technology Park, Pencoed, Bridgend CF35 5AQ
Switchboard Tel: (44) (0)1656 776200, Fax: (44) (0)1656 776201
Customer services Tel: (44) (0)1656 776222, Fax: (44) (0)1656 776220
uk@invacare.com / www.invacare.co.uk