



RESERVEDELSLISTE

Scandi-Puls Luftkompressor

2004



Scandi-Puls Luftkompressor

Position	Varebetegnelse	Varenr.	Antal	Pris
1	Aluminiumslåg	677 03 56	1	75 ,-
2	Jordkabel	677 03 58 *	1	10,-
3	Vibrationsstang inkl. møtrik	677 03 59 *	1 4	20,-
4	Gummimembran med skrue, møtrik og skive	677 03 51 *	1	20,-
5	Lufthammare	677 03 62 *	1	50,-
6	Gummikudde til lufthammare	677 03 52 *	1	10,-
7	Gummipakning til lufthammare	677 03 54	1	10,-
8	Filter til luftkammer	677 03 57	1	10,-
9	Luftbeholder	677 03 55 *	1	60,-

Varenr. mærket med * findes som lagervare på Høng



RESERVEDELSLISTE

Scandi-Puls Pumpe

2004




Scandi-Puls Pumpe

Position	Varebetegnelse	Varenr.	Antal	Pris
1	Lås til synkronmotor	677 03 49 *	1	10,-
2	Sprint till synkronmotor	677 03 48 *	1	5,-
3	Luftfordeler	677 03 44 *	1	80,-
4	Fjeder til synkronmotor	677 03 47 *	4	10,-
5	Metalskive til synkronmotor	677 03 46 *	1	10,-
6	Synkronmotor	677 03 45 *	1	160,-
7	Luftkompressor, komplet	677 03 43 *	1	270,-
8	Bundplade	677 03 39 *	1	100,-
9	Skruer	677 03 38 *	4	10,-
10	Strømafbryder	677 03 42 *	1	30,-
	Reparationskit, madras	677 02 07 *	1	10,-
	PVC slange	677 03 37 *		50,-

Varenr. mærket med * findes som lagervare på Høng



Scandi-Puls Luftkompressor





Scandi-Puls Pump

