



INVACARE® TOP END® PRO™ TENNIS WHEELCHAIR



For new players, adjustability is key. The **Invacare® Top End® Pro™ Tennis Wheelchair** has just that, featuring adjustable front and rear seat heights, center-of-gravity, footrest angle/height, backrest angle/height and anti-tip. Plus, the Top End Pro Tennis Wheelchair is available with either a short or tall frame. The short frame is ideal for children and the tall frame for adults, which makes this chair a great choice for programs that have multiple users. Adjusts in minutes, super easy to use, lightweight and durable... This chair is a winner!

Features

Adjustable front and rear seat heights

Anti-tip/5th wheel: Adjustable height

Fixed back

Adjustable back height

Adjustable tension back upholstery

Adjustable center-of-gravity

Fixed camber 15°, 18° or 20°

High performance Sun Components® wheels:
24", 25", 26"

Anodized pushrims

High pressure clincher tires

Smart fork caster system


72mm casters

Adjustable platform footrest

Powder coat paint

250 lb. weight capacity

Model #: SPTN



1. Front seat height and center-of-gravity adjustments
2. Rear seat height adjustment
3. Quick release lever
4. Single swivel fixed fifth wheel
5. Precision fork caster system
6. Adjustable footrest system



Invacare Corporation
www.invacare.com

USA

One Invacare Way
Elyria, Ohio
44035-4190
(800) 333-6900

Canada

570 Matheson Blvd. E.,
Unit 8
Mississauga, Ontario
L4Z 4G4 Canada
(800) 668-5324

Top End

www.topendwheelchair.com
Sports and Recreation Products
4501 63rd Circle North
Pinellas Park, Florida
33781-5914
(727) 522-8677
(800) 532-8677
Fax: (727) 522-1007

©2014

Invacare Corporation.
All rights reserved.
Trademarks are identified
by the symbols ™ and ®.
All trademarks are owned
by or licensed to Invacare
Corporation unless
otherwise noted. Spinerger
is a registered trademark of
Spinerger Inc. Sun and Equalizer
are registered trademarks of
Sun Components. Kenda is
a registered trademark of
Kenda Rubber Industrial
Co., Ltd. Specifications are
subject to change without
notification.
Form No. 12-507
Rev. 12/15
141157

Options

Clamp-on fixed aluminium side guards

Spinerger® SLX wheels

Sun Components Equalizer® wheels

12 tab pushrims

Kenda® tires

Click straps

O-ring positioning belt

Spoke guards

Frame protectors

Alignment gauge

Telescoping anti-tip bar

Cushion

Folding Back

Fold-in side guards

