


scalamobil / scalacombi

scalamobil S35/S38	244
Components	244
Drive unit	245
E-Box	246
Battery	247
Battery tray	248
Wheel	249
Brakes	250
Handle	251
scalacombi S36/S39	252
scalacombi	252
scalaport X7 für S35	253
scalachair X3	254
Components	254
scalamobil accessories	255
Anti tippers	255
scalacombi eco S34	256
Components	256
Drive unit housing	257
Drive unit housing	258
Drive unit gear box battery	259
Drive unit gear box	260
Drive unit gear box	261
Wheel	262
Brake	263
Seat unit	264
Foot rest	265
Handle	266

scalamobil S35/S38


Components



Order number	Description	Qty	Net in EUR
1456776	Immobiliser key A10/E35/E36/V25 grey	pcs	11,80
1456795	Programming key scalamobil/scalacombi	pcs	1,85
1519809	Charger 24V 1,0A with main cable EU	pcs	114,80
1520182	Mains cable for charger 1.0A EU	pcs	21,60
1520385	Drive unit scalamobil S35	pcs	4.591,45
1520387	Handle complete S35/S38	pcs	704,05
1520389	Battery pack S35/S36/S38/S39	pcs	160,00
1520559	Documentation Kit S35/S38 multilingual	pcs	25,00
1520560	Documentation Kit S35/S38 German	pcs	15,00
1520992	Drive unit scalamobil S38/S39	pcs	5.182,15

scalamobil S35/S38


Drive unit



Order number	Description	Qty	Net in EUR
1456185	Locking ring 14 x 1	pcs	0,15
1456646	Spacer ring 2mm	pcs	0,65
1456648	Spacer ring 10mm	pcs	0,65
1456789	Distance socket scalamobil lower mounting	pcs	1,50
1456854	Set collar galvanized for scalamobil	pcs	10,70
1456873	Screws and bolts for mounting of upper bolts scalamobil (without bolts)	pcs	1,50
1457767	Lower bolt scalamobil (Dummy Code)	pcs	22,35
1457770	Upper bolts scalamobil (Dummy Code)	pcs	32,65
1461655	Side shell cover scalamobil S35	pcs	1,55
1463454	Device feet scalamobil S35	pcs	3,15
1463458	Rubber protector housing scalamobil	pcs	4,10
1489247	Adjusting ring set for scalamobil/viamobil	pcs	16,95

scalamobil S35/S38


E-Box



Order number	Description	Qty	Net in EUR
1456876	Cover for electronic box scalamobil S35	pcs	7,40
1457281	Circuit board scalamobil IQ-Feature with bearing sensor	pcs	379,00
1459071	Cylinder head screw DIN 912 - M 4 x 12 - 8.8	pcs	0,40
1461741	Sticker „Operating cover“ scalamobil S35	pcs	1,30
1519841	Label warning advice scalamobil	pcs	1,00
1520043	Ejot Easyboss 50 x 3,5 x 8	pcs	1,35
1520045	Ejot Delta PT Screw 50 x 18-AZK	pcs	0,30
1520353	Housing E-Box S35/S36	pcs	16,35
1520952	Circuit board scalamobil S38 IQ-Feature with bearing sensor	pcs	383,00

scalamobil S35/S38


Battery



Order number	Description	Qty	Net in EUR
1456389	Fuse 15 A	pcs	1,55
1456799	Batteries scalamobil complete	pcs	112,00
1456875	Locking mechanism battery pack scalamobil	pcs	3,30
1461707	Releasing device battery scalamobil	pcs	3,30
1465938	Fuse PICO II 2.5 A Littelfuse Serie 473	pcs	2,90
1466248	Pressure spring VD-052G	pcs	0,35
1519643	Fuse 4A for circuit board of battery S35	pcs	1,55
1520389	Battery pack S35/S36/S38/S39	pcs	160,00
1592407	#NV	pcs	33,50

scalamobil S35/S38


Battery tray



Order number	Description	Qty	Net in EUR
1456877	Carrying handle scalamobil	pcs	6,40
1463533	Screw Ejot PT K 40 x 16 WN 1451	pcs	0,15
1520342	Sticker „scalamobil“ motor unit S35	pcs	1,30
1520343	Sticker „scalacombi“ motor unit S36	pcs	1,30
1520348	Adapter/Cover for carrying handle S35/S36	pcs	14,75
1520349	Chute for battery pack S35/S36/S38/S39	pcs	20,90

scalamobil S35/S38


Wheel



Order number	Description	Qty	Net in EUR
1456676	Tyre scalamobil S28/S38	pcs	69,10
1456792	Tyre for scalamobil S35/S34	pcs	51,00
1456806	Magnet holder scalamobil S35	pcs	58,15
1463662	Countersunk head screw DIN 7991 - M 4 x 14 - 8.8	pcs	0,40
1463684	Screw for rim - M 4 x 18	pcs	0,20
1463690	Screw M 5 x 16	pcs	0,15
1464450	Outer wheel scalamobil S35 complete	pcs	129,55
1464451	Inner wheel scalamobil S35 left complete	pcs	129,55
1464452	Inner wheel scalamobil S35 right complete	pcs	129,55
1490058	Hub cover S34/S35 beamless	pcs	22,90
1490524	Selftapping countersink screw	pcs	0,10
1490525	Selftapping countersink screw	pcs	0,15
1520957	Rim scalamobil S38	pcs	96,00
1520958	Rim outer wheel scalamobil S38	pcs	97,45
1592719	Magnet holder scalamobil S38/S39	pcs	73,90

scalamobil S35/S38

Brakes





Order number	Description	Qty	Net in EUR
1461721	Washer	pcs	1,15
1465336	Ball bearing compensation disc EPL13	pcs	0,30
1519727	Countersunk head screw DIN 7991 M 5 x 16 - 8.8	pcs	0,25
1520253	Eجت Delta PT-Screw WN 5454-45 x 12-A2K	pcs	0,25
1520391	Brake wheel S35/S36	pcs	10,70
1520394	Brake lining S35/S36	pcs	11,70
1520474	Tension Spring brake scalamobil S35	pcs	2,85
1520769	Brake scalamobil S35/S36 right complete	pcs	168,95
1520770	brake scalamobil S35/S36 left complete	pcs	168,95
1520999	Brake scalamobil S38 left complete	pcs	165,00
1521000	Brake scalamobil S38 right complete	pcs	165,00

scalamobil S35/S38


Handle

1520387
1520388



Order number	Description	Qty	Net in EUR
1456217	Countersunk head screw DIN 7991 - M 6 x 45 - 8.8	pcs	0,50
1456567	Screw Ejot PT K 30 x 8 WN 1451	pcs	0,15
1463679	ISO 2338 - A - 6 M 6 x 30 - St	pcs	0,40
1490801	Alber Logo 17 x 17mm	pcs	1,25
1520345	Adapter for guidance of lever S35/S36	pcs	10,70
1520354	SoftTouch cushion S35/S36	pcs	27,60
1520359	Cover Bridge S35/S36	pcs	20,90
1520362	Sticker Display S35/S36	pcs	4,60
1520363	Silicone switch on/off S35/S36	pcs	3,55
1520387	Handle complete S35/S38	pcs	704,05
1520388	Handle complete S36/S39	pcs	704,05
1520443	Sticker „scalamobil“ for handle S35	pcs	2,90
1520444	Cover	pcs	5,40
1520469	Distance bush	pcs	1,50
1520470	Clamp lever for lateral adjustment S35/S36	pcs	19,70
1520495	Circuit board Display S35/S36	pcs	50,75
1520496	Handle S35 left hand side with selector switch for drive direction up/down	pcs	174,65
1520498	Handle S35 right hand side with drive switch	pcs	174,65
1520499	Hand wheel S35 connection motor unit	pcs	17,50
1520501	Nut	pcs	2,15
1520502	Guidance for lever of lateral adjustment S35	pcs	6,05
1520509	Hand wheel for heightadjustment S35 (kit)	pcs	18,50


scalacombi S36/S39



Order number	Description	Qty	Net in EUR
1456453	Steering wheels scalacombi/scalaport	pcs	29,40
1520346	Clamp for backrest	pcs	1,50
1520386	Drive unit scalacombi S36	pcs	5.007,85
1520455	Headrest S35/S36	pcs	165,00
1520472	Seat plate S36 (padded)	pcs	163,35
1520499	Hand wheel S35 connection motor unit	pcs	17,50
1520503	Arm rest S36	pcs	25,30
1520505	Backrest for scalacombi S36	pcs	174,65
1520506	Hand wheel seat unit scalacombi S36	pcs	13,00
1520508	Footrest S36	pcs	59,50
1520538	Bottom part, assembled	pcs	138,75
1564990	Chest shoulder belt scalacombi S36/S39, consisting of 2 automatic belt retractors and height adjustable buckles on handle	pcs	320,00
1565509	Lap Belt (Automatic) for scalacombi S36	pcs	140,00
1565528	Footrest scalacombi S36 complete (replacement)	pcs	242,00
1566469	Seat Unit scalacombi S36	pcs	368,00

scalaport X7 für S35


scalaport X7 für S35



Order number	Description	Qty	Net in EUR
1456185	Locking ring 14 x 1	pcs	0,15
1456404	Pressure spring for adjustment knob	pcs	1,15
1456453	Steering wheels scalacombi/scalaport	pcs	29,40
1458895	Locking ring 25 x 1,2	pcs	0,15
1463436	Cylinder head screw M 10 x 25 - 8.8	pcs	0,80
1490076	Clamping mechanism wheelchair frame right side	pcs	79,50
1490077	Clamping mechanism wheelchair frame left side	pcs	79,50
1490083	Brake pad X7	pcs	2,95
1490093	Rubber pad X7 small	pcs	2,95
1490098	Rubber pad X7 large	pcs	2,95
1490161	Adaptation Bracket X7 right for S30 incl. snatch	pcs	41,15
1490162	Adaptation Bracket X7 left for S30	pcs	34,30
1490184	Rod X7	pcs	15,00
1490185	Spacer brass X7	pcs	2,35
1490186	Spacer X7	pcs	3,95
1519489	Magnet X7	pcs	7,40
1520507	Clamping mechanism X7 for S35	pcs	336,85

scalachair X3

Components




Order number	Description	Qty	Net in EUR
1456191	Cylinder head screw DIN 912 - M 6 x 20	pcs	0,15
1456251	Armrest for scalachair X3	pcs	57,50
1456355	Cushion seat and back for scalachair X3	pcs	140,90
1456371	Ball bearing 28 x 12 x 8mm	pcs	6,30
1456406	Quickpin quick release bolt	pcs	21,85
1456412	Clamping lever for push handles for scalachair X3	pcs	24,10
1456452	Wheel for scalachair X3	pcs	24,10
1456850	Pair of push handles for scalachair X3 left and right	pcs	110,00
1456862	Armrest for scalachair X3 left and right	pcs	61,25
1456864	Holding device for handles upper part scalachair X3	pcs	53,45
1456866	Release lever for scalachair X3	pcs	210,80
1456867	Footrest for scalachair X3 with anti-slip cover pad	pcs	35,90
1456868	Steering linkage scalachair X3	pcs	55,40
1456870	Holding device for handles lower part scalachair X3	pcs	50,10
1458178	Seat for scalachair X3	pcs	73,85
1458179	Backrest for scalachair X3 grey without cushion RAL 7040	pcs	73,85
1463386	Footrest fixation for scalachair X3	pcs	276,60

scalamobil accessories

Anti tippers


1519747
1519748



Order number	Description	Qty	Net in EUR
1456183	Washer DIN 125 - A 6,4	pcs	0,05
1456233	Cylinder head screw M 6 x 20	pcs	0,25
1456406	Quickpin quick release bolt	pcs	21,85
1456914	Fixation set 2 x 12mm for anti-tippers	pcs	46,45
1456915	Fixation set 2 x 15mm for anti-tippers	pcs	46,45
1456924	Toothed bar for anti-tippers e-fix	pcs	14,40
1459944	Washer DIN 988 - 10 x 16 x 1	pcs	0,05
1467845	Support angle right for anti-tippers	pcs	19,00
1467846	Support angle left for anti-tippers	pcs	19,00
1489365	Wheel for anti-tippers complete	pcs	2,40
1489366	Slide rod anti-tippers mounted	pcs	20,60
1519603	Securing disk DIN 6799-Gr.8 vz	pcs	0,05
1519747	Anti-tippers without jack-up function 12mm (pair)	pcs	262,00
1519748	Anti-tippers without jack-up function 15mm (pair)	pcs	262,00


scalacombi eco S34

Components



Order number	Description	Qty	Net in EUR
1456776	Immobiliser key A10/E35/E36/V25 grey	pcs	11,80
1519564	Footrest S34	pcs	257,40
1519596	Arm rest left S34	pcs	42,85
1519597	Arm rest right S34	pcs	42,85
1519706	Handlebar S34 complete	pcs	469,00
1519707	Drive unit scalacombi S34	pcs	3.105,00
1519719	Handle bar S34 without backrest	pcs	325,80
1519797	Seat unit S34	pcs	349,00
1519799	Documentation Kit S34 multilingual	pcs	25,00
1519809	Charger 24V 1,0A with main cable EU	pcs	114,80


scalacombi eco S34



Order number	Description	Qty	Net in EUR
1463516	Screw Eject PT KA 30 x 5 WN	pcs	0,15
1463520	Screw Eject PT K 30 x 14 WN	pcs	0,15
1467418	Membrane e-fix E25/E26	pcs	6,00
1491037	Front housing motor unit S34	pcs	97,20
1519591	LED of immobiliser	pcs	20,00
1519592	Main circuit board S34	pcs	291,60
1519633	Charging socket with cable S34	pcs	13,20
1519635	Connector S34	pcs	13,80
1519636	Battery cable (negative pole)	pcs	4,95
1519637	Battery cable (positive pole)	pcs	14,30
1519693	On/Off switch S34	pcs	25,70
1519756	Holder for charging socket	pcs	1,25
1520043	Ejot Easyboss 50 x 3,5 x 8	pcs	1,35
1520044	Ejot Easyboss 50 x 2,5 x 11	pcs	2,00

scalacombi eco S34


Drive unit housing



Order number	Description	Qty	Net in EUR
1466990	Screw Remform 4 x 30 WN	pcs	0,15
1491041	Rubber protector for housing S34	pcs	8,40
1519580	Label „scalacombi eco“	pcs	9,85
1519634	Connecting cable circuit board S34	pcs	28,60
1519720	Light conductor status bar display S34	pcs	9,85
1519725	Housing motor unit S34 with light conductor label and rubber protectors	pcs	66,30

scalacombi eco S34


Drive unit gear box battery



Order number	Description	Qty	Net in EUR
1456613	Foam rubber small	pcs	0,20
1456799	Batteries scalamobil complete	pcs	112,00
1461226	Rubber pad for battery S34	pcs	1,50
1489287	Velcro strap 350mm black	pcs	4,50
1489290	Velcro strap 500mm black	pcs	4,50

scalacombi eco S34


Drive unit gear box



Order number	Description	Qty	Net in EUR
1456184	Washer DIN 125	pcs	0,05
1463393	Securing disk DIN 6799	pcs	0,55
1463653	Cylinder head screw DIN 912 - M 8 x 25	pcs	0,95
1519593	Motor unit S34	pcs	668,95
1519605	Frame S34	pcs	451,75
1519646	Punch climbing mechanism S34	pcs	17,25
1519717	Rubber pad for punch S34	pcs	9,85
1519727	Countersunk head screw DIN 7991 M 5 x 16 - 8.8	pcs	0,25
1520777	Wheel S34	pcs	2,80

scalacombi eco S34


Drive unit gear box



Order number	Description	Qty	Net in EUR
1456184	Washer DIN 125	pcs	0,05
1459412	Cylinder head screw DIN 912 - M 8 x 35	pcs	0,35
1461699	Toothed wheel scalamobil S35/S34	pcs	11,70
1463686	Torx screw DIN 7991 - M 5 x 12	pcs	0,20
1467093	Countersunk head screw M 4 x 10	pcs	0,60
1491032	Excenter S34	pcs	75,45
1519568	Chain S34 (30 elements)	pcs	17,25
1519574	Rubber pad for bracket of motor unit S34	pcs	3,95
1519575	Bush for bracket of motor unit S34	pcs	2,95
1519576	Blank sheet for bracket of motor unit S34	pcs	6,90
1519577	Cover for Excenter S34	pcs	15,35
1519588	Brake shaft S34	pcs	20,60
1519718	Bearing	pcs	14,30

scalacombi eco S34


Wheel



Order number	Description	Qty	Net in EUR
1456792	Tyre for scalamobil S35/S34	pcs	51,00
1463683	Screw hubcap - M 3 x 8	pcs	0,15
1463684	Screw for rim - M 4 x 18	pcs	0,20
1490058	Hub cover S34/S35 beamless	pcs	22,90
1490209	Rim outer wheel scalamobil S35/S34	pcs	94,85
1520151	Bearing 61806-2RS	pcs	5,90

scalacombi eco S34


Brake



Order number	Description	Qty	Net in EUR
1461721	Washer	pcs	1,15
1465336	Ball bearing compensation disc EPL13	pcs	0,30
1519727	Countersunk head screw DIN 7991 M 5 x 16 - 8.8	pcs	0,25
1520253	Ejot Delta PT-Screw WN 5454-45 x 12-A2K	pcs	0,25
1520391	Brake wheel S35/S36	pcs	10,70
1520394	Brake lining S35/S36	pcs	11,70
1520474	Tension Spring brake scalamobil S35	pcs	2,85
1520769	Brake scalamobil S35/S36 right complete	pcs	168,95
1520770	brake scalamobil S35/S36 left complete	pcs	168,95

scalacombi eco S34


Seat unit



Order number	Description	Qty	Net in EUR
1456181	Securing disk S - 3 - FSt	pcs	0,05
1456183	Washer DIN 125 - A 6,4	pcs	0,05
1456191	Cylinder head screw DIN 912 - M 6 x 20	pcs	0,15
1456807	Seat plate with cushion scalacombi	pcs	136,95
1519586	Spacer for backrest S34	pcs	1,00
1519587	Clamp for backrest S34	pcs	2,30
1519711	Backrest S34	pcs	137,20


scalacombi eco S34

Foot rest



Order number	Description	Qty	Net in EUR
1456184	Washer DIN 125	pcs	0,05
1456225	Washer DIN 125 - A10,5	pcs	0,15
1460350	Cylinder head screw DIN 7984 - M 8 x 16	pcs	0,55
1519594	Castor wheel S34	pcs	24,60
1521167	Adjustment knob scalacombi eco S34	pcs	12,50

scalacombi eco S34



Order number	Description	Qty	Net in EUR
1456416	Adjustment knob scalamobil scalamobil	pcs	11,75
1459302	Cylinder head screw DIN 912 - M 6 x 25	pcs	0,15
1463520	Screw Eject PT K 30 x 14 WN	pcs	0,15
1489535	Plug for handle S34	pcs	0,55
1489778	Sticker „Alber“ 28 x 28mm	pcs	1,05
1491035	Clamp for handle S34	pcs	17,25
1491040	Profile cover for stanchion handle S34	pcs	4,95
1519573	Reception tube for handle S34	pcs	85,80
1519578	Stanchion for handle S34	pcs	97,20
1519580	Label „scalacombi eco“	pcs	9,85
1519638	Spiral cable with toggle switch S34	pcs	62,95
1519669	Handle S34 with soft cover	pcs	97,20
1519719	Handle bar S34 without backrest	pcs	325,80
1519937	Cover S34	pcs	27,40
1520499	Hand wheel S35 connection motor unit	pcs	17,50