



INVACARE® TOP END®  
**PRO™**  
BASKETBALL WHEELCHAIR




The **Invacare® Top End® Pro™ Basketball Wheelchair** delivers custom-built performance with adjustments that take only minutes. Perfect for competitors still dialing in their game.

## Features

Custom aluminum frame  
Offensive 4 bend wing  
Adjustable front and rear seat heights  
Anti-tip/fifth wheel: Adjustable height  
Fixed back  
Adjustable back height  
Adjustable tension back upholstery  
Adjustable Center-of-Gravity  
Rear seat height adjustment range is 6" and the frame is available in a tall or short variety to accommodate every aspiring athlete

Fixed camber: 15°, 18° or 20°  
High performance Sun® wheels  
Anodized pushrims  
High pressure clincher tires  
Smart fork caster system  
72mm casters  
Adjustable platform footrest  
Powder coat paint  
250 lb. weight capacity  
Quick release lever allows the player to change his or her rear seat height in just minutes

## Model #: SPBB



1. Front seat height and center-of-gravity adjustments
2. Rear seat height adjustment
3. Quick release lever
4. Single swivel fixed fifth wheel
5. Precision fork caster system
6. Adjustable footrest system



**Invacare Corporation**  
www.invacare.com

### USA

One Invacare Way  
Elyria, Ohio  
44035-4190  
(800) 333-6900

### Canada

570 Matheson Blvd. E.,  
Unit 8  
Mississauga, Ontario  
L4Z 4G4 Canada  
(800) 668-5324

### Top End

www.topendwheelchair.com  
Sports and Recreation Products  
4501 63rd Circle North  
Pinellas Park, Florida  
33781-5914  
(727) 522-8677  
(800) 532-8677  
Fax: (727) 522-1007

©2014

Invacare Corporation.  
All rights reserved.  
Trademarks are identified by the symbols ™ and ®. All trademarks are owned by or licensed to Invacare Corporation unless otherwise noted. Spinerger is a registered trademark of Spinerger Inc. Sun and Equalizer are registered trademarks of Sun Components. Kenda is a registered trademark of Kenda Rubber Industrial Co., Ltd. Specifications are subject to change without notification.  
Form No. 04-188  
Rev. 12/14  
141156

## Options

Clamp-on fixed aluminum side guards  
Spinerger® SLX wheels  
Sun® Equalizer™ wheels  
12 tab pushrims  
Kenda® tires  
O-Ring positioning belt  
Close mount/cut pushrim tab

Click straps  
Spoke guards  
Frame protectors  
Alignment gauge  
Double swivel anti-tip bar  
Cushion: soft or hard  
Folding back

