

KÜSCHALL K-SERIES attract



Local Homologation number, ISO, MDA, TUV logos

The product includes all standard features

Maximum user weight: 120kg

We reserve the right to deliver a + or - tolerance of +/- 10mm on all measurements

! We reserve the right to make technical changes without notice



Manual Wheelchair
Local Prescription Form
January 2018

- Please make your selection with a cross in the box
- Standard feature



CODE	DESCRIPTION
<input checked="" type="checkbox"/> DDZ0068	KÜSCHALL K-SERIES attract
SEATING CHOICE	
FRAME	
Material / Angle / Type	
<input checked="" type="checkbox"/> DDZ0293	75° to the ground Aluminium - tapered front frame
<input type="checkbox"/> DDZ0294	90° to the ground Aluminium - tapered front frame
SEAT WIDTH (SB)	
<input type="checkbox"/> DDZ0126	SB 360 mm
<input type="checkbox"/> DDZ0128	SB 380 mm
<input type="checkbox"/> DDZ0130	SB 400 mm
<input type="checkbox"/> DDZ0132	SB 420 mm
<input type="checkbox"/> DDZ0134	SB 440 mm
<input type="checkbox"/> DDZ0136	SB 460 mm
<input type="checkbox"/> DDZ0138	SB 480 mm
SEAT DEPTH (ST)	
<input type="checkbox"/> DDZ0414	ST 375 mm
<input type="checkbox"/> DDZ0417	ST 400 mm
<input type="checkbox"/> DDZ0420	ST 425 mm
<input type="checkbox"/> DDZ0423	ST 450 mm
<input type="checkbox"/> DDZ0426	ST 475 mm
SEAT POSITION (SP)	
<input type="checkbox"/> DDZ0482	SP 125 mm
<input checked="" type="checkbox"/> DDZ0463	SP 110 mm
<input type="checkbox"/> DDZ0464	SP 95 mm
<input type="checkbox"/> DDZ0465	SP 80 mm
BACK HEIGHT (RH)	
<input type="checkbox"/> DDZ0585	Back height depends on the height of the positioning backrest
<input type="checkbox"/> DDZ0503	RH 315 mm
<input type="checkbox"/> DDZ0506	RH 330 mm
<input type="checkbox"/> DDZ0509	RH 345 mm
<input type="checkbox"/> DDZ0512	RH 360 mm
<input checked="" type="checkbox"/> DDZ0514	RH 375 mm
<input type="checkbox"/> DDZ0516	RH 390 mm
<input type="checkbox"/> DDZ0518	RH 405 mm
<input type="checkbox"/> DDZ0520	RH 420 mm
<input type="checkbox"/> DDZ0522	RH 435 mm
<input type="checkbox"/> DDZ0524	RH 450 mm
<input type="checkbox"/> DDZ0526	RH 465 mm
<input type="checkbox"/> DDZ0528	RH 480 mm
FRONT SEAT-TO-FLOOR (SHv)	
<input type="checkbox"/> DDZ0850	SHv 480 mm
<input type="checkbox"/> DDZ0851	SHv 490 mm
<input checked="" type="checkbox"/> DDZ0852	SHv 500 mm
<input type="checkbox"/> DDZ0853	SHv 510 mm
REAR SEAT-TO-FLOOR (SHh)	
<input type="checkbox"/> DDZ0945	SHh 400 mm
<input type="checkbox"/> DDZ0946	SHh 410 mm
<input type="checkbox"/> DDZ0947	SHh 420 mm
<input checked="" type="checkbox"/> DDZ0948	SHh 430 mm
<input type="checkbox"/> DDZ0949	SHh 440 mm
<input type="checkbox"/> DDZ0950	SHh 450 mm
<input type="checkbox"/> DDZ0951	SHh 460 mm
<input type="checkbox"/> DDZ0952	SHh 470 mm
<input type="checkbox"/> DDZ0953	SHh 480 mm
KNEE-TO-HEEL LENGTH (UL)	
75° frame: UL ≤ SHv -30 mm; 90° frame: UL ≤ SHv -40 mm	
<input type="checkbox"/> DDZ1670	UL 290 mm
<input type="checkbox"/> DDZ1671	UL 300 mm
<input type="checkbox"/> DDZ1672	UL 310 mm
<input type="checkbox"/> DDZ1673	UL 320 mm
<input type="checkbox"/> DDZ1674	UL 330 mm
<input type="checkbox"/> DDZ1675	UL 340 mm
<input type="checkbox"/> DDZ1676	UL 350 mm
<input type="checkbox"/> DDZ1677	UL 360 mm
<input type="checkbox"/> DDZ1678	UL 370 mm
<input type="checkbox"/> DDZ1679	UL 380 mm
<input type="checkbox"/> DDZ1680	UL 390 mm
<input type="checkbox"/> DDZ1681	UL 400 mm
<input type="checkbox"/> DDZ1682	UL 410 mm
<input type="checkbox"/> DDZ1683	UL 420 mm
<input checked="" type="checkbox"/> DDZ1684	UL 430 mm
<input type="checkbox"/> DDZ1685	UL 440 mm
<input type="checkbox"/> DDZ1686	UL 450 mm
<input type="checkbox"/> DDZ1687	UL 460 mm
<input type="checkbox"/> DDZ1688	UL 470 mm
<input type="checkbox"/> DDZ1689	UL 480 mm
BACKREST	
BACKREST TYPE	
<input checked="" type="checkbox"/> DDZ0773	Aluminium backrest, with straight backbrace
<input type="checkbox"/> DDZ0772	Backrest with lower rounded bar for fixed backrest mounting
Only without backrest upholstery (0769) and for positioning backrests such as Matrx Only without push handles or with push handles height adjustable integrated	
BACK UPHOLSTERY	
<input checked="" type="checkbox"/> DDZ0725	Tension velcro adjustable back, black cover
<input type="checkbox"/> DDZ0735	Aero back, tension adjustable
<input type="checkbox"/> DDZ0769	Without backrest upholstery
Select if you want to order Positioning backrests like Matrx	
MATRIX BACKRESTS	
Only possible with "without push handles" (DDZ0750) and "push handles integrated" (DDZ0733) when mounting Matrx backrests please make sure that the upper part of the Matrx backrest is positioned not higher than the max RH value indicated on the EPF i.e. for the k-Series attract 480mm	
Height refers to actual height of the Matrx backrest ; Requires additional lead time, always provided with the option "without backrest upholstery"	
<input type="checkbox"/> DDZ0764	Matrx MX2, height 150 mm with quick release hardware
<input type="checkbox"/> DDZ0765	Matrx MX2, height 230 mm with quick release hardware
<input type="checkbox"/> DDZ0766	Matrx MX2, height 305 mm with quick release hardware
<input type="checkbox"/> DDZ0767	Matrx MX2, height 405 mm with quick release hardware
NC for SB 480mm (DDZ0128) NC for SB 480mm (DDZ0128) NC for SB 480mm (DDZ0128) NC for SB 480mm (DDZ0128)	

Frame	SHv	UL / Footrest mounting
	High	Std
75° tapered	48	29-39 40-45
	49	29-40 41-46
	50	29-41 42-47
	51	29-42 43-48
90° tapered	48	29-37 38-44
	49	29-38 39-45
	50	29-40 41-40
	51	29-41 42-47

<input type="checkbox"/>	DDZ0756	Elite Deep, height 310 mm with easy-set mounting hardware	
<input type="checkbox"/>	DDZ0757	Elite Deep, height 410 mm with easy-set mounting hardware	
<input type="checkbox"/>	DDZ0758	Elite Standard, height 260 mm with easy-set mounting hardware	
<input type="checkbox"/>	DDZ0759	Elite Standard, height 310 mm with easy-set mounting hardware	
<input type="checkbox"/>	DDZ0760	Elite Standard, height 360 mm with easy-set mounting hardware	NC for SB 480mm (DDZ0138)
<input type="checkbox"/>	DDZ0761	Elite Standard, height 410 mm with easy-set mounting hardware	
PUSH HANDLES			
<input type="checkbox"/>	DDZ0750	Without push handles	
<input checked="" type="checkbox"/>	DDZ0751	Push handles fixed, short	Possible with straight bar (DDZ0773) only without Matrix backrest
<input type="checkbox"/>	DDZ0702	Push handles fixed, long	Possible with straight bar (DDZ0773) only without Matrix backrest
<input type="checkbox"/>	DDZ0703	Push handles height adjustable, rear set	NC with Matrix backrests
<input type="checkbox"/>	DDZ0704	Foldable push handles	Possible with straight bar (DDZ0773) only without Matrix backrest
<input type="checkbox"/>	DDZ0733	Push handles height adjustable, integrated	Only possible with Matrix backrests
BACKREST ANGLE (RW)			
<input type="checkbox"/>	DDZ0718	78° angle from the seat	
<input checked="" type="checkbox"/>	DDZ0709	82° angle from the seat	
<input type="checkbox"/>	DDZ0719	86° angle from the seat	
<input type="checkbox"/>	DDZ0711	90° angle from the seat	
<input type="checkbox"/>	DDZ0723	94° angle from the seat	
SEATING			
SEAT UPHOLSTERY			
<input checked="" type="checkbox"/>	DDZ0610	Standard	
<input type="checkbox"/>	DDZ0611	Tension adjustable with one central velcro band (soft)	
SEAT BAG			
<input type="checkbox"/>	DDZ0606	Seat bag, click-in and removable	
SEAT CUSHION			
<input type="checkbox"/>	DDZ0600	Standard seat cushion, 60 mm, double foam with visco top	
<input type="checkbox"/>	DDZ0662	Matrix Libra cushion, black smooth cover	NC in ST375mm (DDZ0414); SB 360mm NC in ST450/475mm
<input type="checkbox"/>	DDZ0663	Matrix Vi cushion, black smooth cover	NC in ST375mm (DDZ0414); SB 360mm NC in ST450/475mm
FOOTREST			
FOOTREST TYPE (TUBE / PLATE)			
<input checked="" type="checkbox"/>	DDZ1508	Aluminium footrest tube, rigid	Rigid tube with foot-bow
<input type="checkbox"/>	DDZ1570	Carbon footplate, angle adjustable	
FOOTREST COVER			
<input type="checkbox"/>	DDZ1513	Footrest cover	Only with standard footrest tube (DDZ1510)
SIDE GUARDS			
CLOTHES GUARDS / MUD GUARDS			
<input checked="" type="checkbox"/>	DDZ1424	Composite clothes guards	
<input type="checkbox"/>	DDZ1425	Composite mudguards	
<input type="checkbox"/>	DDZ1401	Carbon clothes guards	
<input type="checkbox"/>	DDZ1408	Carbon mudguards	
ARMREST			
<input type="checkbox"/>	DDZ1416	Armrest, height adjustable, short armpad & composite clothes guards, non-locking	NC (DDZ1424/1425 and DDZ1401/1408) NC with SH<420 mm
<input type="checkbox"/>	DDZ1421	Armrest, height adjustable, long armpad & composite clothes guards, non-locking	NC with SH<420 mm
FRAME			
FRAME COLOURS, STANDARD			
<input checked="" type="checkbox"/>	DDZ0360	Attract black	
<input type="checkbox"/>	DDZ0307	Attract grey	
<input type="checkbox"/>	DDZ0308	Attract blue	
<input type="checkbox"/>	DDZ0311	Attract red	
FRONT WHEELS			
<input checked="" type="checkbox"/>	DDZ0801	4" Everyday wheel with black tyre	
<input type="checkbox"/>	DDZ0860	4" Comfort wheel with black tyre	
<input type="checkbox"/>	DDZ0810	5" Everyday wheel with black tyre	NC with 90° front frame (DDZ0294)
<input type="checkbox"/>	DDZ0863	5" Comfort wheel black tyre	NC with 90° front frame (DDZ0294)
REAR WHEELS			
CAMBER AXLE TYPE			
<input checked="" type="checkbox"/>	DDZ0277	Standard aluminium axle, black	
REAR WHEEL FIXATION			
<input checked="" type="checkbox"/>	DDZ0250	Standard rear wheel position	
<input type="checkbox"/>	DDZ0251	Rear wheel extension (5 positions)	Only with rear wheel extension (DDZ0251)
POSITION OF AXLE			
<input type="checkbox"/>	DDZ0263	1st position from the front	
<input type="checkbox"/>	DDZ0264	2nd position from the front	
<input type="checkbox"/>	DDZ0265	3rd position from the front	
<input type="checkbox"/>	DDZ0266	4th position from the front	
<input type="checkbox"/>	DDZ0267	5th position from the front	
CAMBER ANGLE			
<input type="checkbox"/>	DDZ0270	0° camber	
<input checked="" type="checkbox"/>	DDZ0271	3° camber	
REAR WHEELS TYPE			
<input checked="" type="checkbox"/>	DDZ0921	24" Active wheels	With alu grey handrims with 6 points of fixation (screws) and light tyre 7.5 bar as standard 24 spokes and hub in black ; NC with spoke guards
<input type="checkbox"/>	DDZ0915	24" küschall Starec	Starec hub colourless, rim black - Cross-spoked
<input type="checkbox"/>	DDZ0908	24" High Performance wheels	Starec hub, rim and spokes black
<input type="checkbox"/>	DDZ0906	24" Spinergy LXX	Hub, rim and spokes black
ALBER ADD-ON KIT			
<input type="checkbox"/>	DDZ0916	E-Motion adaptation plate	ALBER antitipper needs to be installed with this option 2 ALBER antitippers included
<input type="checkbox"/>	DDZ0926	Twion adapter plate	2 ALBER antitippers included
REAR WHEEL TYRE			
<input type="checkbox"/>	DDZ1014	24" Schwalbe Rightrun black/grey	
<input type="checkbox"/>	DDZ1010	24" Schwalbe Marathon Plus	SmartGuard Puncture protected (5mm thick layer of highly flexible special rubber)
<input type="checkbox"/>	DDZ1004	24" Solid tyre, profile, grey	Only with küschall Starec wheels (DDZ0915)
<input type="checkbox"/>	DDZ1011	24" Schwalbe ONE	V-Guard puncture resistant (cut-resistant high-tech fiber) Not for active wheels (DDZ0921)
HANDRIMS			
<input type="checkbox"/>	DDZ1102	24" Supergripp, dist. 30 or 40 mm	Long tabs with 2 holes at distance 30 or 40 mm, mounted at distance 30 mm
<input type="checkbox"/>	DDZ1104	24" Steel polished, dist. 30 mm	Short tabs with only one hole at distance 30 mm
<input type="checkbox"/>	DDZ1113	24" Light alu, grey, dist. 30 mm	Short tabs with only one hole at distance 30 mm
<input type="checkbox"/>	DDZ1120	24" The Surge LT 30mm, dist. 30 mm	Short tabs with only one hole at distance 30 mm
SPOKE GUARDS			
<input type="checkbox"/>	DDZ1200	24" Transparent	NC with Active (DDZ0921) and Spinergy LXX (DDZ0906) wheels
SAFETY			
BRAKES			
<input checked="" type="checkbox"/>	DDZ1300	Standard brake, lever bent	
<input type="checkbox"/>	DDZ1302	Performance brake	
<input type="checkbox"/>	DDZ1303	Active brake	
BRAKE EXTENSION LEVER			
<input type="checkbox"/>	DDZ1310	Extension lever for standard brake	Only with standard brakes, levers bent
SAFETY OPTIONS			
ANTITIPPER			
<input type="checkbox"/>	DDZ1702	Antitipper tube, right	NC with ALBER add-on kit NC with cane holder or tipping aid on the same side
<input type="checkbox"/>	DDZ1703	Antitipper tube, left	NC with cane holder or tipping aid on the same side
<input checked="" type="checkbox"/>	DDZ1704	Antitipper tube, both sides	NC with cane holder or tipping aid

TRANSPORTATION KIT		Only available with backrest height (RH) > 400 mm, not with Matrix backrests
<input type="checkbox"/> DDZ1770	Transportation kit incl. posture belt and labels	Without tie down system
POSTURE BELT		
<input type="checkbox"/> DDZ1730	Posture belt	NC with transportation kit (DDZ1770)
PASSIVE LIGHTS		
<input type="checkbox"/> DDZ1710	Passive lights, assembled	
<input type="checkbox"/> DDZ1711	Passive lights, enclosed	
ACCESSORIES		
FRAME PROTECTORS		
<input type="checkbox"/> DDZ1739	Front frame protectors "küsçhall", 2 pieces	NC with UL < 400 mm
TIPPING AID		
<input type="checkbox"/> DDZ1708	Tipping aid, right	NC with cane holder or antitipper on the same side
<input type="checkbox"/> DDZ1709	Tipping aid, left	NC with cane holder or antitipper on the same side
<input type="checkbox"/> DDZ1718	Tipping aid, both sides	NC with cane holder or antitipper
CANE HOLDER		
<input type="checkbox"/> DDZ1706	Cane holder, right	NC with tipping aid or antitipper on the same side
<input type="checkbox"/> DDZ1707	Cane holder, left	NC with tipping aid or antitipper on the same side
PUMP		
<input type="checkbox"/> DDZ1715	Air pump	
KÜSCHALL WATCH		
<input type="checkbox"/> DDZ1712	küsçhall watch, right side	Mounted in castor housing
<input type="checkbox"/> DDZ1713	küsçhall watch, left side	
TOOL KIT		
<input type="checkbox"/> DDZ1717	Tool kit	
WHEELS BAG		
<input type="checkbox"/> DDZ1728	Wheel bag: for 24", 25" and 26" rear wheels	

SUPPORT FOR MEASUREMENTS / MYKÜSCHALL REQUESTS

